

# A HIGH THROUGHPUT LC-MS LIPIDOMICS PLATFORM FOR SCREENING & STRATIFYING CANCER SAMPLES FROM BIOBANKS AND CLINICAL TRIALS

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

Adam M King<sup>1,2</sup>, Robert D Trengove<sup>2</sup>, Lauren G Mullin<sup>3</sup>, Bryan Katzenmeyer<sup>3</sup>, Paul D Rainville<sup>3</sup>, Giorgis Isaac<sup>3</sup>, Robert S Plumb<sup>3</sup>, Lee A Gethings<sup>1</sup>, and Ian D Wilson<sup>4</sup>  
 1 Waters Corporation, Wilmslow, UK. 2 Murdoch University, Perth, AUS. 3 Waters Corporation, Milford, MA, US. 4 Imperial College, London, UK

## INTRODUCTION

Lipids encompass a large group of compounds of various classes, each exhibiting a range of physicochemical properties and biological function. Changes in expression and metabolism of lipids has been linked to numerous diseases including diabetes<sup>1</sup> and various cancers<sup>2</sup>. In order to screen the lipid phenotype of large batches of samples obtained from clinical trials and biobanks requires a high throughput analytical assay in order to analyse the thousands of potential samples.

Conventional discovery LC-MS lipidomic assays have sample acquisition times of >15 minutes per sample. Using these methodologies to analyse larger cohorts of samples from biobanks can lead to weeks of expensive analysis putting pressure on laboratory resources.

Reducing a column's internal diameter, column length and scaling down mobile phase flow rates and gradients can dramatically reduce the overall acquisition time with minimal impact on chromatographic performance. Here we describe the development and application of a high throughput discovery lipidomic assay.




Figure 1. Summary schematic of sample workflow, analysis and data processing.

TO DOWNLOAD A COPY OF THIS POSTER, VISIT [WWW.WATERS.COM/POSTERS](http://WWW.WATERS.COM/POSTERS)

## METHODS

### Sample preparation

- Human plasma samples from breast cancer patients (n=20) and normal control subjects (n=6).
- Pooled QC sample prepared from each study sample (50 µL).
- Lipids extracted from 100 µL of sample with 400 µL of IPA.
- Extracts then incubated at 2–8 °C for 2 hours.
- Extracts centrifuged to remove proteins.
- Supernatant removed for analysis (Fig. 1).

### MS conditions

- Waters Synapt G2-Si with IMS enabled
- Positive ESI, sensitivity mode, 50–1200 m/z
- Capillary voltage: 0.5 kV, Sampling cone: 30 V
- Source temp: 120 °C, desolvation temp: 500 °C
- Desolvation gas: 800 L/hr, cone gas: 50 L/hr.
- IMS wave velocity: 600 m/s, wave height: 40 V

### LC conditions

- Column: Waters BEH C8, 1.0 x 50 mm (1.7 µM)
- Mobile phase:
  - A) H<sub>2</sub>O:IPA:MeCN w/ NH<sub>4</sub>CH<sub>3</sub>CO<sub>2</sub> & CH<sub>3</sub>COOH
  - B) IPA:MeCN w/ NH<sub>4</sub>CH<sub>3</sub>CO<sub>2</sub> & CH<sub>3</sub>COOH

Table 1. Gradient composition and injection volumes for both the conventional and rapid lipid profiling assay.

Conventional method			Rapid lipid method		
Injection	Vol (µL)	%B	Injection	Vol (µL)	%B
Initial	0.6	1.0	Initial	0.25	1.0
0.10	0.6	1.0	0.05	0.25	1.0
2.00	0.6	30.0	0.50	0.25	30.0
11.50	0.6	90.0	2.80	0.25	90.0
12.00	0.6	99.9	3.00	0.25	99.9
12.50	0.6	99.9	3.15	0.25	99.9
12.75	0.6	1.0	3.25	0.25	1.0
13.25	0.6	1.0	3.70	0.25	1.0

## Data processing


- All data was processed using Progenesis Q1 (Nonlinear Dynamics, UK).
- Raw data was aligned, peak picked and then normalized to all compounds.
- Detected ions underwent adduct deconvolution to determine neutral mass values.
- Statistical analysis was performed using EZinfo (Umertrics, SWE).
- Group separation was visualized by principle component analysis (PCA).
- Significant features were determined through orthogonal partial least squared discriminate analysis (OPLS-DA) and S-plots of the features.
- Compound identification was performed through LipidMaps and Waters IROA CCS database searches.

## RESULTS

### Rapid LC methodology

- Scaled chromatography reduced sample acquisition times from 13.25 mins to 3.7 mins and solvent consumption by 75 %.
- Column I.D. Was reduced by a factor of 4 while flow rate reduced by only 2.4 enabling equivalent column volumes and increased linear velocity.
- Chromatographic separation of lipid classes was maintained through the scaling process (Fig. 2).

Figure 2. Comparison of a conventional discovery lipid assay with the scaled down 1.0 mm rapid lipid separation from human plasma.



### Ion mobility spectrometry

- Due to the reduction in chromatographic resolution, enabling ion mobility spectrometry (IMS) gave an additional separation for co-eluting ions.
- The collision cross sectional (CCS) measurements provided increased specificity with lipid classes forming distinct groups by CCS value (Fig. 3).
- With IMS enabled the spectral quality of co-eluting features was improved. (Fig. 4)




Figure 4. A) High collision energy fragment ion spectra without IMS of two co-eluting lipids B) plot of identified features by retention time and CCS with the co-eluting features highlighted and separated by CCS value. The highlighted region shows a co-eluting DG & SM which are clearly separated in the IMS domain. C) Increased specificity provides improved fragmented ion matching and database identifications.

- Using conventional data independent acquisition modes, fragment ions from co-eluting species can be miss assigned.
- The ion mobility separation prior to CID improved the assignment of fragment to precursor ions.
- This in turn improved the spectral matching for database searching.
- 5 features were determined to be up regulated in the breast cancer samples with another 10 lipids shown to be down regulated.
- These features were considered for database searching using lipid maps.
- Additional database searching was conducted against the IROA CCS database.




Figure 3. Plot of identified plasma features by measured m/z and CCS value. Different lipid classes are highlighted.

### Feature Identification

- Ions that had a %CV in the QC samples >30 % were filtered out and the remaining underwent statistical analysis.
- Following OPLS-DA (Fig. 5A), significant features were determined by s-plot (Fig. 5B).




Figure 5. A) OPLS-DA comparison of breast cancer samples and healthy controls. B) Resulting S-plot of significant features.

- Of the lipids identified, a number of PCs and TAGs showed a reduction in expression in breast cancer patients (Table 2).
- The decrease in PCs can be indicative of an increase in phospholipase A2 activity, previously reported by Yunping *et al*<sup>3</sup> to be linked to breast cancer.
- Five PS species were also noted to be increased in breast cancer patients.
- In the literature, PS has been noted as a potential biomarker for cancer, corresponding to the results shown here<sup>4</sup>.

Table 2. List of up regulated (green) and down regulated (red) features in breast cancer patients and potential identifications following database searching

Lipid identifi- cation	Neutral mass (Da)	m/z	Retention time (min)	CCS (Å <sup>2</sup> )	ΔCCS (Å <sup>2</sup> )	Peak width (min)	Anova (p)	q Value	Max Fold Change	Minimum CV%
TG(52:3)	856.75	874.79	3.14	334.1	-	0.20	5.5E-05	0.00166	1.5	6.03
TG(52:4)	-	872.77	3.06	331.4	-	0.23	1.0E-04	0.00201	1.8	4.85
TG(54:5)	880.75	898.78	3.07	337.6	-	0.18	1.5E-04	0.00237	2.6	5.53
PS (40:4)	839.57	822.56	1.62	304.4	-	0.30	2.1E-04	0.00274	2.0	4.38
TG(54:3)	-	902.82	3.22	341.1	-	0.24	3.4E-04	0.00386	2.0	5.54
DG(54:1)	-	577.52	3.2	267.1	-	0.18	3.7E-04	0.00404	1.8	2.97
TG(60:1)	832.75	850.79	3.19	332.5	-	0.16	3.9E-04	0.00408	1.7	5.15
PS(36:2)	773.56	774.56	1.55	295.8	-	0.31	5.0E-04	0.00452	1.9	3.16
PS(38:3)	799.57	782.57	1.62	299.7	-	0.28	6.5E-04	0.00529	1.8	3.09
PS(36:1)	789.55	790.56	1.54	298.7	-	0.41	2.9E-03	0.01119	2.9	2.21
PC(36:4)	781.56	782.57	2.01	306.6	-	0.28	6.6E-03	0.02300	1.5	5.15
PC (38:4)	809.59	810.60	2.19	312.2	7.2	0.37	8.3E-03	0.02592	1.6	4.06
PC (36:2)	785.60	786.60	2.19	304.6	4.6	0.20	1.9E-02	0.04109	1.3	4.69
PS (38:2)	815.57	816.57	1.63	306.3	-	0.41	2.2E-02	0.04428	1.5	4.56
PC(34:2)	757.57	758.57	2.02	296.9	-	0.24	2.7E-02	0.04926	1.3	4.84

## CONCLUSIONS

- A rapid discovery lipidomic method has been developed and applied in the assessment of plasma breast cancer samples.
- The reduction in acquisition time demonstrated potential for larger cohort studies to be acquired in a matter of days rather than weeks.
- The additional of IMS improved the high collision energy spectra.
- Generation of CCS measurements augmented the database searches and increased ID confidence.
- Significant lipids identified corresponded to those documented in the literature.

### References

- Al-Sulaiti, H., I. Diboun, S. Banu, M. Al-Emadi, P. Amani, T. M. Harvey, A. S. Domling, A. Latiff, and M. A. Elrayess. 2018. Triglyceride profiling in adipose tissues from obese insulin sensitive, insulin resistant and type 2 diabetes mellitus individuals. *J Transl Med* 16: 175.
- Yang, L., X. Cui, N. Zhang, M. Li, Y. Bai, X. Han, Y. Shi, and H. Liu. 2015. Comprehensive lipid profiling of plasma in patients with benign breast tumor and breast cancer reveals novel biomarkers. *Anal Bioanal Chem* 407: 5065–507
- Y. Qiu, B. Zhou, M. Su, S. Baxter, X. Zheng, X. Zhao, Y. Yen, W. Jia. 2013. Mass spectrometry-based quantitative metabolomics revealed a distinct lipid profile in breast cancer patients. *Int J Mol Sci*, 14: 8047–8061.
- Sharma, B., and S. S. Kanwar. 2017. Phosphatidylserine: A cancer cell targeting biomarker. *Semi Cancer Biol*.