

# A WORKFLOW DRIVEN PLATFORM SOLUTION FOR MAM-BASED CRITICAL QUALITY ATTRIBUTE MONITORING OF BIOTHERAPEUTICS IN PROCESS DEVELOPMENT AND QC

Waters

THE SCIENCE OF WHAT'S POSSIBLE.®

Nilini S Ranbaduge, Ying Qing, Min Du, Weibin Chen  
Waters Corporation, 34 Maple Street, Milford, MA 01757, USA

## INTRODUCTION

- Biotherapeutics undergo rigorous characterization and monitoring during development and production to maintain product quality and safety.
- Peptide-based monitoring is often used to monitor multiple attributes affecting the product quality.
- A compliance-ready, easy-to-use, high performance LC-MS system with automated methods is highly desired.
- The BioAccord is an SmartMS-enabled LC-MS system purposefully designed with integrated workflow methods for biopharmaceutical analyses including peptide-based quality attribute profiling and monitoring.



- Small footprint benchtop system
- Easy-to-use & maintain
- Workflow-driven
- cGMP compliance-ready

## METHODS

### Sample Preparation:

- The Trastuzumab (Genentech, USA) sample was subjected pH, heat and oxidative stress conditions as given below.

pH 9.0, at 37°C				
1 day	2 days	4 days	6 days	
Heat, 37°C				
4 days			6 days	
H <sub>2</sub> O <sub>2</sub> , room temperature, 1 day				
0.005%			0.05%	

### LC-MS System: BioAccord System

- Column: ACQUITY UPLC CSH C18, 1.7 μm, 2.1 x 100 mm
- Column temperature: 60 °C
- Mobile phase: A. 0.1% FA in water, B. 0.1% FA in Acetonitrile
- Total run Time: 80 min
- Sample temperature: 6 °C
- Injection volume: 5 μL
- Flow rate: 0.2 mL/min
- Gradient: 1%-35% Acetonitrile + 0.1% formic acid over 52 min


- Ionization: ESI+
- Acquisition mode: MS with fragmentation
- Acquisition range: m/z 50-2000
- Capillary voltage: 1.2 kV

## METHODS


- Cone voltage: 30 V
- Collision energy: 60-120 V
- Desolvation energy: 350 °C
- Intelligent data capture: on

- Informatics: UNIFI Scientific Information v1.9.4
  - Peptide mapping method
  - Accurate mass screening method
  - UNIFI scientific library


### End-to-end Workflow




### Separation & detection



### Data processing



### Reporting




**Figure 1:** The end-to-end workflow for peptide attribute profiling and monitoring. The UNIFI processing methods are purposefully designed with automated data acquisition, processing and reporting capabilities.

## RESULTS

### Peptide Mapping Method

#### (A) Chromatographic separation of peptides and detection




#### (B) Identify peptides and their modifications

Protein name	Fragment label	Peptide	Modifiers	Observed RT (min)	Observed m/z	Charge	Observed mass (Da)	Mass error (ppm)	Response	Matched 1st Gen Primary Ions
TnAb	279	FTGADTSK		13.67	485.2478	2	969.4983	-0.5	1420659	
TnAb	2737	GRFSDVAVFVENSNGPENNK		30.11	848.7109	3	2544.1361	1.9	138932	
TnAb	275	GLBNVBAR		17.97	415.7302	2	830.4531	1.4	16312430	
TnAb	2713	GPSVPLAPSSK		22.86	593.8281	2	1186.6489	1.9	20901224	
TnAb	276	IYPTNGYTR		11.57	542.7746	2	1084.5419	-0.3	5447678	
TnAb	2766	IYPTNGYTR	Deamidation N [5]	12.84	543.2662	2	1085.5251	-1.0	1797664	
TnAb	2766	IYPTNGYTR	Deamidation N [5]	13.21	543.2679	2	1085.5285	2.1	421322	
TnAb	175	LLVSAFLYSQVPSR		34.42	591.6587	3	1772.9617	2.0	9212968	
TnAb	2728	LSCAASGNK	Carbamidomethyl C [3]	17.24	584.2964	2	1167.5856	2.5	13685866	
TnAb	2739	LTVDK		5.15	288.1736	2	575.3399	-0.1	3709279	
TnAb	27368	NQVSLTCLVK	Carbamidomethyl C [7]	21.62	581.3192	2	1161.6311	1.3	20449848	
TnAb	27368	NQVSLTCLVK	Carbamidomethyl C [7], Deamidation N [1]	24.34	581.8081	2	1162.6089	-4.1	254887	

- Custom attribute library
- Sequence & modifications
- Mass & retention time
- Detected charge states, etc.

#### (C) Sequence confirmation using fragment ions




**Figure 2:** (A) The peptides and their modifications are chromatographically separated prior to MS analysis. This figure shows the total ion chromatogram of one of the pH stressed samples. (B) The component table provides all identified peptide sequences and modifications with their respective retention times, mass tolerance, m/z, number of primary fragment ions, etc. (C) The figure shows annotated fragment ion spectra for native and modified IYPTNGYTR. The fragmentation information is used in data filtering for high confidence peptide attribute identification.

### Peptide Attribute Monitoring

#### (A) Target peptide attribute list

Component name	Label	Expected RT (min)	Expected neutral mass (Da)	Expected fragment (m/z)	Adducts	Description	Internal standard?
1	YGGFL	22.87	555.2693	171+H <sup>+</sup> - LeuK		IYPTNGYTR	<input checked="" type="checkbox"/>
2	SLSLSPGK - Lysine C-TERM	17.61	659.3490	2742B - Lysine C-TERM [B]+H <sup>+</sup> - TnAb		IYPTNGYTR	<input type="checkbox"/>
3	SLSLSPGK	12.22	787.4440	2742+H <sup>+</sup> - TnAb		IYPTNGYTR	<input type="checkbox"/>
4	DTLMSR	15.44	834.2099	2721+H <sup>+</sup> - TnAb		IYPTNGYTR	<input type="checkbox"/>
5	DTLMSR Oxidation M	12.70	850.4239	2722B Oxidation M [B]+H <sup>+</sup> - TnAb		IYPTNGYTR	<input type="checkbox"/>
6	IYPTNGYTR	11.56	1083.5348	2716+H <sup>+</sup> - TnAb		IYPTNGYTR	<input type="checkbox"/>
7	IYPTNGYTR Deamidation N	12.83	1084.5189	2716B Deamidation N [5]+H <sup>+</sup> - TnAb		IYPTNGYTR	<input type="checkbox"/>
8	IYPTNGYTR Deamidation N_1	13.20	1084.5189	2716B Deamidation N [5]+H <sup>+</sup> - TnAb		IYPTNGYTR	<input type="checkbox"/>
9	NTAYLQMSLR	19.87	1309.6449	2710+H <sup>+</sup> - TnAb		IYPTNGYTR	<input type="checkbox"/>
10	NTAYLQMSLR Deamidation N	21.23	1310.6289	2710B Deamidation N [5]+H <sup>+</sup> - TnAb		IYPTNGYTR	<input type="checkbox"/>
11	NTAYLQMSLR Deamidation N_1	20.41	1310.6289	2710B Deamidation N [5]+H <sup>+</sup> - TnAb		IYPTNGYTR	<input type="checkbox"/>


#### Accurate mass screening-based detection and quantification



**Figure 3:** (A) The target peptide list is imported from custom attribute library created using peptide mapping results. Each library entry contains peptide sequence, neutral mass and retention time information (and fragment ions). Both neutral mass and retention time are used to identify peptide CQAs in the accurate mass screening method. (B) The XICs of target peptides (native and modified) are extracted for the monoisotopic masses or multiple isotopes. The peak area of each XIC is used to generate relative %modifications. (C) The summary plot shows native and deamidated IYPTNGYTR peptide %modification levels across samples. Each relative modification level is determined based on the total MS response for the respective peptide sequence.


#### Limit check capability

Node	Field name	Component	Sample type	Level	Error minimum	Warning minimum	Warning maximum	Error maximum
1	Component	%Peptide Mod-MS (%)	DTLMSR Oxidation M	All	All	25	50	
2	Component	%Peptide Mod-MS (%)	IYPTNGYTR Deamidation N, IYPTNGYTR Deamidation N_1	All	All	6	15	




**Figure 4:** (A) The limits can be applied during data processing by introducing upper/lower thresholds. (B) The limits are used for easy visualization of data that flags both samples and attribute in a preset color. In here, yellow and red display warning and error limits of an attribute respectively.

#### New peak detection



**Figure 5:** Comparison mode feature compares sample to its reference/control. This figure shows a chromatographic comparison of Trastuzumab control and a pH stressed sample. Specific data filters are applied to isolate unique peptide peaks or peptides present at a higher level (>5-fold) in the stress sample compared to the control.

#### Comparable quantification data obtained with and without fragmentation



**Figure 6:** When performed in MS only mode %modification level for an attribute provides comparable results to that of the MS with fragmentation mode. This allows the user to choose either method in routine peptide attribute monitoring.

## CONCLUSION

- The BioAccord is a compact LC-MS system that can be easily deployed across biopharmaceutical organizations from discovery thru QC.
- Integrated workflows (UNIFI Informatics) for intact protein, released-glycan and peptide analysis support automated attribute screening and monitoring workflows.
- Accessibility to non MS experts is facilitated by one-button startup, intelligent diagnostics, and automated tuning/calibration routines.
- Information transfer between peptide characterization and monitoring workflows is enabled by the UNIFI scientific library functionality, allowing lists of attributes to be incrementally updated for later use in targeted monitoring analyses.
- Automated data acquisition, processing, review, reporting, and signoff streamlines operations within regulated laboratories where data integrity is a concern.
- The ability to automatically highlight product attributes exceeding warning or error limits and denote new data features (potential attributes) vs. a reference sample greatly simplifies review and reporting of multiple attribute monitoring (MAM) studies.