

Fully automated chemical imaging with LA-REIMSTM

Gabriel Stefan Horkovics-Kovats^{1,2}, Richard Schäffer¹, Csaba Hajdu¹, Gitta Schlosser², Julia Balog¹; Budapest
¹Waters Research Center Kft., ²Eötvös Loránd University

Introduction

Mass spectrometry chemical imaging is an emerging technology, which can be used for multiple applications in different scientific fields, for example chemical imaging based histological assessment or drug monitoring. The Laser Assisted – Rapid Evaporative Ionization Mass Spectrometry (LA-REIMSTM) is an ambient technique requiring no sample preparation of the target tissue. The technique combined with an imaging setup developed for this purpose is suitable for a micro-scale resolution metabolic profiling through point-by-point laser ablation of the tissue section.

Goal

Creating a functional and fully automated LA-REIMSTM based imaging platform, testing its usability and compare to other imaging modalities (e.g. DESITM).

Instrumentation

- Imaging setup is linked to a Xevo™ G2-QToF-MS with REIMS interface. All experiments and testing procedures were performed in sensitivity and negative ion mode of the mass spectrometer.
- For this purpose, an Opolette™ 2940 from Oportek is used which, due to its applied wavelength of 2940 nm, selectively disrupts hydrogen bonds and thus releases other constituents for MS - analysis.
- To handle microscope slides automatically, the Panoramic™ Scan II - automatic slide sorter from 3DHISTECH has been integrated into the setup.
- A Masslynx™ compatible home made software allows the control of all components of the setup and enables easy to use ROI selection and acquisition handling

REIMS, DESI, Xevo and Masslynx are trademarks of Waters Technologies Corporation. Opolette is a trademark of OPOTEK. Panoramic is a trademark of 3DHISTECH.





Figure 1. Current LA-REIMSTM Imaging setup and schematic hardware structure


Creation and Challenges

- Combine laser, optics and Panoramic Scan II into an integrated system
- Correct laser beam alignment to achieve round spot size
- Testing of different lenses to obtain the best possible spot size and spectra
- Camera adjustment through hardware and software to avoid radial and tangential torsion as much as possible
- Creation of a suitable image preview

Setup components (Figure 1):

- Opolette 2940
- Adjustment plates
- Cage System Iris Diaphragm (Ø0.8 - Ø20 mm)
- CaF₂ Plate Beamsplitter
- Possible mount for energy measurement
- 2x motorized linear translation stage (stepper motor)
- mounted neutral density filter
- ZnSe aspheric lens (f=12.7mm)
- Panoramic Scan II – automatic slide sorter

Results: Proof of Concept



- Successful sampling and image generation (Figure 2)
- First direct comparisons between LA-REIMSTM and DESITM results (Figure 3)
- Correct classification of brain sectors through PCA and LDA (Figure 4 and 5)

Figure 2. Created image on a specific ROI of a cancerous dog testicle section, sampled on 50 microns resolution that visualizes separation of tissue types




Figure 3. LA-REIMSTM VS DESITM imaging comparison on mouse brain sections on 50 microns sampling resolution (500-1000 m/z, single scan comparison from the same cortex region). Showcasing of complementary spectra: LA-REIMSTM spectra featuring more abundant PEs with ammonia loss; DESITM spectra featuring more abundant PIs




Figure 4. PCA 1 of mouse brain using custom built software, AMX

Figure 5. Classification of mouse brain section based on previously calculated LDA classifier

Conclusion

We have created a fully automated, integrated, high throughput (up to 150 slides), LA-REIMSTM based imaging platform and demonstrated its feasibility to achieve 50 um resolution clear chemical images from tissue sections.



PREPARED WITH THE PROFESSIONAL SUPPORT OF THE DOCTORAL STUDENT SCHOLARSHIP PROGRAM OF THE CO-OPERATIVE DOCTORAL PROGRAM OF THE MINISTRY OF INNOVATION AND TECHNOLOGY FINANCED FROM THE NATIONAL RESEARCH, DEVELOPMENT AND INNOVATION FUND.