

Lisa Reid¹, Mark Towers¹, Lee Gethings¹, Royston Goodacre², Ian Wilson² & Paul Richardson²
¹Waters Corporation, Wilmslow, Cheshire, UK
²Department of Biochemistry and Systems Biology, Institute of Systems, Molecular, and Integrative Biology, University of Liverpool, UK

INTRODUCTION

In biomedical research, to profile blood products for large cohort lipidomic studies may consist of long LC/MS analyses requiring several weeks of instrument time. Consequently, increasing the potential for instrument failure, unwanted batch effects and sample degradation.

DESI imaging is traditionally used to investigate compound localization within tissue where, by visualizing the structure and molecular composition of tissues, a greater understanding of organ structure and changes due to disease may be obtained.

Here we have examined the feasibility of using DESI as an alternative high throughput direct analysis method for performing lipid profiling. This technique has been applied to screen sera obtained from a UK human population as part of the Human Serum Metabolome study (HUSERMET).

Each slide of 44 samples required 8 minutes acquisition time giving a total acquisition for all 500 samples of <3 hours per polarity. **Figure 2** shows a typical total ion chromatogram (TIC) from a line of 11 spots. Acquisition speed could be increased if desired, however, on this occasion the scan rate was intentionally limited to ensure excellent signal strength and maximise the information generated per sample.

