

Quantification of Corticosteroids and Androgens in Serum, Utilizing Waters MassTrak™ Steroid Serum Sets 2 and Sets 3 for Clinical Research

Waters™

Laura Moran¹, Dominic Foley², Niall Tobin¹, Norma Breen¹, Leanne Davey¹

¹Waters Technology, IDA Business Park, Drinagh, Wexford, Y35 D431, Ireland, ²Waters Corporation, Stamford Avenue, Altrincham Road, Wilmslow, SK9 4AX UK

Introduction

Steroid hormones encompass a large class of small molecules that play a central role in metabolic processes, such as the regulation of sexual characteristics, blood pressure, and inflammation. Enzymes that form part of the steroid biosynthetic pathway are pivotal in these metabolic processes, and their dysfunction can be examined through the correct measurement of steroid hormones in a clinical research setting. The availability of lyophilized calibrators and QCs reduces sample preparation time, aids in method harmonization, and assists with metrological traceability in accordance with ISO15189:2022, when used alongside an analytically selective chromatographic method.

Outlined is a quantitative clinical research method utilizing Waters MassTrak Steroid Serum Sets 2 and 3, Cals & QCs and Waters Oasis™ PRiME HLB μElution plate technology for the extraction of testosterone, androstenedione, 17-hydroxyprogesterone (17-OHP), dehydroepiandrosterone sulfate (DHEAS), Cortisol, 11-deoxycortisol and 21-deoxycortisol from human serum samples.




Figure 1: Sample preparation workflow.

Chromatographic separation was performed on an ACQUITY™ UPLC™ I Class Plus (FL-I) System, using an ACQUITY UPLC HSS T3 Column, accompanied by a Xevo™ TQ-S micro, mass detector.





Figure 2: Chromatographic selectivity on the ACQUITY UPLC HSS T3 Column for the selection of Steroid Hormones



Figure 3: ACQUITY™ UPLC™ I Class Plus (FL-I) System, accompanied by a Xevo™ TQ-S micro mass detector.


An accurate ($\pm 6\%$) and precise ($\pm 10\%$) research method for the quantification of steroid hormones in serum.




Scan the QR Code to download the full Application Note




Figure 4: MassTrak Steroid Serum Calibrators and QCs



Graph 1: MassTrak Steroid Serum Set 2 & 3 Precision for each QC – two replicates over 5 days



Graph 2: MassTrak Steroid Serum Set 2 & 3 Repeatability for each QC – four replicates over 5 days



Graph 3: MassTrak Steroid Serum Set 2 & 3 Accuracy – Mean % deviation from target concentration at low, medium & high end of the calibration line

Key Takeaways

- Accuracy ($\pm 6\%$) and precision ($\pm 10\%$) have been confirmed through comparison to External Quality Assurance (EQA) schemes, in-house panels and QC material for all seven steroid hormones.
- Set 2 and Set 3 were analysed over a five-day period and were linear with a co-efficient of determination of (r^2) > 0.999 for all analytes.
- Deming and Linear regression analysis were performed. No statistically significant bias was observed for each compound, with a mean method bias of $\pm 1.3\%$ for Set 2 and $\pm 1.0\%$ for Set 3.
- Analytical sensitivity using Signal: Noise (S/N) of the low calibrator (Calibrator 1) of each set, was >10:1 for each analyte across several analytical runs.

For Research Use Only. Not for use in diagnostic procedures.