

Author

Meng Yuan

PerkinElmer, Inc.

Shelton, CT 06484 USA

Analysis of Butylated Hydroxytoluene in Food with Headspace Trap-GC/MS

Introduction

Butylated hydroxytoluene (BHT, 2,6-di-tert-butyl-4-methylphenol) is a common food additive. BHT is found in many types of food including butter, meats, cereals, chewing gum, baked goods, snack foods, dehydrated potatoes and beverages. It is used to preserve food odor, color and flavor. BHT is oxidized preferentially in fats or oils, protecting the foods from spoilage.

Concern exists that long-term human consumption of BHT may have potential health risks. It has undergone the additive application and review process required by the U.S. Food and Drug Administration (FDA); the committee concluded that no evidence in the available information on BHT demonstrates a hazard to the public when it is used at levels that are now current and in the manner now practiced. However, uncertainties exist requiring that additional studies should be conducted.¹ The chemical properties which make BHT an excellent preservative may also be implicated in health effects. The oxidative characteristics and metabolites of BHT may contribute to carcinogenicity. Some people may have difficulty metabolizing BHT, resulting in health and behavioral changes.

This application note will demonstrate a fast and easy analytical technique to determine the amount of BHT in foods. Headspace (HS) sample introduction is used because it provides a means to analyze food without any sample preparation. Headspace eliminates the need for solvents and other sample-preparation steps to reduce cost and complexity of extraction. In this application note, an adsorbent trap is used to concentrate the headspace sample and increase sensitivity, allowing for low-level detection or small sample sizes.

The analysis is carried out with gas chromatography mass spectrometry (GC/MS) – this will allow us to resolve the BHT from other volatile compounds in the food matrices and provide positive identification of the BHT with mass spectral data. Calibration of the system and analysis of food samples will be demonstrated.

Experimental

The instrumental platform for this application is the TurboMatrix™ HS Trap coupled to a Clarus® 680 GC/MS, both platforms from PerkinElmer. The transfer line of the HS was directly connected to the Elite™-17ms column with a universal butt connector. The samples are heated in a sealed vial at 80 °C for 30 minutes to drive the BHT from the food into the headspace. Using automated headspace technology, the gas is extracted from the vial, concentrated on an adsorbent trap (PerkinElmer® Air Toxics), and injected into the GC/MS system. Table 1 shows the detailed instrumental setup parameters for the HS Trap-GC/MS system.

Sample Introduction		PerkinElmer TurboMatrix HS-40 Trap		Gas Chromatograph		PerkinElmer Clarus 680 GC	
Needle Temperature	90 °C	Headspace Connector	Universal Connector	Inlet Temperature	150 °C	Oven Program Initial Temp	50 °C
Transfer Line Temperature	110 °C	Oven Temperature	80 °C	Hold Time 1	1 min	Ramp 1	25 °C/min to 280 °C
Trap Low Temperature	40 °C	Trap Low Temperature	40 °C	Hold Time 2	1.8 min	Vacuum Compensation	On
Trap High Temperature	280 °C	Trap High Temperature	280 °C	Headspace Control	On	Column	Elite-17ms 30 m x 0.25 mm x 0.25 µm
Dry Purge (Helium)	5 min	Dry Purge (Helium)	5 min	Carrier Gas	Helium	Mass Spectrometer	PerkinElmer Clarus 600 MS
Trap Hold Time	6 min	Trap Hold Time	6 min	Mass Range	45-300 u	Solvent Delay Time	0.1 min
Desorb Time	0.5 min	Desorb Time	0.5 min	Scan Time	0.20 sec	Scan Time	0.20 sec
Thermostatting Time	30 min	Thermostatting Time	30 min	InterScan Delay Time	0.02 sec	Transfer Line Temperature	240 °C
Pressurization Time	1 min	Pressurization Time	1 min	Transfer Line Temperature	240 °C	Source Temperature	200 °C
Decay Time	2 min	Decay Time	2 min	Source Temperature	200 °C	Multiplier	500 V
Column Pressure	17 psi	Column Pressure	17 psi	Multiplier	500 V		
Vial Pressure	35 psi	Vial Pressure	35 psi				
Desorb Pressure	10 psi	Desorb Pressure	10 psi				
Universal Capillary Column Connector	Part No. N9302149	Universal Capillary Column Connector	Part No. N9302149				
Transfer Line	Fused Silica 2 m x 320 µm	Transfer Line	Fused Silica 2 m x 320 µm				

Calibration-Standards Preparation

A 10 ng/µL standard stock solution was prepared by diluting 0.1 mL of a 1000 µg/mL BHT standard to 10 mL with methanol. 1 ng/µL, 2 ng/µL and 5 ng/µL standard working solutions were prepared by diluting 0.1 mL, 0.2 mL and 0.5 mL of a 10 ng/µL BHT standard to 1 mL with methanol. 20 ng/µL, 50 ng/µL and 100 ng/µL standard working solution was prepared by diluting 0.02 mL, 0.05 mL and 0.1 mL of a 1000 µg/mL BHT standard to 1 mL with methanol.

The working curve was prepared by injecting 1 µL of each working standard solution into headspace vials. Working calibration standards at 1, 2, 5, 10, 20, 50, and 100 ng were prepared fresh each day.

One gram of each food sample purchased at local Shanghai markets were placed into the headspace vials. All headspace vials were sealed immediately and transferred to the headspace-trap vial tray.




Figure 1. Example chromatogram of a 100 ng standard injection of BHT.

Results and Discussion

The instrument calibration included seven calibration levels in the working curve; the response of this calibration curve was linear (Table 2). Additionally, the method is precise throughout the calibration range, as demonstrated by the relative standard deviation of 3.2% at the calibration limit (1 ng, n=5) and 1.9% at 10 ng (n=5).

Figure 1 is an extracted ion chromatogram, of m/z 205, from the analysis of a 100 ng BHT standard. Figure 2 demonstrates the spectral data of BHT which matches exactly the fragmentation of BHT in the NIST® spectral library.




Figure 2. Background subtracted spectra from the analysis of a BHT standard.

Following the calibration of the system, five food samples were analyzed: a cracker, powdered coffee creamer, instant noodles, sausage, and tea leaves. The BHT concentrations are quantified (Table 3). BHT concentration in the food samples analyzed here was below the quantitation limit of 1 ng/g. It can be seen in Figure 3 (Page 4) that the BHT peak is easily identified in the sample analysis. Each sample was analyzed in triplicate – the area reproducibility achieved (Table 3) demonstrates that the method remains very precise, even below the quantitation limit.

Table 2. Calibration Table for BHT.

Name	Retention Time (min)	Quantifier Ion	Qualifier Ion 1	Qualifier Ion 2	%RSD (n=5 at 1 ng)	%RSD (n=5 at 10 ng)	r^2
BHT	7.60	205	220	57	3.2	1.9	0.9980

Table 3. %RSD of BHT in Food Samples.

Sample	BHT (ng/g) in 1 g of Sample	BHT (ng/g) in 1 g of Duplicate	BHT (ng/g) in 1 g of Triplicate	Values Mean (ng/g)	%RSD
Crackers	0.65	0.55	0.54	0.58	10.5
Coffee Creamer	0.66	0.73	0.69	0.69	5.1
Instant Noodles	0.67	0.67	0.70	0.68	2.5
Sausage	0.67	0.53	0.56	0.59	12.6
Tea Leaves	0.62	0.54	0.53	0.56	8.8




Figure 3. Resultant chromatogram from the analysis of instant noodles for BHT.

Conclusion

BHT is a common food additive used to prevent spoilage. Analysis of BHT is needed for both food quality and safety reasons. Food is often a complicated sample matrix which is time consuming to prepare and analyze. This method uses headspace technology to virtually eliminate sample preparation and reduce the cost and labor of the analysis. In addition to eliminating sample preparation, the method is both sensitive and precise as demonstrated by the analysis of standard reference materials and a variety of food samples. The throughput of the system is further improved by the Clarus 680 GC/MS with a fast-cooling GC oven, further improving throughput and productivity. The MS data provides positive confirmation of BHT in sample matrices.

References

1. Database of Select Committee on GRAS Substances (SCOGS) Reviews-Butylated Hydroxytoluene (BHT), available from <http://www.accessdata.fda.gov/scripts/fcn/fcnDetailNavigation.cfm?rpt=scogsListing&id=41>