

Repeatability in High Temperature Polyethylene Analysis using Agilent PLgel MIXED-B

Technical Overview

Introduction

A gel permeation chromatography system comprising Agilent PLgel MIXED-B columns and the Agilent PL-GPC 220 integrated, high temperature instrument is ideally suited to the analysis of polyethylene. The ability of the system is demonstrated in a repeatability study using eight polyethylene injections.

A commercial sample of high density polyethylene (HDPE) was prepared at 2 mg/mL using the Agilent PL-SP 260 sample preparation system with a dissolution temperature of 160 °C and a dissolution time of two hours. Eight aliquots of the master batch solution were dispensed into the Agilent PL-GPC 220 autosampler vials and placed in the autosampler carousel of the PL-GPC 220 where the hot zone temperature was 160 °C and the warm zone was 80 °C.

Figure 1 shows an overlay of the raw data chromatograms obtained for eight consecutive injections of the HDPE sample.




Figure 1. An overlay of the raw data chromatograms obtained for eight consecutive polyethylene injections.



The data were analyzed against a polystyrene standards calibration using the following Mark-Houwink parameters to obtain the polyethylene equivalent molecular weight averages. The molecular weight averages, peaks areas and a statistical analysis of the results for the eight injections are shown in Table 1.

Polystyrene in TCB¹ $K = 12.1 \times 10^{-5}$ $\alpha = 0.707$

Polyethylene in TCB² $K = 40.6 \times 10^{-5}$ $\alpha = 0.725$

Table 1. Calculated Molecular Weights for Eight Injections of Polyethylene and Calculated % Variation

Injection number	Mp	Mn	Mw
1	17,289	76,818	333,851
2	16,988	77,434	335,496
3	17,428	77,514	332,616
4	17,521	77,052	335,635
5	17,348	76,520	334,212
6	17,487	77,728	333,511
7	16,898	77,578	335,642
8	17,457	77,288	334,923
Mean	17,302	77,241	334,485
Std dev	220	387	1,048
% Variation	1.3	0.5	0.3

Conditions

Sample	High density polyethylene
Columns	3 × Agilent PLgel 10 μm MIXED-B, 7.5 × 300 mm (p/n PL1110-6100)
Eluent	TCB + 0.0125% BHT
Flow rate	1.0 mL/min
Inj vol	200 μL
Temp	160 °C
System	Agilent PL-GPC 220

Figure 2 is an overlay of the molecular weight distribution calculated for the eight consecutive injections of the HDPE sample, and illustrates the excellent repeatability obtained with the PL-GPC 220 using Agilent PLgel 10 μm MIXED-B columns.




Figure 2. Overlay of the molecular weight distribution calculated for eight consecutive injections of polyethylene that illustrates the excellent repeatability obtained with the Agilent PL-GPC 220 using Agilent PLgel 10 μ m MIXED-B columns.

References

1. H. Col and D. K. Giddings, J. Polym. Sci., (A2) 8 (1970) 89.
2. T. G. Scholte et al., J. Appl. Polym. Sci., 29 (1984) 3763.

GPC/SEC columns and calibrants from Agilent

Agilent offers a comprehensive portfolio of GPC/SEC columns and calibrants for high-performance separations based on molecular size in solution. Agilent delivers leading solutions for characterizing and separating polymers by GPC/SEC, and manufactures all components for accurate polymer analysis.

Look at the Agilent Literature Library on www.agilent.com/chem/gpc-sec for a comprehensive range of application notes and technical overviews to help you get the best from your Agilent GPC/SEC columns and instruments.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2015
Printed in the USA
April 30, 2015
5990-8496EN



Agilent Technologies