



# A Discovery Proteomics Workflow for the Elucidation of Prostate Cancer Biomarkers

## Application Note

### Authors

Claire Tonry and Stephen Pennington  
School of Medicine,  
UCD Conway Institute  
University College Dublin

Mark Sartain and Christine Miller  
Agilent Technologies, Inc.  
Santa Clara, California, USA

### Introduction

In biological systems, the host microenvironment is profoundly altered during tumor growth, and this includes becoming hypoxic due to insufficient blood supply. The hypoxic microenvironment correlates with increased tumor aggressiveness, invasiveness, and resistance to clearance [1]. In addition to other effects, the deprivation of molecular oxygen under hypoxic conditions may modulate the tumor cell proteome, leading to alterations in cell proliferation, and dynamics of the cell cycle [2]. Given the importance of androgen-regulated proteins, it is anticipated that further characterization of the role of hypoxia and androgen sensitivity may provide further insight into the mechanisms that drive aggressive, cellular hypertrophy/tumor growth [3].

This application note specifically demonstrates the LC/MS-based discovery portion of a clinical research workflow to comprehensively compare the proteome of androgen-independent and androgen-sensitive cell lines under both hypoxic and normoxic conditions, to identify potential protein biomarkers that may be indicative of important changes within the cellular microenvironment that drive abnormal cellular growth.



**Agilent Technologies**

## Methods

### Cell culture

The PCa cell lines LNCaP, LNCaP-abl (abl), and LNCaP-abl-Hof (Hof) were gifted to the Irish Prostate Cancer Research Consortium, Dublin, Ireland from the laboratory of Professor Helmut Klocker (Department of Urology, University of Innsbruck, Austria). Culturing the above cell lines was conducted in a Class II laminar flow cabinet. Cells were maintained in T175 cm<sup>2</sup> flasks with ventilation (Sarstedt) in a 5 % CO<sub>2</sub>-humidified atmosphere at 37 °C. LNCaP cells were maintained in Advance RPMI 1640 media (GIBCO Life Technologies), and supplemented with 10 % fetal calf serum (FCS) (Sigma-Aldrich), 2 μM/mL L-glutamine (GIBCO Life Technologies), 50 unit/mL penicillin, and 50 μg/mL streptomycin (GIBCO Life Technologies). Abl and Hof cells were maintained in Advance RPMI 1640 media supplemented with 10 % charcoal stripped FCS (Sigma-Aldrich), 2 μM/mL L-glutamine (GIBCO Life Technologies), 50 unit/mL penicillin, and 50 μg/mL streptomycin (GIBCO Life Technologies). Cell media were changed every 3–4 days.

### Simulation of hypoxia in PCa cell lines

LNCaP, Abl, and Hof cell lines were seeded into 10 cm<sup>2</sup> culture dishes and grown to 70–80 % confluence. For each cell line, medium was removed and replaced with appropriate medium supplemented with 1 mM dimethylloxaloglycine (DMOG; Cambridge Bioscience) or 1 mM dimethyl sulfoxide (DMSO; Sigma Aldrich) as a control. Cells were incubated in treated media for either 8 or 24 hours prior to cell lysis. For each cell line, three biological replicates were performed for each time point.

### Sample preparation for nano LC-MS/MS analysis

Whole cell lysates were prepared for nano LC-MS/MS analysis according to the filter-aided sample preparation (FASP) method described by Wisniewski; *et al.* [4]. Briefly, 50 μg of cell lysate protein were reduced through boiling (95 °C for 5 minutes) with DTT in a final concentration of 0.1 M. Then, 200 μL of UA buffer (8 M urea, 0.1 M Tris-HCL, pH 8.5) was added to each sample, and samples were transferred to 30,000 MWCO Vivacon 500 spin filters (Sartorius) and centrifuged at 14,000 g for 40 minutes at 21 °C. Bound proteins were alkylated by incubating the spin filters for 5 minutes in 0.05 M iodoacetamide (IAA), followed by centrifugation at 14,000 g for 30 minutes at 21 °C. Spin filter membranes were then washed three times by adding UB (8 M urea, 0.1 M Tris/HCL, pH 8.0) and centrifuging at 14,000 g for 40 minutes at 21 °C. For maximum protein identification, sample protein was digested with both Lys-C (Wako Chemicals GmbH) and trypsin (Promega) enzymes. Proteins were initially digested overnight with Lys-C (enzyme:substrate 1:50) in a wet chamber. Digestion was completed by a 3-hour incubation with trypsin (enzyme:substrate 1:100) in a thermomixer set to 37 °C and 600 rpm. Digestion was stopped by acidifying the samples through the addition of trifluoroacetic acid (TFA) to a final concentration of 1 %. Peptide material from digested cell lysates were purified using C18 resin ZipTips (Millipore). Each ZipTip contains C18 resin packed into a 10-μL pipette tip with a loading capacity of 5 μg protein/peptide per tip. This allows for purification of peptide material of molecular weight between 0–50 kDa. For purification of cellular peptides, the C18 resin was activated with 10 μL of acetonitrile (x10). The resin was then equilibrated by pipetting 10 μL of 0.5 % trifluoroacetic acid (TFA) (x10). Peptides were then bound to the resin by pipetting 15 μL of digested sample through the resin (x10). Bound peptides were eluted into fresh Eppendorf tubes in 25 μL of Elution Buffer (70 % acetonitrile, 0.1 % TFA) (x2). To ensure maximum yield of purified peptide for nano LC-MS/MS analysis, this process was repeated four times for each sample. Eluted peptides were dried down under vacuum for approximately 1 hour at 30 °C, and resuspended in 30 μL of Buffer A (3 % ACN, 0.1 % formic acid) to allow for ~3 μg of peptide per 5-μL injection on the Agilent 6550 Q-TOF mass spectrometer.

## Chromatography

### Chromatography parameters

LC System		
Analytical pump	Agilent Nano Pump (G2226A)	
Loading pump	Agilent Capillary Pump (G1376A)	
Autosampler	Agilent Micro Well-plate Sampler (G1377A)	
Column	Agilent Polaris-HR-3C18 3 $\mu$ m High Performance Chip (G4240-62030) 150 mm $\times$ 75 $\mu$ m separation column w/ 360 nL enrichment column	
Injection volume	5.00 $\mu$ L	
Autosampler temperature	4 $^{\circ}$ C	
Needle wash	10 seconds in wash port (1.0 % formic acid in 50 % methanol)	
Mobile phase	A) 0.1 % formic acid in H <sub>2</sub> O B) 0.1 % formic acid in 90 % acetonitrile	
Flow rate	300 nL/min (analytical pump) 2.5 $\mu$ L/min (loading pump)	
Loading pump solvent composition	3 % Solvent B	
Analytical pump gradient program	Time (min)	%B
	0.00	3
	60.00	25
	90.00	40
	95.00	90
	100.00	90
	100.10	3
HPLC-Chip Cube timetable	Time (min)	Inner valve position
	110.00	Enrichment
Stop time	120 minutes (analytical pump)	
Post time	3 minutes (analytical pump)	

## Software

The Agilent MassHunter data analysis suite (B.06.00) was used for qualitative analysis. The Q-TOF data were first processed using Agilent Spectrum Mill software (B.05.00.181 SP1). MS/MS spectra were searched against the SwissProt human database (downloaded on 5/18/2015) and validated at the spectral level using 1.2 % FDR as the criteria. Data were searched with trypsin specificity, with parameters shown in Table 1. Validation of peptide identifications was based on a maximum FDR of 1.2 % across each LC run, minimum peptide length of six amino acids, and precursor charge range of 2–6. These data were also searched with  $\pm 20$  ppm precursor and  $\pm 50$  ppm fragment ion tolerance. Spectrum Mill determines the area for each precursor and, following peptide spectral matching, assigns the peak area to the corresponding peptide. Protein intensity was calculated as the median intensity of the peptides from the protein.

## Mass spectrometry

### MS Parameters

Instrument	Agilent G6550A Q-TOF mass spectrometer		
Ion source	Agilent HPLC-Chip Cube (G4240A)		
Instrument mode	2 GHz, extended dynamic range, $m/z$ 1,700		
Acquisition mode	Auto MS/MS		
Polarity	Positive		
Gas temperature	280 $^{\circ}$ C		
Drying gas	11 L/min		
Capillary voltage	1,900 V		
Fragmentor	360 V		
Skimmer	45 V		
Octapole RF	750 V		
Acquisition speed	3 spectra/second (MS), 3 spectra/s (MS/MS)		
Isolation width	Narrow ( $\sim 1.3$ amu)		
Precursor selection			
Max precursor per cycle	20		
Threshold (Abs)	1,000		
Precursor abundance base scan speed	Yes		
Targets (counts/spectrum)	25,000		
MS/MS accumulation time limit	Yes		
Isotope model	Peptides		
Active exclusion	Exclude after 1 spectrum, release 0.50 minutes		
Charge state preference	2, 3, >3		
Ramped collision energy			
	Charge	Slope	Offset
	2	3.1	1
	3	3.6	-4.8
	>3	3.6	-4.8
Reference correction	$m/z$ 1,221.99063700		

Table 1. Data Search Parameters

Parameter	Value
Maximum missed cleavages	2
Fixed modification	carbamidomethylation (C)
Variable modification	oxidized methionine (M)
Minimum matched peak intensity	50 %
Precursor mass tolerance	$\pm 20$ ppm
Product mass tolerance	$\pm 50$ ppm
Maximum ambiguous precursor charge	3
Minimum detected peaks	4

Protein intensity results were imported into Agilent Mass Profiler Professional software (MPP) version B.14.8 for differential analysis. MPP parameters were as follows: Normalization: none; Baseline: z-transform; Filter by Abundance: Greater than or equal to 2 and present in 66 % of 1 in 6 conditions.

Differential results were mapped to pathways using the Pathway Architect module in MPP. The source for pathways was KEGG, and averaged results were used for mapping.

## Experimental

Figure 1 illustrates the sample preparation strategy, and is summarized as follows:

- As a model platform, the androgen-sensitive LNCaP cell line as well as two androgen-independent sub-lines, LNCaP-abl and LNCaP-abl-Hof, were treated with dimethylxalylglycine (DMOG) to mimic hypoxic conditions.
- After 8-hour and 24-hour incubations, cells were lysed in 1 % SDS buffer and digested with trypsin and Lys C according to the FASP protocol.
- Three biological replicates were prepared for each condition.
- The 8-hour incubation replicates and 24-hour replicates were pooled before trypsin digestion to create two sample quality controls (Sample QC).
- Peptides were purified using C18-packed stage tips (prepared in house).
- A technical quality control sample (Technical QC) was prepared by pooling the sample QCs following trypsin digestion and peptide cleanup.
- For LC/MS worklists, all biological replicates generated for each time point were analyzed in a randomized order along with Sample QCs and the Technical QCs at the beginning, middle, and end.




Figure 1. Experimental design for LC/MS-based identification of potential prostate cancer biomarkers. Samples were used for both discovery and targeted analyses. -T and +T denote treatment with and without dimethylxalylglycine, respectively.

## Workflow overview

Figure 2 depicts the discovery and targeted-based proteomics workflow planned for this long-term study. Key steps are summarized as follows:

- Discovery proteomics raw data files were processed through Spectrum Mill software, and a file was generated for import into MPP.
- Following differential analysis, protein expression changes were visualized within the Pathway Architect module of MPP.
- For the verification of selected proteins of interest, MRM assays will be designed using Skyline software on an Agilent 6495 triple quadrupole mass spectrometer.
- MassHunter Quantitative Analysis and MPP software will be used to quantitate and visualize MRM data, respectively.




Figure 2. Data acquisition and analysis workflow to characterize and confirm potential prostate cancer biomarkers.

## Results and Discussion

### Protein identification

LC-MS/MS analysis with the 6550 Q-TOF mass spectrometer on cell lysate generated after 8 hours of incubation led to the identification of 1,225 entities common to all sample groups, while analysis of the lysates generated after 24 hours of incubation led to the identification of 1,324 entities common to all samples groups. Comparison of the datasets revealed that 1,080 entities were commonly identified at both time points (Figure 3).





Figure 3. Venn diagram, comparing protein identifications from 8 and 24 hours of incubation, showing total proteins identified from the Sample QC replicates.



### Differential analysis

Spectrum Mill protein identification results and abundances were exported to MPP for differential analysis. MPP can normalize data against a housekeeping protein when there is significant biological variation in protein expression. Since these data are from cell lines under well controlled conditions, normalization was found to have no effect. One-way ANOVA analysis was performed across all samples with a p-value cut-off set to 0.05, fold change >2.0, and no multiple testing correction used. Principal component analysis was performed with the biological replicates from each cell line (treated and untreated) at 8 hours and 24 hours (Figure 4A). This analysis shows clear separation between many of the treated and untreated cell lines. Clustering analysis of proteins identified as significant following ANOVA analysis at 24 hours (34 differentially regulated proteins) showed clear separation between the androgen-sensitive (LNCaP) and androgen-independent (Abl, Hof) cell lines (Figure 4B). Taken together, these results suggest that differences in protein expression from the androgen-dependent versus androgen-independent cell lines are more pronounced than the expression differences induced by hypoxia.




Figure 4. Differential analysis of the three cell lines treated with DMOG and control (no drug) at 8 hours and 24 hours. A) Principal component analysis; B) Hierarchical clustering results of significant features at 24 hours.

To identify specific proteins induced by hypoxic conditions, an unpaired t-test analysis was performed between control and +DMOG replicates for each cell line. P-value cut-off was set to 0.05, fold change >2.0. For all three cell lines, a greater number of significantly changing proteins were identified after 8 hours of incubation with DMOG (Figure 5A), as opposed to 24-hour incubation (Figure 5B).




Figure 5. Volcano plots showing significant features in hypoxic (DMOG treated) versus normoxic (control) conditions for each cell line. Comparisons were made for the 8-hour (A) and 24-hour (B) timepoints.



## Pathway analysis

Significant differential proteins identified by ANOVA analysis were mapped by Pathway Architect Software to curated KEGG pathways ([www.genome.jp/kegg/](http://www.genome.jp/kegg/)). Mapping the significant proteins and the abundance results provided the ability to view pathways impacted by cell line and treatment. For instance, the oxidative phosphorylation pathway was found to be impacted following 8-hour incubation with DMOG

(Figure 6). It has been shown that oxidative phosphorylation is altered by hypoxia [6]. Cancer cells shift from oxidative phosphorylation to glycolysis as the primary means of producing ATP in conditions of metabolic stress, such as hypoxia, to maintain cellular energetics for survival [5]. This phenomenon would be worthy of further investigation to understand cancer progression under hypoxic conditions.




Figure 6. The KEGG oxidative phosphorylation pathway overlaid with averaged protein abundances from significant proteins identified in the 8-hour incubation.



## Conclusions

This application note introduces a discovery-to-targeted workflow to identify potential cancer biomarkers with Q-TOF-based label-free shotgun proteomics.

The results shown focused on the discovery phase, where Agilent Mass Profiler Professional software identified significant differential proteins between relevant cancer cell lines and growth conditions. To add biological context, these differential proteins were mapped onto biological pathways and visualized with Agilent Pathway Architect software. Targeted experiments to confirm these differential proteins through peptide MRM-based quantitation are currently underway.

## References

1. A. Fraga, *et al.* "Hypoxia and Prostate Cancer Aggressiveness: A Tale with Many Endings" *Clin. Genitourin Cancer* [Internet]. Elsevier Inc.; **13**, 295–301 (2015). doi: 10.1016/j.clgc.2015.03.006
2. L. Harrison, K. Blackwell. "Hypoxia and anemia: factors in decreased sensitivity to radiation therapy and chemotherapy?" *Oncologist* [Internet], **9 Suppl 5**, 31–40 (2004). doi: 10.1634/theoncologist.9-90005-31
3. L. Yin, Q. Hu, R. W. Hartmann. "Recent progress in pharmaceutical therapies for castration-resistant prostate cancer" *Int. J. Mol. Sci.* [Internet]. [cited 2014 Jun 3]; **14**, 13958–78 (2013). doi: 10.3390/ijms140713958
4. J. R. Wisniewski, *et al.* "Universal sample preparation method for proteome analysis" *Nature Methods*, **6(5)**, 359-62 (2009). doi: 10.1038/nmeth.1322
5. K. E. Allison, B. L. Coomber, B. W. Bridle. "Metabolic Reprogramming in the Tumour Microenvironment: a Hallmark Shared by Cancer Cells and T Lymphocytes" *Immunology*, 1–10 (2017). <https://doi.org/10.1111/imm.12777>
6. P. Sontakke, *et al.* "Hypoxia-like signatures induced by BCR-ABL potentially alter the glutamine uptake for maintaining oxidative phosphorylation" *PLoS ONE*, **11(4)**, 1–16 (2016). <https://doi.org/10.1371/journal.pone.0153226>

[www.agilent.com/chem](http://www.agilent.com/chem)

For Research Use Only. Not for use in diagnostic procedures.

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2017  
Printed in the USA  
November 20, 2017  
5991-8448EN



**Agilent Technologies**