

ISO 9001  
REGISTERED



# Capillary GC Troubleshooting: a Practical Approach

Katherine Stenerson

Gas Separations R&D

Supelco, Supelco Park, Bellefonte, PA 16823

# Outline

- ✓ Basic Troubleshooting Strategy
- ✓ Preventing Problems
- ✓ Identifying Common Problems
- ✓ Recommended Reading
- ✓ Discussion

# Troubleshooting Strategy

- ✓ Have appropriate equipment and supplies on hand.
- ✓ Establish a systematic approach.
- ✓ Know what to look for first.
- ✓ Record what you did to correct the problem.

# Suggested Equipment to Have on Hand for Troubleshooting:

---

- ✓ Electronic Leak Detector
- ✓ Flow Meter
- ✓ “Test” Column
- ✓ Replacement Accessories (Syringes, Ferrules, Septa)
- ✓ Replacement Purifiers

# Five Major Sources of Chromatographic Problems:


---

- ① Operator Error
- ② The Sample
- ③ The Column
- ④ The GC Electrical System
- ⑤ The Gas Flow System (both internal and external to the GC)

# Approaching the Problem

- ✓ Check first to see if a “fix” for the problem is already known.
- ✓ Check the Supelco Capillary GC Troubleshooting Guide (Bulletin 853.)
- ✓ Check back in the instrument maintenance record.
- ✓ Talk to others in your lab.

# Isolate the Source of the Problem



## When reviewing operating method parameters consider the following:

- ? Is my starting temp. low enough to allow sufficient sample focusing?
- ? For splitless injections, is my splitter opening at the appropriate time?
- ? Is my column flow set to give me maximum efficiency at the most critical point?
- ? Are heated zones (injectors, detectors, interfaces) set appropriately?



# The Best Way to Solve Problems is to Prevent Them!

---

- ✓ Install and maintain proper purification for *all* gases in the GC system.
- ✓ Maintain the injector by periodically inspecting and changing the liner, septa, and seal (H/P™.)
- ✓ Use the proper injection technique-this includes using the right liner for the job.
- ✓ When necessary, use a guard column to protect the analytical column.

# Gas Purification

## ✓ Carrier Gas

- At minimum, remove hydrocarbons, water, and oxygen.

## ✓ Hydrogen (FID)

- At minimum, remove hydrocarbons.

## ✓ Air (FID)

- At minimum, remove water and hydrocarbons.

## ✓ Nitrogen make-up (FID, ECD)

- At minimum, remove hydrocarbons.

## ✓ P-5 make-up (ECD)

- At minimum remove hydrocarbons (especially halogen-containing), oxygen.

# Acceptable Purity Levels for Chromatography Grade Gases

## Impurity / Maximum Concentration

Gas	O <sub>2</sub>	H <sub>2</sub> O	CO <sub>2</sub>	CO	Total Hydrocarbons
Helium	<1.0 ppm	<1.0 ppm	<1.0 ppm	<1.0 ppm	<1.0 ppm
Nitrogen	<1.0 ppm	<1.0 ppm	<1.0 ppm	<1.0 ppm	<1.0 ppm
Air	20-22%	<1.0 ppm	<1.0 ppm	<1.0 ppm	<1.0 ppm
Hydrogen	<1.0 ppm	<1.0 ppm	<1.0 ppm	<1.0 ppm	<1.0 ppm
Argon/ Methane	<1.0 ppm	<1.0 ppm	<1.0 ppm	<1.0 ppm	<1.0 ppm

# Suggested purifiers

	Hydrocarbons	Water	Oxygen
Carrier	Supelcarb HC Supelpure HC	Mole Sieve 5A	OMI -2
H <sub>2</sub>			
Air		Mole Sieve 5A	
N <sub>2</sub> makeup			
P-5	OMI -2		OMI -2

# What are some signs that my purifiers need to be changed?

## Hydrocarbon Traps

Noise in the baseline (FID)

Increase in background peaks on tune (MSD)

Higher than normal baseline reading on FID

Extra peaks visible in run

## Molecular Sieve 5A

Increase in column bleed

Water visible in MS background

Poor peak shapes for gaseous VOCs (purge and trap)

Extra peaks visible in run

OMI™-2 color change

# Injector Maintenance

✓ Change (as needed):

1. Liner and O-ring\*
2. Seal and washer\*

\*H/PTM GCs

✓ Inspect the inlet periodically

- Look for contamination in the liner
- Look for residue on the seal

# Using the right liner and injection technique can also prevent problems:

## ✓ Split Injection


- used for concentrated samples
- high flow of carrier gas through liner during injection
- should use liner designed for split injection

## ✓ Splitless Injection


- used for trace analysis
- low flow of carrier gas through liner during injection
- inertness and internal volume of liner used are critical

# Some liners used for split injection


Cup (unpacked)



Cup (wool packed)



Split/splitless- wool packed





# Some liners used for splitless injection

2 mm ID, straight



Dual-tapered



Single-tapered



## If you *must* clean a liner....

- ✓ Handle liners with gloves or forceps.
- ✓ Use *clean* compressed gas and/or a fine brush to remove particles.
- ✓ Rinse liner in an appropriate solvent and dry with *clean* compressed gas.
- ✓ Use mineral acid and/or detergent only if absolutely necessary.

# Using a Guard Column

- ✓ Choose a guard column which has been deactivated.
- ✓ Usually, the ID of the guard matches the analytical column.
- ✓ A 5-10 meter length is normally used.
- ✓ Connect with either a GlasSeal™ or butt connector.


# Common Problems

- 1 Poor Peak Shapes (either tailing, fronting, or just generally ugly.)
- 2 Nonlinearity
- 3 Baseline Noise and /or Drift
- 4 Ghost Peaks
- 5 Missing Peaks / Poor Response
- 6 Insufficient Resolution

# 1. Poor Peak Shapes

---

✓ Fronting can indicate column overload.




✓ Tailing can indicate activity in the system or improper column installation.



# 1. Poor Peak Shapes (cont.)

---

- ✓ Generally ugly peaks, such as  $\alpha,\alpha$ -dimethylphenethylamine, can be caused by a variety of problems.



## 2. Nonlinearity

---

The most common causes are:

- ➔ Column overload
- ➔ Detector overload
- ➔ Standards preparation
- ➔ Poor peak shape resulting in improper integration

# An Example of Overload:





## Preventing column overload for splitless injections:


---

- ✓ Inject a smaller amount / use a 1 ul syringe.
- ✓ Use a thicker film column.
- ✓ Use a column with a wider ID.
- ✓ Decrease upper limit of calibration range.
- ✓ Use a column of slightly different polarity.

# An example of poor peak shape affecting linearity:

---

- ✓ Benzoic acid is typically of poor shape when doing splitless injections.





### 3. Common causes of baseline noise / drift.

---

- ✓ Column bleed
- ✓ Dirty detector
- ✓ Contaminants in carrier gas / carrier gas purity

# Effect of carrier gas purity on baseline noise:






Column bleed results from the normal degradation of the stationary phase.

---

- ✓ All columns bleed to some extent.
- ✓ Bleed increases with temperature.
- ✓ The amount of bleed will increase in the presence of oxygen.


# A Typical Bleed Profile:



# Column Bleed on an MSD


- ✓ Visible as baseline rise in the TIC.
- ✓ Check spectra for key bleed ions:
  - PTE™-5: 207, 281
  - SPB™-1: 73, 207, 281
  - SPB™-624: 207, 269
- ✓ Make sure interface temp. is < column max. temp.

Ion 207 corresponds to a fragment known as D3:






# MS spectra of bleed from a PTE™-5 Column




# MS spectra of bleed from an SPB™-1 Column




# MS spectra of bleed from an SPB™-624 Column

Abundance



m/z-->

# MS Spectra of Septa Bleed




# Prevent column bleed!

- ✓ Sufficiently purge column with carrier before ramping it up in temperature.
- ✓ Make sure carrier gas is filtered for water and oxygen.
- ✓ Check integrity of all fittings leading to the column.
- ✓ Do not heat the column above its maximum temp.
- ✓ Precondition the column prior to use.


## 4. Ghost Peaks

- ✓ Residue in the inlet liner and at the head of the column
- ✓ Contaminated syringe / and or wash solutions on an autosampler.
- ✓ Sample carryover
- ✓ Contaminated carrier gas

# If pieces of septa get into an inlet liner...



...even a simple analysis can be ruined.






## 5. Missing Peaks / Poor Response

- ✓ Sample decomposition
  - Activity in the inlet or column
  - Injection port temperature too high
  - Sample not stable enough for GC
  - Standards not stable
- ✓ Coelution
- ✓ Insufficient run time / final temperature
- ✓ Sample not volatile enough for GC
- ✓ Improper column installation


# Nasty samples can damage a column by creating active sites:



# Here, the responses of some acids were affected:


Abundance

After sample injection




Time-->


Response can also be affected by the position of the column in the inlet:



Here, the column was not inserted far enough:



Here, the column was inserted too far:



## 6. Insufficient Resolution

### ✓ Wrong column

- Longer columns increase resolution.
- Smaller ID columns increase resolution.
- A different phase altogether may be needed.

### ✓ Wrong Conditions

- Carrier gas flow too fast or slow .
- Oven ramp rate too fast.

# Recommended Reading Supelco Bulletins

1. **741**: The Supelco Guide to Leak-Free Connections
2. **783**: Cleaning Flame Ionization Detectors
3. **853**: Capillary Troubleshooting Guide
4. **875**: Supelco Capillary GC Selection Guide
5. **895**: Installation and Maintenance Instructions for .25 mm and .32 mm ID Fused Silica Capillary Columns
6. **897**: Installation and Maintenance Instructions for .53 mm ID Fused Silica Capillary Columns
7. **898**: Gas Management Systems for GC
8. **899**: Capillary GC Inlet Sleeve Selection Guide
9. **916**: Purge and Trap System Guide
10. **918**: Selecting Purifiers for Gas Chromatography



Help is just a phone call or mouse click away!

✓ **Supelco Technical Service**

phone: 1-800-359-3041

email: [techservice@supelco.sial.com](mailto:techservice@supelco.sial.com)

✓ **Supelco Customer Service**

phone: 1-800-247-6628

email: [supelco@sial.com](mailto:supelco@sial.com)

✓ **Sigma-Aldrich Website**

[www.sigma-aldrich.com](http://www.sigma-aldrich.com)