

DETERMINATION OF QUANTITATIVE PROTEIN SIGNATURES FOR DUCTAL CARCINOMA (BREAST CANCER) BY LC/MS PROTEOME ANALYSIS

Hans Vissers¹, Marc Kipping¹, Toralf Reimer², Annika Kasten², Cornelia Koy², Jim Langridge¹, Michael Glocker²

¹ Waters Corporation, MS Technologies Centre, Manchester, United Kingdom

² Proteome Center Rostock, Department for Proteome Research, University of Rostock, Germany

INTRODUCTION

With the availability of the human genome sequence, data-driven research for tackling the molecular grounds of multi-factorial, polygenic diseases can be considered a realistic challenge to the scientific community. In most recent research projects, protein expression profiles are obtained using sophisticated MS-based equipment, producing read-outs termed protein signatures rather than single protein markers.

In this study, a comprehensive MS-based discovery strategy is applied for a polygenic disease. The method employs the separation and detection of non-labelled tryptic fragments by means of an LC/MS acquisition. During the acquisition, the collision energy within the gas cell is continuously switched from low to elevated energy and no precursor isolation is applied.

The low-energy functions contain all detectable peptide pseudo molecular ions. In a complementary fashion, the resulting high-energy data provides extensive multiplexed fragmentation information. The high-energy fragment ions are aligned to their related precursor ions in chromatographic space by time and profile. Relative quantification is achieved via normalization of the MS datasets and comparison of the peptide intensities across injections and between samples. Identification of peptides exhibiting a change in expression level is made using the peptide exact mass and the fragment ion information from the high-energy dataset.

An initial study was conducted on a small patient group. Quantitative multi-variance analysis was performed. Initial results on samples from patients who suffer from ductal carcinoma, breast cancer, indicate that expression levels of the newly-found potential protein signatures might become useful in diagnosis and possibly prognosis.

EXPERIMENTAL

Sample preparation

75 μL of breast cancer and healthy tissue protein extract samples were taken up in 50 mM NH_4HCO_3 , 0.1% RapiGest™ SF, pH 8.5 to a final concentration of $\sim 1 \mu\text{g}/\mu\text{L}$. Reduction and alkylation was with 2.5 μL 100 mM DTT and 2.5 μL 300 mM IAA, respectively. The proteins were digested with 1:25 (w/w) sequence grade trypsin overnight (16 hr). Trypsin was added immediately after the addition of DTT and IAA to limit endogenous protease activity.

RapiGest was removed by the addition of 2 μL conc. HCl, followed by centrifugation, and the supernatant collected. Samples were diluted with 0.1% formic acid to an appropriate final working concentration prior to analysis, corresponding to an 0.7 μg of protein digest on-column load.




Figure 1. X-ray pictures (mammograms) illustrating the location of the cancer and healthy tissue; blue = cancer tissue; yellow = healthy, unaffected control tissue.

LC/MS conditions

LC/MS quantification experiments were conducted using a 1.5 hr reversed-phase gradient at 250 nL/min (5 to 40% acetonitrile over 90 min) on the Waters® Expression^E High Definition Proteomics™ System, using as an inlet the nanoACQUITY UPLC® System and an Atlantis® 3µm C₁₈ NanoEase™ 75 µm x 15 cm nanoscale LC column. Samples were run in triplicate.

The Expression^E System also included the Q-ToF Premier™ Mass Spectrometer, which was programmed to step between normal (5 eV) and elevated (25 to 40 eV) collision energies on the gas cell, using a scan time of 1.5 s per function over 50 to 1990 m/z. Protein identifications and quantitative information were generated by the use of dedicated algorithms, part of the Expression^E System informatics, and searching human-specific databases.

RESULTS

Protein and peptide replication

Figure 2 displays the peptide replication rate across both healthy and tumor tissue for patient A.




Figure 2. Accurate mass retention time pair replication rate (peptide count) vs. injection replicate # (tumor and unaffected conditions) for patient A; Grey = replication rate <math>< 2/3</math> injections. Blue = replication rate >math>> 2/3</math> injections (unique to tumor, 6675; unique to unaffected, 3204; common to both conditions, 7858).

Figures 3 and 4 show the peptide and protein replication rates for the healthy tissue from a single patient.




Figure 3. Peptide replication rate for condition G (unaffected control tissue), patient A; from 2 of 3 injections = 57%; from all 3 injections = 34%.




Figure 4. Protein replication rate condition G (unaffected control tissue), patient A; from 2 of 3 injections = 50%; from all 3 injections = 34%.

Clustering analysis

Figures 5, 6, and 7 display the PCA plots obtained from analysis of the low-energy precursor ion (peptide) information from the LC/MS experiments; retention time, mass, and intensity.




Figure 5. PCA scores plot of accurate mass retention time pairs (peptides) from all patients (A, B, and C) and conditions (tumor, T and unaffected, G); $t[1]$ = 1st PCA component; $t[2]$ = 2nd PCA component. The plot shows clear separation between unaffected and tumor tissue from each patient.




Figure 6. Loadings PCA plot of accurate mass retention time pairs (as in Figure 5); $M1. p[1]$ = 1st PCA component; $M1. p[2]$ = 2nd PCA component.




Figure 7. Log-log intensity visualization of accurate mass intensity pairs from Figure 5 and 6 (condition G (x-axis) vs. condition T (y-axis) for patient A). Cluster annotations as for PCA in Figure 6. Interestingly, the peptides (unique to patient and condition shown above) visualized in this plot are of relative high abundance.

Relative quantification

Presented in Figure 8 is a log-log plot of the peptide precursor ion intensity between the unaffected and tumor tissue for a single patient. Displayed are those ions that are statistically up- or down-regulated ($p < 0.05$ and $p > 0.95$). The peptides were subsequently searched utilizing both the peptide mass and fragment ion information, of which an example is shown in Figure 9.




Figure 8. Significantly up- and down-regulated peptides shown here were selected for protein identification by database search. Log-log intensity condition T vs. G of patient C.

Cluster	Intensity	Average MS	Average RT	OK	P	Score	Locus	CT2/CO2
12591	933 5134	54.0	33.5	...	342 (1.23+1.21) [0.98]
1651	2494 1516	66.3	33.5	...	536 (1.17+1.22) [1.00]
46	1622 9698	70.1	8.7	...	0.07 (2.68+0.14) [0.00]
953	1598 8516	47.2	56.9	...	0.34 (1.24+0.24) [0.00]
1259	1598 8502	29.5	56.9	...	0.31 (1.17+0.47) [0.00]
1259	1352 7148	47.1	56.9	...	0.34 (1.09+0.31) [0.00]
664	1309 7226	68.2	56.9	...	0.44 (0.93+0.29) [0.00]
1350	1283 7883	58.7	42.4	...	0.08 (2.88+1.13) [0.00]
894	1395 5547	34.7	42.4	...	0.05 (2.88+1.05) [0.00]
2068	1625 7468	68.8	114.2	...	20.20 (2.31+0.20) [1.00]
6642	1213 6138	48.5	114.2	...	14.73 (2.69+0.36) [1.00]
1412	1643 9408	64.8	114.2	...	10.26 (2.24+0.25) [1.00]
13330	2035 8420	49.1	40.2	...	4.76 (1.58+0.43) [0.99]
13838	1789 9930	61.6	40.2	...	19.20 (2.88+0.78) [1.00]
13895	1634 7858	48.5	40.2	...	11.25 (2.42+0.72) [1.00]
1871	1956 8348	59.1	40.2	...	8.08 (2.09+0.25) [1.00]
2432	2406 1824	48.1	47.7	...	9.45 (2.25+0.24) [1.00]
4701	1622 8668	57.4	57.7	...	11.47 (2.44+0.35) [1.00]
5274	1107 6047	45.1	57.7	...	10.70 (2.37+0.47) [1.00]
2228	1211 6884	47.6	57.7	...	5.21 (1.89+0.25) [1.00]
1598	1598 7457	55.8	229.7	...	0.96 (1.84+0.17) [1.00]
1968	1154 5728	52.4	229.7	...	7.85 (2.38+0.22) [1.00]
2491	1055 5377	56.7	229.7	...	7.46 (2.01+0.44) [1.00]
2493	1278 5812	54.2	229.7	...	8.23 (2.12+0.24) [1.00]
8402	1247 8182	43.4	229.7	...	9.43 (2.28+0.54) [1.00]
8488	848 4983	38.2	229.7	...	3.71 (1.31+0.52) [1.00]
555	1646 8927	78.0	229.7	...	3.23 (1.28+0.31) [1.00]
153	2663 6022	71.4	154.2	...	0.23 (0.48+0.13) [1.00]
158	2813 3844	74.3	154.2	...	0.23 (0.48+0.17) [0.98]
620	2238 1022	61.8	154.2	...	0.12 (2.12+0.28) [0.66]
767	1038 6138	43.7	154.2	...	0.15 (1.68+0.31) [0.68]
758	2170 5778	70.3	154.2	...	0.11 (2.28+0.48) [0.66]
3201	872 4738	36.3	154.2	...	0.20 (0.62+0.71) [0.68]
71	1813 3338	74.8	144.2	...	0.29 (0.33+0.12) [0.86]
2980	1778 8394	45.6	272.2	...	1.60 (0.52+0.29) [0.99]
4152	2984 2878	78.0	272.2	...	0.30 (2.24+1.93) [1.00]
4971	1400 7782	68.1	272.2	...	60.91 (0.93+0.38) [1.00]
11589	1823 8584	56.5	272.2	...	34.47 (3.54+0.94) [1.00]
13837	2292 1975	71.0	272.2	...	49.40 (2.89+0.65) [1.00]
14094	1174 6497	52.4	272.2	...	14.01 (2.88+0.72) [1.00]
2698	2354 3315	68.2	272.2	...	47.47 (3.84+0.30) [1.00]
10895	1195 8181	48.4	10.1	...	10.52 (2.38+0.71) [1.00]
3743	2227 1802	68.3	10.1	...	16.44 (2.89+0.38) [1.00]
7269	1588 7251	24.3	10.1	...	4.90 (1.59+0.59) [1.00]
1386	1819 8834	67.7	10.1	...	2.61 (0.97+0.25) [1.00]
2578	587 3722	42.1	21.0	...	0.38 (1.03+0.59) [0.80]
340	3658 3888	73.4	9.2	...	3.85 (1.03+0.81) [0.80]
5913	1385 7846	48.8	9.2	...	0.55 (1.88+0.83) [1.00]
351	1821 8272	51.8	15.9	...	0.88 (0.22+0.11) [1.00]

Figure 9. Protein identification results with peptides from each protein listed (accurate mass retention time pair search results of high-energy fragmentation data). Peptides identified from down-regulated osteoinductive factor protein in condition T of patient C are highlighted.

[APPLICATION NOTE]

In Figure 10, peptide-level information is displayed in further detail for the osteoinductive factor protein highlighted in Figure 9. Of particular note is the consistency of the intensity profile across the peptides (blue bars). These peptides are then annotated on a log-log precursor ion distribution plot, which is shown in Figure 11.




Figure 10. Measure of the consistency in relative abundance (blue), mass measurement error (yellow), and protein identification probability (green) for the highlighted peptides (shown here by cluster number) in Figure 9. The results display a high consistency for all peptides matched to osteoinductive factor protein.




Figure 11. Log-log intensity condition T vs. G of patient C. Annotations of peptides from osteoinductive factor (identified in Figures 8 and 9) is shown to be down-regulated in tumor tissue of all patients.

Protein signatures

Absolute protein amounts were estimated and expressed as the $^2\log$ ratio vs. a protein spike at the 150 fmol level, providing both an instrument-specific absolute concentration response factor and condition signatures. These condition signatures are displayed in Figure 12, where it can be clearly seen that these profiles are consistent across patients A and C. These signatures do not require comparative analysis and therefore could easily be extended for larger scale studies.




Figure 12. Protein signatures from unaffected and tumor tissue for patients A and C. Colors represent different proteins identified across all patients and conditions and their absolute concentration.

CONCLUSIONS

The presented approach takes into consideration that both tumor and control materials are from the same persons.

Inter-individual protein abundance differences are excluded by a stringent selection procedure in which:

- Only those proteins are considered as disease-related that occur in all three investigated patients
- Only those proteins are considered that are regulated synchronously in all three comparisons

Future work

Future works will focus on the complementary analysis of invasive ductal carcinoma samples by means of 2D PAGE analysis followed by MALDI-TOF PMF analysis. The same selection criteria will be applied as within the approach presented in this study, allowing us to identify proteomics analysis technique independent signatures.

Acknowledgements

Kieran Neeson and Scott Geromanos are kindly acknowledged for their assistance and contributions.

References

Silva JC, Denny R, Dorschel CA, Gorenstein M, Kass IJ, Li GZ, McKenna T, Nold MJ, Richardson K, Young P, Geromanos S. Quantitative Proteomic Analysis by Accurate Mass Retention Time Pairs. *Anal Chem*. 2005 Apr 1; 77 (7): 2187-200.

Silva JC, Gorenstein MV, Li GZ, Vissers JP, Geromanos SJ. Absolute Quantification of Proteins by LC/MS^E; a Virtue of Parallel MS Acquisition. *Mol Cell Proteomics*. 2006 Jan; 5(1): 144-56.

Hughes MA, Silva JC, Geromanos SJ, Townsend CA. Quantitative Proteomic Analysis of Drug-Induced Changes in Mycobacteria. *J Proteome Res*. 2006 Jan; 5 (1): 54-63.

Silva JC, Denny R, Dorschel C, Gorenstein MV, Li GZ, Richardson K, Wall D, Geromanos SJ. Simultaneous Qualitative and Quantitative Analysis of the *Escherichia coli* Proteome: A Sweet Tale. 2006 Apr; 5 (4): 589-607.

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™

Waters, Atlantis, and nanoACQUITY UPLC are registered trademarks of Waters Corporation. Expression[®] High Definition Proteomics, Q-ToF Premier, RapiGest, and NanoEase are trademarks of Waters Corporation. All other trademarks are the property of their respective owners.

©2006-2007 Waters Corporation Produced in the U.S.A.
Nov. 2007 720001831EN AG-PDF

Waters Corporation
34 Maple Street
Milford, MA 01757 U.S.A.
T: 1 508 478 2000
F: 1 508 872 1990
www.waters.com

