

Advantages of Oasis PRiME HLB for the LC-MS/MS Analysis of Vitamin D Metabolites in Serum for Clinical Research

Robert Wardle and Lisa Calton
Waters Corporation, Wilmslow, UK



GOAL

To successfully develop an LC-MS/MS method for the analysis of 25-hydroxyvitamin D₃ (25OHD₃), 25-hydroxyvitamin D₂ (25OHD₂), 24,25-dihydroxyvitamin D₃ (24,25(OH)₂D₃), and C3-epi-25-hydroxyvitamin D₃ (C3-epi-25OHD₃), using the Xevo® TQ-S micro Mass Spectrometer and utilizing Waters® Oasis® PRiME HLB to minimize phospholipid interferences and improve analytical sensitivity.

BACKGROUND

Matrix interferences can be challenging when measuring serum levels of vitamin D metabolites by LC-MS/MS. In particular, lysophosphatidylcholines (LysoPCs 16:0, 18:1, and 18:0) – which have similar hydrophobicity to 25-hydroxyvitamin D (25OHD) – have been shown to cause ion suppression in mass spectrometric methods. Despite being structurally very different, LysoPCs have proved challenging to remove during sample preparation and are difficult to separate chromatographically from vitamin D metabolites.

For LC-MS/MS vitamin D metabolite analysis, Oasis PRiME HLB minimizes phospholipid interferences and improves analytical sensitivity, without the need for lengthy derivitization.



Figure 1. Oasis PRiME HLB μElution solid phase extraction (SPE) Plate.

THE SOLUTION

Here we describe an approach for method optimization using Oasis PRiME HLB μElution solid phase extraction (SPE) Plates (Figure 1) that assesses the reduction of LysoPCs compared to Oasis HLB SPE and protein precipitation. Using a Waters ACQUITY UPLC® HSS PFP Column, chromatographic separation of 25OHD₂, 25OHD₃, 24,25(OH)₂D₃, and C3-epi-25OHD₃ was achieved and peak area profiles were compared with targeted LysoPCs (16:0, 18:1, and 18:0, having precursor ions of m/z 496, m/z 522, and m/z 524 and product ions of m/z 184) for the three sample extraction methods. Performance of the optimized extraction method was assessed using a Waters ACQUITY UPLC I-Class (FTN)/ Xevo TQ-S micro System.

SAMPLE PREPARATION AND ANALYSIS

Stable labeled internal standards were added to 100 μL of calibrator, QC, or sample and protein precipitation performed using a methanol/zinc sulfate_(aq) solution.

Following centrifugation, the supernatant was transferred to an Oasis PRiME HLB μ Elution Plate. An elution profile was generated by analyzing the eluate from varying concentrations (0–100%) of acetonitrile_(aq) with triplicate preparations for each level passed through the SPE sorbent.

An ACQUITY UPLC I-Class System was used to separate diluted eluate using a 2.1 x 100 mm ACQUITY UPLC HSS PFP Column with a water, methanol, ammonium acetate, and formic acid gradient and analyzed on a Xevo TQD using the MRM transitions and precursor scan parameters listed in Table 1.

RESULTS

Mean peak areas from the triplicate preparations for each vitamin D metabolite and LysoPCs tested were plotted on scatter charts to view their elution profiles (Figure 2).

The elution profile demonstrates that LysoPCs begin to elute from the Oasis PRiME HLB sorbent at a similar organic concentration to the vitamin D metabolites, as they have very similar hydrophobic properties. However, at concentrations >90% acetonitrile_(aq) almost all LysoPCs are retained on the Oasis PRiME HLB SPE sorbent, while the vitamin D metabolites are released. Therefore, 25% acetonitrile_(aq) wash and 100% acetonitrile elution conditions were selected for the optimized method.

LysoPCs and vitamin D metabolites mean peak areas obtained from the optimized Oasis PRiME HLB protocol were compared to a simple protein precipitation extraction using methanol/zinc sulfate_(aq) and an optimized Oasis PRiME HLB protocol, which used methanol as the elution solvent. The results are summarized in Figure 3.

Using Oasis PRiME HLB with acetonitrile, >99% of all targeted LysoPCs were removed when compared to Oasis HLB with methanol and protein precipitation. This resulted in an increase of 5x in vitamin D metabolites peak areas when compared to Oasis HLB and an increase of 10x when compared to protein precipitation.

Analyte	Transition (m/z)	Cone voltage (V)	Collision energy (eV)
25OHD ₃ and C3-epi-25OHD ₃	401.3 > 159.1 (Quan)	24	24
	401.3 > 365.3 (Qual)	24	10
[² H ₃]-25OHD ₃ and [² H ₃]-C3-epi-25OHD ₃	404.3 > 162.1 (ISTD)	24	24
25OHD ₂	413.3 > 355.3 (Quan)	24	10
	413.3 > 83.1 (Qual)	24	24
[² H ₃]-25OHD ₂	416.3 > 358.3 (ISTD)	24	10
	417.3 > 159.1 (Quan)	20	24
24,25(OH) ₂ D ₃	417.3 > 159.1 (Quan)	20	24
	417.3 > 381.3 (Qual)	20	10
[² H ₆]-24,25(OH) ₂ D ₃	423.3 > 159.1 (ISTD)	20	24
LysoPCs	Precursors of m/z 184	20	20

Table 1. MRM transitions and precursor scan parameters.




Figure 2. Elution profile plots for vitamin D metabolites and LysoPCs with different concentrations (0–100%) of acetonitrile being passed through the Oasis PRiME HLB SPE Plate.




Figure 3. Comparison of vitamin D metabolites and LysoPC peak areas from samples prepared by protein precipitation (PPT), Oasis HLB (HLB MeOH), and Oasis PRiME HLB (PRiME MeCN).

Performance of the optimized extraction method was assessed using an ACQUITY UPLC I-Class/Xevo TQ-S micro System. Over five days, all calibration curve correlation coefficients (r^2) were >0.99 across the following concentration ranges (Table 2) for 25OHD₃, 25OHD₂, C3-epi-25OHD₃ and 24,25(OH)₂D₃, with conversion factors from nmol/L to ng/mL being shown.

Separation of the C3-epi-25OHD from 25OHD was achieved for both 25OHD₂ and 25OHD₃ using the ACQUITY UPLC HSS PFP Column. An additional isobaric peak present in both the quantifier and qualifier transitions for 24,25(OH)₂D₃ was also separated and confirmed to be 25,26(OH)₂D₃. An example chromatogram is shown in Figure 4, demonstrating a runtime of <8 minutes injection-to-injection.

Total precision and repeatability of the method was assessed by extracting and quantifying serum samples using five replicates at low, mid and high concentrations across five days (n=25). All results were ≤8.2%CV as shown in Table 3.

Accuracy was assessed by analyzing DEQAS (Vitamin D External Quality Assessment Scheme) and NIST SRM972a (National Institute of Standards and Technologies Standard Reference Materials 972a) samples and calculated

Analyte	Calibration curve range (nmol/L)	Conversion factor (nmol/L to ng/mL)
25OHD ₃	2.5 – 374	0.401
25OHD ₂	2.4 – 363	0.413
C3-epi-25OHD ₃	0.50 – 74.9	0.401
24,25(OH) ₂ D ₃	0.48 – 72.0	0.417

Table 2. Calibration curve ranges for each vitamin D metabolite.




Figure 4. Chromatogram of an extracted serum sample containing endogenous 25OHD₂, 25OHD₃, 24,25(OH)₂D₃, C3-epi-25OHD₃ and 25,26(OH)₂D₃.

Compound	Low			Mid			High		
	Conc. (nmol/L)	Total (%CV)	Repeatability (%CV)	Conc. (nmol/L)	Total (%CV)	Repeatability (%CV)	Conc. (nmol/L)	Total (%CV)	Repeatability (%CV)
25OHD ₃	4.9	5.9	4.6	24.5	4.0	3.8	243	3.7	2.3
25OHD ₂	4.7	5.4	3.7	23.8	3.4	3.4	235	3.2	2.3
C3-epi-25OHD ₃	1.0	8.2	7.0	4.6	4.6	4.6	49.0	3.8	3.0
24,25(OH) ₂ D ₃	0.97	4.8	4.3	4.4	4.2	3.4	47.3	3.1	2.6

Table 3. Total precision and repeatability assessment for the analysis of vitamin D metabolites in serum.

concentrations were compared to the NIST assigned values. For 25OHD₃, a Deming regression equation of $y=1.02x-1.09$ and Altman Bland analysis demonstrated good agreement with minimal bias (-0.1%). Only seven data points were available for 25OHD₂, with an overall mean percentage difference of -0.9% (range -9.7%–10.0%) obtained when compared to the NIST assigned values. A large scatter was observed for C3-epi-25OHD₃ DEQAS samples, but there was good agreement when compared to NIST SRM972a material, having a bias of within ±6.7% for all samples tested (triplicate preparations of three samples containing C3-epi-25OHD₃). A large scatter was also observed in 15 samples assessed for 24,25(OH)₂D₃, when compared to LC-MS/MS ALTM values of a DEQAS pilot study (NIST assigned values are not available). The LC-MS/MS ALTM was derived from only six laboratories, with returned results indicating a large standard deviation between the participating laboratories.

Analytical sensitivity was assessed by extracting ten replicates of stripped serum samples spiked at low to high concentrations over three days (n=30). A precision of <20%CV and S:N (ptp) of >10:1 were obtained at 1nmol/L for 25OHD₃, 25OHD₂ and 0.5nmol/L for 24,25(OH)₂D₃ and C3-epi-25OHD₃.

SUMMARY

A clinical research method has successfully been developed for the LC-MS/MS analysis of 25OHD₃, 25OHD₂, C3-epi-25OHD₃ and 24,25(OH)₂D₃ in serum. The advantages of the phospholipid removal properties of Oasis PRiME HLB have been demonstrated for the analysis of these vitamin D metabolites and include the following benefits:

- Removal of phospholipids from serum samples using a simplified extraction workflow
- Excellent precision performance for all analytes at all concentration levels tested (≤8.2% CV)
- Good agreement was observed to the DEQAS NIST assigned values for 25OHD₃
- Increased analytical sensitivity for 25OHD₃, 25OHD₂, C3-epi-25OHD₃, and 24,25(OH)₂D₃ negating the need for derivatization

For Research Only. Not for use in diagnostic procedures.

Waters

THE SCIENCE OF WHAT'S POSSIBLE.®

Waters, The Science of What's Possible, ACQUITY UPLC, Xevo, and Oasis are registered trademarks of Waters Corporation. All other trademarks are the property of their respective owners.

©2017 Waters Corporation. Produced in the U.S.A. May 2017 720005965EN TC-PDF

Waters Corporation
34 Maple Street
Milford, MA 01757 U.S.A.
T: 1 508 478 2000
F: 1 508 872 1990
www.waters.com