

Application Note # CA-270384

Fast Analysis of Paraffins, iso-Paraffins, Olefins, iso-Olefins, Naphthenes and Aromatics in Hydrocarbon Streams

Introduction

The Bruker PIONA+™ Analyzer is a development to characterize the complete hydrocarbon sample composition, including individual oxygenates, in a single analysis of spark ignition engine fuels by multi-dimensional gas chromatography. The sample is separated in the component groups per carbon number and in individual components through the use of multiple columns and traps. In the PIONA+ system, paraffins, iso-paraffins and iso-olefins, olefins, naphthenes and aromatics are identified. However, analysis time in the PIONA mode is about 180 minutes, which severely limits the number of samples that can be analyzed per day. A unique aspect of the design of the Bruker PIONA+ Analyzer is the ability to independently heat the individual traps (concurrent heating). This application note describes work to determine if concurrent heating could be used to significantly reduce the analysis time and improve sample throughput.

Instrumentation:

Bruker PIONA+ Analyzer
CompassCDS Chromatography software
PIONA+ plug-in software

Results and discussion

As shown in Figure 1, when the PIONA+ system is operated in the conventional mode, a total analysis time of 180 minutes is required to elute all of the component groups. However, as shown in Figure 2, through the use of concurrent heating of both Bruker Molsieve 5A and 13X traps, a reduction in total analysis time from 180 minutes to 95 minutes is obtained, i.e. ~50 % less.




Figure 1: Chromatogram of a test sample 1 in conventional PIONA mode.

The reduction in analysis time is achieved by carefully optimizing the temperature settings of the various traps and columns and invoking simultaneous heating of the Molsieve traps. This approach, referred to as the "FastPIONA" mode, results in the elution of the paraffins immediately after their naphthene and iso-paraffin counterparts. Yet the elution integrity of the component groups remains intact with no negative influence on either naphthene or iso-paraffin groups. Figure 3 shows a close-up of the resulting elution sequence. The same results occur for the olefin group separations. The cyclic and iso-olefins normally elute from 100–130 minutes and the n-olefins, in conventional PIONA mode, from 150–180 minutes.

However, through the use of concurrent heating they completely elute between 55 and 70 minutes (Figure 4). In the example of Figure 5, a commercial standard was analyzed for O-PIONA (oxygenates, paraffins, iso-paraffins, olefins, iso-olefins, naphthenes and aromatics) utilizing the concurrent heating approach described above. Note that the oxygenates elute separately from the other component groups and are easily identified and calculated by the PIONA+ software.




Figure 2: Chromatogram of the test sample in fast PIONA mode.




Figure 3: Close-up view of C8 and C9 components. Separation in both conventional and fast PIONA mode.

Conclusion

Temperature programming a Molsieve 5A trap (concurrent heating) together with a Molsieve 13X trap greatly reduces the total analysis time of PIONA+ type analysis and, in some cases, can improve the performance. Total analysis time is reduced by one third in the case of O-PIONA and by as much as half for PIONA. The reduction in analysis time is possible without any hardware changes and can be achieved without any negative effect on the quality of the component class separation. In addition, the analyses remain compliant with the method EN ISO 22854:2008.




Figure 4: Close-up view of olefins in conventional and fast PIONA mode.




Figure 5: Chromatogram of a commercial standard in fast O-PIONA mode.

References

EN ISO 22854:2008 Liquid petroleum products. Determination of hydrocarbon types and oxygenates in automotive-motor gasoline. Multidimensional gas chromatography method.

Keywords	Instrumentation & Software
EN ISO 22854	Bruker PIONA+ Analyzer
spark ignition engine fuels	CompassCDS Chromatography Software
hydrocarbons	PIONA+ plug-in software
oxygenates	

For research use only. Not for use in diagnostic procedures.

● Bruker Daltonik GmbH

Bremen · Germany
Phone +49 (0)421-2205-0
Fax +49 (0)421-2205-103
sales@bdal.de

www.bruker.com

Bruker Daltonics Inc.

Billerica, MA · USA
Phone +1 (978) 663-3660
Fax +1 (978) 667-5993
ms-sales@bdal.com

Fremont, CA · USA
Phone +1 (510) 683-4300
Fax +1 (510) 490-6586
cam-sales@bruker.com