
Metrosep RP 2 Guard – 6.1011.030

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Guard column for anion and cation columns; comprising:

- 2 × **A** Thread stopper
- 1 × **B** Filter screw
- 1 × **C** Frit
- 1 × **D** RP 2 filter
- 1 × **E** Filter housing
- 1 × **F** Connector

Connection to separation column

1. Remove thread stoppers **A** from filter screw **B** and filter housing **E**.
2. Connect the separation column to the filter screw **B** using the connector **F**.
3. Connect eluent supply to the filter housing **E**.

Exchange of filters

The filters **C** and **D** built into the guard column have to be replaced at a considerable pressure rise or no later than after 4 weeks. The 6.1011.130 spare filters available as an option (10 pieces each) have to be replaced as follows:

1. Unscrew filter screw **B** from filter housing **E** using two wrenches.
2. Remove old RP 2 filter **D** from filter housing **E**.
3. Remove old frit **C** from filter screw **B** using a pointed tool.
4. Hold filter screw **B** vertically and put new frit **C** into the depression on the screw.
Note: To avoid leakage, press frit **C** completely into the depression (e. g. by using reverse side of 6.2617.010 tool for piston seal).
5. Place new RP 2 filter **D** onto frit **C**.
6. Screw filter housing **E** onto filter screw **B** from above and tighten by hand.
7. Tighten filter housing **E** for an additional ¼ turn using two wrenches.