6.01051.2X0 Metrosep C 6 - XXX/2.0

6.01051.210	Metrosep C 6 - 100/2.0
6.01051.220	Metrosep C 6 - 150/2.0
6.01051.230	Metrosep C 6 - 250/2.0

DE**Säulenmaterial**

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EN

Column material

Silica gel with carboxyl groups, particle size 5 µm

Dimensions

6.01051.210	100 x 2.0 mm
6.01051.220	150 x 2.0 mm
6.01051.230	250 x 2.0 mm

pH range

2 - 7 (eluent)

2 - 10 (sample)

Maximum pressure

20.0 MPa (200 bar)

Maximum flow rate

6.01051.210	1.0 mL/min
6.01051.220	0.7 mL/min
6.01051.230	0.4 mL/min

Recommended flow rate

0.25 mL/min (for all lengths)

Application

Determination of monovalent and divalent cations. Optimized separation of sodium and ammonium with sizable differences in concentrations. Large-volume injections.

Eluent

- Standard eluent:
Mixture made up of 1.7 mmol/L dipicolinic acid and 1.7 mmol/L nitric acid.

- Permitted eluent additives:
0 - 100% acetone and acetonitrile, no alcohol.
- Permitted organic additives for the sample matrix:
0 - 100% acetone, acetonitrile and alcohol.

Preparation

- Rinse the column with eluent for 2 h.

Guard column

- Metrosep C 6 Guard/2.0 (6.01051.600)
- Metrosep C 6 S-Guard/2.0 (6.01051.610)

Storage

Store the column in standard eluent at ambient temperature.

Regeneration

NOTE

Ensure that the maximum pressure is never exceeded during regeneration.
If the pressure becomes too high, reduce the flow rate.

1. Disconnect the column outlet from the detector inlet.
2. Rinse the column with ultrapure water before and after regeneration.
3. Depending on the type of contamination, regenerate the column as follows:
 - *Organic contamination*
Regenerate the column in the direction opposite to the flow at a flow rate of 0.25 mL/min with the following solutions in succession:
 - 1 h with ultrapure water.
 - 1 h with acetonitrile-water mixture (40:60).
 - 1 h with ultrapure water.
 - *Inorganic contamination*
Regenerate the column in the direction opposite to the flow at a flow rate of 0.25 mL/min for 1 h with 10 mmol/L nitric acid + 4 mmol/L dipicolinic acid.

General notes

- Sample solutions must be microfiltered (0.45 µm).
- The column is not resistant to alcohol additives in the eluent.
- The eluent has only a limited shelf life due to bacterial and fungal growth. We therefore recommend preparing fresh eluent regularly.

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