Adjusting the sample rack and needle / Probenrack und Nadel justieren

EN

Adjusting the sample rack

It could happen that the rack position openings of the 874 Oven Sample Processor do not precisely match the opening of the oven. The sample vial could become jammed as a result.

This lack of exactitude can be rectified by adjusting the rack. The rack can be adjusted by correcting the rack offset in *tiamo*[™].

Moving to a sample position

Proceed as follows:

- Open the manual control in *tiamo*[™].
- Under **Sample changers ▶ 874_1**, click on **Tower 1**.
- On the tab **Move**, enter the **Target position 1** and click on **[Start]**.

If the opening of rack position 1 is now located considerably to the left or right of the middle of the oven, then the value of the rack offset must be corrected in the rack configuration.

Note


Carry out the correction of the rotation angle of the rack (**Rack offset -5.0° ... $+5.0^{\circ}$**) in *tiamo*.

1 mm deviation of the rack to the left requires a correction of the rotation angle by -0.3° .

1 mm deviation of the rack to the right requires a correction of the rotation angle by $+0.3^{\circ}$.

Correcting the rack offset

- In *tiamo*[™], select the **Configuration**.
- Double-click on the device name **874_1**.
- Click the tab **Rack**.
- Click **[Rackdata]** to open the settings of the rack parameters.

Rack data

Rack name:

Rack code:

Number of positions:

Rack parameters Lift positions Special beakers

Beaker radius samples: mm

Beaker sensor:

Rack offset: °

- Correct the value for **Rack offset** according to the observed deviation.
- Close the rack data configuration with **[OK]**.
- Click on **[Initialize rack]**.

Aligning the needle position

The needle might be in an off-center position in relation to the sample vial. Please correct this setting by aligning the needle holder (1-2).

Note

Please adjust the rack prior to aligning the needle.


Figure 1: Sample vial and needle in alignment position

- 1 Screws
- 2 Needle holder
- 3 Needle
- 4 Sample vial

Moving to sample position

Proceed as follows:

- Place a sample vial (1-4) in the first rack position.
- Move to the first rack position (see description under "Moving to a sample position", page 1).
- Use the lift to position the needle (1-3) approx. 3 mm above the sample vial (see Figure 1, page 2).

Adjusting the holder

Caution

Take care when tightening the screws, as they run in a plastic thread (needle holder).

- Loosen the screws (1-1) until you can move the needle holder.
- Move the needle holder (1-2) until the needle is centered as precisely as possible above the sample vial.
- Fasten the screws again.

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