

## Application Note # CA-270379

# Paraffins, Olefins, iso-Paraffins, iso-Olefins, Naphthenes and Aromatics (PIONA) in Hydrocarbon Streams

### Introduction

This application note describes the quantitative analysis of n- and iso-saturates, n- and iso-olefins, and aromatics in spark-ignition engine fuels by multi-dimensional gas chromatography using the Bruker PIONA+™ Analyzer. Through the use of this system, hydrocarbon types (paraffins, iso-paraffins, olefins, iso-olefins, naphthenes and aromatics) are analyzed and reported, based on carbon number or as a total.

#### Instrumentation:

Bruker PIONA+ Analyzer  
Bruker 450-GC Gas Chromatograph  
PIONA+ multi column module

#### Software:

compassCDS Chromatography Software from Bruker with PIONA+ plug-in software

### Conditions

All conditions (temperatures, valve/switch timings) for the different columns and traps are set and tuned at the Bruker factory per method/mode to obtain an optimized chromatographic separation. In this application, the PIONA mode for the analyzer has been selected. In the PIONA mode of operation, the MoISieve 5A trap is not only used to separate the n-paraffins from the iso-paraffins but also the n-olefins from the iso-olefins as well. An example is shown in Figure 1.

Table 1: Elution scheme for PIONA.

From	To (min)	Components	Column Route
0	30.0	C1 to C12 N + iP	1 <sup>st</sup> OV-275 fraction via 5A and olefin trap to 13x
30.0	35.0	C6 to C8 A and pN	2 <sup>nd</sup> OV-275 fraction via arom/eth to non-polar column
35.0	40.0	>200 °C fraction	Back flush non-polar column of 2 <sup>nd</sup> OV-275 fraction
50.0	80.0	C3 to C12 nP	5A in flow to 13x
80.0	86.5	C8 to C10 A	3 <sup>rd</sup> back flush OV-275 fraction via arom/eth to non-polar column
86.5	95.0	>200 °C fraction	Back flush non-polar column of 3 <sup>rd</sup> OV-275 fraction
100.0	130.0	C3 to C12 cO+iO	Olefin trap in flow via 5A trap to 13x
130.0	150.0	None	Cooling 13x column
150.0	180.0	C3 to C12 nO	5A in flow to 13x




Figure 1: PIONA+ chromatogram of a test sample analysis in PIONA mode.

## Results and discussion

An example of a PIONA analysis is shown in the chromatogram below. A test sample is used to calibrate the PIONA+ analyzer as shown in Figure 1. Clearly visible is the group type separation per carbon number. The high resolution allows for easy identification and thus accurate and precise quantification.

compassCDS software, together with the PIONA+ plug-in, generate weight% and volume% profile reports as shown in Tables 2 and 3. In one view, the amounts of the different groups as well as the totals per group and per carbon number can be seen.

Another example is shown in Figure 2. A second test sample containing additional C5 components and Unsaturation, was analyzed using the PIONA+ system. It is noteworthy that in spite of the fact that additional components are present (compared to the example shown in Figure 1), excellent peak resolution and group separation are still produced. This enables the analyst to achieve precise identification and subsequent quantification, as shown in Table 4 (weight% report) and Table 5 (volume% report). Thus, in one overview, the saturated and unsaturated groups are reported as a total and per carbon number. In addition, the totals per carbon number are revealed.

## Conclusion

The data and results clearly demonstrate that the Bruker PIONA+ Analyzer is capable of providing both mass% and volume% in accordance with the DIN 51448-2 standard method.

Table 2: Weight% report of calibration sample 1.

Carbon	Saturates			Unsaturation			Aromatics	Total
	Cyclic	Iso	Normal	Cyclic	Iso	Normal		
3	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
4	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
5	0.00	0.00	0.00	0.00	0.78	0.00	0.00	0.00
6	1.20	0.87	0.74	1.16	0.85	0.35	2.83	8.00
7	3.56	2.14	2.13	1.80	0.40	0.41	2.90	13.35
8	5.76	4.19	4.04	0.00	0.63	0.00	10.71	25.33
9	4.65	0.70	3.02	0.00	0.20	0.00	6.22	14.59
10	2.45	0.45	5.20	0.00	0.00	0.00	5.30	13.40
11	0.00	0.00	4.32	0.00	0.00	0.00	0.00	4.32
Total	17.62	8.35	19.45	2.97	1.89	0.77	27.95	78.98
Fraction >200 °C	14.22							
Polynaphthenes	6.80							

Table 3: Volume% report of calibration sample 1.

Carbon	Saturates			Unsaturation			Aromatics	Total
	Cyclic	Iso	Normal	Cyclic	Iso	Normal		
3	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
4	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
5	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
6	1.22	1.02	0.87	1.15	0.96	0.40	2.49	8.13
7	3.64	2.42	2.41	1.77	0.45	0.46	2.59	13.73
8	5.74	4.61	4.44	0.00	0.68	0.00	9.56	25.03
9	4.57	0.75	3.23	0.00	0.00	0.00	5.51	14.06
10	2.35	0.47	5.51	0.00	0.00	0.00	4.62	12.95
11	0.00	0.00	4.00	0.00	0.00	0.00	0.00	4.00
Total	17.52	9.26	20.46	2.93	2.09	0.86	24.77	77.90
Fraction >200 °C	14.95							
Polynaphthenes	7.15							




Figure 2: Piona+ chromatogram of test sample 2 in PIONA mode.

Table 4: Weight% report of calibration sample 2.

**Normalized Weight Percent Profile**

Carbon	Saturates			Unsaturates			Aromatics	Total
	Cyclic	Iso	Normal	Cyclic	Iso	Normal		
3	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
4	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
5	1.42	0.00	0.95	0.00	0.78	0.00	0.00	3.15
6	4.05	4.49	3.49	1.47	0.00	1.36	9.15	14.01
7	4.07	1.02	3.88	0.00	0.09	3.81	2.02	14.90
8	5.12	1.08	4.81	0.00	0.18	4.08	4.47	19.74
9	5.05	1.08	4.99	0.00	0.20	3.92	5.19	20.43
10	4.10	1.04	4.05	0.00	0.00	2.68	2.20	14.07
11	0.00	3.70	0.00	0.00	0.00	0.00	0.00	3.70
Total	23.81	12.40	22.17	1.47	1.26	15.85	23.03	100.00
Fraction >200 °C	0.00							
Polynaphthenes	0.00							

Table 5: Volume% report of calibration sample 2.

**Normalized Volume Percent Profile**

Carbon	Saturates			Unsaturates			Aromatics	Total
	Cyclic	Iso	Normal	Cyclic	Iso	Normal		
3	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
4	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
5	1.45	0.00	1.16	0.00	0.90	0.00	0.00	3.51
6	4.07	5.17	4.02	1.43	0.00	1.52	7.89	24.10
7	4.07	1.13	4.30	0.00	0.10	4.13	1.77	15.51
8	5.00	1.17	5.18	0.00	0.19	4.33	3.91	19.78
9	4.86	1.13	5.23	0.00	0.21	4.05	4.51	20.00
10	3.85	1.08	4.21	0.00	0.00	2.73	1.88	13.75
11	0.00	3.36	0.00	0.00	0.00	0.00	0.00	3.36
Total	23.30	13.04	24.10	1.43	1.41	16.76	19.97	100.00
Fraction >200 °C	0.00							
Polynaphthenes	0.00							

**References**

DIN (Deutsches Institut für Normung e. V) 51448-2,  
 "Testing of liquid petroleum hydrocarbons - Determination of  
 hydrocarbon types – Part 1: Gas chromatographic analysis by  
 column switching procedure", Berlin, Germany. [www2.din.de](http://www2.din.de).

**Keywords**

DIN 51448-2  
 Quantitation  
 engine fuels  
 PIONA analysis  
 Hydrocarbon streams

**Instrumentation & Software**

Bruker PIONA+ Analyzer  
 PIONA plug in software  
 compassCDS Chromatography  
 Software

For research use only. Not for use in diagnostic procedures.

**● Bruker Daltonik GmbH**

Bremen · Germany  
 Phone +49 (0)421-2205-0  
 Fax +49 (0)421-2205-103  
[sales@bdal.de](mailto:sales@bdal.de)

**Bruker Daltonics Inc.**

Billerica, MA · USA  
 Phone +1 (978) 663-3660  
 Fax +1 (978) 667-5993  
[ms-sales@bdal.com](mailto:ms-sales@bdal.com)

[www.bruker.com/chemicalanalysis](http://www.bruker.com/chemicalanalysis)