



Semivolatile organic compounds

Analysis of semivolatile compounds in indoor air

Application Note

Environmental

Authors

Agilent Technologies, Inc.

Introduction

Gas chromatography with an Agilent CP-Sil 8 CB column separates 27 semivolatiles in indoor air in 40 minutes.




Agilent Technologies

Conditions

Technique : GC
Column : Agilent CP-Sil 8 CB, 0.32 mm x 50 m fused silica (df = 0.12 µm) (Part no. CP7751)
Temperature : 35 °C → 160 °C, 25 °C/min → 280 °C, 5 °C/min, (26 min), → 300 °C, 5 °C/min, (22 min)
Carrier Gas : Helium, 2.7 mL/min
Injector : TDAS 5000 thermal desorber, using Tenax tubes
Detector : MS
Sample Size : 48 - 300 L of indoor air
Concentration Range : low ppb level
Courtesy : Axel Clausen, National Institute of Occupational Health, Copenhagen, Denmark

Peak identification

1. decane	14. o-terphenyl
2. undecane	15. dibutyl phthalate
3. naphthalene	16. eicosane
4. tridecane	17. pyrene
5. isomer 1 of 2,2,4-trimethyl-1,3-pentanediol monoisobutyrate (texanol)	18. stearic acid
6. isomer 2 of texanol	19. docosane
7. biphenyl	20. tricosane
8. 1-dodecanol	21. tetracosane
9. 2,6-di-tert-butyl-1-hydroxy-toluene (BHT)	22. di-(2-ethylhexyl)-phthalate (DEHP)
10. dodecanoic acid	23. hexacosane
11. di-isobutyrate	24. octacosane
12. heptadecane	25. 1,3,5-triphenylbenzene
13. octadecane	26. triacontane
	27. dotriacontane



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01901



Agilent Technologies