

Hardness Determination in Decarburized Steel

LECO Corporation; Saint Joseph, Michigan USA

Instrument: AMH43 Vickers Microindentation Package with 50X and 10X Objectives



Decarburization

Decarburization can be defined as the loss of carbon from the surface of iron-based alloys as the result of heating in a medium that reacts with the carbon.

While "Case Hardening" refers to a process by which the surface layer of iron-based alloys is substantially harder than the interior, the opposite occurs during decarburization. The surface of decarburized steel is SOFTER than the core. At times it is necessary to quantify the degree of decarburization.

Sampling and Sample Preparation

Sample Identification: Decarburized Steel

Sectioning

| | |
|------------------------|--|
| Saw | MSX300R1 |
| Vises | 808-957, 809-388 |
| Blade | 14" x 0.098", Resin Bonded, Al ₂ O ₃ , 809-346 |
| Speed | 2410 RPM |
| Coolant/Rust Inhibitor | 812-469 |

Mounting

| | |
|----------------|------------------------|
| Press | PR-32 |
| Media | Bakelite (Red) 811-118 |
| Other Comments | "Normal Cycle" |

Manual Preparation

Grinding - GPX200, 10" Wheel

| | Time (Minutes:Seconds) | Hand Pressure (Manual) | Wheel Direction | Wheel Speed (FPM) |
|--------------------|---------------------------|---------------------------|--------------------|----------------------|
| 180 grit SiC | 1:00 | Medium | CCW | 250 |
| 812-206-PRM, water | | | | |
| 320 grit SiC | 1:00 | Medium | CCW | 250 |
| 812-210-PRM, water | | | | |
| 600 grit SiC | 1:00 | Medium | CCW | 250 |
| 812-214-PRM, water | | | | |

Polishing

| | Time (Minutes:Seconds) | Hand Pressure (Manual) | Wheel Direction | Wheel Speed (FPM) |
|--|---------------------------|---------------------------|--------------------|----------------------|
| 3μ Diamond Compound, Ultra Silk, Microid Extender, 810-872, 812-438, 811-003 | 1:30 | Medium | CCW | 200 |
| 1μ Diamond Compound, Red Felt, Microid Extender, 810-870, 812-225, 811-003 | 0:30 | Light | CCW | 200 |

Etching

| | Time (Minutes:Seconds) |
|------------------|---------------------------|
| 2% Nital Reagent | 00:20 to 00:30 |

Automatic Prep

Sample Identification: Decarburized Steel

Grinding (Fixed Holder) - GPX200, 10" Wheel

| | Time (Min:Sec) | Head Dir. | Head Pressure (lb.) | Head Speed (RPM) | Wheel Dir. | Wheel Speed (FPM) |
|-------------------------------|-------------------|--------------|------------------------|---------------------|---------------|----------------------|
| Platinum #1 812-337, Water | 2:00 | CW | 35 | 75 | CCW | 200 |

Pre-Polishing - FAS Magnetic System, 10", 812-382

| | Time (Min:Sec) | Head Dir. | Head Pressure (lb.) | Head Speed (RPM) | Wheel Dir. | Wheel Speed (FPM) |
|---|-------------------|--------------|------------------------|---------------------|---------------|----------------------|
| Silver Disk, 6μ Cameo Suspension, Microid Extender, 812-340, 812-356, 811-003 | 2:00 | CW | 35 | 75 | CCW | 200 |

Polishing

| | Time (Min:Sec) | Head Dir. | Head Pressure (lb.) | Head Speed (RPM) | Wheel Dir. | Wheel Speed (FPM) |
|---|-------------------|--------------|------------------------|---------------------|---------------|----------------------|
| 3μ Premium Suspension, Ultra Silk, Microid Extender, 810-997-016, 812-438, 811-003 | 3:00 | CW | 40 | 100 | CCW | 200 |
| 1μ Premium Suspension, Red Felt, Microid Extender, 810-998-016, 812-225, 811-003 | 0:30 | CW | 40 | 100 | CCW | 200 |

Etching

| | Time (Minutes:Seconds) |
|------------------|---------------------------|
| 2% Nital Reagent | 00:20 to 00:30 |

Calibration Samples


Calibration is performed using a Certified Micro-indentation Hardness Standard from LECO, NIST, or other suitable reference material. The Hardness Standard should be certified with the same load and indenter type as the sample analysis, and have a hardness close to the material being analyzed. A multi-point calibration is available if needed. This would allow multiple hardness points to be calibrated within a given load.

Accessories

863-586-B/O AMH43 Advanced Analysis Software Module

Application Parameters—Individual Patterns

Load: 300 g
 Indenter Type: Vickers (Diamond Pyramid)
 Dwell: 13 seconds
 Patterns: 10 evenly-spaced patterns along surface edge (expected Decarb surface)



Pattern Configuration


5 impressions spaced 100 microns apart; 6th impression @ 1200 microns total depth

Decarb Transition Determination

88% of Core Average (Core Average = Avg. of 5th & 6th impression); This formula is input into the Case Depth Hardness area in each pattern, and is automatically calculated per pattern

Application Parameters—Complete Sample Surface

Load: 300 g
 Indenter Type: Vickers (Diamond Pyramid)
 Dwell: 13 seconds
 Patterns: 5 Segment Patterns based off the sample surface profile



Pattern/Segment Configuration

5 Segment Patterns with 250 micron spacing between impressions; 75 micron spacing between segment lines; 5th segment @ 1200 microns from surface

Decarb Transition Determination

Color Hardness Contour Plot is used to visually display Hardness/Location relationships along the Decarb layer

Analysis Procedure

Follow the macro buttons at the top of the AMH43 Application page from left to right in order to run the analysis.




New—Input sample data

Trace Part—Shows sample outline

Load—Select pre-defined pattern


Verify—Position each pattern



Push the "Analyze" button and the hardness impressions will be made and measured as instructed by the pattern's detail.

| Vector | Indents | Units | | | | | | | |
|--------|----------|---------|--------|------|-------------|---------------|---------------|-------|-----------|
| Indent | Distance | Spacing | Offset | Load | Est. Length | Est. Hardness | Circle Radius | Dwell | Objective |
| 1 | 100.00 | 100.00 | 0.00 | 300 | 43.06 | 300.00 | 0.00 | 13 | 50x |
| 2 | 200.00 | 100.00 | 0.00 | 300 | 43.06 | 300.00 | 0.00 | 13 | 50x |
| 3 | 300.00 | 100.00 | 0.00 | 300 | 43.06 | 300.00 | 0.00 | 13 | 50x |
| 4 | 400.00 | 100.00 | 0.00 | 300 | 43.06 | 300.00 | 0.00 | 13 | 50x |
| 5 | 500.00 | 100.00 | 0.00 | 300 | 43.06 | 300.00 | 0.00 | 13 | 50x |
| 6 | 1200.00 | 700.00 | 0.00 | 300 | 43.06 | 300.00 | 0.00 | 13 | 50x |
| * | | | | | | | | | |


Typical Results—Individual Patterns



Push the "Excel Report" macro button and the test data will automatically be inputted into an Excel template.

| | | | | | | | | | | | | | | |
|-----------|----------|------|-----------|----------|------|-----------|----------|------|-----------|----------|------|-----------|----------|------|
| Pattern 1 | | | Pattern 2 | | | Pattern 3 | | | Pattern 4 | | | Pattern 5 | | |
| Hardness | Distance | HRC | Hardness | Distance | HRC | Hardness | Distance | HRC | Hardness | Distance | HRC | Hardness | Distance | HRC |
| 319 | 100 | 32.1 | 313 | 74.6 | 31.4 | 341 | 99.8 | 34.5 | 300 | 75 | 29.8 | 295 | 101 | 29.1 |
| 352 | 200 | 35.8 | 341 | 174 | 34.5 | 325 | 199 | 32.8 | 319 | 175 | 32.2 | 330 | 189 | 33.3 |
| 352 | 300 | 36.9 | 357 | 274 | 36.3 | 343 | 299 | 34.8 | 341 | 274 | 34.5 | 326 | 305 | 32.8 |
| 359 | 400 | 36.6 | 380 | 373 | 36.6 | 336 | 399 | 34 | 350 | 374 | 35.5 | 352 | 401 | 35.8 |
| 350 | 501 | 35.5 | 357 | 474 | 36.3 | 382 | 499 | 36.9 | 357 | 475 | 36.3 | 359 | 501 | 36.6 |
| 380 | 1200 | 38.8 | 328 | 1200 | 33.1 | 377 | 1200 | 36.5 | 383 | 1200 | 39.1 | 383 | 1200 | 39.1 |
| 321 | 106 | | 301 | 0 | | 325 | 0 | | 325 | 203 | | 327 | 181 | |

| | | | | | | | | | | | | | | |
|-----------|----------|------|-----------|----------|------|-----------|----------|------|-----------|----------|------|------------|----------|------|
| Pattern 6 | | | Pattern 7 | | | Pattern 8 | | | Pattern 9 | | | Pattern 10 | | |
| Hardness | Distance | HRC | Hardness | Distance | HRC | Hardness | Distance | HRC | Hardness | Distance | HRC | Hardness | Distance | HRC |
| 288 | 59.9 | 25.4 | 278 | 95.7 | 26.5 | 286 | 74.2 | 28 | 288 | 74 | 28.2 | 297 | 81.6 | 29.3 |
| 313 | 195 | 31.4 | 293 | 210 | 28.9 | 300 | 190 | 29.8 | 311 | 214 | 31.2 | 323 | 212 | 32.6 |
| 311 | 291 | 31.2 | 313 | 299 | 31.4 | 336 | 275 | 34 | 348 | 300 | 35.3 | 321 | 290 | 32.4 |
| 345 | 375 | 35 | 309 | 400 | 30.9 | 321 | 374 | 32.4 | 345 | 400 | 35 | 338 | 375 | 34.3 |
| 341 | 475 | 34.5 | 343 | 498 | 34.8 | 348 | 474 | 35.3 | 345 | 500 | 35 | 336 | 474 | 34 |
| 380 | 1200 | 40.7 | 389 | 1200 | 37.7 | 377 | 1200 | 36.5 | 380 | 1200 | 36.8 | 391 | 1200 | 39.9 |
| 326 | 327 | | 313 | 411 | | 319 | 234 | | 319 | 232 | | 320 | 193 | |



Typical Results—More Complicated Patterns

Line Number 1

| Hardness | Distance | HRC |
|----------|----------|------|
| 952 | 253 | 35.8 |
| 948 | 504 | 35.3 |
| 948 | 753 | 35.3 |
| 980 | 1000 | 38.8 |
| 941 | 1250 | 34.5 |
| 923 | 1500 | 32.6 |
| 952 | 1750 | 35.8 |
| 928 | 2000 | 33.1 |
| 964 | 2250 | 37.2 |
| 934 | 2500 | 33.8 |
| 936 | 2750 | 34.0 |
| 936 | 3000 | 34.0 |
| 948 | 3250 | 35.3 |
| 945 | 3500 | 35.0 |
| 948 | 3750 | 35.3 |
| 952 | 4000 | 35.8 |
| 932 | 4260 | 33.6 |
| 930 | 4500 | 33.3 |
| 936 | 4750 | 34.0 |
| 932 | 5010 | 33.6 |
| 926 | 5250 | 32.8 |
| 938 | 5500 | 34.3 |
| 928 | 5760 | 33.1 |
| 941 | 6000 | 34.5 |
| 941 | 6250 | 34.5 |
| 934 | 6500 | 33.8 |
| 930 | 6750 | 33.3 |
| 932 | 7000 | 33.5 |
| 941 | 7250 | 34.5 |
| 919 | 7500 | 32.1 |
| 938 | 7750 | 34.3 |
| 930 | 8000 | 33.3 |
| 917 | 8250 | 31.9 |
| 917 | 8500 | 31.9 |
| 921 | 8750 | 32.4 |

Line Number 2

| Hardness | Distance | HRC |
|----------|----------|------|
| 952 | 253 | 35.8 |
| 948 | 504 | 35.3 |
| 948 | 753 | 35.3 |
| 980 | 1000 | 38.8 |
| 941 | 1250 | 34.5 |
| 923 | 1500 | 32.6 |
| 952 | 1750 | 35.8 |
| 928 | 2000 | 33.1 |
| 964 | 2250 | 37.2 |
| 934 | 2500 | 33.8 |
| 936 | 2750 | 34.0 |
| 936 | 3000 | 34.0 |
| 948 | 3250 | 35.3 |
| 945 | 3500 | 35.0 |
| 948 | 3750 | 35.3 |
| 952 | 4000 | 35.8 |
| 932 | 4260 | 33.6 |
| 930 | 4500 | 33.3 |
| 936 | 4750 | 34.0 |
| 932 | 5010 | 33.6 |
| 926 | 5250 | 32.8 |
| 938 | 5500 | 34.3 |
| 928 | 5760 | 33.1 |
| 941 | 6000 | 34.5 |
| 941 | 6250 | 34.5 |
| 934 | 6500 | 33.8 |
| 930 | 6750 | 33.3 |
| 932 | 7000 | 33.5 |
| 941 | 7250 | 34.5 |
| 919 | 7500 | 32.1 |
| 938 | 7750 | 34.3 |
| 930 | 8000 | 33.3 |
| 917 | 8250 | 31.9 |
| 917 | 8500 | 31.9 |
| 921 | 8750 | 32.4 |

Line Number 3

| Hardness | Distance | HRC |
|----------|----------|------|
| 930 | 251 | 33.3 |
| 948 | 501 | 35.3 |
| 941 | 752 | 34.5 |
| 955 | 1000 | 36.1 |
| 945 | 1250 | 35.0 |
| 945 | 1500 | 35.0 |
| 941 | 1750 | 34.5 |
| 957 | 2000 | 36.3 |
| 943 | 2250 | 34.8 |
| 932 | 2500 | 33.5 |
| 943 | 2750 | 34.8 |
| 948 | 3000 | 35.3 |
| 969 | 3250 | 37.7 |
| 936 | 3500 | 34.0 |
| 934 | 3750 | 33.8 |
| 962 | 4000 | 36.9 |
| 941 | 4250 | 34.5 |
| 943 | 4500 | 34.8 |
| 943 | 4750 | 34.8 |
| 945 | 5000 | 35.0 |
| 947 | 5250 | 35.3 |
| 947 | 5500 | 35.3 |
| 948 | 5750 | 35.3 |
| 964 | 6000 | 37.2 |
| 945 | 6250 | 35.0 |
| 950 | 6500 | 35.5 |
| 945 | 6750 | 35.0 |
| 943 | 7000 | 34.8 |
| 945 | 7250 | 35.0 |
| 936 | 7500 | 34.0 |
| 943 | 7750 | 34.8 |
| 936 | 8000 | 34.0 |
| 928 | 8250 | 33.1 |
| 932 | 8500 | 33.6 |
| 930 | 8750 | 33.3 |

Individual hardness data from complicated patterns (involving hundreds of impressions) does not lend itself to easy viewing and interpretation (see example at left).

When this is the case, it is useful to utilize the Color Hardness Contour mapping feature of the AMH system (see image below).

