

Thermo. Titr. Application Note No. H-009

Title: Determination of Sulfate in Brines

Scope: Determination of the sulfate content of brines

Principle: An aliquot of brine is acidified with nitric acid and titrated with standard barium chloride solution to a single thermometric endpoint. In concentrated brines, the endpoint is subject to some rounding. For accurate results, it is necessary to determine the method blank on aliquots of a sample brine solution.

Reagents: Standard 1mol/L BaCl₂ solution
5 mol/L HNO₃
Anhydrous Na₂SO₄ A.R.

Method: Basic Experimental Parameters:

Data rate (per second)	10
Titrant delivery rate (mL/min.)	2
No. of endothermic endpoints	1
Data smoothing factor	55

Procedure:

Pipette a 25.00 mL or 50mL aliquot of brine into a titration vessel. Add 1 mL 5 mol/L HNO₃, and titrate with 1mol/L BaCl₂ solution to an exothermic endpoint.

Determination of method blank:

Titrate aliquots of 20, 25, 30, 40 and 50mL of a selected typical brine sample according to Section 4.2. Subject the results to regression analysis, plotting aliquot volume on the x-axis and BaCl₂ titre on the y-axis. The y-intercept is the method blank in mL, and must be subtracted from all titres. It will be noted in the example given here; the intercept is negative, meaning that this amount must be effectively added to the titre.

Standardization of BaCl₂ titrant:

Dry anhydrous A.R. Na₂SO₄ for 2 hours at 200°C. Cool in a dessicator. Weigh accurately 5 amounts ranging from approximately 0.13g to 0.65g in roughly equal increments directly into titration vessels. Add 30mL D.I. water and 1 mL 5 mol/L HNO₃ and titrate. Convert masses of Na₂SO₄ titrated to mmole, and plot on the x-axis, with

corresponding titres of BaCl₂ on the y-axis. Perform a regression analysis, and compute the gradient of the regression curve. The molarity of the BaCl₂ is the reciprocal of the gradient. In this instance, the y-intercept is not used as the method blank, due to the need to *match the sample matrix*.

Results (Example):	Analysis of brines:		
	Sample No.	Sample Aliquot, mL	Sulfate as SO ₄ ²⁻ , g/L
	1	50	4.73, 4.75 φ = 4.74
	2	20 - 50	12.16, 12.13, 12.17, 12.14, 12.08 φ = 12.14, S.D. = 0.033
	3	50	2.80, 2.77 φ = 2.79
	4	25	7.70, 7.70 φ = 7.70
	5	25 & 50	3.10, 3.14, 3.11 φ = 3.12

Determination of Method Blank (Example based on sample #2)		
(see Fig. 1)	Sample aliquot volume, mL	Titre BaCl ₂ , mL
	20	2.533
	25	3.155
	30	3.792
	40	5.032
	50	6.330

Standardization of BaCl ₂ Titrant			
(see Fig. 2)	Mass of Na ₂ SO ₄ , g	*Equiv. mmole Na ₂ SO ₄	BaCl ₂ titre, mL
	0.1320	0.925	1.012
	0.2690	1.884	1.971
	0.5444	3.813	3.909
	0.4119	2.885	2.986
	0.6662	4.667	4.751

* Assumes reagent purity of 99.5%




Fig.1. Regression analysis to determine method blank
y-intercept = method blank = -0.0017mL




Fig. 2. Regression analysis to standardize BaCl₂
 Molarity = 1/gradient = 1/1.0004 = 0.9996 mol/L


Calculation:

$$\text{SO}_4^{2-} \text{ g/L} = \frac{((\text{titre, mL} - \text{blank, mL}) \times M \text{ BaCl}_2 \times \text{FW SO}_4^{2-})}{\text{sample volume, mL}}$$

Example:

$$\begin{aligned} \text{SO}_4^{2-} \text{ g/L} &= \frac{((2.004 - (-0.0017)) \times 0.9996 \times 96.058)}{25.00} \\ &= 7.70 \end{aligned}$$

Thermometric Titration Plot:



Legend:

Red = solution temperature curve

Black = second derivative curve