

Thermo. Titr. Application Note No. H-017

Title:	Determination of Bromide and Chloride in Photographic Developer Solutions
---------------	--

Scope:	Determination of bromide and chloride in photographic developer solutions
---------------	---

Principle:	An aliquot of photographic developer solution is titrated with 0.1M AgNO ₃ . Endpoints are determined using thermometric and potentiometric sensors simultaneously. A combination silver billet/reference electrode was used as the potentiometric sensor
-------------------	--

Reagents:	0.1M AgNO ₃ solution Concentrated HNO ₃
------------------	--

Method:	<p>Basic Experimental Parameters:</p> <table> <tr> <td>Data rate (per second)</td> <td>10</td> </tr> <tr> <td>Titrant delivery rate (mL/min.)</td> <td>1</td> </tr> <tr> <td>No. of exothermic endpoints</td> <td>2</td> </tr> <tr> <td>Data smoothing factor (thermo)</td> <td>55</td> </tr> <tr> <td>No. of potentiometric endpoints</td> <td>2</td> </tr> <tr> <td>Data smoothing factor (potentiometric)</td> <td>110</td> </tr> </table> <p>Procedure: (Fresh Developer or Replenisher). Pipette a 10mL aliquot of undiluted solution into a titration beaker equipped with a spin ring. Add 15 mL deionized water and neutralize with 2mL concentrated nitric acid. Mix, and allow temperature to equilibrate in a room temperature water bath. Titrate to the first endpoint for bromide assay only, or to the second endpoint for combined bromide and chloride determination.</p>	Data rate (per second)	10	Titrant delivery rate (mL/min.)	1	No. of exothermic endpoints	2	Data smoothing factor (thermo)	55	No. of potentiometric endpoints	2	Data smoothing factor (potentiometric)	110
Data rate (per second)	10												
Titrant delivery rate (mL/min.)	1												
No. of exothermic endpoints	2												
Data smoothing factor (thermo)	55												
No. of potentiometric endpoints	2												
Data smoothing factor (potentiometric)	110												

Results:		Bromide g/L		Chloride g/L	
		Ag billet	Thermo	Ag billet	Thermo
	Mean (n=8)	1.62	1.61	1.39	1.40
	Std. Dev.	0.005	0.007	0.004	0.003

Calculation:

$$\text{halide } g / L = \frac{((\text{titre} - \text{blank}) \times M \text{ AgNO}_3 \times \text{FW halide})}{\text{aliquot, mL}}$$

Thermometric/Potentiometric Titration Plot:



Legend:

Potentiometric:

Red = mV curve

Green = first derivative curve

Thermometric:

Plum = solution temperature curve

Brown = second derivative curve