


Fully automated water analysis by OMNIS



The determination of the physical and chemical parameters as electrical conductivity, pH value, alkalinity, the calcium and magnesium hardness as well as the total hardness are necessary for evaluating the water quality.

A fast and accurate determination in tap water is realized using an automated OMNIS System working in parallel on different workstations. An 856 Conductivity Module with Dosinos extends the system.

Method description

Sample

Tap water Herisau

Five-ring conductivity measuring cell with Pt1000, $c = 1.0 \text{ cm}^{-1}$

6.0915.130

Sample preparation

No sample preparation is required.

Configuration

| | |
|--|-------------|
| Main module Pick&Place L | 2.1012.0010 |
| Module base M/L | 6.02003.010 |
| Rack base, 2x | 6.02001.010 |
| Pick&Place module, 4x | 2.1014.0010 |
| "Peristaltic" (4-channel) pump module, 2x | 2.1016.0110 |
| Gripper fingers 50 - 72 mm | 6.02601.020 |
| OMNIS sample rack 9 x 200 mL, 7x | 6.02041.020 |
| Retainer for 200 mL sample beaker | 6.02064.020 |
| Beaker adapter for 200 mL sample beaker, 4x | 6.01404.020 |
| Titration head 6xNS14 / 3xNS9 (P&P), 4x | 6.01403.000 |
| Stirring propeller 20 mm ETFE, 3x | 6.01900.030 |
| OMNIS Professional Titrator without stirrer, 2x | 2.1001.0010 |
| OMNIS Dosing Module without stirrer, 2x | 2.1003.0010 |
| OMNIS Titration Module without stirrer, 2x | 2.1002.0010 |
| Digital measuring module, 6x | 6.02100.010 |
| Electrode cable plug-in head Q / plug P, 1.5 m, 6x | 6.02104.310 |
| OMNIS 10 mL cylinder unit, 3x | 6.03001.210 |
| OMNIS 5 mL cylinder unit, 3x | 6.03001.150 |
| 807 Dosing unit 20 mL, 3x | 6.3032.220 |
| OMNIS Rod stirrer "Sample Robot", 3x | 2.1006.0010 |
| OMNIS | 6.0600.0000 |
| OMNIS Stand-alone license (including one instrument license) | 6.06003.010 |
| OMNIS instrument license, 2x | 6.06002.010 |
| 800 Dosino, 3x | 2.800.0010 |
| 856 Conductivity module | 2.856.0010 |
| Combined dCa ISE | 6.00502.300 |
| dAquatrode with Pt1000 | 6.00202.300 |

Solutions

| | |
|--------------------|--|
| Titration 1 | $c(\text{HCl}) = 0.1 \text{ mol/L}$, if possible this solution should be bought from a supplier. |
| Titration 2 | $c(\text{EDTA}) = 0.1 \text{ mol/L}$, if possible this solution should be bought from a supplier. |
| Auxiliary solution | 24.2 g of TRIS is weighed in and transferred into a 1 L volumetric flask and dissolved in approx. 500 mL dist. H_2O . 12 mL acetyl acetone is added and the solution is made up to the mark with dist. H_2O . This solution can only be used for a few days. It masks Fe^{3+} and Al^{3+} for a better differentiation of Ca^{2+} and Mg^{2+} . |

Analysis

100 mL sample is pipetted into a 200 mL disposable beaker. First, the conductivity cell is dipped 3 times for preconditioning into the sample. Then, the conductivity is measured. After pH measurement, the alkalinity was determined by a DET titration using $c(\text{HCl}) = 0.1 \text{ mol/L}$. Directly after the alkalinity determination, 15 mL of the auxiliary solution is added. The hardness of the sample solution is determined by titration with $c(\text{EDTA}) = 0.1 \text{ mol/L}$ until after the second equivalence point.

Parameters

Conductivity

| Mode | MEAS Cond |
|--------------------|------------------------------------|
| Drift | Time controlled (no drift control) |
| Measuring time | 60 s |
| Measuring interval | 2 s |
| Volume increment | 50 μL |

pH measurement

| Mode | MEAS pH |
|-------------------|----------|
| Drift | 2 mV/min |
| Min. waiting time | 10 s |
| Max. waiting time | 110 s |

Method description

Alkalinity

| | |
|---------------------|-----------|
| Mode | DET pH |
| Drift | 50 mV/min |
| Max. waiting time | 26 s |
| Meas. Point density | 4 |
| Min. increment | 10 µL |
| Stop volume | 7 mL |
| EP criterion | 5 |
| EP recognition | Greatest |

Hardness

| | |
|---------------------|-----------|
| Mode | DET U |
| Drift | 50 mV/min |
| Max. waiting time | 26 s |
| Meas. Point density | 4 |
| Min. increment | 10 µL |
| Max. increment | 100 µL |
| Stop volume | 4 mL |
| EP criterion | 5 |
| EP recognition | All |

Results

| Sample (n = 63) | Mean value | s(abs) | s(rel) / % |
|---------------------------|------------|--------|------------|
| Conductivity / (mS/cm) | 0.546 | 0.001 | 0.27 |
| pH | 8.079 | 0.046 | 0.58 |
| Alkalinity / (mmol/L) | 5.674 | 0.015 | 0.27 |
| Ca / (mmol/L) | 2.133 | 0.015 | 0.71 |
| Mg / (mmol/L) | 0.814 | 0.014 | 1.77 |
| Total hardness / (mmol/L) | 2.947 | 0.013 | 0.44 |

| | |
|-----------------|----------|
| Sample (n = 63) | Time |
| Total duration | 04:31:46 |