

Agilent CrossLab Start Up Services

Agilent Cary WinUV 5.3 Software Site Preparation Checklist

Thank you for purchasing software from **Agilent Technologies**. CrossLab Start Up is focused on helping customers shorten the time it takes to start realizing the full value of their instrument system investment.

Correct site preparation is the key first step in ensuring that your instruments and software systems operate reliably over an extended lifetime. This document is an **information guide and checklist** prepared for you that outlines the software requirements for the system set up in your lab.

Introduction

Customer Information

- If you have questions or problems in providing anything described as part of Customer Responsibilities below, please contact your local Agilent or partner support / service organization for assistance prior to delivery. In addition, Agilent and/or its partners reserve the right to reschedule the installation dependent upon the readiness of your laboratory.
- Should your site not be ready for whatever reasons, please contact Agilent as soon as possible to re-schedule any services that have been purchased.
- Other optional services such as additional training, operational qualification (OQ) and consultation for user-specific applications may also be provided at the time of installation when ordered with the system but should be contracted separately.

Customer Responsibilities

Ensure that your site meets the following specifications before the installation date. For details, see specific sections within this checklist, including

- The **computing environment** and the necessary space is made available.
- The number and location of **electrical outlets** for your computer systems and peripherals are planned.
- Your site meets the **software, hardware, and networking specifications** described later in this document.
- Locate your **sales order information**, software authorization codes, and/or software licenses/certificates.
- Check for Hotfix, Microsoft updates, and patch compatibility appropriate for the customer's system. **Critical Note: Current Microsoft updates and Windows version must be installed in the customer provided PC before installation can begin.**
- The necessary **software media**, disks etc. are available including upgrade/update disks.
- That a suitable **backup solution** is identified for your software and **data**.
- If Agilent is delivering **Installation and Introduction services**, users of the instrument system should be present throughout these services. Otherwise, they will miss important operational, maintenance, and safety information.
- Availability of a **system/network administrator** as needed to connect to your intranet.

- Please consult the **Special Requirements** section for other product-specific information.

Important Customer Web Links

- To access **Agilent University**, visit <http://www.agilent.com/crosslab/university/> to learn about training options, which include online, classroom and onsite delivery. A training specialist can work directly with you to help determine your best options.
- To access the **Agilent Resource Center** web page, visit <https://www.agilent.com/en-us/agilentresources>. The following information topics are available:
 - Sample Prep and Containment
 - Chemical Standards
 - Analysis
 - Service and Support
 - Application Workflows
- The **Agilent Community** is an excellent place to get answers, collaborate with others about applications and Agilent products, and find in-depth documents and videos relevant to Agilent technologies. Visit <https://community.agilent.com/welcome>
- Videos about specific preparation requirements for your instrument can be found by searching the **Agilent YouTube** channel at <https://www.youtube.com/user/agilent>
- **Need to place a service call?** [Flexible Repair Options | Agilent](#)

Important Version and Hardware Information

- Please refer to the following table for supported Cary WinUV Pharma software, firmware and hardware configurations

Cary WinUV Version (Firmware)	Supported Instrument Hardware
5.3.0 (FW 6.4.0.141)	Agilent Cary 60 UV-Vis (Cary 50 is not tested or supported)

Site Preparation

Software Specifications for Workstations

Special notes

- If you have purchased a system including hardware instrumentation, refer to the hardware and peripheral device Site Preparation Checklist for requirements regarding laboratory bench space requirements.
- **Implement a backup strategy or backup all your data prior to making any changes.**
- No other instrument software can be installed.
- **Upgrade information:**
 - **Direct upgrades from versions earlier than Cary WinUV 5.1.3 are not supported.**
 - **To upgrade from Cary WinUV 5.1.3 and earlier, the old version must be uninstalled, and the new version installed.**
 - **Data from versions Cary WinUV 5.1.3 and earlier are compatible with Cary WinUV 5.3.**
- Upgrade services (if available) only include:
 - upgrade on the original PC.
 - Installing to a new PC.
- Upgrade services **do not** include:
 - Backing up or restoring data.
 - Transferring data.
 - Installing or upgrading operating systems.
 - Verifying data
 - IQ/OQ
 - Please contact your local Agilent representative for more information.

Workstation Requirements

Software Specification Description	Supported	Comments
Operating system name, version	Microsoft Windows 10 revision 21H2 or higher Microsoft Windows 11 Revision 21H2 or higher	64-bit only Windows 10: Evaluated on Enterprise (LTSC) and Pro. Windows 11: Evaluated on Enterprise and Pro.
O/S .NET and other add-ons	NET 4.8 Advanced Services .NET 3.5.1	Installed by the System Preparation Tool if required.
Language settings/compatibility	English Japanese Chinese Spanish Portuguese	English, Japanese, Simplified Chinese, Spanish and Portuguese Microsoft Windows operating systems are supported.
Account settings/privileges	Full Local Administrator privileges	Required for installation, configuration and upgrading
Microsoft Windows Updates	Outstanding updates must be completed or must not start before installing software	Disconnect the PC from the network while performing the installation and configuration
Web Browser	Google Chrome 98 or higher Microsoft Edge	Cary UV Help now supports HTML5
PDF Reader	Adobe Acrobat Reader DC	Required to view manuals and documentation
Antivirus software	Symantec Endpoint Protection	
Power Saving	None - disable sleep or power saving modes. If enabled, this may result in data loss or unstable instrument control.	Any PC power saving functions for hardware such as hard disks, monitor, LAN etc. should be disabled.

Computer Hardware Specifications for Workstations

Special notes

- Cary UV systems are evaluated under Agilent's ISO9001 procedures using Cary WinUV software and a desktop computer. We strongly recommend a computer configuration that matches the system used during evaluations.

Hardware Specification Description	Evaluated Configuration	Evaluated Operating System	Comments
Desktop Computer	PC: HP z2G9	Windows 10 Pro 21H2 Windows 10 Ent (LTSC) Windows 11 Pro 21H2 Windows 11 Pro 22H2 Windows 11 Ent 22H2	
Laptop Computer	Elitebook 655 G10	Windows 10 Pro 22H2 Windows 11 Pro 22H2	
Processor type and speed	Intel® i5 or equivalent 3.0 GHz or greater 6 Core		
Memory	8Gb		
Internal storage/devices/media	500Gb solid-state-drive		
Graphics resolution	1600 x 900		
Pointing devices	Computer Mouse		
USB Port	USB 2.0		Required for installation
LAN	100/1000 MBIT		Required for network connectivity

Special Requirements

- For further details on supported hardware, software, instruments, and firmware requirements please refer to the **Agilent Cary WinUV 5.3.0 Release Notes and SSB** documents which can be found on the installation media or by contacting your local Agilent representative.
- For further details on installation and upgrades, please consult the **Cary WinUV 5.3.0 Installation Instructions** which can be found on the installation media or by contacting your local Agilent representative.

Service Engineer Review (Optional)

If the Service Engineer completed a review of the Site Preparation requirements with the customer, the Service Engineer should complete the following Comments section.

Service Engineer Comments

If there are any specific points that should be noted as part of performing the site preparation review or other items of interest for the customer, please write in this box.

Site Preparation Verification

Service Request Number:

Date of Review:

Service Engineer Name:

Customer Name:

Service Engineer Signature:

Total number of pages in this document: