

Reaction Cell Frontier: Selenium Isotope Dilution Analysis using Innovative Reaction Cell ICP-MS

Naoki Sugiyama Agilent Technologies, Inc., 9-1 Takakura-cho, Hachioji-shi, Tokyo 192-0033 JAPAN.



FP25

Introduction

Selenium Isotope Dilution Analysis using Innovative Reaction cell ICP-MS

Selenium (Se) is one of major elements required to be quantified in various types of samples due to its essentiality in animal/human and toxicity when ingested in excess. However it is a difficult element to accurately quantify. There would be three reasons. The first is that required Limit of Quantification (LOQ) for Se is commonly as low as ppb to sub ppb range. The second is Se has high ionization potential (IP=9.75eV), causing low sensitivity due to poor ionization in plasma, susceptibility to matrix suppression and difficulty in finding suitable Internal Standard (ISTD) element. The third is that all Se isotopes suffer severe interference as summarized in table 1. Isotope Dilution (ID) method is well known technique for its absolute accuracy, however it is also a method that is least susceptible to signal drift and matrix suppression, which are common problems in ICP-MS, deteriorating accuracy of quantification. However application of ID method to Se quantification has been limited because of the third reason stated above; interference.

Table 1: Interference on major Se isotopes

m/z	Abundance %	isobar	argide	oxide	hydride	chloride	doubly charged	dimer
77	7.83	⁷⁷ Co ⁺ , ⁷⁷ Ar ⁺ , ⁷⁷ Ca ⁺ , ⁷⁷ Al ⁺	⁷⁷ Se ⁺ , ⁷⁷ Br ⁺ , ⁷⁷ K ⁺	⁷⁷ SeO ⁺ , ⁷⁷ SeO _{2⁺}	⁷⁷ SeH ⁺ , ⁷⁷ SeH ₂ ⁺	⁷⁷ SeCl ⁺ , ⁷⁷ SeCl ₂ ⁺	⁷⁷ Se ²⁺ , ⁷⁷ SeO ²⁺	⁷⁷ Se ₂ ⁺
78	23.77	⁷⁸ K ⁺	⁷⁸ Ar ⁺ , ⁷⁸ Ca ⁺ , ⁷⁸ Al ⁺	⁷⁸ Se ⁺ , ⁷⁸ Br ⁺ , ⁷⁸ K ⁺	⁷⁸ SeH ⁺ , ⁷⁸ SeH ₂ ⁺	⁷⁸ SeCl ⁺ , ⁷⁸ SeCl ₂ ⁺	⁷⁸ Se ²⁺ , ⁷⁸ SeO ²⁺	⁷⁸ Se ₂ ⁺
80	49.61	⁸⁰ K ⁺	⁸⁰ Ar ⁺ , ⁸⁰ Ca ⁺ , ⁸⁰ Al ⁺	⁸⁰ Se ⁺ , ⁸⁰ Br ⁺ , ⁸⁰ K ⁺	⁸⁰ SeH ⁺ , ⁸⁰ SeH ₂ ⁺	⁸⁰ SeCl ⁺ , ⁸⁰ SeCl ₂ ⁺	⁸⁰ Se ²⁺ , ⁸⁰ SeO ²⁺	⁸⁰ Se ₂ ⁺
82	8.73	⁸² K ⁺	⁸² Ar ⁺ , ⁸² Ca ⁺ , ⁸² Al ⁺	⁸² Se ⁺ , ⁸² Br ⁺ , ⁸² K ⁺	⁸² SeH ⁺ , ⁸² SeH ₂ ⁺	⁸² SeCl ⁺ , ⁸² SeCl ₂ ⁺	⁸² Se ²⁺ , ⁸² SeO ²⁺	⁸² Se ₂ ⁺

The innovative reaction cell of ICP-QQQ allows us to measure more than one Se isotopes free from interference, making it possible to quantify Se using ID method. Recently Giuseppe et al.¹⁾ applied ID method for multi element analysis focusing on less susceptibility to matrix suppression rather than the absolute accuracy. They proposed Online Isotope Dilution Analysis (OIDA) that removes time consuming step for spiking enriched-isotope standard into samples.

Experimental




Figure 1: Configuration of Agilent ICP-QQQ

Agilent ICP-QQQ

Agilent Technologies developed a new high-end ICP-MS, Triple Quadrupole ICP-MS (ICP-QQQ). As shown in Fig.1, it has two Quadrupoles before and after Collision/Reaction Cell. The 1st Quadrupole selects ions entering the cell, providing consistent reaction condition to changing sample composition. It solves problem of current reaction cell technologies, allowing analyst to use reaction cell for more elements/applications, more effectively. The Agilent ICP-QQQ has following features.

- The best performance cell for solving interference problem.
 - > Advanced reaction cell overcoming weak point of existing reaction cell due to unique QQQ configuration.
 - > Versatile He collision cell. The performance is better than Agilent 7700 due to MS/MS scan.
- Sensitivity is much higher than Agilent 7700 series ICP-MS. A lab-prototype achieved 810 Mcps/ppm with oxide ratio (CaO/Ce⁺) of 1.47% by x-axis and 1.2 Gcps/ppm with oxide ratio of 2.47% by s-axis for Yttrium.
- Random back ground noise is as low as <0.2cps.
- Abundance sensitivity is immeasurably low; <10-9
- It keeps Agilent 7700 compatible robustness and durability, being applicable to developed Agilent 7700 applications.

Results and Discussion

1/ Cell gas mode investigation

Interference removal capability was studied to determine optimum cell gas condition and the best isotopes of Se for ID method. An experimental result is given in Fig.2. The Figure shows Background Equivalent Concentration (BEC) of four reaction gas modes in various synthetic matrix samples. All modes used MSMS scan that operated the 1st Quad as 1 amu-window band pass mass filter, selecting ions entering reaction cell. While Se isotopes were detected at the nature mass in the three cell gas modes, He, H₂ and NH₃+H₂ modes, they were detected as SeO⁺ in O₂+H₂ mode. In this mode Se⁺ reacts with O₂, converting to SeO⁺. We call it mass shift method. Fig.3 illustrates the mass shift method. Reaction between Se⁺ and O₂ is endothermic though (ΔH=0.769eV), the reaction was promoted at octopole bias of -18V effectively forming SeO⁺. Addition of small amount of H₂ into cell was found effective for further reducing BEC.




Fig.2 Result of interference removal capability study

From Fig.2 we can see O₂+H₂ mass shift method gives the best interference removal capability, allowing us to measure three Se isotopes, ⁷⁷Se, ⁷⁸Se and ⁸²Se nearly free from interference. O₂ cell gas can be used with existing cell-ICP-MS though, the method doesn't work fine when sample contains either Zr, Mo or Ru since their isotopes overlaps SeO⁺. With ICP-QQQ, due to ion selection by Q1, the method can be applied to samples regardless of the composition. It is clear advantage of ICP-QQQ over existing cell-ICP-MS, allowing us to apply the developed method to various types of samples.




Fig.3 Illustration of O₂ mass shift method of ICP-QQQ

2/ Selenium Online Isotope Dilution Analysis (OIDA)

⁸²Se enriched standard (⁸²Se 97.43%, ⁸⁰Se 1.65% and ⁷⁸Se 0.51%) was purchased from Oak Ridge national Laboratory. The standard was dissolved in ultrapure HNO₃, and diluted to be appropriate concentration. It was spiked to all tested samples online using Agilent online-ISTD mixing kit (IPA was added to the spike solution to enhance Se⁺ signal). Please refer to reference [1] for detail of OIDA. Se concentration of sample is calculated by formula [1] using measured isotope ratio of two Se isotopes, R_m. I determined Se concentration in samples using two isotope ratios, ratio of ⁸²Se to ⁸⁰Se and ratio of ⁸²Se to ⁷⁸Se for comparison.

$$C_x = C_y \left(\frac{R_m - R_n}{R_n - R_s} \right) \left(\frac{R_s - R_n}{R_x - R_n} \right) \quad [1]$$

R_m : measured isotope ratio of mixed reference standard and spike solution
 R_n : measured isotope ratio of mixed unknown sample and spike solution
 C_x : concentration of unknown sample
 C_y : concentration of natural reference standard
 R_s : natural isotope ratio
 R_x : isotope ratio of spike

3/ Quantification Error study in ID method

For accurate quantification, ID method normally needs correction of measured isotope ratio such as mass bias correction and dead time correction. However we can expect good accuracy with OIDA method without corrections when C_x = C_r is met. Theoretically Error = 0 is achieved when C_x = C_r without correction. I investigated the accuracy of measurement without correction. I analyzed 0.1, 0.2, 2 and 20 ppb natural Se solution using 2ppb solution as natural reference standard. From Fig.4 it is found that quantification error becomes significant when C_x << C_r and C_x >> C_r and accuracy of 95-105% will be achievable to concentration ranging from 100ppt to 5ppb without correction. Since accuracy of measured isotope ratio of ⁸²Se/⁷⁸Se suffer severe mass bias effect, it gave larger uncertainty in determined concentration.




Fig.4 Quantification error of ID method

4/ Se analysis by OIDA using O₂+H₂ mass shift method of ICP-QQQ

Se concentration was determined in various Certified Reference Material (CRM) using developed method. Analyzed CRM contain environmental water, rock, sediment, food, bio-sample and plant. Those CRM are all prepared in appropriate way, i.e. some of them were microwave digested with HNO₃, HCl and H₂O₂. The tested samples were all diluted so that expected Se concentration falls between 100ppt and 5ppb. Integration time of 1 sec. was applied to three Se isotopes measurement and ten replicates were done for averaging. Se concentration is calculated by formula [1] using measured isotope ratio, R_m.

Fig.5 shows Se quantification result in recovery % to certified value. The developed method determined Se concentration of all CRM with good accuracy except two milk samples. At this moment, I have no clear idea to explain the result. It may be simply because of human error during sample prep, or due to unknown interference or signal enhancement. One of plan I have for next step is to utilize ID method spiking isotope-enriched standard before sample preparation.




Fig.5 Result of Se analysis in CRM by OIDA using O₂+H₂ mass shift method of ICP-QQQ

Conclusions

- 1/ We found that reaction cell using O₂+H₂ can effectively solve the interference by REE⁺ on Se isotopes. The developed method on ICP-QQQ can be applied to samples containing Zr, Mo and Ru, since they are removed by the first Quad.
- 2/ We demonstrated the OIDA using O₂+H₂ mass shift method of ICP-QQQ can be used for Se quantification in wide range of samples.
- 3/ Next our plan is to investigate the high recovery of Se in milk samples.

Reference

- 1) Giuseppe Centineo, Jose Angel, Rodriguez Castrillon and Esther Munoz Agudo, Agilent Technologies application note 5990-917EN
- 2) A. Henriot, Fresenius J Anal Chem, 1994, 350, 657-658