



Mestrelab Research

# Gears Chrom Calibration 1.0

---

## STARTING GUIDE



Document Number

P/N 345 R1.1

## **COPYRIGHT**

**©2023 MESTRELAB RESEARCH S.L.**

All rights reserved. No parts of this work may be reproduced in any form or by any means - graphic, electronic, or mechanical, including photocopying, recording, taping, or information storage and retrieval systems - without the written permission of the publisher.

Products that are referred to in this document may be either trademarks and/or registered trademarks of the respective owners. The publisher and the author make no claim to these trademarks.

While every precaution has been taken in the preparation of this document, the publisher and the author assume no responsibility for errors or omissions, or for damages resulting from the use of information contained in this document or from the use of programs and source code that may accompany it. In no event shall the publisher and the author be liable for any loss of profit or any other commercial damage caused or alleged to have been caused directly or indirectly by this document.

## CONTENTS


<b>1. THE WORKFLOW .....</b>	<b>3</b>
1.1. INPUT .....	3
1.2. PROCESSING.....	8
1.3. PLUGINS.....	8
1.3.1. <i>The Chromatogram tab</i> .....	9
1.3.2. <i>The CSV Input tab</i> .....	10
1.3.3. <i>The Controls tab</i> .....	13
1.3.4. <i>The Report tab</i> .....	14
1.4. DESIGN .....	14
1.5. OUTPUT .....	14
<b>2. THE OUTPUT FOLDER.....</b>	<b>15</b>
2.1. THE HTML FILE .....	15
2.2. THE CSV REPORTS.....	16
2.3. THE MNOVA AND PDF REPORTS .....	18
2.4. OTHER OUTPUT.....	19
<b>3. MNOVA GEARS RESULTS VIEWER .....</b>	<b>19</b>
3.1. THE MGEARS VIEWER.....	20
3.1.1. <i>The Well plate view</i> .....	20
3.1.2. <i>The Results section</i> .....	21
3.1.3. <i>Assigning peak/Modifying peak assignment</i> .....	21
3.2. THE CHROM CAL VIEWER.....	23
3.2.1. <i>The Calibration and Cross-validation points</i> .....	23
3.2.2. <i>The Calibration parameters</i> .....	24
3.2.3. <i>The Plot Area</i> .....	24
3.2.4. <i>Updating the calibration/cross-validation results</i> .....	25

Generating a calibration curve is a standard and routine process when performing quantitative chromatography. With Gears Chrom Cal, this process has been fully automated to allow you work on different analytes at the same time and calculate calibration curves quickly and efficiently in a single, simple run.

In this document, we will go through the different configuration steps right to the output and result revision to help you get started with Chrom Cal.

## 1. The workflow


Launch Mgears from the Mnova **Automation** ribbon. The dialog with the usual six tabs will open.



### 1.1. Input


Chrom Cal analysis follows the general Mnova Gears workflow based on its five-step setup, the first step being the choice of input data. Mgears can read data from your **Disk**, **Database**, or from **Real-Time** acquisition. In this guide, we will work with data from disk directories (*please refer to the [Mnova Gears manual](#) for further detail about other input types*).

Click on the  button and select the data files to upload as your **Main Directory**.




Use the **Automatic Detection** mode to detect you input files.



When using **Automatic Detection**, if your data folder contains different types of data files, it is recommended that **Experiment selection** is used in order to restrict detection to GC/LC/MS and avoid analysis of other undesired files. In this case, you must select the **GC/LC/MS** experiment type, as shown below.




**Important!** When using **Automatic Detection**, make sure your data file extensions are included in the **Mask Manager** under the **Settings** tab so that Mgears can recognize the appropriate files correctly. Otherwise, manually add your standard files extension or custom file formats and press **OK**.






With Chrom Cal, the use of a **Mapping file** is mandatory to provide the necessary input information (i.e., compound information and concentrations) that allows the grouping and evaluation of the calibration curve of each compound in the samples. It also serves to provide the compound and concentration information for the reference samples used for cross-validation of the calibration curves.


To add and use a **Mapping file**, click on the **Filter** icon at the bottom left-hand side of the Mgears dialog. The two available filtering options (**Filtering String** and **Mapping File**) will appear. Click on  and upload your **Mapping File** (in either **.txt** or **.csv** format) with the data and parameters you want to work with. Then, click on  to set up the parameters to read the mapping file. The **Advanced Filtering Options** dialog will open.



In the **Advanced Filtering Options** dialog, you must set the column with the datafile names for matching.


 In our example, the **Well plate** position included in the file names (and not the whole file name) will be used for matching. For this to work correctly, the **Substring Matching** option must be enabled.

To configure matching columns, you can either indicate the number of the column directly by clicking on this box , or open the assistant  to visualize the CSV and select the corresponding column, as shown in the image below. Click on **OK** to save the selection.




**Note.** In the **CSV Column Detection** dialog, only the first 10 rows of the filtering file will be loaded.


**Important!** Remember that in order to be correctly read by Mgears, the CSV separator (comma, semicolon, tab, or vertical line) should be configured in the **Settings** tab (in this case, the CSV is using a comma as a separator, therefore we select the comma in the **Settings**).




If needed, you can set the **Starting Row** at which matching should start and decide whether to **Read Until** a specific row in the CSV. You can also match the **well plate position** so that the results will be displayed in the [Result Viewer](#) according to their actual well plate position.



Once you've completed configuring your filtering options, click on **OK**.

You can now perform an automatic inspection of the selected directory by clicking on  to verify that your input files will be found and filtered correctly. As you can see, the MS data files and sample positions are correctly detected.



**Watch out!** Depending on the operating system you have on your computer, Mnova may support different dataset format types. You can check the [MS supported formats](#) on our website.




## 1.2. Processing


In the **Processing** tab, you can upload a processing script to apply your own customized processing settings (please refer to the [Mnova Gears manual](#) for more details about this option).

## 1.3. Plugins

In the **Plugins** section, select and add the Chrom Cal plugin.




Then, click on **Chrom Cal Plugin Settings** to configure your analysis and report parameters.



A dialog with four tabs should appear.


**Note.** At this point, if a **Mapping File** is not configured in the **Input** tab, the following message will appear to remind you to add it.




### 1.3.1. The Chromatogram tab

In the **Chromatogram** tab, click on **Add** to configure the chromatogram extraction preferences. A dialog will open.


Select the analytical trace type (DAD, DAD range, PDA, or other trace) for chromatogram extraction. Then, choose whether to align the chromatograms automatically or manually; we recommend the use of a manual alignment. In such case, the value introduced will be added to the chromatogram in order to align it with the TIC. (Please refer to this [article](#) to see how you can calculate the time-shift required to align chromatograms manually with Mnova.)



In this tab, it is also possible to **Apply Trace Baseline Correction**. The **Folder** and **Smooth Factor** can be set manually by: (1) typing the value, (2) using the cursor, or automatically by (3) autodetection.



Once the **Chromatogram tab** is configured, the chosen parameters are displayed in the main **Settings** window, and the **Add** button is replaced by an **Edit** button.



### 1.3.2. The CSV Input tab

In the **CSV Input tab**, Mgears will automatically detect the **Mapping file** and will allow you to configure matching columns with the input information required for the analysis. You must match:

- The **Compound Id** with the column containing the name/number/ID of the analyte. This will allow grouping and evaluation of the datapoints for each analyte separately to generate the corresponding calibration curve. If the column is set to “zero”, all datafiles will be evaluated together and a single calibration curve will be obtained.

**Column 7**

D	E	F	G	H	I	J
Pos	Tube BC	ELN ID	Compound ID	Vol (uL)	Conc ( mM)	Parent MW
1	1.1E+09	1	Compound_357	30	3.333333333	259
2	1.2E+09	2	Compound_357	30	1.111111111	259
3	1.3E+09	3	Compound_357	30	0.37037037	259
4	1.4E+09	4	Compound_357	30	0.12345679	259
13	1.5E+09	13	Compound_148	30	3.333333333	608
14	1.6E+09	14	Compound_148	30	1.111111111	608
15	1.7E+09	15	Compound_148	30	0.37037037	608
16	1.8E+09	16	Compound_148	30	0.12345679	608
25	1.9E+09	25	Compound_505	30	3.333333333	343
26	2E+09	26	Compound_505	30	1.111111111	

Compound Id:

- The **Molecular Formula** and/or **Monoisotopic Mass** to allow correct peak assignment to the different analytes.

**Column 11                      13**

I	J	K	L	M	N	O
Conc ( mM)	Parent MW	Full MW	Parent MIW	Full Formula	Salt	Proje
3.333333333	259	295.81	259	C16H21NO2	HCl	
1.111111111	259	295.81	259	C16H21NO2	HCl	
0.37037037	259	295.81	259	C16H21NO2	HCl	
0.12345679	259	295.81	259	C16H21NO2	HCl	
3.333333333	608	608.69	608	C33H40N2O9	Parent	
1.111111111	608	608.69	608	C33H40N2O9	Parent	
0.37037037	608	608.69	608	C33H40N2O9	Parent	
0.12345679	608	608.69	608	C33H40N2O9	Parent	
3.333333333	343	343	343	C20H29		

Molecular Formula:   
Monoisotopic Mass:

- The samples' **Concentrations**.

**Column 20**

Q	R	S	T
ndor	ghed Mass	Validation	Conc(mg/mL)
YYY	4.13		0.986033333
YYY	4.13		0.328677778
YYY	4.13		0.109559259
YYY	4.13		0.036519753
YYY	5.84		2.028966666
YYY	5.84		0.676322222
YYY	5.84		0.225440741
YYY	5.84		0.075146914
YYY	3.59		1.112222222

Concentration: 20


- **Validation Point Marker** for result validation, if cross-validation is required. The column matched to this parameter must contain a "Y" or "YES" in the rows with the validation samples.

**Column 19**

G	H	I	J	K	L	M	N	O	P	Q	R	S	T
Compound ID	Vol (uL)	Conc ( mM)	Parent MW	Full MW	arent MIW	Full Formula	Salt	Project	Chemist	Vendor	ghed Mass	Validation	onc(mg/mL)
Compound_505	30	0.37037037	343	343	343	C20H29N3O2	Parent		XXXX	YYYY	3.59	Y	0.127037
Compound_357	30	3.333333333	259	295.81	259	C16H21NO2	HCl		XXXX	YYYY	4.13	Y	0.986033
Compound_148	30	1.111111111	608	608.69	608	C33H40N2O9	Parent					Y	0.676322

Validation Point Marker: 19


To configure matching columns, you can either indicate the number of the column directly by clicking on this box , or click on **Inspect CSV** to define input columns using the CSV assistant.



When using the **Inspect CSV** option, a **CSV File Input Definition** dialog will open with a table display of the **Parameters**, their assigned **Column**, column **Header**, and the **First Value** found in the file.

Parameter	Column	Header	First Value
1 Compound ID	7	Compund ID	Compound_357

Double-click on a “Column” cell to open the **CSV Column Detection** dialog and visualize the CSV, then select the column corresponding to the parameter and click **OK** to save the selection.



In the **CSV input** tab, you have also an option to **Evaluate All Compounds Together**. This option is interesting if a universal detector such as CAD is used, and when enabled, a calibration curve with all the compounds is generated in addition to the separate analyte calibration curves.

ID	Compound ID	Vol (uL)	Conc (mM)	Parent MW	Full MW	Area
1	Compound_357	30	3.333333333	259	295.81	2
2	Compound_357	30	1.111111111	259	295.81	2
3	Compound_357	30	0.37037037	259	295.81	2
4	Compound_357	30	0.12345679	259	295.81	2
5	Compound_148	30	3.333333333	608	608.69	6
6	Compound_148	30	1.111111111	608	608.69	6
7	Compound_148	30	0.37037037	608	608.69	6
8	Compound_148	30	0.12345679	608	608.69	6
9	Compound_505	30	3.333333333	343	343	3
10	Compound_505	30	1.111111111	343	343	3
11	Compound_505	30	0.37037037	343	343	3
12	Compound_505	30	0.12345679	343	343	3
13	Compound_357	30	1.111111111	259	295.81	2
14	Compound_357	30	0.37037037	259	295.81	2
15	Compound_357	30	0.12345679	259	295.81	2
16	Compound_148	30	3.333333333	608	608.69	6
17	Compound_148	30	1.111111111	608	608.69	6
18	Compound_148	30	0.37037037	608	608.69	6
19	Compound_148	30	0.12345679	608	608.69	6
20	Compound_505	30	3.333333333	343	343	3
21	Compound_505	30	1.111111111	343	343	3
22	Compound_505	30	0.37037037	343	343	3
23	Compound_505	30	0.12345679	343	343	3

Compound Id: 7

Analyzed together

Analyzed together

Analyzed together

Analyzed together

Evaluate All Compounds Together

## 1.3.3. The Controls tab


In the **Controls** tab, result acceptance criteria and cross-validation rules can be defined, and control flags can be configured to receive a color warning when a result does not fit the configured parameters. The available controls include:

- **Accepted R2:** If enabled, a threshold must be defined. A **Pass** ✓/**Fail** ✗ message will be displayed with the calibration result for the **R2** value.
- **Accepted Slope Range:** If enabled, a slope range must be defined. A **Pass** ✓/**Fail** ✗ message will be displayed with the calibration result for the **Slope** value.
- **Accepted Diff Conc (%):** If enabled, a threshold must be defined. A flag is obtained when the variation of the concentration (C) (calculated in percentage as  $[C_{\text{Calculated}} - C_{\text{Expected}}] / C_{\text{Expected}}$ ) exceeds the accepted threshold. The flag color and text can be customized.

Note that for this control to work, a **Validation Point Marker** must be duly defined in the **CSV Input** tab.

- **Peak assigned:** If enabled, a flag will be obtained when no peaks are detected for a component. The flag color and text can be customized.


Enable the option of interest, then click on the color selection button and choose the flag color. Click on **OK** to save your preferences.



Finalize your plugin settings setup by pressing **OK** and move to the next step.

### 1.3.4. The Report tab

In the **Report** tab, choose the report formats you wish to obtain for your analysis results. **Mnova**, **PDF**, **CSV** reports can all be generated.




## 1.4. Design

In the **Design** tab, it is possible to load an Mnova layout template or use a custom script to produce a final customized report.


In this analysis, we are not going to use custom designs. The results will be laid out using the default configuration.

## 1.5. Output

Choose a directory to which Mgears can save the results of the analysis.




Optionally, enable the **Add Nickname to the Results Folder** and type the nickname of your choice. You can also **Add Incremental Numbering** after the chosen nickname.




Enable the Mnova document creation to be able to view/review the spectra along with the results in the **Mgears Result Viewer** (if not enabled, when loading the results in the Mgears Result Viewer, only the data and Calibration curves will be visible, with no spectra).

You can also choose to create PDF reports with the individual results or save the results to a database (*please refer to the [Mnova Gears manual](#) for more details on the output configuration options*).

Now that everything is set, click on  **Run** to start the analysis.

## 2. The Output Folder

The analysis results obtained are saved in the directory selected in the **Output** tab and is named according to the date and time stamp at which the execution was carried out, in addition to the nickname and number added.



### 2.1. The HTML file

In the HTML report, an overview of the results is available. Each sample is reported in a row with **Controls charts**, a preview of the **MS spectrum**, and hyperlinks to the **Mnova** and **PDF** files generated.

### Mgears ChromCal Results

Parameters

Parameter	Value
Results Directory	C:/Users/Usuario/Desktop/Results/Chrom Cal/2022-05-05T12.30.18_Mytest-3
Started On	2022-05-05T12:30:18
Completed On	2022-05-05T12:31:22

Detailed Results

Show  entries

Copy
CSV
Columns
PDF
Print


Search:

#	fileName	name	Controls Chart	MS	Mnova File	Pdf File
1	STD01-03-B03.lcd	B03			<a href="#">B03.mnova</a>	<a href="#">B03.pdf</a>
2	STD01-02-B02.lcd	B02			<a href="#">B02.mnova</a>	<a href="#">B02.pdf</a>
3	STD01-01-B01.lcd	B01			<a href="#">B01.mnova</a>	<a href="#">B01.pdf</a>
4	STD31-04-C04.lcd	C04			<a href="#">C04.mnova</a>	<a href="#">C04.pdf</a>
5	STD31-03-C03.lcd	C03			<a href="#">C03.mnova</a>	<a href="#">C03.pdf</a>



You can change the display of the columns and adapt them to your preferences. Click on **Columns** and select/unselect the parameters you wish to display/hide.


You can also **Copy**, **Print**, or save those results into **CSV** or **PDF** formats.



## 2.2. The CSV reports

Three types of CSV reports are obtained for each calibration analyte when the **CSV** option is checked in the [Reports](#) tab of the settings. These reports are saved under a separate sub-folder for each analyte under the “ChromCalReports” folder:


- The **LinearFitResults** reports the calibration curve parameters.



### LinearFitResults

	A	B	C	D
1	Slope	Intercept	R-Square	Validation R-Square
2	640565.1	41360.54923	0.984563711	NO


- The **PointWiseResults** reports the results for each data file included in the calibration of the compound.



### PointWiseResults

	A	B	C	D	E
1	File Name	Peak Area	Concentration	Estimated Peak Area	Area Diff ^2
2	STD01-03-B03.lcd	569114.2342	0.676322222	474588.9467	8935029969
3	STD01-02-B02.lcd	-	0.225440741	-	-
4	STD01-01-B01.lcd	55016	0.075146914	89497.03813	1188941990
5	STD01-04-B08.lcd	1348666	2.028966666	1341045.742	58068336.66
6	STD01-03-B07.lcd	569114.2342	0.676322222	474588.9467	8935029969
7	STD01-02-B06.lcd	-	0.225440741	-	-
8	STD01-01-B05.lcd	55016	0.075146914	89497.03813	1188941990
9	STD01-04-G04.lcd	1280416	2.028966666	1341045.742	3675965578
10	STD01-02-G02.lcd	118691	0.225440741	185770.0153	4499594290
11	STD01-01-G01.lcd	-	0.075146914	-	-

- The **CrossValidation** reports are obtained only when [Cross-validation points](#) are configured.



### CrossValidation

	A	B	C	D	E	F	G
1	File Name	Peak Area	Estimated Concentration	Expected Concentration	% Conc Diff	Validation	Flag
2	STD01-03-G03.lcd	579811.2	0.840586958	0.676322222	24.28794001	NO	Cross-validation Failed

## 2.3. The Mnova and PDF reports

Both Mnova and PDF reports are obtained for each calibration analyte and saved under the corresponding sub-folder when the **Mnova and PDF** option is checked in the [Reports](#) tab of the settings:

**Mnova**

**PDF**

**POINT-WISE RESULTS**

File Name	Peak Area	Concentration	Estimated Peak Area	Area Diff %2
STD01-03-803.kci	569114.23	0.6783	474588.95	8.94e+9
STD01-02-802.kci	-	0.2254	-	-
STD01-01-801.kci	55016.00	0.0751	89497.04	1.19e+9
STD01-04-804.kci	1348666.00	2.0290	1341045.74	5.81e+7
STD01-03-807.kci	569114.23	0.6783	474588.95	8.94e+9
STD01-02-806.kci	-	0.2254	-	-
STD01-01-805.kci	55016.00	0.0751	89497.04	1.19e+9
STD01-04-004.kci	1280416.00	2.0290	1341045.74	3.88e+9
STD01-02-002.kci	118891.00	0.2254	185770.02	4.50e+9
STD01-01-001.kci	-	0.0751	-	-

**LINEAR FIT RESULTS**

Slope	Intercept	R-Square	Validation R-Square
640565.08	41360.55	0.9848	X

Validation R-Square  $R^2 = 0.98$

**POINT-WISE RESULTS**


File Name	Peak Area	Concentration	Estimated Peak Area	Area Diff %2
STD01-03-803.kci	569114.23	0.6783	474588.95	8.94e+9
STD01-02-802.kci	-	0.2254	-	-
STD01-01-801.kci	55016.00	0.0751	89497.04	1.19e+9
STD01-04-804.kci	1348666.00	2.0290	1341045.74	5.81e+7
STD01-03-807.kci	569114.23	0.6783	474588.95	8.94e+9
STD01-02-806.kci	-	0.2254	-	-
STD01-01-805.kci	55016.00	0.0751	89497.04	1.19e+9
STD01-04-004.kci	1280416.00	2.0290	1341045.74	3.88e+9
STD01-02-002.kci	118891.00	0.2254	185770.02	4.50e+9
STD01-01-001.kci	-	0.0751	-	-

**LINEAR FIT RESULTS**

Slope	Intercept	R-Square	Validation R-Square
640565.08	41360.55	0.9848	X

Validation R-Square  $R^2 = 0.98$

**Note.** When the [Evaluate All Compounds Together](#) option is enabled, **CSV, Mnova** and **PDF** reports for the evaluation, with all the available data points together, are saved in a sub-folder named “allComponents”.



## 2.4. Other output

- A “Chrom Cal” directory is created containing the analysis output (JSON files) with data and calibration curves that will be loaded and reviewed in the **Mgears Result Viewer**. These files can also be used as input in downstream applications that need to import the calibration parameters, e.g., [Chrom Reaction Optimization](#).
- A “documents” directory, containing the output individual Mnova files for each data file (unless Mgears is configured to save Mnova files in another location).
- A log file of the execution.
- A copy of the settings used in the current evaluation.
- A resume file of the steps followed in the execution.
- A CSS folder, a data folder, a JS folder, and an images folder.

## 3. Mnova Gears Results Viewer

Chrom Calibration has two different – but interconnected – result viewers, the **Mgears Viewer** and the **Chrom Cal Viewer**, to display individual calibration curve results.


Open the **Mgears Viewer** from the Mnova **Automation** tab, then click on  to open your analysis results. When the experiment is open, both **Mgears** and **Chrom Cal Viewers** are launched automatically.






### 3.1. The Mgears viewer

In the **Mgears Viewer** you can see the names of the dataset files, their positions in the well plate, and any associated numerical results. Click on a specific dataset to view the corresponding results and spectrum in Mnova.



#### 3.1.1. The Well plate view

The well plate can show either data file **Names** or **Controls charts**. You can change the display by selecting the corresponding option for the **Result in Well Plate** under the **Results** section.



#### File names

Well Plate		1	2	3	4	5	6	7	8
A		undefined	undefined	undefined	undefined	A05	A06	A07	A08
B		B01	B02	B03	undefined	B05	B06	B07	B08
C		C01	C02	C03	C04	C05	C06	C07	C08
D									
E									
F		F01	F02	F03	F04				
G		G01	G02	G03	G04				
H		H01	H02	H03	H04				

Failed Cross-Validation

#### Controls charts


Well Plate		1	2	3	4	5	6	7	8
A		undefined	undefined	undefined	undefined				
B					undefined				
C									
D									
E									
F									
G									
H									

Failed Peak assignment

### 3.1.2. The Results section

The **Results** section can include up to three tabs with information on:

- Peak-specific results relative to the datafile in question under the **Calibration** tabs; when the [Evaluate All Compounds Together](#) option is enabled, two **Calibration** tabs are available, one for the analyte evaluation (**Compound\_148**) and another for the overall evaluation (**All**).
- **Controls**, as defined in the [Controls tab](#) of the settings. This tab is only visible if at least one analysis control is enabled.



The screenshot displays three overlapping windows from the software interface. The top window is titled 'Calibration (Compound\_148)' and shows a table with 7 rows of parameters and their values. The middle window is titled 'Calibration (All)' and also shows a table with 7 rows of parameters and their values. The bottom window is titled 'Controls' and shows a table with 1 row of control information, including a 'Peak Assignment' control that has passed.

Parameter	Value
1 Calibration	Compound_148
2 Concentration	0.68
3 Area	569114.23
4 Area Calc	474588.95
5 Area Diff	94525.29
6 Area Diff^2	8935029968.9
7 Peak RT	2.33


Parameter	Value
1 Calibration	allComponents
2 Concentration	0.68
3 Area	569114.23
4 Area Calc	492150.06
5 Area Diff	76964.17
6 Area Diff^2	5923483732.6
7 Peak RT	2.33

Control	Value	Passed	Flag
1 Peak Assignment		✓	

### 3.1.3. Assigning peak/Modifying peak assignment

Not only you can visualize analysis results with the **Mgears viewer**, but also you can assign peaks to samples that failed peak assignment, or modify peak assignment manually and update your results.

To do so, open your sample of interest (in this example, we are going to assign a peak to a sample that failed automatic assignment), then click on the  button. A **Manual Peak Selection** dialog will open to allow you choose a peak in the spectrum and assign it to your compound. Select **peak** and press **OK**. Mgears will automatically recalculate and update the results (including output files and reports) for the selected well, as can be seen in the image below.




The screenshot shows the Mgears Viewer interface. In the 'Results' section, the 'Chrom Cal' viewer displays a 'Controls Chart' with a table of controls. The first control, 'Peak Assignment', has a red 'X' in the 'Passed' column and the text 'Peak Assignment Failed'. A red circle '1' highlights a magnifying glass icon. A 'Manual Peak Selection' dialog box is open, showing a table of peaks from 'Chromatogram 0 Da Peaks':


RT	Area
1 0.582	1888.000
2 1.499	2670.000
3 2.340	59784.000
4 2.582	8904.000
5 2.840	11760.000
6 3.165	143503.000

A red circle '2' points to the peak at RT 2.340. A 'Recalculating...' dialog box is also visible with a red circle '4' and a 'Cancel' button. A red circle '3' points to the 'OK' button in the 'Manual Peak Selection' dialog.

The screenshot shows the Mgears Viewer interface after manual peak selection. The 'Controls Chart' table now shows the 'Peak Assignment' control with a green checkmark in the 'Passed' column. A red circle '5' highlights the 'Peak Assignment' row. The 'Chrom Cal' viewer shows the 'Chromatogram 0 Da Peaks' table with the peak at RT 2.340 highlighted in green. A red circle '5' points to this peak. The 'Recalculating...' dialog box is still present.

You can also manually edit the spectrum peaks using the Mnova tools. In this case, it will be necessary to press the **Analyze again**  button in the **Mgears Viewer** to update the results according to your changes.

**Note.** When modifying peak assignment manually, the selection of the peak in the different calibration curves is not affected. Therefore, it may be necessary to review the calibration result in the **Chrom Cal** viewer and update peak selection (to include or exclude it), as explained in [section 3.2.4](#).

If at any point you close the Chrom Cal Viewer, it is possible to open it again by clicking on the  icon in the **Results** section.

## 3.2. The Chrom Cal viewer

In the **Chrom Cal Viewer**, you can view the calibration results (calibration and cross-validation points, calibration parameters, and plot) for each analyte.

The screenshot shows the Chrom Cal Viewer interface for 'Compound\_505'. It is divided into four main sections:

- Calibration points:** A table with columns 'Sample', 'C', 'Area', and 'Area Diff^2'. It lists three samples: H04, H03, and H02.
- Cross-validation points:** A table with columns 'Validation', 'Area', 'C', 'C (Ref)', and 'C (% Dev)'. It lists one validation point: C02.
- Calibration parameters:** A table with columns 'Parameter', 'Value', and 'Passed'. It lists R2, Slope, and Intercept.
- Plot Area:** A scatter plot with 'Concentration' on the x-axis and 'Peak Area' on the y-axis. It shows four data points and a linear regression line.

You can select the results to display from the **Calibration** options list.

This screenshot shows the 'Calibration' dropdown menu in the Chrom Cal Viewer. The menu is open, showing a list of options: 'allComponents', 'Compound\_148', 'allComponents', 'Compound\_505', and 'Compound\_357'. The 'Compound\_505' option is currently selected and highlighted in blue. The background table shows calibration points for samples B03, B02, and B01.

### 3.2.1. The Calibration and Cross-validation points

For each calibration, the individual samples are listed in the calibration points table with their **Concentration**, **Peak Area**, and **Area Diff^2** (indicating the level of deviation of a point from the calibration curve).

	Sample	C	Area	Area Diff^2
8	<input checked="" type="checkbox"/> H04	1.143	9.32e+5	4.49e+8
9	<input checked="" type="checkbox"/> H03	0.381	3.39e+5	4.55e+8
10	<input checked="" type="checkbox"/> H02	0.127	1.12e+5	6.46e+7



The cross-validation samples are also listed in a separate table with their **Peak Area**, **Concentration**, **Reference Concentration**, and the **Percentage of Deviation** in the concentration value.

	Validation	Area	C	C (Ref)	C (% Dev)
1	<input checked="" type="checkbox"/> C02	1.19e+5	0.125	0.127	1.6

Selected samples  are included in the global calibration result and plot, whereas unselected samples  are excluded from the global result and plot.

### 3.2.2. The Calibration parameters

For each calibration, the **R2**, **Slope**, and **Intercept** values are displayed in a table with the test result (if configured in the [Controls tab](#)), as can be seen below:




#### Compound\_505


Parameter	Value	Passed
<b>R2</b>	0.995	✓
<b>Slope</b>	777367.2	
<b>Intercept</b>	21654.7	

#### Compound\_148

Parameter	Value	Passed
<b>R2</b>	0.986	✗
<b>Slope</b>	645493.3	
<b>Intercept</b>	33139.0	


### 3.2.3. The Plot Area

For each calibration, the different samples' **Concentrations** and **Peak areas** are plotted, and the fitted calibration curve is represented as can be seen in the example below. **Zoom in**  and **Zoom out**  tools are made available to be able to inspect the graph with more flexibility. This button  will reset the calibration graph to its original size fitting to the plot area.




When hovering the cursor over either a calibration or cross-validation point in the plot:


- a tool tip with the name of the sample or validation point, the concentration, and peak area will be displayed: “Name (Concentration, Area)”; and,
- the corresponding row in the results table will be highlighted.



### 3.2.4. Updating the calibration/cross-validation results

It is possible to update the calibration results by selecting/unselecting samples to include/exclude them from the analysis. Chrom Cal will automatically update the results in the viewer as well as in the output files and reports.



If manual changes are introduced to the spectra using the Mnova tools the **Analyze again**  button in the **Chrom Cal Viewer** should be used to save the updated results to the output folder and reports.

**For more details on the Mnova Gears options, please refer to the [Mnova Gears manual](#).**