

A Novel Benchtop Time-of-Flight GC-MS System for High Performance Analysis of Human Urine

David E. Alonso, Joseph E. Binkley, Lorne M. Fell | LECO Corporation, St. Joseph, MI

Introduction

The development of improved profiling workflows are critical for effective investigation of environmental exposure and its relationship to chronic diseases. Metabolomic methods are particularly useful in this endeavor due to the proximity of metabolites to system phenotype, and their relatively quick insight into system perturbations from diet, drugs, gut microbiota, etc. (Figure 1). Unfortunately, a major problem with these methods continues to be incomplete sample characterization due to metabolite diversity, wide concentration range, and the overall complexity of biological matrices.

Urine is an ideal bio-fluid for studying the effects of exposure, since it is easy to obtain in large volumes, and is relatively free from interfering proteins and lipids. In addition, high concentrations of drugs and other xenobiotic materials can be detected over extended periods of time. In this study, we used a high performance, benchtop gas chromatography time-of-flight mass spectrometer for comprehensive characterization of derivatized urine samples (Figure 2).




Figure 1. Human Exposure (Schantz and coworkers, *Anal. Bioanal. Chem.*, 2015, 407, 2945-2954).



Figure 2. Pegasus® BT GC-TOFMS.

Objectives

- Develop and implement an effective workflow to facilitate the comparison of urine samples
- Use time-proven chemical protocol to expand metabolite coverage by producing stable, GC-MS amenable compounds
- Use high performance GC-TOFMS and powerful software tools to quickly and confidently annotate metabolites in urine

Sample Preparation

Urine (200 µL) was treated with urease (10 mg; 37°C for 15 minutes) to remove excess urea. This was followed by the addition of methanol (800 µL) to produce a solution that was vortexed for 10 minutes, and centrifuged at 12,000 rpm for 10 minutes. The supernatant was transferred to a GC vial with a 400 µL insert, and evaporated to dryness using a Speed Vac. 100 µL of MTBSTFA + 1% TBMCS was added to the dry sample which was heated at 100°C for 1 hour, and then injected into a high performance GC-TOFMS using a split ratio of 20:1. The total acquisition time was 12 minutes (Figure 3).

Instrumentation Parameters and Methods

Gas Chromatograph	Agilent 7890 with MPS2 Autosampler
Injection	1 µL, Split 20:1
Carrier Gas	He @ 0.8 ml/min, Constant Flow
Column	Rxi-5ms, 20 m x 0.18 mm i.d. x 0.18 µm (Restek, Bellefonte, PA, USA)
Temperature Program	60°C (0.5min), ramped 36°C/min to 330°C, held 4 min
Mass Spectrometer	LECO Pegasus® BT
Ion Source Temperature	250°C
Ionization Mode	EI
Mass Range (m/z)	45-600
Acquisition Rate	20 spectra/s




Figure 3. TICs for NIST SRMs 3673 (top) and 3672 (bottom).

Results (Non-Smokers' Urine)




Figure 4. (A) Analytical Ion Chromatogram (AIC) and (B) TIC of Non-Smokers' Urine (NIST SRM 3673). (C) Table of Selected Compounds in SRM 3673.




Figure 5. Peak True (Deconvoluted) and Library Spectra for BHT (A,B) and Ibuprofen (C,D).




Figure 6. Peak True and Library Spectra for Salicylic Acid (A,B) and Hippuric Acid (C,D).

Results (Smokers' Urine)




Figure 7. (A) Analytical Ion Chromatogram (AIC) and (B) TIC of Smokers' Urine (NIST SRM 3673).




Figure 8. Peak True and Library Spectra for Theobromine (A,B) and Naproxen (C,D).




Figure 9. Peak True and Library Spectra for p-Cresol (A,B) and Santonox (C,D).

Target Analyte Finding (TAF) – Quick Processing for Sample Comparison




Figure 10. (A) ChromaTOF 5.0 TAF Method and (B) Target Ion(s).




Figure 11. TAF AICs (A,B) and Results (C,D). Showing Greater Quantities of Target Compounds in Smokers' versus Non-Smokers Urine.

Conclusions

A workflow was developed and implemented for the effective comparison of urine. It involved the preparation of stable compounds through derivatization, data acquisition using a high performance benchtop TOFMS, comprehensive data processing (NonTarget Deconvolution™), and quick retrospective analysis of rich datasets using Target Analyte Finding.

