


New Synthetics and New Innovations:


Synthetic fentanyl analogues and the new Agilent Ultivo Tandem Mass Spectrometer



Julie Cichelli, Ph.D.
Agilent Technologies,
Wilmington, DE.



Introduction




- Opioid crisis is becoming more severe, with reportedly approx. 100 deaths from OD per day across the USA.
- Intelligence suggesting that historically abused opiates such as Heroin, are being 'cut' with many new synthetic, more potent opioids such as Fentanyls.
- Keeping up with the new synthetic illicit and analyzing for these is difficult.
- Herein we describe a Tandem LC/MS/MS method for a list of 12x new synthetic fentanyl opioids and 4-ANPP, the precursor scheduled by the DEA.
- Details of the methodology will be presented together with results obtained from Serum and Urine matrices using the new Innovative Agilent Ultivo Tandem Mass Spectrometer.

Coming together of New Innovations, some good, some bad!

- The recent emergence of synthetic fentanyl-type opioids has led to the development and continued development of new, innovative certified standards for the forensic laboratories;
- Thankfully, we are constantly and increasingly able to accurately confirm and quantify more synthetics using such certified standards;
- Modern LC/MS Instrumentation for such analyses is also being introduced;
- Agilent Ultivo Tandem MS is such instrumentation, and therefore, it was put through this test to see how well this new innovation would prove to the new and emerging situations of the forensic toxicology community;
- The following summary will outline typical results obtained from such an application using the Agilent Ultivo TMS and will detail its own innovations.

List of Synthetic Fentanyl Opioids for this study

- 4-ANPP* (Synthetic precursor)
- 3-methylfentanyl
- Acetylfentanyl*
- Acetylnorfentanyl*
- Acrylfentanyl*
- Butyrylfentanyl*
- Carfentanil
- Furanylfentanyl
- N-desmethyl U-47700
- Norcarfentanyl
- Para-fluorobutyrylfentanyl
- Valerylfentanyl*
- U-47700*
- **W-18 RM**



*Asterix denotes ISTD also obtained

Sample Preparation

(Standards and ISTDs obtained from Cerilliant & Cayman Chemicals)

(Serum and negative urine obtained from Golden West Biologicals)

Urine Matrix:

1. Negative Urine diluted 1/10 using deionized water;
2. 1/10 diluted negative Urine Centrifuged at 10,000 RPM through nanosep 3K filters then supernatant combined as a stock solution;
3. ISTDs Spiked into diluted negative urine stock at a concentration of 25ng/ml;
4. Serial Dilutions made for calibrators from 500ng/ml down to 1pg/ml using ISTD-spiked filtered negative Urine stock. (Equivalence dilution = 1/10)

Serum Matrix:

1. Gold Serum (10mL) mixed with 20mL of cold acetonitrile;
2. Vortexed for 2 minutes and centrifuged at 5,000 RPM for 15 minutes;
3. Supernatant (20mL) spiked with ISTD at a concentration of 50ng/ml;
4. Supernatant & ISTDs diluted with 20mL deionized water. ISTD concn. 25ng/ml.
5. Serial Dilutions made for calibrators 500ng/ml to 1pg/ml.
6. (Equivalence dilution = 1/6)

LC Conditions - 1290 UHPLC

HPLC: Agilent 1290 Binary pump, well plate autosampler, thermostatted column compartment

| LC Conditions | | | |
|-------------------------|--|-------|--------------------|
| Column | Poroshell 120 EC-C18, 2.7um, 2.1x100 mm | | |
| Column temperature | 55 °C | | |
| Injection volume | 5.0 µL (injector loop 20uL) | | |
| Draw speed (uL/min) | 100 | | |
| Eject speed (uL/min) | 100 | | |
| Autosampler temperature | 6 °C | | |
| Needle wash | 5 seconds in wash port (100% methanol) | | |
| Mobile phase | A = 5mM Ammonium Formate/0.01% Formic acid in water B = 0.01% formic acid in Methanol | | |
| Gradient program | Time | B (%) | Flow rate (mL/min) |
| | 0.00 | 10 | 0.500 |
| | 0.50 | 15 | 0.500 |
| | 3.50 | 50 | 0.500 |
| | 5.50 | 95 | 0.500 |
| | 6.00 | 95 | 0.500 |
| Stop time | 6.0 min | | |
| Post time | 1,0 min | | |




MS Conditions- Agilent Ultivo TMS

Instrument and Source conditions:

| Parameters | |
|------------------------|-------------------|
| Ion mode | AJS ESI, positive |
| Gas temperature | 325 °C |
| Drying gas flow | 9 L/min |
| Nebulizer gas | 35 psi |
| Sheath gas temperature | 350°C |
| Sheath gas flow | 11 L/min |
| Capillary voltage | 3500 V |
| Nozzle voltage | 0 V |

| MRM Parameters | |
|--------------------------------|-------------------------|
| Delta EMV | N/A |
| Q1 and Q3 resolution (amu) | 0.7/0.7 (unit/unit) |
| Fragmentor Voltage | Various – See MRM Table |
| Cell Accelerator voltage | 9 V |
| Cycle time | 400 msec |
| Total MRMs | 38 |
| Dwell Time each MRM Transition | 10 msec |



MS Conditions (Contd.) - Agilent Ultivo TMS

MRM Instrument Settings:

| | Compound | Compound Name | ISTD | Precursor Ion | MS1 Resolution | Product Ion | MS2 Resolution | Dwell | Fragmentor | CAV | Collision Energy | Polarity |
|----|----------|----------------------------|------|---------------|----------------|-------------|----------------|-------|------------|-----|------------------|----------|
| 1 | | W-18 RM | No | 422.1 m/z | Unit | 174.9 m/z | Unit | 10 ms | 175 V | 9 V | 32 V | + |
| 2 | | W-18 RM | No | 422.1 m/z | Unit | 110.9 m/z | Unit | 10 ms | 175 V | 9 V | 56 V | + |
| 3 | | Furanyl fentanyl | No | 375.2 m/z | Unit | 188.1 m/z | Unit | 10 ms | 170 V | 9 V | 24 V | + |
| 4 | | Furanyl fentanyl | No | 375.2 m/z | Unit | 105 m/z | Unit | 10 ms | 170 V | 9 V | 48 V | + |
| 5 | | Valeryl fentanyl-D5 | Yes | 370.3 m/z | Unit | 105 m/z | Unit | 10 ms | 180 V | 9 V | 48 V | + |
| 6 | | para-fluorobutyrylfentanyl | No | 369.2 m/z | Unit | 188.1 m/z | Unit | 10 ms | 180 V | 9 V | 24 V | + |
| 7 | | para-fluorobutyrylfentanyl | No | 369.2 m/z | Unit | 105 m/z | Unit | 10 ms | 180 V | 9 V | 52 V | + |
| 8 | | Valeryl fentanyl | No | 365.3 m/z | Unit | 188.1 m/z | Unit | 10 ms | 180 V | 9 V | 24 V | + |
| 9 | | Valeryl fentanyl | No | 365.3 m/z | Unit | 105 m/z | Unit | 10 ms | 180 V | 9 V | 48 V | + |
| 10 | | Butyrylfentanyl-D5 | Yes | 356.3 m/z | Unit | 188.1 m/z | Unit | 10 ms | 180 V | 9 V | 24 V | + |
| 11 | | Butyrylfentanyl-D5 | Yes | 356.3 m/z | Unit | 105.1 m/z | Unit | 10 ms | 180 V | 9 V | 48 V | + |
| 12 | | cis-3-methyl fentanyl | No | 351.2 m/z | Unit | 202.1 m/z | Unit | 10 ms | 180 V | 9 V | 24 V | + |
| 13 | | Butyrylfentanyl | No | 351.2 m/z | Unit | 188.1 m/z | Unit | 10 ms | 180 V | 9 V | 24 V | + |
| 14 | | Butyrylfentanyl | No | 351.2 m/z | Unit | 105 m/z | Unit | 10 ms | 180 V | 9 V | 48 V | + |
| 15 | | cis-3-methyl fentanyl | No | 351.2 m/z | Unit | 105 m/z | Unit | 10 ms | 180 V | 9 V | 48 V | + |
| 16 | | Acrylfentanyl-D5 | Yes | 340.2 m/z | Unit | 188.1 m/z | Unit | 10 ms | 165 V | 9 V | 24 V | + |
| 17 | | Acrylfentanyl-D5 | Yes | 340.2 m/z | Unit | 105 m/z | Unit | 10 ms | 165 V | 9 V | 44 V | + |
| 18 | | U-47700-D6 | Yes | 335.2 m/z | Unit | 284 m/z | Unit | 10 ms | 120 V | 9 V | 16 V | + |
| 19 | | Acrylfentanyl | No | 335.2 m/z | Unit | 188.1 m/z | Unit | 10 ms | 165 V | 9 V | 24 V | + |
| 20 | | Acrylfentanyl | No | 335.2 m/z | Unit | 105 m/z | Unit | 10 ms | 165 V | 9 V | 44 V | + |


MS Conditions (Contd. 2) - Agilent Ultivo TMS

MRM Instrument Settings:


| | | | | | | | | | | | | |
|----|--|-------------------------|-----|-----------|------|-----------|------|-------|-------|-----|------|---|
| 21 | | Acetylfentanyl-13C6 | Yes | 329.2 m/z | Unit | 105 m/z | Unit | 10 ms | 170 V | 9 V | 44 V | + |
| 22 | | U-47700 | No | 329.1 m/z | Unit | 284 m/z | Unit | 10 ms | 120 V | 9 V | 16 V | + |
| 23 | | U-47700 | No | 329.1 m/z | Unit | 172.9 m/z | Unit | 10 ms | 120 V | 9 V | 36 V | + |
| 24 | | Acetylfentanyl | No | 323.2 m/z | Unit | 188.1 m/z | Unit | 10 ms | 170 V | 9 V | 24 V | + |
| 25 | | Acetylfentanyl | No | 323.2 m/z | Unit | 105.1 m/z | Unit | 10 ms | 170 V | 9 V | 44 V | + |
| 26 | | N-desmethyl U-47700 | No | 315.1 m/z | Unit | 284 m/z | Unit | 10 ms | 120 V | 9 V | 16 V | + |
| 27 | | N-desmethyl U-47700 | No | 315.1 m/z | Unit | 172.9 m/z | Unit | 10 ms | 120 V | 9 V | 36 V | + |
| 28 | | Norcarfentanil | No | 291.2 m/z | Unit | 231.1 m/z | Unit | 10 ms | 95 V | 9 V | 12 V | + |
| 29 | | Norcarfentanil | No | 291.2 m/z | Unit | 142 m/z | Unit | 10 ms | 95 V | 9 V | 16 V | + |
| 30 | | 4-ANPP-D5 | Yes | 286.2 m/z | Unit | 105.1 m/z | Unit | 10 ms | 140 V | 9 V | 36 V | + |
| 31 | | 4-ANPP | No | 281.2 m/z | Unit | 188.1 m/z | Unit | 10 ms | 140 V | 9 V | 16 V | + |
| 32 | | 4-ANPP | No | 281.2 m/z | Unit | 105 m/z | Unit | 10 ms | 140 V | 9 V | 36 V | + |
| 33 | | Acetyl norfentanyl-13C6 | Yes | 225.2 m/z | Unit | 84.1 m/z | Unit | 10 ms | 120 V | 9 V | 20 V | + |
| 34 | | Acetyl norfentanyl | No | 219.1 m/z | Unit | 84.1 m/z | Unit | 10 ms | 120 V | 9 V | 20 V | + |
| 35 | | Acetyl norfentanyl | No | 219.1 m/z | Unit | 55.2 m/z | Unit | 10 ms | 120 V | 9 V | 44 V | + |
| 36 | | Carfentanil-D5 | Yes | 400.3 m/z | Unit | 113 m/z | Unit | 10 ms | 130 V | 9 V | 36 V | + |
| 37 | | Carfentanil | No | 395.2 m/z | Unit | 335.2 m/z | Unit | 10 ms | 135 V | 9 V | 16 V | + |
| 38 | | Carfentanil | No | 395.2m/z | Unit | 113 m/z | Unit | 10 ms | 135 V | 9 V | 36 V | + |

Overlaid EIC of each Fentanyl Opioid Transition

(Serum Mix Calibrant 100ng/ml)




Typical Results – acetyl fentanyl in Urine (1pg/ml→500ng/ml)



LLOQ of 50pg/mL

$R^2=0.9999$, $N=3$


Typical Results – acetyl fentanyl in Serum (1pg/ml→500ng/ml)



LLOQ of 10pg/mL

$R^2=0.9997, N=3$


Typical Results – acrylfentanyl in Urine (1pg/ml→500ng/ml)



LLOQ of 50pg/mL

$R^2=0.9998$, $N=3$


Typical Results – acrylfentanyl in Serum (1pg/ml→500ng/ml)



LLOQ of 10pg/mL

$R^2=0.9994$, $N=3$


Typical Results – 4-ANPP in Urine (1pg/ml→500ng/ml)



LLOQ of 100pg/mL

$R^2=0.9998, N=3$


Typical Results – 4-ANPP in Serum (1pg/ml→500ng/ml)



LLOQ of 10pg/mL

$R^2=0.9979$, $N=3$


Typical Results – Carfentanyl in Urine (1pg/ml→500ng/ml)



LLOQ of 50pg/mL

$R^2=0.9992$, $N=3$

Typical Results – Carfentanyl in Serum (1pg/ml→500ng/ml)



LLOQ of 10pg/mL

$R^2=0.9995$, $N=3$

LLOQ Results Table for Fentanyl synthetics – Agilent Ultivo TMS

| Analyte | Actual LLOQ (Urine) pg/ml | Actual LLOQ (Serum) pg/ml | Equivalence LLOQ (Urine) pg/ml (1/10 dilution) | Equivalence LLOQ (Serum) pg/ml (1/6 dilution) |
|----------------------------|---------------------------|---------------------------|--|---|
| 4-ANPP | 100 | 10 | 1000 | 60 |
| 3-methylfentanyl | 100 | 50 | 1000 | 300 |
| Acetylfentanyl | 50 | 10 | 500 | 60 |
| Acetylnorfentanyl | 100 | 50 | 1000 | 300 |
| Acrylfentanyl | 50 | 10 | 500 | 60 |
| Butyrylfentanyl | 100 | 50 | 1000 | 300 |
| Carfentanil | 50 | 10 | 500 | 60 |
| Furanylfentanyl | 50 | 50 | 500 | 300 |
| Para-fluorobutyrylfentanyl | 50 | 50 | 500 | 300 |
| Norcarfentanyl | 100 | 50 | 1000 | 300 |
| N-desmethyl U-47700 | 100 | 50 | 1000 | 300 |
| Valerylfentanyl | 100 | 50 | 1000 | 300 |
| U-47700 | 100 | 50 | 1000 | 300 |
| W-18 RM | 50 | 50 | 500 | 300 |

*Equivalence refers to the actual concentration of the Urine or Serum prior to dilution

Innovations: Ultivo Triple Quad LC/MS

AGILENT IS RESHAPING THE
FUTURE OF
MASS SPECTROMETRY.



JOIN US AND WELCOME IN A NEW
ERA OF "FIT FOR PURPOSE" MASS
SPECTROMETRY.

Typical LC/TQ Lab Layout




6 feet
1.83 meters

LC/TQ Instrumentation

Vacuum $\sim 10^{-5}$ torr

Source & Ion Transfer Optics

Quadrupoles




U.S. Patent Nov. 18, 1980 Sheet 2 of 3 4,234,791

Collision Cell


Detector

Electronics

Ultivo Ion Guide, Mass Analyzer and Detector





ULTIVO: Designed for Robust Performance



Ultivo: VacShield

- Fast, ion injector exchange without system venting
- Easy-to-use, increased throughput, less downtime for maintenance





Agilent Technologies

Pre MS1 Ion Beam Compressor – “Cyclone Ion Guide”




Shorter than octopole.

Transport ions through multiple Vacuum stages


Transport Ions with minimal ion losses

Twisted and tapered


Reduce ion beam size



Cyclone Transmission vs m/z – outer rods off




Cyclone Transmission vs m/z – outer rods on



Benefits of Cyclone Ion Guide

- Twisted, tapered (2) hexapole set
 - Convergence section – greater ion acceptance
 - Outer rods – increased high mass transmission efficiency (additional confinement)
 - Resistive rods – DC bias to gently pull ions towards the exit
 - Parallel section – superior collimation of the ion beam
- ***More ions, increased signals, better results***









Agilent Technologies

Traditional Collision Cell

- Consists of **multiple** parallel rods held at alternating positive or negative voltage
- Controls orientation before electron transfer
- Uses gas (N₂, argon) to stimulate collisions
- Design criteria:
 - Simple and **reliable** design
 - High ion transmission - **Improved**
 - Broad mass range transmission - **Improved**
 - Short ion transmission time - **Improved**
 - Exiting beam conditions indistinguishable between MS and MS/MS
 - **Compress ion beam for improved Q2 acceptance**








Collision Cell Results

Significantly shorter (60 mm less) collision cell affords reduced bench space


Compresses ion beam to improve transmission into smaller diameter quads

Operates at multifrequency for broad range transmission

Nearly equivalent MS/MS spectra compared to 6420/60



Benefits of Vortex Collision Cell



- Twisted, tapered hexapole
 - Taper – maximum ion collection from Q1 and collimates ion beam for transmission to Q2
 - Twist – optimal pressure and fragmentation of ions at front, higher pressure for cooling and compression of ion in the middle
 - Crosstalk is prevent by applying a DC-bias to resistive rods (more kinetic energy) to prevent stalling
 - Openings at end enable a fast pressure drop, while ions are confined (prevents scattering before entry into Q2)
- ***Optimum Ion Transmission and better MS/MS performance***



Agilent Technologies




Virtual Pre/Post Filters

20x thinner than traditional pre/post filters


Reversed polarity U+/U- DC only, no RF required.

Noding is eliminated


Consistent tuning with improved instrument performance



10%



Benefits of Virtual Pre/Post-Filters



- Pre-filters are necessary to reduce fringe fields
- Virtual filters cancels the DC component in the fringe field region by applying an opposite sign DC voltage to the filter
- ***Maximum ion transmission at entrance and exit of resolving quadrupoles***
- ***Quicker switching → Faster MRMs***

Redesigned Hyperbolic Quadrupole




A quadrupole's performance are defined by the following:

1. Electrode shape
2. Radius
3. Length
4. Frequency
5. Voltage



Agilent's Next Gen hyperbolic Quadrupoles operate at a higher frequency and voltage which allows a smaller sized quadrupole to perform with exceeding performance.

Easy Change Detector Assembly



Summary from Initial Ultivo TMS Synthetic Fentanyl Analyses.

- Quick 6 minute MRM method has been developed;
- Minimal sample preparation techniques, 1/10 dilution and filtering for Urine Matrix & protein precipitation with 1/6 dilution for Serum Matrix;
- Sensitivity measured as 100pg/ml or less for every analyte in both matrices (5uL injection volume);
- Precision data less than 10%RSD at the LLOQ for each analyte;
- Completely new and innovative miniature Tandem Mass Spectrometer performance demonstrated herein.

Ultivo Triple Quad LC/MS

AGILENT IS RESHAPING THE
FUTURE OF
MASS SPECTROMETRY.

Questions?



JOIN US AND WELCOME IN A NEW
ERA OF "FIT FOR PURPOSE" MASS
SPECTROMETRY.