

Instrument: Pegasus® BT and ChromaTOF® Sync**Characterization and Comparison of Whiskey* Aroma Profiles with GC-TOFMS and ChromaTOF Sync**

LECO Corporation; Saint Joseph, Michigan USA

Key Words: Whiskey, GC, TOFMS, ChromaTOF Sync, Alignment, PCA, Aroma Profile, Deconvolution

Introduction

Characterizing and comparing the aroma profiles of food and beverage products, like whiskeys, can be very interesting for better understanding the products and the specific chemicals that may contribute to the sensory observations of the products. Gas chromatography with mass spectrometry (GC-MS) is an excellent analytical technique for probing the aroma profile of these types of samples. The aroma contributing analytes tend to be volatile and semi-volatile and are well-suited to GC analysis. Chromatography effectively separates individual analytes in these complex samples, and MS detection then provides spectral information and good tentative identifications. With full m/z range acquisition and sensitive detection, LECO's Pegasus BT Time-of-Flight (TOF) MS produces rich data that describes these complex samples. Analytical software tools, like ChromaTOF Sync, are powerful for then probing the data and uncovering similarities, differences, and trends between the samples. ChromaTOF Sync performs sample set peak finding and incorporates deconvolution to produce a combined peak table that compiles and aligns analyte information across the entire sample set. Combining sample set information facilitates data review and comparing features through the set of samples. ChromaTOF Sync also includes additional comparative tools, like Principal Component Analysis (PCA) for general characterization and exploring trends.

In this application note, six different whiskeys were analyzed with LECO's Pegasus BT GC-TOFMS. The associated data was compared with ChromaTOF Sync to uncover interesting analytes and trends in the whiskey samples. Several examples are highlighted and discussed here.




Figure 1. Six different whiskey samples were analyzed (in duplicate). Overlaid Total Ion Chromatograms (TIC) are shown.

*Please note that all samples are referred to as "whiskey" and "whiskeys," throughout this application note for simplicity and consistency, because the samples were sourced from three geographic origins—US, Scotland, and Ireland. It is important to acknowledge that in Scotland (and Canada), the product "Scotch," is referred to as "Whisky" or "Whiskies," whereas in the US and Ireland the products are referred to as "Whiskey," and "Whiskeys."



Figure 3. Many ethyl esters are observed in the samples and details about this compound class are compiled.

Reviewing this compound class provides interesting context about the whiskey samples in terms of their aroma attributes, and it can also provide additional context for the identifications of the other analytes in the samples. While this series is tentatively identified, the consistency across the series and the very good similarity scores makes it a good option for retention index calculations. Essentially, the ethyl ester series was used as a naturally occurring series to calculate retention index values for the sample overall. Thus, the observed retention times and known retention index values for this series were used with *ChromaTOF Sync* for calculating retention index and supporting identifications of other analytes in the samples.

For example, 2-methyl 1-butanol is shown in Figure 4. This analyte had a very good similarity to the library with a score of 897. Retention index further supports this identification with an observed RI value of 738 compared to the library RI value of 739. The trends across the sample set can also be observed in Figure 4 with highest levels of this analyte in whiskey 4. This analyte has an aroma type of ethereal with descriptors of ethereal, whiskey, fusel, alcoholic, fatty, greasy, winey, leathery, and cocoa.^[1] This identification and the relative trends can provide additional insight to these whiskey samples.




Figure 4. 2-methyl 1-butanol. Spectral information, overlaid chromatographic profiles, and relative trends (as a bar graph) are shown.

ChromaTOF Sync also incorporates deconvolution into the peak finding, allowing for additional separation and mathematical resolution in instances of chromatographic overlap. An example is shown in Figure 5. A TIC chromatogram is shown in the top left corner of this figure. Two peak markers are indicated, but it is not readily apparent that two peaks are eluting in this section of the chromatogram when viewing the TIC. Deconvolution effectively resolves these coeluting features from each other, and from other analytes and background, and provides pure spectra for each as well as pure chromatographic profiles for each by indicating unique masses per feature. In this case, benzaldehyde and 1-heptanol coelute and are resolved with deconvolution. Benzaldehyde has a similarity score of 917 while 1-heptanol has a similarity score of 838. Both identifications are supported with retention index with observed RI values of 973 and 975 and library RI values of 970 and 975 for benzaldehyde and 1-heptanol, respectively. Neither of these analytes are very clear in the TIC, but their good identifications and relative trends across the sample set are readily determined with deconvolution. Both of these vary between the whiskey types and have potentially interesting aroma contributions. The alcohol has a green aroma type with descriptors of musty, leafy, violet, herbal, green, sweet, woody, and peony, and benzaldehyde has a fruity aroma type with descriptors of strong, sharp, sweet, bitter, almond, and cherry.^[1]




Figure 5. Deconvolution example. Benzaldehyde and 1-heptanol coelute and are mathematically isolated with ChromaTOF Sync deconvolution.

Benzaldehyde has fruity aroma notes and interesting trends across the sample set. Fruity tends to be an important descriptor of many whiskeys and reviewing other analytes that have a fruity aroma type may also be interesting. Several analytes with "fruity" aroma types are tabulated in Figure 6. Identification metrics (similarity scores, retention index, etc.) and aroma notes can be observed in the table and trends across sample set can be observed in the heat map.

Name	Formula	Similarity	CAS	Quant mass	RI. calc	RI. lib	aroma	descriptor	Med RT	1	2	3	4	5	6	7	8	9	10	11	12
Propanoic acid, ethyl ester	C ₅ H ₁₀ O ₂	924	105-37-3	57.05	710	710	fruity	sweet fruity rum juicy fruit grape pineapple	389.1	5.69	1.77	17.19	0.86	8.14	1.90	0.09	2.93	36.00	47.78	30.76	11.21
Propanoic acid, 2-methyl-, ethyl ester	C ₆ H ₁₂ O ₂	911	97-62-1	43.08	758	756	fruity	sweet ethereal fruity alcoholic fusel rummy	444.6	4.97	6.92	19.76	4.43	19.12	4.40	16.98	12.01	36.15	59.34	35.00	15.19
Isobutyl acetate	C ₆ H ₁₂ O ₂	897	110-19-0	43.04	774	772	fruity	sweet fruity ethereal banana tropical	462.8	8.94	3.93	19.66	19.18	13.84	4.58	18.33	15.91	17.52	10.51	40.10	10.16
Butanoic acid, ethyl ester	C ₆ H ₁₂ O ₂	927	105-54-4	71.06	802	802	fruity	fruity juicy fruit pineapple cognac	495.5	4.52	2.35	14.60	17.63	13.33	3.25	12.11	8.75	37.94	24.41	37.56	14.10
Butanoic acid, 2-methyl-, ethyl ester	C ₇ H ₁₄ O ₂	914	7452-79-1	57.09	851	849	fruity	sharp sweet green apple fruity	554.9	4.35	9.25	11.05	14.28	10.28	3.63	19.02	0.76	35.05	19.18	12.70	18.94
Butanoic acid, 3-methyl-, ethyl ester	C ₇ H ₁₄ O ₂	887	108-64-5	88.07	855	853	fruity	fruity sweet apple pineapple tutti frutti	559.9	1.93	2.07	11.04	19.03	12.64	0.98	13.26	0.54	10.26	54.27	21.16	15.96
1-Butanol, 3-methyl-, acetate	C ₇ H ₁₄ O ₂	902	123-92-2	70.10	877	876	fruity	sweet fruity banana solvent	586.3	0.66	1.02	16.22	16.74	12.26	0.09	18.65	15.44	14.14	41.31	25.54	17.21
Pentanoic acid, ethyl ester	C ₇ H ₁₄ O ₂	907	539-82-2	85.09	900	900	fruity	sweet fruity apple pineapple green tropical	614.5	1.74	3.22	10.57	10.72	18.29	5.91	17.80	1.67	38.38	74.77	23.05	11.94
Ethyl tiglate	C ₇ H ₁₂ O ₂	950	5837-78-5	83.07	941	939	fruity	sweet fruity tutti frutti tropical berry floral caramel	659.1	8.84	8.99	12.79	13.97	10.72	6.92	8.82	4.12	35.00	47.74	40.37	12.61
Benzaldehyde	C ₇ H ₆ O	917	100-52-7	106.07	975	962	fruity	strong sharp sweet bitter almond cherry	695.1	8.97	3.25	16.53	10.28	18.91	5.43	10.66	15.84	40.39	72.52	44.39	19.52
Furan, 2-pentyl-	C ₇ H ₁₀ O	807	3777-69-3	81.06	995	993	fruity	fruity green earthy beany vegetable metallic	717.2	3.44	6.00	11.69	18.82	13.03	3.90	14.31	13.44	11.68	75.59	37.86	16.78
Hexanoic acid, ethyl ester	C ₈ H ₁₆ O ₂	925	123-66-0	88.07	999	999	fruity	sweet fruity pineapple waxy green banana	721.5	0.13	5.39	12.03	9.68	13.03	8.36	15.19	7.95	51.99	54.06	59.65	18.38
Acetic acid, hexyl ester	C ₈ H ₁₆ O ₂	875	142-92-7	43.05	1013	1011	fruity	fruit green apple banana sweet	735.9	3.44	0.38	0.00	0.00	13.51	9.36	17.01	2.58	74.75	74.47	35.63	12.01
2-Nonanone	C ₉ H ₁₈ O	845	821-55-6	58.06	1095	1092	fruity	fruity tutti frutti tropical berry floral herbal	817.1	9.68	1.71	17.75	11.24	19.94	3.70	18.39	10.54	37.58	18.96	75.47	13.59
Heptanoic acid, ethyl ester	C ₉ H ₁₈ O ₂	886	106-30-9	88.08	1098	1098	fruity	fruity pineapple cognac rum wine	820.2	2.65	9.09	13.62	16.35	17.66	0.16	13.88	2.02	30.05	50.12	49.23	12.21
Hexanoic acid, 2-methylpropyl ester	C ₉ H ₁₈ O ₂	877	105-79-3	99.11	1152	1149	fruity	fruity pineapple green apple skin green apple sour tropical peach	870.6	3.48	0.12	0.00	0.00	15.51	9.28	18.89	18.69	34.96	45.10	26.61	18.97
Butanedioic acid, diethyl ester	C ₈ H ₁₆ O ₄	887	123-25-1	101.05	1178	1181	fruity	mild fruity cooked apple ylang	894.6	1.07	5.96	15.97	15.21	12.30	8.18	19.15	8.56	14.91	63.86	71.27	16.20
Ethyl (E)-2-octenoate	C ₁₀ H ₁₈ O ₂	877	7367-82-0	125.14	1248	1249	fruity	fruity pear skin green waxy tropical plum skin fatty	956.3	4.16	3.25	0.00	0.00	12.37	0.98	11.21	13.81	37.49	76.62	10.87	17.20
Isopentyl hexanoate	C ₁₁ H ₂₀ O ₂	853	2198-61-0	43.08	1251	1250	fruity	fruity banana apple pineapple green	958.8	1.16	9.82	11.09	12.52	16.90	8.40	17.83	12.53	32.02	86.61	35.48	12.07
n-Caprylic acid isobutyl ester	C ₁₁ H ₂₀ O ₂	885	5461-06-3	57.08	1348	1348	fruity	fruity green oily floral	1038.6	7.53	6.56	16.66	19.78	11.85	7.55	17.05	16.82	45.79	79.85	33.79	18.65
Ethyl 9-decanoate	C ₁₁ H ₂₀ O ₂	852	67233-91-4	88.07	1388	1388	fruity	fruity fatty	1069.9	2.81	9.98	0.00	0.00	19.91	5.80	0.00	0.00	39.24	94.63	36.25	15.50
Octanoic acid, 3-methylbutyl ester	C ₁₁ H ₂₀ O ₂	903	2035-99-6	70.09	1449	1446	fruity	sweet oily fruity green soapy pineapple coconut	1116.1	4.08	6.50	11.07	17.79	14.35	0.70	22.89	17.04	10.01	107.28	32.28	14.40
Limonene	C ₁₀ H ₁₆	887	138-86-3	68.08	1040	1030	fruity	citrus herbal terpene camphor	762.8	3.24	8.48	17.66	11.42	17.79	1.44	14.08	5.87	16.84	63.76	35.96	11.28
Propyl octanoate	C ₁₁ H ₂₀ O ₂	852	624-13-5	145.15	1292	1290	fruity	coconut cacao gin	993.4	4.52	7.66	18.29	10.70	19.89	4.77	19.47	1.60	39.06	100.38	32.97	13.86

Figure 6. Analytes with "fruity" aroma types.

While it is interesting to look at specific compound classes (like the ethyl esters), individual analytes (like benzaldehyde), and analytes with similar aroma attributes (like those with fruity notes), *ChromaTOF Sync* also supports overall unsupervised sample comparisons such as Principal Component Analysis (PCA). Features in the compiled peak table can be used as variables to explore general trends in the samples. In this case, tentatively identified analytes (similarity > 800 and library RI within 30 units) were used for PCA. The scores and loadings plots are shown in Figure 7. The scores suggest that the whiskeys 1, 3, 5, and 6 were more similar to each other while whiskey 2 (red) and whiskey 4 (pink) were most distinct. It is potentially interesting to note that sample 2 is from a distillery in Ireland, sample 4 is from a distillery in the United States, and the other whiskeys are all from Scotland.




Figure 7. PCA Scores (left) and Loadings (right) for the 6 whiskeys.

The scores plot can indicate which samples are most similar or different from each other and the associated loadings can help provide insight to those differences. For example, whiskey 2 (red) has the lowest PC1 scores, and 4 analytes with the lowest loadings on PC1 are also indicated in Figure 7. These analytes are distinct to whiskey 2 and are listed in Figure 8. These analytes are particularly interesting as they all have aroma types of herbal with additional descriptors of camphor, eucalyptol, and other spicy notes.^[1] Relative to the other whiskeys, whiskey 2 also had some distinct sensory descriptors with more spicy and clove notes than the other whiskeys. It is likely that some of these distinct analytes are connected to those distinct sensory notes.

Name	Formula	Similarity	CAS	Quant mass	R.I. calc	R.I. lib	aroma	descriptor	Med RT	1	2	3	4	5	6	7	8	9	10	11	12	
cis-Dihydrocarvone	C ₁₀ H ₁₆ O	864	3792-53-8	95.08	1213	1195	herbal	herbal warm	926.4	0.00	0.00	5.95	2.11	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
Isocineole	C ₁₀ H ₁₆ O	826	470-67-7	111.11	1025	1016	herbal	cooling pine minty camphor terpene green	747.2	0.00	0.00	4.55	7.71	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
Bicyclo[3.1.0]hexan-3-one, 4-methy	C ₁₀ H ₁₆ O	861	471-15-8	110.13	1132	1114	herbal	herbal warm	852.2	0.00	0.00	2.31	8.11	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00
Eucalyptol	C ₁₀ H ₁₆ O	858	470-82-6	111.11	1046	1032	herbal	eucalyptus herbal camphor medicinal	768.0	0.00	0.00	3.05	4.82	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00	0.00

Figure 8. Four analytes distinct in whiskey 2.

Whiskey samples are very complex with many chemical components. *ChromaTOF Sync* helps to uncover specific analytes and facilitates exploring these interesting analyte trends in the data.

Conclusion

In this application note, LECO's *Pegasus BT* and *ChromaTOF Sync* were used to compare and characterize six different whiskey samples. Aroma profile information was determined, and the software tools helped to uncover interesting analytes and trends within the data.

References

^[1] Good Scents Company database (www.thegoodscentscompany.com)

ChromaTOF[®] | SYNC



LECO Corporation | 3000 Lakeview Avenue | St. Joseph, MI 49085 | Phone: 800-292-6141 | 269-985-5496
info@leco.com • www.leco.com | ISO-9001:2015 Certified | LECO is a registered trademark of LECO Corporation.
Pegasus, ChromaTOF are registered trademarks of LECO Corporation.