

High Throughput Peptide Mapping with the Vanquish UHPLC System and the Q Exactive HF Mass Spectrometer

Martin Samonig,¹ Kai Scheffler,² Remco Swart,¹ Jonathan Josephs³ and Stephane Houel³

¹Thermo Fisher Scientific, Germering, Germany; ²Thermo Fisher Scientific, Dreieich, Germany;

³Thermo Fisher Scientific, San Jose, USA.

Overview

Purpose: Fast peptide mapping of biopharmaceuticals.

Methods: Efficient chromatographic peptide separations with the Vanquish UHPLC system, in combination with fast and high resolution quadrupole-Orbitrap mass spectrometry.

Results: Information about correct sequence, glycosylation and post-translational or artificial modification of recombinant monoclonal antibodies within 5 minutes.

Introduction

Monoclonal antibodies, or mAbs, (Figures 1 and 2) are the major element in the fastest growing sector of biopharmaceuticals within the pharma industry. By 2016, eight of the top ten drugs will be therapeutic proteins. Their manufacture is accomplished in bacterial or eukaryotic expression systems, requiring extensive purification of the target product. During drug development and production, the quality of biotherapeutics needs to be closely monitored.

Various analytical methods have been used to study quality attributes such as structural integrity, aggregation, glycosylation pattern or amino acid degradation. Because of their high information content and versatility, characterization methods based on high-performance liquid chromatography and mass spectrometry are among the most powerful protein characterization techniques. Proteins can be enzymatically digested to obtain peptides enabling their analysis by means of peptide mapping experiments.

Here, we report a fast and sensitive approach that combines enzymatic digestion, fast chromatographic separation, high-resolution mass spectrometry, and rapid data processing to handle the large amount of samples in diverse biopharma workflows. In this study we have analyzed two commercially available drug products: rituximab (trade names MabThera and Rituxan®) and denosumab (trade names Prolia® and XGEVA®).

FIGURE 1. General structure of mAbs and their biological and physico-chemical characteristics.

Physicochemical Characteristics Biological Characteristics

N-terminal heterogeneity

Pyroglutamate formation
Other modifications

Amino acid modifications

Deamidation, oxidation,
glycosylation, isomerization

Fragmentation

Cleavage in hinge region

Oligosaccharides

Fucosylation, sialylation,
galactosylation,...

Disulfide bonds

Free thiols, disulfide shuffling, thioether

C-terminal heterogeneity

Lysine processing, Proline amidation





FIGURE 2. Nomenclature of carbohydrate structures commonly observed on antibodies.



Methods

Sample Preparation

The two drug products Rituximab and Denosumab were denatured for 30 min in 7 M Urea and 50 mM Tris HCL at pH 8.0. The samples were reduced with 5 mM DTT for 30 min at 37°C, alkylation was performed with 10 mM IAA for 30 min at room temperature and the reaction was quenched by addition of 10 mM DTT. Thermo Scientific™ Pierce™ Trypsin Protease (MS Grade) was added and digestion allowed to proceed over night at 37°C. Digests were stopped by addition of TFA to approximately pH 3.0.

Liquid Chromatography and Mass Spectrometry

A Thermo Scientific™ Vanquish™ UHPLC system with a 2.1 x 250 mm i.d. Thermo Scientific™ Acclaim™ RSLC 120, C18, 2.2 μm column and gradients of water and acetonitrile (ACN) with 0.1% formic acid (FA) each were used to separate the peptide mixtures. Five different separation times were applied and compared: 5, 8, 13, 20, and 30 min for the gradient ramping from 4% to 55% eluent B (0.1% FA in 8:2 acetonitrile/water (v/v)). Flow rates were adapted accordingly using 1.1 (5 min), 1.0 (8 min), 0.6 (13 min), 0.4 (20 min), and 0.4 mL/min (30 min). The Thermo Scientific™ Q Exactive™ HF mass spectrometer (MS) equipped with a HESI-II probe was used for mass spectrometric detection using a Full MS / dd-MS2 (Top5) experiment.

Data Analysis

The data were acquired with Thermo Scientific™ Xcalibur™ 3.0 software in combination with Thermo Scientific™ SII for Xcalibur 1.1 software. Data analysis was performed using Thermo Scientific™ PepFinder™ 2.0 and FreeStyle™ 1.0 software packages.

Results

Peptide mapping experiments were performed using the rituximab digest for assessing the sequence coverage for light and heavy chain, as well as for identification and (relative) quantification of a specific set of modifications: a) oxidation, b) glycosylation and c) deamidation.^{1,2} For all five gradient times from 30 min down to 5 min, a very good separation was achieved (Figure 3) and resulting sequence coverages of 100% were obtained from all separation times both for light and heavy chain, even for the very short gradient of 5 min. The sequence coverage map (Figure 4) shows the overlap of the different peptides identified in different intensities and in different lengths due to missed cleavages.

FIGURE 3. Total ion chromatograms obtained from peptide mapping experiments of Rituximab applying gradient lengths of 30, 20, 13, 8 and 5 min. Flow rates and resulting pressures are indicated in the individual traces.




FIGURE 4. Sequence coverage map of rituximab, obtained using the 5 min gradient for peptide separation. The colored bars show the identified peptides, with the numbers in the bars reflecting the retention time. The different colors indicate the intensity of the peptide in the MS1 scan: red = high abundant $>1.6e^{+007}$, yellow $>2.0e^{+006}$, green $>2.3e^{+005}$, light blue $>2.8e^{+004}$, cyan = low abundant $>3.3e^{+003}$.





TABLE 1. Comparison of the oxidation, deamidation and glycosylation modifications identified in runs obtained from the different gradient times.

Modification	Gradient Time					σ	
	30 min	20 min	13 min	8 min	5 min		
Light Chain	Q1+NH3 loss	91.95%	91.17%	89.69%	90.93%	26.57%	28.80%
Heavy Chain	~Q1+NH3 loss	99.62%	99.67%	99.61%	99.68%	99.69%	0.04%
Heavy Chain	N3+Deamidation	0.52%	0.51%	0.58%	-	0.51%	0.03%
Heavy Chain	M34+Oxidation	1.64%	1.54%	1.73%	1.42%	1.45%	0.13%
Heavy Chain	N301+A1G0F	4.32%	4.42%	3.83%	3.52%	3.38%	0.46%
Heavy Chain	N301+A1G1F	1.87%	1.91%	1.72%	3.32%	1.46%	0.73%
Heavy Chain	N301+A2G0	1.09%	1.02%	1.02%	-	0.98%	0.05%
Heavy Chain	N301+A2G0F	37.88%	37.11%	38.59%	40.48%	43.12%	2.41%
Heavy Chain	N301+A2G1F	42.06%	41.89%	43.42%	43.20%	43.35%	0.75%
Heavy Chain	N301+A2G2F	10.23%	10.17%	9.81%	10.36%	10.05%	0.21%
Heavy Chain	N301+A2S1G0F	0.83%	0.86%	-	-	-	0.02%
Heavy Chain	N301+A2S1G1F	2.14%	-	-	-	-	-
Heavy Chain	N301+A3Sg1G0	1.30%	-	-	-	-	-
Heavy Chain	N301+M5	1.61%	1.59%	1.66%	1.87%	1.86%	0.14%
Heavy Chain	N301+Unglycos.	0.54%	0.90%	0.76%	0.83%	0.97%	0.16%
Heavy Chain	G450+Lys	3.57%	3.56%	3.92%	3.40%	3.15%	0.28%
median						0.19%	

Table 1 shows the identification and comparison of a subset of monitored modifications across the different separation times applied. A tilde (~) before the modification indicates the modification was found on the tryptic peptide, but could not be localized on a specific amino acid with MS/MS spectra. The modification is labeled with recovery "Good" when the total peak area, including modified and unmodified forms of the peptide, is at least 10% of the most abundant peptide from the same protein. The recovery "Fair" means it is at least 1%. The relative abundance of the detected modifications in the five different methods has a standard deviation of 0.19% and shows that important information regarding post-translational modifications (PTMs) can be obtained equally and accurately at all separation times.


Since the quantification of modified peptides performed in PepFinder 2.0 software requires their identification based on MS/MS spectra, special care has to be taken in choosing the appropriate ion injection times in the method setup enabling the acquisition of high quality MS/MS spectra required for a positive identification. This is especially true for glycopeptides analyzed using HCD fragment ion spectra, which contain exclusively ions representing sequential loss of glycan residues and no fragment ions representing the decomposition of the peptide as shown in Figure 5. The identification of peptides and modified peptides using PepFinder 2.0 software is based on the comparison of a simulated and the measured spectrum. The strength of the implemented algorithm for spectra and especially fragment ion intensity prediction is displayed in Figure 5 showing the MS/MS spectrum of the low abundant glycosylated peptide TKPREEQFN*STFR (*=A2G0) identified in the 5 min run with the typical fragmentation pattern: the two oxonium ions 204 (HexNAC), 366 (Hex-HexNAC), and the sequence ladder of the fragmented glycan attached to the intact peptide. The precursor ion of this glycopeptide with 735.58 m/z and a +4 charge state is indicated with the dotted blue line in the measured MS/MS spectrum.

FIGURE 5. Simulated (top) and measured (bottom) MS/MS spectra of the glycopeptide aa 290-302 (TKPREEQFN*STFR, *=A2G0) in the 5 min gradient run.



Even with ultra short gradients down to 5 min, as shown in Figure 6A, spectacular separation efficiency and peak widths of less than 1 s were obtained. Figure 6B is highlighting the number of scans obtained across one chromatographic peak. Typically 6 Full MS spectra and 25 (5xTop 5) MS/MS spectra were acquired. The achieved scan speed is key to the success in obtaining full sequence coverage.

Denosumab and B) data point distribution for a Full MS / ddMS² Top5 method during a representative chromatographic peak.



Conclusion

- The applied hardware setup chosen for the experiments consisting of the Vanquish UHPLC system with an Acclaim RSLC 120 2.1 x 250 mm column installed with Thermo Scientific™ Viper™ Fingertight Fitting connections, attached online to the Q Exactive HF mass spectrometer equipped with the HESI-II source, provides a very robust system allowing for very high reproducibility and long term stability.
- Flow rates between 0.4 and 1.1 mL/min, depending on the chosen gradient lengths between 5 and 30 min, delivered chromatograms with peak widths at half maximum of less than 1 second.
- The accelerated scan speed of the Q Exactive HF mass spectrometer delivered sufficient data points over a chromatographic peak and clearly illustrates the compatibility with fast UHPLC separations.
- For all gradient times ranging from 5 to 30 min, 100% sequence coverage for the light and heavy chains for both rituximab and denosumab was obtained.
- The analysis of the most commonly targeted modifications such as the presence/absence of a C-terminal Lys, the N-glycosylation of asparagine on the heavy chain, the N-terminal pyro-glutamine formation on heavy and light chains, oxidations and deamidations, were successfully identified and relatively quantified. However, the obtained results suggest using the slightly longer gradient times of 20 to 30 min for in-depth analyses to also capture the very low abundant modifications.
- The data presented in this study clearly demonstrate the capability of the applied LC-MS setup to significantly speed up peptide mapping experiments enabling high throughput analyses as required e.g. during clone selection in the development phase of biopharmaceuticals.

www.thermoscientific.com

©2015 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. MabThera and Rituxan are registered trademarks of F.Hoffmann-LaRoche Ltd. Prolixa and Xgeva are registered trademarks of Amgen Inc. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
 400 650 5118

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 10 3292 200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 45 453 9100
Korea +82 2 3420 8600
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00

Russia/CIS +43 1 333 50 34 0
Singapore +65 6289 1190
Spain +34 914 845 965
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
UK +44 1442 233555
USA +1 800 532 4752

PN64413-EN 0615S

Abbreviations

ACN – acetonitrile, DTT - dithiothreitol, FA – formic acid, IAA –iodoacetamide, mAb – monoclonal Antibody, PTMs – post translational modifications, PWHM – peak with at half maximum, TFA – trifluoroacetic acid.

References

1. Z. Zhang, Large-scale identification and quantification of covalent modifications in therapeutic proteins. *Anal. Chem.* (2009) 81, 8354-64.
2. Shah, X. G. Jiang, L. Chen, Z. Zhang, LC-MS/MS peptide mapping with automated data processing for routine profiling of N-glycans in immunoglobulins. *J. Am. Soc. Mass. Spectrom.* (2014) 25, 999-1011.

Thermo
 SCIENTIFIC

A Thermo Fisher Scientific Brand