

Influence of Ion Storage Location and Multiple Fill/Laser Pulse Sequence on Top-down UVPD Using an Orbitrap Fusion Mass Spectrometer

Chad R. Weisbrod,¹ Dustin D. Holden,² John E. P. Syka,¹ Jean-Jacques Dunyach,¹ Nick Izgarian,¹ Eugene Zhuk,¹ Jae C. Schwartz¹

¹Thermo Fisher Scientific, San Jose, CA; ²Department of Chemistry & Biochemistry, The University of Texas Austin, TX

Overview

Purpose: Evaluation of performance characteristics of UVPD on the Thermo Scientific™ Orbitrap Fusion™ Tribrid™ Mass Spectrometer.

Methods: Top-down protein characterization with UVPD capable Orbitrap Fusion MS

Results: Multiple laser pulses and multiple ion fills coupled with multiple laser pulses results in higher top-down sequence coverage using UVPD.

Introduction

Ultraviolet photo-dissociation (UVPD) has been demonstrated as an effective means for generating sequence informative polypeptide product ions. Improvement in sequence coverage and decreased charge state dependence on sequence coverage and product ion yield are chief reasons for choosing UVPD^[1] over any other single dissociation technique available for intact protein interrogation. Here we utilize the flexible configuration of the Orbitrap Fusion to enable UVPD reactions within any of its three ion storage regions, the low and high pressure cells (LPC and HPC respectively) of the quadrupole linear ion trap analyzer and the collision cell (HCD cell or Ion Routing Multipole, IRM). Multiple ion packet fill/laser pulse sequences are evaluated for improvement in S/N and protein sequence coverage.

Methods

Sample Preparation

Apomyoglobin was used as received (Sigma Aldrich). Myoglobin was diluted in ultrapure water to a stock concentration of 100 μ M. A final dilution step brought the working solution to a concentration of 0.5 μ M in spray solution (H₂O:Methanol:Acetic Acid; 49:49:2).


Mass Spectrometry

A 193-nm ArF excimer laser (Coherent ExciStar XS) was coupled to an Orbitrap Fusion™ Mass Spectrometer to enable UVPD. Laser optics utilized for coupling instruments resembles the configuration described by Brodbelt *et al.* A blank flange behind the (LPC) was replaced with a UV transparent window flange for laser access. This laser has an energy feedback circuit to enable fixed energy laser pulses within a certain range (0.4-60 mJ/pulse). Pulse width is 5 ns (FWHM). The site of irradiation has been explored among the three ion storage devices bore sighted to the beam. A variety of combinations of laser pulse sequences and (multiple) ion injection events were explored for performance optimization.

Data Analysis

All .raw files were analyzed using Hardklor^[2] for peak deconvolution. Database searches were carried out using ProSight Lite^[3]. Some of the data presented was processed using embedded Lua directly on the Orbitrap Fusion MS instrument computer.

FIGURE 1. UVPD interface with Orbitrap Fusion MS. UVPD can be conducted in all three locations: HPC, LPC, and IRM (HCD).



Results

UVPD reaction location

The site of UVPD was explored to discern whether any performance advantage is conferred via UVPD location. Sequence coverage and number of matching fragment ions was surveyed in the LPC and HPC. Laser beam alignment proved critical for optimized results in each location.

Multiple fill/ laser pulse sequence

Various instrument operation modes using combinations of ion injection events combined with laser pulse sequences have been explored. Sequence coverage was ultimately utilized to determine which configuration offered a benefit to top-down analysis. In addition, preliminary results illustrating proof of principal concept behind RF based laser irradiation protection are shown.

FIGURE 2. Diagram representing the experimental sequence of events regarding ion injection events and laser pulse events: **A) Single ion injection/ single laser pulse experiment, B) Single ion injection followed by N laser pulses, C) N ion injection events (single fill target/N) which equate to a single ion injection event followed by N laser pulses.**




FIGURE 3. Estimated energy deposition as a function of the number of pulses and energy per pulse. The red dashed line represents the threshold for damage for UV-grade fused silica components along the optical path. The yellow area represents the regime which was used to collect all data presented here. The energy was sufficient to induce widespread fragmentation in all intact proteins studied.




Figure 4. UVPD fragmentation pattern produced through irradiation in the HPC with 5 laser pulses of 0.6 mJ each. Inset spectra reveal the product ion peak density produced through this process. Sequence coverage map shows extensive coverage throughout the protein sequence.




Figure 5. (Left) Apomyoglobin [M+19H]¹⁹⁺ was used to study sequence coverage as a function of pulse energy in both the LPC and the HPC. (Right) Number of matching fragment ions are also reported as a function laser pulse energy. These data are preliminary results in an effort to characterize the best location for UVPD on the Orbitrap Fusion MS platform. Fine adjustment to laser alignment was performed for each device without breaking vacuum through monitoring the depletion of precursor.




Figure 6. Product ion coverage maps from two selected experimental data points from Figure 5 (0.2 and 1.0 mJ). The HPC exhibits a much slower roll-off in sequence coverage as a function of laser energy than the LPC. Also, sequence coverage overall is relatively stable across a wide energy range in the HPC. However, smaller fragment ions are observed as the energy is increased (both HPC and LPC) as illustrated by the absence of product ion coverage throughout the middle of the protein sequence at 1.0 mJ. This is as expected if multiple product ion generations are produced during a single laser pulse.




Figure 7. A) Laser energy was varied in 0.1 mJ/pulse steps and the precursor intensity monitored. The insets illustrate the extent of fragmentation achieved at each energy level. As expected for a first-order process, depletion of precursor is an exponential decay with energy delivery. B) Semi-log plot with exponential fit indicating linear increase in total energy with increasing energy.





FIGURE 8. Apomyoglobin sequence coverage and number of matching fragment ions were monitored as a function of laser pulse energy and pulse/fill sequence (see Figure 2)



www.thermoscientific.com

©2015 Thermo Fisher Scientific Inc. All rights reserved. ISO is a trademark of the International Standards Organization. All other trademarks are the property of Thermo Fisher Scientific and its subsidiaries. This information is presented as an example of the capabilities of Thermo Fisher Scientific products. It is not intended to encourage use of these products in any manner that might infringe the intellectual property rights of others. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.




Africa +43 1 333 50 34 0
Australia +61 3 9757 4300
Austria +43 810 282 206
Belgium +32 53 73 42 41
Canada +1 800 530 8447
China 800 810 5118 (free call domestic)
 400 650 5118
 PNG4484-EN 0615S

Denmark +45 70 23 62 60
Europe-Other +43 1 333 50 34 0
Finland +358 10 3292 200
France +33 1 60 92 48 00
Germany +49 6103 408 1014
India +91 22 6742 9494
Italy +39 02 950 591

Japan +81 45 453 9100
Korea +82 2 3420 8600
Latin America +1 561 688 8700
Middle East +43 1 333 50 34 0
Netherlands +31 76 579 55 55
New Zealand +64 9 980 6700
Norway +46 8 556 468 00

Russia/CIS +43 1 333 50 34 0
Singapore +65 6289 1190
Spain +34 914 845 965
Sweden +46 8 556 468 00
Switzerland +41 61 716 77 00
UK +44 1442 233555
USA +1 800 532 4752

FIGURE 9. Experiment indicates proof of principal for protection of ions from irradiation of the laser. This will enable a direct mechanism for over-fragmentation control in UVPD experiments. A single frequency resonant dipolar excitation was applied varying the amplitude to bring the precursor ion out of the path of the laser. As the excitation amplitude is increased the onset of CID is observed. A) UVPD specific and CID specific product ions are monitored from apomyoglobin [M+24H]²⁺.



Conclusion

- UVPD capable Orbitrap Fusion MS system can efficiently fragment polypeptide ions to produce informative product ions
- UVPD experiments were successfully conducted in the HPC and LPC
- Proper control over laser pulse energy will be required to ensure reproducibility and robust product ion generation
- Protection of product ions generated early during multiple laser pulse experiments will be required for control of over-fragmentation

References

- Shaw, JB et al., *J. Am. Chem. Soc.*, **2013**, 135 (34), pp 12646–12651
- Hoopmann, MR et al., *Anal. Chem.*, **79**, 5630-5632 (2007)
- Fellers RT et al., *Proteomics*, **2015**, 15(7), pp 1235-1238
- Holden DD et al., *Poster presented at 62nd ASMS Conference*, Baltimore, MD (2014)

Acknowledgements

We would like to thank the Jennifer Brodbelt and her laboratory for their assistance with initial setup and helpful discussion.

Thermo
SCIENTIFIC

A Thermo Fisher Scientific Brand