


# Optimizing GCxGC Parameters for Petroleum Analysis Using a Free Web-Based Tool

Christina N. Kelly, Joseph E. Binkley, Lorne M. Fell | LECO Corporation, Saint Joseph, MI, USA

## Introduction

Comprehensive two-dimensional gas chromatography (GCxGC) combined with time-of-flight mass spectrometry (TOFMS) is a well-established technique for characterization of petroleum samples because of the ability to separate structural classes of hydrocarbons into distinct bands of analytes that can then be identified with library-matched mass spectra. However, transferring traditional one-dimensional gas chromatography methods to optimized two-dimensional methods can be daunting because of the sheer number of options available for setting up commercial GCxGC systems—choosing orthogonal column phases, secondary oven temperature offsets, modulator temperatures, and second dimension separation period lengths. Using a new online tool, Simply GCxGC™, the process of translating an existing analysis into an optimized GCxGC method is experimentally determined for a sample of diesel fuel. Calculations that maximize use of two-dimensional space and minimize analyte wrap-around are applied to user-controlled parameters to make the most of peak-capacity increases and further enhance the chromatographic resolution gains from the additional dimension of separation.

## Instrumentation



### Pegasus® BT 4D High-Performance TOFMS

- Dual-Stage Thermal Modulator
- Acquisition rates up to 500 spectra/s
- Mass precision within 0.1 Da
- Benchtop size

## Simply GCxGC

Guide Me

- Convert 1D to GCxGC
- Stationary Phase Combinations
- Evaluate Sample Loading
- Determine Secondary Oven Offset
- Evaluate Stationary Phase
- Evaluate Peak Capacity
- Increase Peak Capacity
- Increase Secondary Column Peak Capacity
- Decrease Run Time
- Units / Settings
- Glossary
- Contact Us

### Introduction to Guide Me

The purpose of the Guide Me mode is to provide logical, step-by-step instructions that guide you through the process of developing an optimized GCxGC method.

**More Info**

The Guide Me mode has two components:

- 1) A computational engine which calculates column dimensions and operating parameters for optimum peak capacity.
- 2) A logical experimental scheme for evaluating and modifying the separation.

A GCxGC separation depends on column dimensions and operating parameters (peak capacity) which can be calculated; however, the actual distribution of peaks in the two-dimensional separation plane depends largely on the stationary phase chemistry and the particular sample components. Therefore, the development scheme guides you through the experimentation necessary to adjust the retention in the second dimension, optimize the stationary phases, and modify (increase, decrease, redistribute) the peak capacity as required by the goals of your application.

Previous Next  
Step 1 of 19

### Summary

Detailed

**Overview**

Peak Capacity Gain: 6.2  
 Run Time: 26 min  
 Acquisition Rate: 298 sps

**Columns**

Primary Column  
 Stationary Phase: Rxi-5SiI MS  
 30 m x 0.25 mm x 0.25 µm

Secondary Column  
 Stationary Phase: Rxi-17SiI MS  
 0.30 m x 0.25 mm x 0.25 µm

Transfer Line  
 0.21 m x 0.25 mm x 0.25 µm

Void Time  
 Secondary Column: 0.21 sec

**Inlet**

Carrier Gas: Helium  
 Flow: 1.4 mL/min  
 Inlet Pressure at 300 °C: 28 psi  
 Primary Column Velocity: 40 cm/sec

**Temperature**





Oven Temperature Program  
 Initial Temperature 40 °C for 0.20 min  
 10.2 °C/min to 300 °C  
 Secondary Oven Offset: 5 °C  
 Transfer Line Temperature: 340 °C

**GCxGC**

Optimal Modulation Period: 1.1 sec  
 Method Modulation Period: 1.1 sec  
 Method Hot Pulse Time: 0.34 sec


#	Run Time	PC 1	PC 2	PC Net
1	26 min	305	13	2566

## Workflow




## Normal Phase

Rxi-5MS 30 m x 0.25 mm x 0.25 µm coupled to Rxi-17SiI MS 2 m x 0.25 mm x 0.25 µm



Optimal GCxGC parameters of 0.3 m in secondary oven, 10.3 °C/min ramp rate from 50 °C to 300 °C, carrier gas flow rate of 1.4 mL/min, +16 °C secondary oven offset, 2<sup>nd</sup> dimension separation time of 1.1 s with 0.33 s hot pulse, and spectral acquisition rate of 309 spectra/s. Total peak capacity calculated at 2523.


Slow Down Ramp Rate → More Peak Capacity



Optimal GCxGC parameters of 0.3 m in secondary oven, 4.2 °C/min ramp rate from 50 °C to 300 °C, carrier gas flow rate of 1.4 mL/min, +9 °C secondary oven offset, 2<sup>nd</sup> dimension separation time of 2.4 s with 0.73 s hot pulse, and spectral acquisition rate of 192 spectra/s. Total peak capacity calculated at 3858.


Change 2D Column ID → More 2D Peak Capacity

Rxi-5MS 30 m x 0.25 mm x 0.25 µm coupled to Rxi-17SiI MS 2 m x 0.18 mm x 0.18 µm



Optimal GCxGC parameters of 0.3 m in secondary oven, 4.7 °C/min ramp rate from 50 °C to 300 °C, carrier gas flow rate of 1.4 mL/min, +10 °C secondary oven offset, 2<sup>nd</sup> dimension separation time of 2.2 s with 0.60 s hot pulse, and spectral acquisition rate of 234 spectra/s. Total peak capacity calculated at 4145.


Slow Down Ramp Rate → More Peak Capacity  
Add Loop in 2D Oven



Optimal GCxGC parameters of 0.45 m in secondary oven, 2.7 °C/min ramp rate from 50 °C to 300 °C, carrier gas flow rate of 1.4 mL/min, +15 °C secondary oven offset, 2<sup>nd</sup> dimension separation time of 3.6 s with 1.10 s hot pulse, and spectral acquisition rate of 178 spectra/s. Total peak capacity calculated at 5620.


## Reverse Phase

Rxi-17SiI MS 30 m x 0.25 mm x 0.25 µm coupled to Rxi-5MS 2 m x 0.25 mm x 0.25 µm



Optimal GCxGC parameters of 0.3 m in secondary oven, 10.3 °C/min ramp rate from 50 °C to 300 °C, carrier gas flow rate of 1.4 mL/min, +32 °C secondary oven offset, 2<sup>nd</sup> dimension separation time of 1.1 s with 0.33 s hot pulse, and spectral acquisition rate of 309 spectra/s. Total peak capacity calculated at 2523.


Slow Down Ramp Rate → More Peak Capacity



Optimal GCxGC parameters of 0.3 m in secondary oven, 4.1 °C/min ramp rate from 50 °C to 300 °C, carrier gas flow rate of 1.4 mL/min, +9 °C secondary oven offset, 2<sup>nd</sup> dimension separation time of 2.4 s with 0.73 s hot pulse, and spectral acquisition rate of 193 spectra/s. Total peak capacity calculated at 3903.


Change 2D Column ID → More 2D Peak Capacity  
Add Loop in 2D Oven

Rxi-17SiI MS 30 m x 0.25 mm x 0.25 µm coupled to Rxi-5MS 2 m x 0.18 mm x 0.18 µm



Optimal GCxGC parameters of 0.45 m in secondary oven, 9.6 °C/min ramp rate from 50 °C to 300 °C, carrier gas flow rate of 1.4 mL/min, +9 °C secondary oven offset, 2<sup>nd</sup> dimension separation time of 2.2 s with 0.30 s hot pulse, and spectral acquisition rate of 349 spectra/s. Total peak capacity calculated at 3055.

Slow Down Ramp Rate → More Peak Capacity



Optimal GCxGC parameters of 0.45 m in secondary oven, 4.7 °C/min ramp rate from 50 °C to 300 °C, carrier gas flow rate of 1.4 mL/min, +9 °C secondary oven offset, 2<sup>nd</sup> dimension separation time of 2.2 s with 0.65 s hot pulse, and spectral acquisition rate of 248 spectra/s. Total peak capacity calculated at 4485.

## Conclusions

Using the Simply GCxGC tool to set up experiments yields excellent two-dimensional separations. While initial parameters may not result in the highest peak capacity, Simply GCxGC provides a helpful guide for the next steps that lead to greater peak capacity results and better chromatographic resolution for any given sample, regardless of the initial column phase used.

