

Analysis of Tequila Carbohydrates

Application Note

Food and Beverage

Authors

Stephen Ball, Linda Lloyd
Agilent Technologies, Inc.

Introduction

Tequila is a spirit made from the blue agave plant, which is known to contain raffinose, a trisaccharide commonly found in a variety of vegetables and whole grains. The Agilent Hi-Plex Ca ligand-exchange chromatography column is able to separate raffinose from the other sugars and ethanol that make up this alcoholic beverage.



Agilent Technologies

Conditions

Column	Agilent Hi-Plex Ca, 7.7 × 300 mm, 8 μm (p/n PL1170-6810)
Sample	Sugars in tequila
Sample size	20 mg/mL
Mobile phase	100% DI H ₂ O
Flow rate	0.6 mL/min
Injection volume	20 μL
Temperature	85 °C
Detector	RI




Figure 1. Separation of sugars in tequila using an Agilent Hi-Plex Ca 8 μm column. See Table 1 for peak identification.

Table 1. Peak Identification for Figure 1

Peak	Name	Time (min)	Height (μV)	Area (%)	Width 50% (min)	As. USP	10% Asymmetry	Res. HW	Plate counts	Plates/m
1	Raffinose	8.27	338184.7	12.365	0.21	1.01	1.01	0.00	8972	29906
2	Sucrose, maltose	9.13	608399.4	29.454	0.29	0.98	1.06	2.08	5657	18858
3	Lactose	9.45	273004.4	14.285	0.30	3.67	3.19	0.65	5607	18690
4	Glucose	10.92	184788.9	14.740	0.44	1.04	1.04	2.35	3401	11337
5	Fructose	13.43	161682.3	14.199	0.49	1.01	1.00	3.19	4217	14055
6	Mannitol	17.24	219822.6	14.958	0.35	0.96	0.96	5.34	13125	43750
Total			1785882.3	100.000						

For More Information

These data represent typical results. For more information on our products and services, visit our Web site at www.agilent.com/chem.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice. This publication was originally published in 2008.

© Agilent Technologies, Inc., 2011
Published in USA, June 30, 2011
SI-01677



Agilent Technologies