

## Instrument: TGA801

### Determination of Moisture and Ash in Plant Tissue

LECO Corporation; Saint Joseph, Michigan USA

#### Introduction

The moisture content of plant tissue is often used to monitor the drying and curing process for commercial plant materials. Monitoring moisture values in plant tissue allows producers to ensure that the drying process meets the required moisture content criteria for a particular plant product. By monitoring moisture content, it is possible to prevent mold and fungal growth from high moisture levels, as was over-drying of plant materials. Both high moisture levels and over-drying lead to plant material deterioration and a reduction in product quality. The ash concentration in plant tissues can be used to evaluate the quality and effectiveness of cleaning and/or separation processing of the plant material. Additionally, the determination of a variety of analytically important constituents within Plant Tissue (Carbon, Nitrogen, Sulfur, etc.) require moisture correction utilizing an accurate moisture value.

Thermogravimetric analysis (TGA) is an analytical technique in which changes in sample mass, due to changes in physical and chemical properties of materials, is measured as a function of temperature and/or time. TGA is commonly used to determine selected characteristics of materials that exhibit either mass loss, or gain, due to decomposition, oxidation, or loss of volatile materials such as moisture.

The LECO TGA801 is a macro thermogravimetric analyzer designed to determine moisture and ash content of materials by measuring the change in mass of the sample as a function of the oven temperature while controlling the atmosphere and ventilation rate. The TGA801 allows up to 19 samples to be analyzed simultaneously.

#### Sample Preparation

Samples must be of a uniform consistency to produce suitable results.

#### Accessories

621-331 Ceramic Crucibles, 621-011-507 Double Sided Spoon.

**Sample Mass** ~1.0 g

**Analysis Time** ~4.5 h

#### General Method Parameters

|                         |         |
|-------------------------|---------|
| Crucible Type           | Ceramic |
| Minimum Crucible Weight | 19.0000 |
| Maximum Crucible Weight | 30.0000 |
| Crucible Density        | 3.0     |
| Lid Density             | 3.0     |
| Sample Type             | Leaf    |
| Sample Density          | 1.5     |
| Minimum Sample Weight   | 0.8000  |
| Maximum Sample Weight   | 1.2000  |

#### Method Step Parameters - Moisture

|                    |              |
|--------------------|--------------|
| Step Type          | Preset       |
| Preset Method Step | Moisture     |
| Cooling Option     | Active       |
| Crucible Lids      | No           |
| Start Temperature  | 25.0 °C      |
| End Temperature    | 80.0 °C      |
| Ramp Rate          | 6.0 °C/min   |
| Hold Time          | 120 min      |
| Maximum Time       | 240 min      |
| Atmosphere         | Air          |
| Flow Rate          | 10.0 LPM     |
| Final Weight       | At Constancy |
| Constancy Window   | 9 min        |
| Constancy Level    | 0.0005 g     |

#### Method Step Parameters - Ash

|                    |            |
|--------------------|------------|
| Step Type          | Preset     |
| Preset Method Step | Ash        |
| Cooling Option     | Active     |
| Crucible Lids      | No         |
| Start Temperature  | 80.0 °C    |
| End Temperature    | 500.0 °C   |
| Ramp Rate          | 6.0 °C/min |
| Hold Time          | 60 min     |
| Maximum Time       | 240 min    |
| Atmosphere         | Oxygen     |
| Flow Rate          | 10.0 LPM   |
| Final Weight       | Constancy  |
| Constancy Window   | 9 min      |
| Constancy Level    | 0.0005 g   |

### Method Step Calculations - Moisture

|                      |   |
|----------------------|---|
| Calculation Type     | Preset  |
| Preset Method Step   | Moisture  |
| Measurement Type     | Mass Ratio  |
| Enable Calibration   | Disabled  |
| Moisture Calculation | $((\text{Initial Mass} - \text{Moisture Mass}) \div \text{Initial Mass})$ |

### Method Step Calculations - Ash

|                      |  |
|----------------------|--|
| Calculation Type     | Preset                                       |
| Preset Method Step   | Ash  |
| Measurement Type     | Mass Ratio                                   |
| Enable Calibration   | Disabled                                     |
| Moisture Calculation | $(\text{Ash Mass} \div \text{Initial Mass})$ |

### Method Step Calculations - Ash Dry

|                      |   |
|----------------------|---|
| Calculation Type     | Preset                                      |
| Preset Method Step   | Ash Dry                                     |
| Measurement Type     | Mass Ratio                                  |
| Enable Calibration   | Disabled                                    |
| Moisture Calculation | $(\text{Ash} * ((1 \div \text{Moisture})))$ |

### Procedure

1. Create and/or select a method, using the Method Step Parameters listed above, following the procedure outlined in the LECO TGA801 Instruction Manual.
2. Login and load samples following the procedure outlined in the LECO TGA801 Instruction Manual.

### Typical Results

|                | Initial Mass (g) | % Moisture  | % Ash       | % Ash Dry    |
|----------------|------------------|-------------|-------------|--------------|
| Alfalfa        | 1.0347           | 6.45        | 9.95        | 10.63        |
| 502-273 LRM®   | 1.0039           | 6.46        | 9.99        | 10.68        |
| Lot: 1026      | 1.0209           | 6.46        | 9.98        | 10.66        |
|                | 1.0110           | 6.46        | 9.97        | 10.66        |
|                | 1.0167           | 6.46        | 10.01       | 10.70        |
| <b>Avg =</b>   |                  | <b>6.46</b> | <b>9.98</b> | <b>10.67</b> |
| <b>s =</b>     |                  | <b>0.01</b> | <b>0.02</b> | <b>0.03</b>  |
| Tobacco        | 1.0242           | 2.93        | 9.38        | 9.67         |
| 502-082 LRM    | 1.0090           | 2.93        | 9.43        | 9.72         |
| Lot: 1018      | 1.0109           | 2.91        | 9.47        | 9.76         |
|                | 1.0389           | 2.91        | 9.43        | 9.72         |
|                | 1.0108           | 2.93        | 9.43        | 9.71         |
| <b>Avg =</b>   |                  | <b>2.92</b> | <b>9.43</b> | <b>9.71</b>  |
| <b>s =</b>     |                  | <b>0.01</b> | <b>0.03</b> | <b>0.03</b>  |
| Orchard Leaves | 1.0038           | 4.13        | 5.67        | 5.91         |
| 502-931 LCRM®  | 1.0524           | 4.20        | 5.65        | 5.90         |
| Lot: 1000      | 1.0092           | 4.20        | 5.64        | 5.89         |
|                | 1.0040           | 4.14        | 5.68        | 5.92         |
|                | 1.0776           | 4.17        | 5.67        | 5.92         |
| <b>Avg =</b>   |                  | <b>4.17</b> | <b>5.66</b> | <b>5.91</b>  |
| <b>s =</b>     |                  | <b>0.03</b> | <b>0.02</b> | <b>0.01</b>  |

