

Sulfur in Biomass and Biofuel

LECO Corporation; Saint Joseph, Michigan USA

Instrument: TruSpec[®] S Module

Sample Preparation

A representative, uniform sample is required.

Accessories

528-203 Crucibles, 502-321 Com-Cat[™]

Calibration Samples

LECO Sulfur in Oil samples; 502-420, 502-419, 502-418, 502-417, or other suitable reference materials.

Method Parameters*

External Sulfur Analysis Parameters

Furnace Temperature	1350°C
Analysis Stabilize Comparator	2
Manual Load Baseline Delay Time (seconds)	5
Baseline Time (seconds)	1
Auto Detect Data Missed Time (seconds)	4
Endline Time seconds	5
Minimum Analysis Time (seconds)	60
Comparator Level	0.30 %
Conversion Factor	1.00
Significant Digits	5

System Configuration

Gas Conservation Timeout (minutes)	15
Auto Increment Sample Name	Disable
Nominal Mass	1.0000
Lance Delay Time (seconds)	20
Lance Limits	50000

*Refer to TruSpec/Sulfur Module Operator's Instruction Manual for Method Parameter definitions.

Procedure

1. Prepare instrument for operation as outlined in the operator's instruction manual.
2. Condition the system by analyzing a three to five ~0.25 gram coal samples.
3. Determine blank.
 - a. Enter 1.0000 g mass into Sample Login (F3) using Blank as the sample name.
 - b. Add ~1 gram of 502-321 Com-Cat into a 528-203 Crucible.
 - c. Initiate the analysis sequence (F5), when "Load Sample into Furnace" message appears on the display, click OK, and then slide crucible into the combustion tube until it reaches crucible stop. *Alternately, place crucible onto the appropriate position of carousel if equipped with autoloader system.*
 - d. Repeat steps 3a through 3c a minimum of three times.
 - e. Set the blank following the procedure outlined in the operator's instruction manual.



4. Calibrate/Drift Correct.
 - a. Weigh ~0.1 gram calibration sample into a 528-203 Crucible containing ~0.5 gram of 502-321 Com-Cat, enter mass and sample identification into Sample Login (F3).
 - b. Cover with ~0.5 gram 502-321 Com-Cat.
 - c. Initiate the analysis sequence (F5), when "Load Sample into Furnace" message appears on screen, click OK, then slide crucible into the combustion tube until it reaches crucible stop. *Alternately, place crucible onto the appropriate position of carousel if equipped with autoloader system.*
 - d. Repeat steps 4a through 4c a minimum of three times for each calibration/drift sample used.
 - e. Calibrate or Drift Correct the instrument following the procedure outlined in the operator's instruction manual.
5. Analyze Samples/Biofuels.
 - a. Weigh ~0.1 gram sample into 528-203 Crucible containing ~ 0.5gram of 502-321 Com-Cat, enter mass and sample identification into Sample Login (F3).
 - b. Cover with ~0.5 gram of 502-321 Com-Cat.
 - c. Initiate the analysis sequence (F5), when "Load Sample into Furnace" message appears on screen, click OK, and then slide crucible into the combustion tube until it reaches crucible stop. *Alternately, place crucible onto the appropriate position of carousel if equipped with autoloader system.*
6. Analyze Samples/Biomass.
 - a. Weigh ~0.15 gram sample into 528-203 Crucible and enter mass and sample identification into Login (F3).
 - b. Add ~1 gram of 502-321 Com-Cat to the crucible and thoroughly mix with sample.
 - c. Initiate the analysis sequence (F5). When "Load Sample into Furnace" message appears on screen, click OK and then slide crucible into the combustion tube until it reaches crucible stop. *Alternately, place crucible onto the appropriate position of carousel if equipped with autoloader system.*

LECO[®]

Organic Application Note

Typical Results

Sample	Mass	% S
Wood Chips #1	0.1545	0.138
	0.1555	0.134
	0.1551	0.153
	0.1548	0.137
	0.152	0.134
	X =	0.139
s =	0.008	
Wood Chips #2	0.1553	0.061
	0.1539	0.053
	0.1531	0.058
	0.1523	0.054
	0.1555	0.057
	X =	0.057
s =	0.003	
Biofuel	0.0985	0.01
	0.0971	0.009
	0.1017	0.012
	0.0989	0.01
	0.0995	0.011
	X =	0.010
s =	0.001	