

Screening for Xenobiotics in Postmortem Biological Samples using High Resolution Liquid Chromatography Mass Spectrometry

Aria McCall (1), Robert D. Johnson (1), Sarah R. Olive (2), Neil Loftus (3), Rachel Lieberman (2)

(1) Tarrant County Medical Examiner's Office, Fort Worth, TX (2) Shimadzu Scientific Instruments, Columbia, MD (3) Shimadzu MS/BU, Manchester, United Kingdom

1. Overview

- High resolution LCMS for more selective, more sensitive sample screening in forensic toxicology laboratories
- Method package complete with instrument parameters and optional spectral library containing >900 compounds
- Easy spectral library creation and expansion capabilities
- Faster data processing using Insight Explore software

2. Introduction

- The ability to screen biological samples quickly and accurately is of utmost importance in a Medical Examiner's laboratory.
- To keep up with the workflow, a sensitive and robust analytical technique combined with spectral library matching and user-friendly software is required
- With the constant discovery of novel drugs, forensic laboratories need the ability to retrospectively search through a data file to find a substance not yet added to their library

3. Methods

- Shimadzu Nexera LC40 UHPLC coupled with 9030 QToF Mass Spectrometer (Figure 1)
- Chromatographic separation using Shimadzu Forensic Toxicology Database® parameters
- The acquisition method design was set-up with one MS scan and several DIA-MS/MS events. HRAM data acquired using sequential mass scans with a cycle time less than 1 second (MS and subsequent DIA-MS/MS mass scans with an isolation width of m/z 20 Da). In each MS/MS event, a collision energy (CE) spread (5-55V) was applied to generate precursor MS and merged product ion MS/MS spectra.



Figure 1



4. Results/Workflow

4.1 Library Creation

- An in-house library containing 425 entries was created using certified reference material and the analytical conditions listed in Section 3. Each library entry included compound metadata information and retention time. (Figure 2)




Figure 2

4.2 Sample Extraction

- Drugs from Postmortem blood and tissue samples were extracted using a protein precipitation method, centrifuged, dried down, and reconstituted in mobile phase prior to injection.
- End results were compared with the laboratory's previous methodology.
- This dry down technique requires a much lower sample volume, a quicker sample extraction, and 5x lower injection volume.

4.3 Data Analysis and Reporting

- Data was analyzed using Shimadzu Insight Explore Screening software
- An analyte was deemed a match based on the following criteria:
 - Spectral match >60% similarity index score
 - Mass error < 10 ppm of the calculated mass
 - Retention time +/- 0.2 min from the target
- A 1-point calibrator was included as a pass/fail concentration threshold
- Insight Flagging criteria used to accelerate data review process by allowing review by exception (Figure 3)





Figure 3



5. Conclusion

- Faster, more-encompassing screening methodology for forensic casework
- Less sample volume, lower injection volume required
- Easily expandable spectral library for novel compounds
- Data-independent acquisition methodology provides retrospective data analysis capability

Acknowledgements

Shimadzu Scientific Instruments would like to thank Tarrant County Medical Examiner's Office for their collaborative efforts in making this poster.

For Research Use Only. Not for diagnostic purposes.