

Enhancement of Undergraduate Chemistry Laboratory Experiences Through a Collaboration of Industry and Academia

Alonso, D.E.¹, Ahlberg L.², Clear K.³, Muna G.³, and Binkley J.¹; 1) LECO Corporation, St. Joseph, MI, 2) Andrews University, Berrien Springs, MI, 3) Indiana University at South Bend, South Bend, IN



Background

- 1) Analytical instrumentation in undergraduate laboratories is critical for the scientific training of students at colleges and universities.
- 2) LECO Corporation is currently involved in various projects:
 - a) Cooperative education
 - b) Technology presentations
 - c) Mobile bus visits – Showcase technology to researchers
 - d) **Laboratory experiments** – “hands-on” instrument activities for students and faculty

Objective

To enhance education through exposure to modern analytical instrumentation (GC-TOFMS) and software

Collaboration Institutions





Figure 1: Project Interactive Strategy



Indiana University at South Bend

Public University, 3rd largest IU System Campus
 Location: South Bend, IN
 Enrollment :
 • Undergraduate 4,375
 • Graduate 567



Andrews University

Private Institution, Seventh-Day Adventists
 Location: Berrien Springs, MI
 Enrollment:
 • Undergraduate 1,434
 • Graduate 1,955



Figure 2: LECO Corporation Markets

LECO Corporation

Instrument Manufacturer
 Location: St. Joseph, MI
 Instruments:
 • Analytical Science
 • Metallographic Science
 • Separation Science

Project Tools



Figure 3: LECO Mobile Laboratories



Figure 4: Pegasus BT 4D

Mobile Lab Instruments

- GC-TOFMS and GCxGC-TOFMS
- Low fg sensitivity
- 10⁹ Linear dynamic range
- L-PAL3 Autosampler
- Liquid
- Headspace
- SPME

Results and Discussion

1) Organic Chemistry




Figure 5: Acid-catalyzed dehydration of 2-methylcyclohexanol

Experimental Concepts

- a) Dehydration Reactions
 - Mechanism of Elimination Reactions
 - Stability of Intermediates
 - Hydride/Alkyl Shifts
 - Thermodynamic Stability of Isomeric Alkenes
- b) Gas Chromatography – Mass Spectrometry
 - Total Ion Current (Separation science)
 - Molecular Ion (Structural information)
 - Fragment Ions (Carbocation/radical)




Figure 6: Simple Distillation Glassware for Dehydration Reaction




Figure 7: Product Mixture Total Ion Chromatogram




Figure 8: Mobile Lab Analysis Report

Results and Discussion

2) Instrumental Analysis: Quantitative Analysis of Acetaminophen and Caffeine




Figure 9: Sample Preparation – Extraction

Experimental Concepts

- a) Quantitative Analysis
 - Gas Chromatography
 - Chromatographic Parameters
 - Instrument Parameters
- b) Mass Spectrometry
 - Ionization Methods
 - Mass Analyzers
 - Detectors



Figure 10: Analytical Ion Chromatogram Showing Selected Excedrin Constituents




Figure 11: Acetaminophen and Caffeine Calibration Curves

Summary

- 1) The collaboration between academia and industry enhanced laboratory experiences through the utilization of modern-day instrumentation and software
- 2) Sample analyses were conducted using GC-TOFMS instruments on a mobile laboratory
- 3) New experiments will be developed for future mobile laboratory university projects




Figure 12: SPME, GC-TOFMS AIC for a Citrus Beverage