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Rob Lee
Principal Scientist
Robert_lee@waters.com

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Látky a metabolity
z potu

Cévy (z krve)

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
Figure taken from
Pragst et al., Clin Chim Acta 370 (2007)

Society of Hair Testing - doporučení

Forensic Science International xxx (2011) xxx-xxx

Contents lists available at SciVerse ScienceDirect

Forensic Science International

journal homepage: www.elsevier.com/locate/forensiint

Society of Hair Testing guidelines for drug testing in hair[☆]

Gail A.A. Cooper^{a,*}, Robert Kronstrand^{b,c}, Pascal Kintz^d

^aForensic Medicine and Science, University of Glasgow, G12 8QQ, Glasgow, Scotland
^bDepartment of Forensic Genetics and Forensic Toxicology, National Board of Forensic Medicine, SE-58758 Linköping, Sweden
^cDepartment of Medical and Health Sciences, Faculty of Health Sciences, Linköping University, SE-58185 Linköping, Sweden
^dX-Perfise Consulting, 84 route de Saverne, 67205 Oberhausbergen, France

ARTICLE INFO

Article history:
 Received 24 January 2011
 Accepted 24 January 2011
 Available online xxx

Keywords:
 Hair testing
 Guidelines

ABSTRACT

The Society of Hair Testing (SoHT) Guidelines for Drug Testing in Hair provide labor recommended best practice guidelines whether they are currently offering drug testing in hair or plan to offer a hair testing service in the future. The guidelines include reference to recommend collection and storage procedures, through sample preparation, pre-treatment and analysis of cut-offs.

Crown Copyright © 2011 Published by Elsevier Ireland Ltd. All rights reserved.

Table 1

Recommended cut-offs for substances and metabolites in hair to identify use.

Screening		Confirmation	
Group	Cut-off (ng/mg)	Target analyte	Cut-off (ng/mg)
Amphetamines	0.2	Amphetamine	0.2
		Methamphetamine	0.2
		MDA	0.2
		MDMA	0.2
Cannabinoids	0.1	THC	0.05
		THC-COOH	0.0002
Cocaine	0.5	Cocaine	0.5
		BZE, EME, CE, NC	0.05
Opiates	0.2	Morphine	0.2
		Codeine	0.2
		6-acetylmorphine	0.2
Methadone	0.2	Methadone	0.2
		EDDP	0.2
Buprenorphine	0.01	Buprenorphine	0.01
		Norbuprenorphine	0.01

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
[TECHNOLOGY BRIEF]

Waters
THE SCIENCE OF WHAT'S POSSIBLE™

Using UPLC-MS/MS for the Determination of Carboxy-THC (cTHC) in Hair for Forensic Toxicology

Massimo Gottardi¹, Francisco Ferron², Rob Lee², Andrea Gardumi¹, Simone Donzelli² and Michelle Wood³
¹Comedical, Trento, Italy; ²Waters Corporation, Milano, Italy; ³Waters Corporation, Wilmslow, UK

A robust and sensitive UPLC-MS/MS method for the determination of cTHC in hair.

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





Poděkování

- Francisco Ferron (Waters, Italy)
- Massimo Gottardi and Andrea Gardumi (Comedical, Italy)

Waters

THE SCIENCE OF WHAT'S POSSIBLE.™