



# Toxtyper 3.0

VITATOX 2020, 7. - 9. 9. 2020

Spolehlivý, rychlý a jednoduchý screening

**RADANAL, RECIPE, BRUKER**


**Ing. Daniel Vláčil**  
Country Sales Manager Bruker s.r.o.

# Toxtyper workflow for toxicology screening




Push-Button Sample Submission

UHPLC-MS<sup>n</sup> analysis



Result:  
List of identified compounds



Description	Retention Time [min]	Peak Number	Cmp Name	d RT [min]	m/z [Da]	d m/z [Da]	Purity	Intensity	ID	Comment
AutoMSn(468.32)	4.27	8	Buprenorphine	-0.14	468.32	-0.01	947	6.6E7	MS2	
AutoMSn(304.04)	3.60	7	Cocaine	-0.19	304.04	0.11	997	3.2E7	MS2/MS3	
AutoMSn(234.96)	3.16	5	Lidocaine	-0.12	234.96	0.22	998	2.7E7	MS2	
AutoMSn(289.94)	5.69	10	D5-Diazepam	-0.20	289.94	0.17	957	1.8E7	MS2	
AutoMSn(289.93)	3.25	6	Benzoylcegonine	-0.12	289.93	0.21	988	1.7E7	MS2/MS3	
AutoMSn(199.90)	0.57	1	Ecgoninemethylester	-0.07	199.90	0.23	988	1.0E7	MS2/MS3	MS2 unspecified
AutoMSn(150.00)	2.92	3	Methamphetamine	-0.15	150.00	0.13	995	6.9E5	MS2	
AutoMSn(219.89)	3.14	4	Ritalinic acid	-0.07	219.89	0.24	993	5.5E5	MS2	
AutoMSn(315.93)	4.40	9	Bromazepam	-0.13	315.93	0.08	977	1.6E5	MS2	
AutoMSn(295.92)	6.12	11	Diclofenac	-0.22	295.92	0.10	961	1.3E5	MS2	MS2 unspecified
AutoMSn(180.88)	2.47	2	Theobromine	-0.12	180.88	0.19	919	8.9E5	MS2	
AutoMSn(314.96)	6.40	12	Progesterone	-0.23	314.96	0.27	838	5.8E5	MS2	
AutoMSn(343.11)	6.99	13	THC-COOH	-0.26	343.11	0.08	970	3.4E5	tentative	

- Fully automated data processing and reporting
- MS<sup>n</sup> spectra extraction
  - Spectral library matching
  - Result report generation

Simple and fully automated workflow from sample injection to final result in 12 min

# Choice of sample preparation




- The Toxtyper works with various proven sample preparations, depending on sample specimen and main application.
- A selection of SOPs is part of the Toxtyper Tutorial, e.g.:
  - **Protein precipitation (PP)** for fast and simple sample preparation and preservation of glucuronides
  - **Liquid-liquid extraction (LLE)** or
  - **Solid phase extraction (SPE)**

# Drug identification in urine after LLE




Detection of Fentanyl and Buprenorphine at concentrations < 1 ng/mL after liquid-liquid extraction (LLE, 1 ml urine)


ng/mL	Compound
20	Oxycodone
20	Scopolamine
20	Atropine
4	Zopiclone
20	Tramadol
20	Nortilidine
4	Zolpidem
ISTD	D3-LSD
<b>0.8</b>	<b>Fentanyl</b>
20	Opipramol
<b>0.8</b>	<b>Buprenorphine</b>
4	Zaleplon




# Easy to use software



- Toxtyper software designed for non-MS experts
- Intuitive user interface for easy sample management



# HPLC-MS<sup>n</sup> data acquisition



## Bruker Column Kit

- Acclaim® RSLC 120 C18  
2.2  $\mu\text{m}$ , 120Å 2.1 x 100 mm



## Gradient

- H<sub>2</sub>O/ACN, 2 mM ammonium formate, 0.1% HCOOH
- 11 minutes total run time

## AutoMS<sup>2</sup>/MS<sup>3</sup> Data Acquisition

- Electrospray Ionization
- Zero delay alternating polarity (ZDA™)
- UltraScan 70-800  $m/z$  @ at 32.000  $m/z$  s<sup>-1</sup>
- **Scheduled Precursor List (SPL)** containing precursor mass and retention time information of all library compounds to trigger data dependent acquisition of MS<sup>n</sup> spectra

# Compound identification using library search



Compound	CAS	Monoisotopic neutral mass	Formula	Polarity	RT [min]
2-Amino-5-chloropyridine	1072-98-6	128.01	C5H5ClN2	POSITIVE	2.05
2-Hydroxyethylflurazepam	20971-53-3	332.07	C17H14ClFN2O2	POSITIVE	10.11
3-Hydroxybromazepam	13132-73-5	331.00	C14H10BrN3O2	POSITIVE	6.07
6-O-Acetylnorpine	2784-73-8	327.15	C18H21NO4	POSITIVE	3.56
7-Aminoflurazepam	34084-50-9	283.11	C16H14FN3O	POSITIVE	5.17
Alpha-hydroxyalprazolam	37115-43-8	324.08	C17H13ClN4O	POSITIVE	9.48
Alprazolam	28981-97-7	308.08	C17H13ClN4	POSITIVE	10.18

## Toxtyper & DOA Drugs of Abuse

butabral	11-26-9	224.12	C11H16N2O3	NEGATIVE	7.28
Clobazam	22316-47-8	300.07	C16H13ClN2O2	POSITIVE	10.78
Clonazepam	1622-61-3	315.04	C15H10ClN3O3	POSITIVE	9.98
Cocacethylene	529-38-4	317.16	C18H23NO4	POSITIVE	6.13
Cocaine	50-36-2	303.15	C17H21NO4	POSITIVE	5.16
Codine	76-57-3	299.15	C18H21NO3	POSITIVE	3.32
Creatinine	60-27-5	113.06	C4H7N3O	POSITIVE	0.50
D3-Clomipramine	136765-29-2	317.17	C19H20ClN2O3	POSITIVE	10.42
D4-Haloperidol	136765-35-0	379.17	C21H19ClF4NO2	POSITIVE	8.33
D5-Diazepam	65854-76-4	289.10	C16H8ClN2O5	POSITIVE	11.43
Desallyflurazepam	2895-65-9	288.05	C15H10ClFN2O	POSITIVE	10.48

## Psychotropic Drugs

Application Note LCMS-06  
**Adopting Forensic Analyses to Specific Lab Needs – Identifying Psychotropics in Serum Using the Toxtyper Open Library Concept**

**Abstract**  
 Various instrument types and LC-MS/MS approaches have been adopted in an effort to identify a comprehensive analysis method for detecting and identifying an array of compounds as possible in a single run. However, existing analyses often require the injection of a defined set of substances that are specific to a case of drug. To address this need, the Toxtyper™ Open Library Concept™ (OLC) was developed. The aim of this project was to develop an analysis method for psychotropic substances based on the Toxtyper™ Open Library Concept™. Psychotropic drug data from the Toxtyper library (TL) and opening standard from literature, especially from the Toxtyper™ Open Library Concept™ (OLC) were used to compile a scheduled open library (OL) for the Toxtyper. The OLC was used to identify a comprehensive list of psychotropic substances. A conventional GC-MS and an ion mobility™ analysis were used to identify all other parameters were adapted from the original Toxtyper approach.

**Keywords:** Non-scheduled and scheduled psychotropic substances, open library, GC-MS, Ion Mobility™, LC-MS/MS, Toxtyper™, Library search.

**Non-scheduled and scheduled psychotropic substances**  
 open library  
 GC-MS  
 Ion Mobility™  
 LC-MS/MS  
 Toxtyper™  
 Library search

**Maurer/Wissenbach/Weber LC-MS<sup>n</sup> Library of Drugs, Poisons, and Their Metabolites**

Hans H. Maurer  
 Dirk K. Wissenbach  
 Armin A. Weber

## Maurer/Wissenbach/Weber

**Various libraries available**

**Journal of MASS SPECTROMETRY**

**TRAPPING SPICE**

**WILEY**

**EXPRESS PUBLICATION**

## Synthetic Cannabinoids

# Toxtyper 3.0 library: 1187 compounds from various classes

## Antidepressants

Amitriptyline  
Amitriptylinoxide  
Clomipramine  
Desipramine  
Dibenzepin  
Dosulepine  
Doxepin  
Duloxetine  
Imipramine  
Melitracen  
Nortriptyline  
Oxypertine  
Protriptyline  
Reboxetine

## Analgesics

Acetylsalicylic acid  
Amidopyrine  
Fentanyl  
Nifenazone  
Salicylamide  
Salicylic acid

## Benzodiazepines

7-Aminoflunitrazepam  
Alprazolam  
Bromazepam  
Brotizolam  
Chlordiazepoxide  
Clobazam  
Clobenzepam  
Clonazepam  
Clotiazepam  
Delorazepam  
Desalkylflurazepam  
Diazepam  
Estazolam  
Flunitrazepam  
Lorazepam  
Lormetazepam  
Medazepam  
Metaclazepam  
Midazolam  
Nitrazepam  
Nordazepam  
Oxazepam  
Prazepam  
Temazepam  
Tetraepam

## Opiates

Codeine  
Heroin  
Hydrocodone  
Morphin-3-beta-D-glucuronide  
Morphine  
Morphine-6-beta-D-glucuronide  
Oxycodone  
Oxymorphone  
Papaverin

## Antifungal, Antibacterial ...

Griseofulvin  
Isoconazole  
Naftifine  
Ornidazole  
Phthalylsulfathiazole  
Sulfadiazine  
Sulfadoxine  
Sulfaethidole  
Sulfalene  
Sulfamerazine  
Sulfamethizole  
Sulfapyridine  
Terconazole  
Tolnaftate

## Amphetamines

Amphetamine  
MBDB  
MCPP  
MDA  
MDDMA  
MDEA  
MDMA  
Methamphetamine  
PMMA  
Pseudoephedrine  
Sibutramin

## Metabolites

Desmethyl-chlordiazepoxide  
Desmethyl-Citalopram  
Desmethylclobazam  
Desmethylclomipramine  
Desmethylclozapine  
Desmethyl-Mirtazapine  
Desmethylvenlafaxine  
Nordoxepine  
THC-COOH  
THC-OH

## Antihistaminics

Bamipine  
Cyclizine  
Isothipendyl  
Mequitazine  
Methaphenilene  
Oxomemazine  
Pyribenzamine  
Pyrilamine  
Tritoqualine

## Pesticides

Alachlor  
Flurochloridone  
Isoproturon  
Metamitron  
Methabenzthiazuron  
Methoprotryne  
Metsulfuron-methyl  
Monolinuron  
Napropamide  
Sebuthylazine  
Terbutylazine  
Terbutryn  
Triasulfuron

## Synthetic Cannabinoids

AM-1220  
AM-2233  
JWH-007  
JWH-015  
JWH-018  
JWH-019  
JWH-020  
JWH-072  
JWH-073  
JWH-081  
JWH-122  
JWH-122-5-fluoropentyl-derivate  
JWH-200  
JWH-210  
JWH-250  
JWH-307  
JWH-387  
JWH-398  
JWH-412  
Methanandamide  
RCS-4 ortho-isomer


**And many more ...**



# Identification by MS<sup>2</sup> and MS<sup>3</sup>: Nitrazepam





Extracted ion chromatogram



MS Spectra

experimental vs. library


Compound Spectra



# Identification by MS<sup>2</sup> and MS<sup>3</sup>: Nitrazepam





Extracted ion chromatogram



MS Spectra



experimental vs. library

Compound Spectra




MS<sup>2</sup> Spectra

experimental vs. library



# Identification by MS<sup>2</sup> and MS<sup>3</sup>: Nitrazepam



**Extracted ion chromatogram**



**MS Spectra**



experimental vs. library

Compound Spectra





**MS<sup>2</sup> Spectra**

experimental vs. library



**MS<sup>3</sup> Spectra**


experimental vs. library




**Result fully confirmed  
by MS<sup>2</sup> AND MS<sup>3</sup>**

# MS<sup>n</sup> reproducibility independent of concentration


Example: Amitriptyline



**Highly reproducible MS/MS spectra,  
 independent of concentration**



# Automated result reporting



- Toxtyper provides a **simple list with identified compounds** and few tentative hits which require manual validation.
- **Color coding** can be applied in order to facilitate the review process.
- The Toxtyper software allows for fast **interactive result evaluation** using chromatogram and MS<sup>n</sup> spectra views.
- All results are available as a **PDF report** that can be directly printed or sent via e-mail.

# Toxtyper review perspective: All at a glance



**Toxtyper - Review**

Home Submit Queue Review Admin

**Batches**

Name	Process Status	Creation
20190828	RUNNING	Wed Aug
Master Synth HT	DONE	Fri Aug
Master DVD HT	DONE	Thu Aug
Master SnthCann	DONE	Thu Aug
TT Lib pos	DONE	Thu Aug
bsc	DONE	Wed Aug

**Analyses**

Batch Name: Master DVD HT  
Batch Description: Test HT methods

Sample ...	Status	Method Name	Position	Description	Cmp IDs	Max Cmp	Creatinine [n
Blank	DONE	Psychotropics HT	1:1:1		0	-	0
QC	DONE	Psychotropics HT	1:1:2		3	Mirtazapine	0
QC	DONE	Psychotropics HT...	1:1:2		4	Mirtazapine	0
Blank	DONE	MWW HT	1:1:1		0	-	0
QC	DONE	MWW HT	1:1:2		6	Mirtazapine	0
QC	DONE	MWW HT_custom	1:1:2		7	Mirtazapine	0

**Identified Compounds**

Peak Number	Cmp Name	Retention Time [min]	m/z [Da]	d m/z [Da]	Intensity	Purity'	S/N	ID	Comment
1	Paracetamol	2.92	151.99	0.08	1.8E7	997	223	MS2/MS3	Mass Spectrum Comparison
3	Paracetamol-dimer	2.93	303.07	-0.07	3.9E6	996	67	MS2/MS3	
4	Mirtazapine	4.58	266.05	0.12	3.7E7	980	522	MS2/MS3	
5	Risperidone	5.51	411.22	-0.00	7.8E6	988	160	MS2/MS3	
6	D4-Haloperidol	7.91	380.23	-0.06	1.3E7	924	235	MS2	
7	Desalkylflurazepam	10.59	289.04	0.01	2.1E7	950	324	MS2	

**Toxtyper PDF Report**

**MS<sup>n</sup> Spectra Comparison**

# Batch submission for semi-quantification



1. Injection  
of calibration standards

2. Sample analysis including  
semi-quantitation

Batch Submission Wizard

**Sample Assignment**  
Add more samples to the batch or click "Submit" to finish the batch submission.

Create new Samples:

Sample Identifier Prefix:  Injection Volume (µL):

#	Sample Identifier	Description	Position	Method	Injection Volume	New Calibration	Calibration Table Folder	Calibration Table
1	Blank	MeOH	RA1	Toxtyper_Custom	5.0			
2	Cal Mix 1	10 cpds 2 ISTD	RA2	Toxtyper_One-point calibration	5.0	<input checked="" type="checkbox"/>	D:\Methods\Toxtyper2...	Toxtyper_Cal Mix 1.csv
				Toxtyper_One-point calibration	5.0	<input type="checkbox"/>	D:\Methods\Toxtyper2...	Toxtyper_Cal Mix 1.csv
3	Serum Sample 170420-1			Toxtyper_SemiQuant	5.0		D:\Methods\Toxtyper2...	
4	Serum Sample 170420-2			Toxtyper_SemiQuant	5.0		D:\Methods\Toxtyper2...	
5	Serum Sample 170420-3			Toxtyper_SemiQuant	5.0		D:\Methods\Toxtyper2...	
6	Serum Sample 170420-4			Toxtyper_SemiQuant	5.0		D:\Methods\Toxtyper2...	
7	Serum Sample 170420-5			Toxtyper_SemiQuant	5.0		D:\Methods\Toxtyper2...	
8	Serum Sample 170420-5			Toxtyper_SemiQuant	5.0		D:\Methods\Toxtyper2...	


< Back    Next >    **Submit**    Cancel

# Semi-quantification results



Sample-ID SQM 50 Station Toxtyper2  
 Submitter Tony Technician Method Toxtyper\_SemiQuant (1.1.1)  
 Analysis Name SQM 50\_RB2\_01\_507.d Acquisition Date 11/9/2016 3:10:50 PM  
 Sample Description

## Base Peak Chromatogram



## Library Search Results



Cmp Name	Cmp #	Purity'	RT [min]	d RT	m/z [Da]	d m/z	Intensity	Semi-quant Conc.	ID	Comment
Buprenorphine	6	917	4.33	-0.08	468.31	0.00	1.0 E8	90 ng/mL	MS2	
Methadone	8	982	4.76	-0.02	310.11	0.11	4.3 E7	54 ng/mL	MS2/MS3	
Diazepam	10	972	5.81	-0.11	284.89	0.19	1.5 E7	164 ng/mL	MS2	
D5-Diazepam	9	986	5.79	-0.10	289.95	0.16	1.5 E7	109 ng/mL	MS2	
Fentanyl	5	994	4.15	-0.13	337.21	0.02	1.5 E7	12 ng/mL	MS2	
D3-6-O-Acetylmorphine	3	970	2.90	-0.08	331.17	0.00	1.3 E7	44 ng/mL	MS2	
Amphetamine	2	993	2.85	0.07	136.07	0.03	2.6 E6	262 ng/mL	MS2/MS3	
Morphine	1	887	2.18	-0.02	285.99	0.15	2.6 E6	80 ng/mL	MS2	
Diphenhydramine	7	866	4.55	0.16	256.02	0.15	4.9 E5		MS2/MS3	
Progesterone	11	822	6.57	-0.26	314.99	0.24	4.8 E5		MS2	
Venlafaxine	4	835	4.08	0.06	278.01	0.20	3.4 E5		tentative	MS2 unspecific

Concentration values are listed for all calibrated target compounds



# Maurer/Wissenbach/Weber library for urine drug screening

Identification of **parent drugs** and **metabolites** including glucuronides



# Reporting of results MWW library in alphabetical order

Extract of Toxtyper MWW report. Grouping of drugs and their metabolites.

Butylscopolamine and metabolites

Butylscopolamine	38	983	5.59	360.28	-0.18	1.1 E9	MS2/MS3
Butylscopolamine artifact (-H2O)	53	934	8.10	342.22	-0.12	1.7 E8	MS2/MS3
Butylscopolamine-M (330)	44	956	7.03	330.20	-0.10	1.7 E8	MS2/MS3
Butylscopolamine-M (406)	28	931	4.26	406.25	-0.15	1.3 E8	MS2/MS3
Butylscopolamine-M (HO-aryl)		831	4.26	406.25	-0.15	1.3 E8	MS2/MS3
Butylscopolamine-M (HO-aryl-sulfate)	27	873	4.25	456.23	-0.13	8.0 E7	MS2/MS3
Butylscopolamine-M (sulfate)	39	958	5.78	440.19	-0.09	1.1 E9	MS2/MS3

Citalopram and metabolites

Butyrylcarnitine	13	845	2.87	232.06	0.04	6.1 E6	MS2/MS3
Chavicine-M (demethylenyl-)	62	721	10.03	274.09	0.01	1.3 E7	MS2/MS3
Chavicine-M (demethylenyl-methyl)	56	929	8.43	288.15	-0.05	3.0 E7	MS2/MS3
Chavicine-M (demethylenyl-methyl-glucuronide)		875	8.43	288.15	n/a	3.0 E7	MS2/MS3
Chavicine-M (demethylenyl-methyl-glucuronide) isomer 1	42	897	6.82	464.28	-0.18	9.0 E7	MS2/MS3
Chavicine-M (demethylenyl-methyl-glucuronide) isomer 2		880	6.82	464.28	-0.18	9.0 E7	MS2/MS3
Chavicine-M (demethylenyl-methyl-glucuronide) isomer 2	52	995	8.05	464.25	-0.15	8.2 E7	MS2/MS3
Chavicine-M (demethylenyl-methyl-glucuronide) isomer 1		993	8.05	464.25	-0.15	8.2 E7	MS2/MS3


Citalopram	51	981	7.90	325.18	-0.07	1.2 E9	MS2/MS3
Citalopram-M (N-oxide)	57	912	8.45	341.07	0.03	1.5 E7	MS2/MS3
Citalopram-M (bis-nor-)	48	857	7.44	297.16	-0.06	7.1 E7	MS2/MS3
Citalopram-M (glucuronide)	43	953	7.01	501.20	-0.10	4.1 E7	MS2/MS3
Clobazam-M/artifact (HO-glucuronide +NH4)	64	810	11.54	510.35	-0.25	5.2 E6	MS2
Corticosterone	47	792	7.40	347.26	-0.21	8.9 E6	MS2/MS3
Creatinine-M/artifact (227)	4	997	0.53	227.01	0.09	2.9 E8	MS2/MS3
Cropropamide	7	782	0.78	241.04	0.06	3.1 E7	MS2/MS3
Desmethyl-Citalopram	50	991	7.70	311.25	-0.19	8.0 E8	MS2/MS3
Diisooctylphthalate	72	978	16.70	391.27	-0.17	3.0 E7	MS2/MS3
Diphenhydramine	49	910	7.56	256.11	-0.28	2.8 E8	MS2/MS3
Diphenhydramine-M/artifact (N-oxide dimer)	55	946	8.29	543.34	-0.34	1.9 E9	MS2/MS3
Isovalerylcarnitine	17	796	3.56	245.91	0.19	4.9 E7	MS2/MS3
Lactose	2	841	0.50	365.12	-0.02	2.0 E7	MS2/MS3
Meptazinol	30	959	4.67	234.12	-0.13	3.8 E8	MS2/MS3
Meptazinol-M (HO-ring-glucuronide)	25	890	4.15	426.26	-0.16	8.9 E8	MS2/MS3
Meptazinol-M (glucuronide)	21	979	3.88	410.26	-0.16	3.8 E9	MS2/MS3
Meptazinol-M (nor-glucuronide)	22	974	4.00	396.27	-0.17	1.2 E9	MS2/MS3
Meptazinol-M (sulfate)	32	972	4.83	314.16	-0.06	4.1 E8	MS2/MS3

Meptazinol and metabolites

# Toxtyper versus Toxtyper MWW method Summary



	Toxtyper	MWW (2018)
Total no. of library compounds	1187	> 4500
Metabolites	Some	3000
Polarity	Positive and negative ion mode MS and MS <sup>n</sup>	Positive ion mode MS <sup>n</sup> spectra only
Retention time information	Yes, for all compounds	No
Use of SPL* for MS <sup>n</sup> acquisition	Yes	No, completely untargeted, data dependent acquisition
Identification based on	Precursor <i>m/z</i> MS <sup>n</sup> spectral information Retention time	Precursor <i>m/z</i> MS <sup>n</sup> spectral information Related Metabolites



- **Toxtyper users**  
and their main applications

# Selected Toxyper References



Forensics/Clinical – Forensics – Clinical – Police

- **Dr. Michael Böttcher**, MVZ Labor Dessau GmbH, Germany
- **Dr. Eric Franssen**, Onze Lieve Vrouwe Gasthuis OLVG Amsterdam, Netherlands
- **Dr. George W. Hime**, Miami Dade County Medical Examiner, USA
- **Prof. Dr. Thomas Krämer**, Forensic & Pharmacology, Zürich, Switzerland
- **Dr. Natalia Krupina**, Moscow Regional Bureau of Forensic Medicine, Russia
- **Prof. Dr. Volker Auwärter**, Jürgen Kempf, Laura Huppertz; Institute of Forensic Medicine Freiburg; Germany
- **Dr. Ivana Černá, Ústřední vojenská nemocnice, Praha, Česká republika**
- **Prof. Markus R. Meyer, (Prof. Dr. Dr. Hans Maurer)**, Experimental & Clinical Toxicology, Homburg, Germany
- **Prof. Dr. Wieland**, Katharinenhospital, Stuttgart, Germany → SYNLAB Leinfelden
- **Prof. Dr. Denis Hochstrasser, Dr. Pierre Lescuyer**; Hôpitaux Universitaires de Genève, Switzerland

# Selected Toxyper References

Forensics/Clinical – Forensics – Clinical – Police

- **Dr. Christian Vidal**, State Police (LKA) Lower Saxony, Germany
- **Dr. Michael Pütz, Dr. Lars Müller**, Federal Police (BKA) Wiesbaden, Germany
- **Colonel Dr. Theerin Sinchai**, Royal Thai Police, Thailand
- **Dr. Kristian Wittke**, State Police (LKA) Berlin, Germany

## **User Profile**

- Biggest private drug testing lab in Germany (Limbach group)
- Therapeutic drug monitoring, drugs of abuse testing, workplace drug testing, intoxication cases, clinical drugs testing, especially for addiction medicine.
- Use of mainly GC-MS for the General Unknown Screening before Toxyper

## **Toxyper value for the general unknown screening (GUS)**

- Rapid sample preparation for urine samples: no hydrolysis, no derivatization
- Comprehensive MWW library including metabolites (glucuronides!!!)
- Ease of training and simple data evaluation
  - saves time and money by increasing through-put
  - facilitates overnight and weekend service, e.g. in case of intoxications

Since 2016 a second Toxyper system is in use for saliva screening.

## Bruker "Customer Insight"



### LC-MS for Toxicology: Pushing the Limits of Speed and Sensitivity in Drug Screening

*"We were extremely impressed with the sensitivity, ease of use and rapid sample preparation."*

Ground-breaking LC-MS solutions enable general unknown analysis of drugs in biological samples, at the specialized laboratory of Dr. Michael Böttcher, MVZ Labor-Dessau GmbH



Dr. Michael Böttcher's specialist medical laboratory uses cutting-edge toxicology instrumentation to offer screening services for clinics and organizations across Europe.

*"We rely on Bruker's robust LC-MS/MS solution, the Toxtyper™, to rapidly analyze biological samples for drugs and drugs of abuse, with unprecedented accuracy."*

*"Post-run time is also greatly reduced on the Toxtyper™ compared to GC-MS. Data-mining is more efficient, and you don't need as much experience."*



## Streamlining sample preparation: GC-MS vs. Toxyper

### Timetable for GC/MS (1 urine sample)

- Enzymatic hydrolysis:	2.25 h
- Extraction:	0.60 h
- Derivatization:	0.30 h
- Chromatographic run + cooling:	0.30 h (17 min +2 min)
	<b>3.45 h = 3 h 25 min</b>

- All sample preparation steps with different recovery
- Different „on board“ stability for the different substances
- **Urine samples to be analysed batchwise**
- Acidic substances need separate extraction and derivatization

### Timetable for Toxyper (1 sample)

- Enzymatic hydrolysis:	-
- Extraction:	0.50 h
- Derivatization:	-
- Chromatographic run:	0.50 h (TT-)
	<b>1.00 h</b>

**Allows continuous analysis for all kind of samples**

- No extraction losses
- No additional stability problems
- Glucuronides detectable (MWW 2014 library)

*“It became clear that the hydrolyzing process for urine analysis is really limiting in GC-MS. If you can avoid it, it's a big step forward. Oral fluid is such a clean matrix; the Toxyper™ can achieve extreme sensitivity, which was unexpected at the beginning.”*

- Rapid sample preparation enables random access for urgent samples.
- For GC-MS this is uneconomical as sample prep is cumbersome and thus done batch-wise.

# Prof. Dr. Wieland / Dr. Shipkova



## User Profile

- Central Institute for Laboratory Medicine and Clinical Chemistry, Klinikum Stuttgart
- Toxtyper was purchased in 2016 by Prof. Wieland and Dr. Shipkova



## Toxtyper value

- Use of Toxtyper for screening of drugs and toxins in blood, saliva and urine.
- Replacement of GC-MS and a couple of targeted LC-MS/MS methods, e.g. for benzodiazepines, pregabalin and as well as synthetic opioids
- *"The Toxtyper is user-friendly, provides a reasonable solution for a 24 h/365 day emergency setting, and will not overburden the technical skills of laboratory technicians with minimal LC-MS experience"<sup>1)</sup>*

In 2018 Prof. Wieland and Dr. Shipkova changed to SYNLAB (Leinfelden) and purchased a Toxtyper to offer the same analytics at the new site.

<sup>1)</sup> *Clinical Mass Spectrometry* 4–5 (2017) 11–18

# Identification and detection limits in urine



Analyte	LLOI/LLOD Toxyper “MWW” [µg/L]	LLOI/LLOD Toxyper “DOAL” [µg/L]	Cut Off Immunoassay [µg/L]	LLOD of GC- or LC with MS Detection [µg/L]
Acetylcodeine	11/40	10/20	300	10
Amphetamine	n.d.	100/200	500	150
Benzoylcegonine	22/50	9/10	300	75 <sup>1</sup>
Methadone	3/5	2/5	300	40
Nordiazepam	80/200	50/100	200	125
Sufentanil	2/10	1/1	n.a.	1
Pregabalin	430/2000	430/500	n.a.	200

Urinary creatinine concentrations were between 50 and 350 mg/dl. n.d. = not determined; n.a. = not applicable. LLOI = lower limit of identification; LLOD = lower limit of detection.

<sup>1</sup> confirmation by GC-MS was based on cocaine detection.

\* LLOI was calculated, LLOD was experimentally established.

Simplified automated sample preparation for urine with Tecan robot: Hydrolysis, precipitation (incl. addition IS), 10fold dilution

- MWW:  
Maurer/Wissenbach/Weber method and library
- DOAL  
Toxyper drugs of abuse method and library (83 compounds)

# Publications on urine (2017) and oral fluid screening (2018)



Contents lists available at [ScienceDirect](#)

Clinical Mass Spectrometry

journal homepage: [www.elsevier.com/locate/clinms](http://www.elsevier.com/locate/clinms)



## Detection of drugs of abuse in urine using the Bruker Toxtyper™: Experiences in a routine clinical laboratory setting

M. Ott, K. Berbalk, T. Plecko, E. Wieland, M. Shipkova\*

Zentralinstitut für Klinische Chemie und Laboratoriumsmedizin, Klinikum Stuttgart, Stuttgart, Germany

### ARTICLE INFO

#### Keywords:

Drugs of abuse screening  
Assay performance  
Ion trap mass spectrometry  
Routine laboratory

### ABSTRACT

Urine screening can be used to detect misuse of illicit drugs and validate opioid replacement therapy. It is common that immunochemical assays are combined with GC-MS for these applications. Bruker has released an ion trap mass spectrometer, called Toxtyper™, with the potential to replace current screening algorithms to detect drug misuse.

Here, we compare our current strategy of urine screening for misuse of cannabis, amphetamines, opiates, benzodiazepine, methadone, sufentanil, and pregabalin to the Toxtyper protocols provided by the manufacturer.

The analytical performance of the instrument was determined on a selected drug panel and with 188 samples being compared to establish concordance between our currently established approach and the Toxtyper.

The lower limits of detection and identification for acetylcodeine, amphetamine, benzoylecgonine, fentanyl, and nordiazepam were below the common cut-offs for immunological screening assays and comparable to GC-MS. Imprecision and accuracy, both within- and between-series, were consistently < 25%. The screening for pregabalin and sufentanil was less sensitive than a targeted LC-MS/MS assay. Concordance in the screening of the Toxtyper assay protocols used and the inherent imprecision of the assay.

Our study has revealed that a considerable portion of our current time-consuming protocol for screening drugs of abuse in urine, based on the combination of multiple analytical methods, could be consolidated to the Toxtyper for a majority of the most-relevant substances in our patient population.

*Ther Drug Monit.* 2018 Jun 6. doi: 10.1097/FTD.0000000000000544. [Epub ahead of print]

## Evaluation of an Ion Trap LC-IT/MS Instrument (Toxtyper) for Drug of Abuse Screening in Oral Fluid.

Plecko T<sup>1</sup>, Berbalk K, Wieland E, Shipkova M.



### Author information

1 Zentralinstitut für Klinische Chemie und Laboratoriumsmedizin, Klinikum Stuttgart, Stuttgart, Germany.

### Abstract

**BACKGROUND:** Oral fluid (OF) is increasingly used as an alternative sample matrix in drug of abuse (DOA) screening. Screening is commonly performed by immunoassays and results confirmed using laborious GC-MS based methods. Therefore, an easy to operate ion trap mass spectrometric (IT-MS) commercial screening method (Toxtyper, Bruker Daltronik, Bremen, Germany) combined with a laboratory-developed sample preparation procedure have been evaluated for their application to OF.

**METHODS:** Oral fluid samples were subjected to protein precipitation followed by HybridSPE-phospholipid extraction. Chromatographic separation was achieved by ultra-high performance liquid chromatography (UHPLC); MS2/MS3 spectra were recorded by IT-MS and analyzed using a library provided by the manufacturer (Bruker Daltronik). The lower limit of detection (LLOD), linearity, imprecision, inaccuracy, and specificity (interferences, matrix effects) were investigated for methadone, buprenorphine, pregabalin, fentanyl, amphetamine, MDMA (3,4-Methylenedioxy-N-methylamphetamin), cocaine, acetylcodeine and nordiazepam, after

# Prof. Maurer, Experimental and Clinical Toxicology, Saarland University



Prof. Hans Maurer

RESEARCH ARTICLE

WILEY

Blood plasma level determination using an automated LC-MS<sup>n</sup> screening system and electronically stored calibrations exemplified for 22 drugs and two active metabolites often requested in emergency toxicology

Achim T. Caspar | Markus R. Meyer | Hans H. Maurer

Department of Experimental and Clinical Toxicology, Institute of Experimental and Clinical, Pharmacology and Toxicology, Saarland University, Homburg, Germany

**Correspondence**

Hans H. Maurer, Department of Experimental and Clinical Toxicology, Saarland University, D-66421 Homburg, Germany  
Email: hans.maurer@uks.eu

**Abstract**

Fast and comprehensive qualitative and quantitative methods preferably by gas chromatography-mass spectrometry (GC-MS) and/or liquid chromatography-mass spectrometry (LC-MS) are needed to support the (differential) diagnosis of acute poisonings in emergency toxicology. One option is a commercially available qualitative screening solution based on LC-MS<sup>n</sup> (Bruker Daltonik Toxtyper™, TT). Identified and toxicologically relevant compounds should be quantified to assess severity of poisonings. The aim of the present study was to test the TT system for quantification simultaneous with the screening process in blood plasma exemplified for 22 relevant drugs and two active metabolites. A standard liquid-liquid extraction was used for sample work-up followed by 1:5 dilution of the final extracts. They were analyzed using the TT system consisting of a Bruker amaZon speed ion trap and a Thermo Fisher Dionex Ultimate 3000 LC system. Plasma levels were assessed using full-scan data and an electronically stored five-point calibration. The calibration model was linear for the studied ranges and could be used for at least two months. The method was validated according to international guidelines. The acceptance criteria recommended for emergency toxicology for accuracy and precision were fulfilled for all tested compounds, but bromazepam, lorazepam, oxycodone, and prothipendyl could reliably be determined only above the therapeutic range. In conclusion, the presented procedure allowed the combination of a comprehensive LC-MS<sup>n</sup> screening with fast automated assessment of plasma levels for emergency toxicology of tested compounds.

Received: 30 May 2018 | Revised: 6 July 2018 | Accepted: 9 July 2018

DOI: 10.1002/dta.2466

# Calibrator and QC concentrations of 24 target compounds



**TABLE 1** Final plasma concentrations in µg/L of the analytes in calibrators (Cal 1-5) and quality control samples (QC) LOW, therapeutic (TH), HIGH, and dilution (DILU) as well as the used weightings of linear calibration model and the therapeutic ranges according to Schulz et al.<sup>32</sup>

Analyte	Weighting	Cal 1	Cal 2	Cal 3	Cal 4	Cal 5	QC LOW	QC TH	QC DILU	QC HIGH	Therapeutic Range
9-Hydroxyrisperidone	1/x <sup>2</sup>	100	500	1,000	1,500	2,000	120	100	800	1,600	10-100
Amitriptyline	1/x <sup>2</sup>	100	500	1,000	1,500	2,000	120	200	800	1,600	50-200
Bromazepam	equal	1,000	1,500	2,000	2,500	3,000	1,200	200	1,200	2,400	50-200
Carbamazepine	1/x <sup>2</sup>	2,500	5,000	10,000	15,000	20,000	3,000	12,000	-	16,000	2,000-12,000
Citalopram	1/x <sup>2</sup>	100	500	1,000	1,500	2,000	120	200	800	1,600	20-200
Clozapine	1/x <sup>2</sup>	100	500	1,000	1,500	2,000	120	600	800	1,600	100-600
Codeine	1/x <sup>2</sup>	250	500	1,000	1,500	2,000	300	250	800	1,600	30-250
Diazepam	equal	500	1,000	1,500	2,000	2,500	600	2,000	1,000	2,000	200-2,000
Diphenhydramine	equal	500	1,000	1,500	2,000	2,500	600	1,000	1,000	2,000	100-1,000
Doxepin	1/x <sup>2</sup>	100	500	1,000	1,500	2,000	120	150	800	1,600	20-150
Fluoxetine	1/x	250	500	1,000	1,500	2,000	300	500	800	1,600	120-500
Lorazepam	equal	500	1,000	1,500	2,000	2,500	600	250	1,000	2,000	20-250
Mirtazapine	1/x <sup>2</sup>	150	500	1,000	1,500	2,000	180	300	800	1,600	50-300
Nordiazepam	equal	750	1,000	1,500	2,000	2,500	900	800	1,000	2,000	200-800
O-Desmethyltramadol	1/x <sup>2</sup>	100	500	1,000	1,500	2,000	120	300	800	1,600	100-300
Oxazepam	equal	1,000	1,500	2,000	2,500	3,000	1,200	2,000	1,200	2,400	500-2,000
Oxycodone	1/x <sup>2</sup>	250	500	1,000	1,500	2,000	300	100	800	1,600	5-100
Paracetamol (acetaminophen)	1/x <sup>2</sup>	2,500	30,000	60,000	90,000	120,000	3,000	20,000	-	96,000	10,000-20,000
Prothipendyl	1/x <sup>2</sup>	50	125	250	375	500	60	10	200	400	5-10
Quetiapine	1/x <sup>2</sup>	100	500	1,000	1,500	2,000	120	500	800	1,600	100-500
Risperidone	1/x	50	125	250	375	500	60	60	200	400	20-60
Tramadol	1/x	250	500	1,000	1,500	2,000	300	300	800	1,600	100-300
Venlafaxine	1/x <sup>2</sup>	100	500	1,000	1,500	2,000	120	400	800	1,600	100-400
Verapamil	1/x <sup>2</sup>	100	500	1,000	1,500	2,000	120	350	800	1,600	50-350

Calibration covers the therapeutic to toxic range


# Toxtyper Results and Conclusions



## Toxtyper Analysis Report

Sample-ID 114773\_P\_LLE\_TTQ Station  
 Submitter Method toxtyper\_semiquant (2.0)  
 Analysis Name 114773\_P\_LLE\_TTQ\_BC3\_01\_5294.d Acquisition Date 5/17/2016 5:54:52 AM  
 Sample Description

### Base Peak Chromatogram



### Library Search Results

Cmp Name	Cmp #	Purity'	RT [min]	d RT	m/z [Da]	d m/z	Intensity	Semi-quant	Conc.	ID	Comment
3-Hydroxybromazepam	6	985	4.05	-0.04	331.94	0.06	6.9 E5			MS2/MS3	MS2 unspecific
Bromazepam	9	976	4.44	-0.09	315.91	0.10	3.6 E6	<1000 ng/mL		MS2	
Caffeine	2	971	3.00	-0.09	194.86	0.23	4.2 E6			MS2/MS3	
Carbamazepine	14	990	4.73	-0.15	236.89	0.21	3.3 E7	3089 ng/mL		MS2	
Carbamazepine10,11-epoxide	7	918	4.21	-0.08	252.94	0.16	5.0 E6			MS2	
D3-Trimipramine	15	922	4.75	-0.22	298.14	0.10	1.8 E8			MS2	
Diazepam	20	936	5.76	-0.16	285.00	0.08	4.1 E7	2299 ng/mL		MS2	
EDDP	8	970	4.44	-0.20	278.10	0.09	5.9 E6			MS2	
Lorazepam	17	899	5.02	-0.15	320.97	0.05	6.4 E5	<500 ng/mL		MS2/MS3	MS2 unspecific
Methadone	13	956	4.72	-0.06	310.14	0.08	5.8 E7			MS2/MS3	
Nordiazepam	18	977	5.27	-0.10	270.92	0.14	6.6 E7	>2500 ng/mL		MS2	
Oxazepam	16	977	4.92	-0.15	286.88	0.18	4.0 E6	<1000 ng/mL		MS2/MS3	MS2 unspecific
Temazepam	19	977	5.34	-0.17	300.96	0.11	4.3 E6			MS2/MS3	MS2 unspecific
Theobromine	1	900	2.49	-0.10	180.84	0.23	1.5 E6			MS2	

## 4 | CONCLUSIONS

The current method allowed a reliable and fast blood plasma screening and level assessment for 24 analytes in 24/7 emergency toxicology. The method was successfully integrated into the Toxtyper™ screening solution and validated according to international guidelines and recommendations. The blood levels of the analytes could be assessed automatically and quickly, based on a stored five-point calibration model and were given together with the blood screening results in a pdf report. However, the method was limited to blood level assessment in emergency toxicology to identify acute overdosing or poisoning.

# George W. Hime, Elisa Shoff Miami Dade Medical Examiner



*George W. Hime, M.S., Assistant Laboratory Director and Elisa Shoff, Toxicologist II, Miami-Dade Medical Examiner Dept. Toxicology Laboratory, Miami, USA*



*"The Toxtyper has been our go-to instrument for postmortem casework containing novel substances, substances that are otherwise undetectable via GCMS, and unknown substances we have not yet seen in the lab. This instrument has been a blessing to our laboratory over the last several years of dealing with the flood of new illegal drugs in Miami. From synthetic cathinones to synthetic fentanyls we have detected and identified them all reliably and with high confidence."*

**Benefits of using Toxtyper:**

*"Reliable, rugged, sensitivity beyond expectation, and ease of maintenance and operation are all features of this instrument. Full spectral MS3 data at pg/mL sensitivities in postmortem whole blood or tissues cannot be beat in our work by any other instrument."*

Toxtyper in the "opioid crisis":  
Use of smaller sub-libraries and SPE sample preparation for urine, blood and tissue samples to enhance sensitivity for the detection of highly potent NPS such as designer opioids.



# George W. Hime, Elisa Shoff Miami Dade Medical Examiner



## Qualitative Identification of Fentanyl Analogs and Other Opioids in Postmortem Cases by UHPLC-Ion Trap-MS<sup>n</sup>

Elisa N. Shoff\*, M. Elizabeth Zaney, Joseph H. Kahl, George W. Hime, and Diane M. Boland

Miami-Dade County Medical Examiner Department, 1851 NW 10th Avenue, Miami, FL 33136, USA

*Journal of Analytical Toxicology*, 2017;41:484–492

doi: 10.1093/jat/bkx041

Advance Access Publication Date: 13 June 2017

Article

OXFORD

\*Author to whom correspondence should be addressed

**Table III.** Comparison of number of positive cases for synthetic opioids on LC-Ion Trap MS<sup>n</sup> and GC-MS with respective limits of detection

Analog	No. of cases detected by LC-Ion Trap MS <sup>n</sup>	No. of cases detected by GC-MS	Limits of detection on LC-Ion Trap MS <sup>n</sup> (ng/mL)	Limits of detection on GC-MS (ng/mL)
Acetyl fentanyl	13	4	0.2	5
Beta-hydroxythiofentanyl	9	0	0.1	<sup>a</sup>
Butyryl fentanyl	3	1	0.2	0.5
Carfentanil	134	30	0.1	1
Furanyl fentanyl	37	9	0.5	2.5
4-Fluoroisobutyryl fentanyl	22	16	0.5	1
U-47700	4	1	0.5	–

<sup>a</sup>Beta-hydroxythiofentanyl cannot be detected via GC-MS.

### Abstract

Since 2013, the Miami-Dade County Medical Examiner Department has seen an increase in the number of opioid-related deaths. The majority of these deaths are attributed to fentanyl and its analogs. Fentanyl has infiltrated the local heroin supply. From 2014 to 2016, there was a significant increase in fentanyl-related deaths, followed by the identification of several novel fentanyl analogs. In 2016, four additional fentanyl analogs were identified: acetyl fentanyl, butyryl fentanyl, furanyl fentanyl and carfentanil. To address this epidemic, a method was developed for the qualitative identification of analgesic compounds in postmortem samples using UHPLC-Ion Trap mass spectrometry with MS<sup>n</sup> capabilities. The method was validated from 0.1 to 5 ng/mL, with a majority having limits of detection from ~500 postmortem cases were submitted for analysis based on case history and/or initial screening results. Of those cases, 375 were positive for illicit fentanyl and/or one or more fentanyl analogs. Due to the potency of these compounds, they were almost always included in the cause of death. Worth emphasizing and extremely alarming is the detection of carfentanil in 134 cases, 104 of which were initially missed by gas chromatography mass spectrometry. By incorpo-

ration of UHPLC-Ion Trap mass spectrometry, all 134 cases were identified. This method provides a sensitive and specific means for the identification of fentanyl and its analogs in postmortem samples. The identification of carfentanil in 134 cases, 104 of which were initially missed by gas chromatography mass spectrometry, highlights the importance of using sensitive and specific methods for the identification of synthetic opioids in postmortem cases.

# Prof. Volker Auwärter, Dr. Jürgen Kempf

## Inst. of Forensic Medicine, Freiburg

### User Profile

- Leading national institute in forensic toxicology
- World-renowned expertise in the analysis of New Psychoactive Substances (NPS)
- Drug testing (forensic and clinical), intoxication and death case analysis, full post-mortem toxicology, drug market monitoring

### Toxtyper value

- Comprehensive and cost-effective pre-screening
- Quick and easy update of Toxtyper libraries to keep pace with emerging NPS
- Use for all kind of biological samples as well as for drug preparations

*"It is often the case that, if they could, customers would always like to have the whole picture of xenobiotics present in the sample. But money is a restriction for them: if you want to conduct a full toxicological analysis, the price will be very high. This is why it's valuable to us to be able to offer broad screening with the Toxtyper™, covering almost all relevant analytes, at a lower price. It's a benefit to us and to them." explains Prof. Dr. Auwärter.*

*From: Bruker Customer Insights*

# Example: Semi-quantification of designer benzodiazepines using Toxtyper



Dr. Jürgen Kempf

- Use of a benzodiazepine sub-library
- LODs in serum between 1-25 ng/mL
- Two exceptions: 3-OH-Bromazepam and Demoxepam (50 ng/mL)

Name	LOD [ng/mL]
2-OH-Ethylflurazepam	25
3-OH-Bromazepam	50
3-OH-Flubromazepam	25
3-OH-Phenazepam	25
7-Aminoclonazepam	10
7-Aminoflunitrazepam	1
7-Aminonitrazepam	1
Adinazolam	5
$\alpha$ -OH-Alprazolam	5
$\alpha$ -OH-Midazolam	5
$\alpha$ -OH-Triazolam	5
Alprazolam	1
Bromazepam	10
Chlordiazepoxide	5
Clobazam	1
Clonazepam	5
Clonazolam	5
Cloniprazepam	1
Clotiazepam	1
Delorazepam	5

Name	LOD [ng/mL]
Demoxepam	50
Desalkylflurazepam	10
Deschloroetizolam	1
Diazepam	5
Diclazepam	1
Etizolam	1
Flubromazepam	5
Flubromazolam	5
Fludiazepam	10
Flunitrazepam	1
Flunitrazolam	5
Flurazepam	1
Fonazepam	5
Loprazolam	1
Lorazepam	5
Lormetazepam	5
Meclonazepam	5
Medazepam	5
Metizolam	1
Midazolam	1


Name	LOD [ng/mL]
Norflunitrazepam	5
Nifoxipam	-
Nimetazepam	10
Nitrazepam	25
Nitrazolam	5
Norclobazam	10
Nordazepam	5
Oxazepam	25
Phenazepam	5
Prazepam	1
Pyrazolam	5
RO-5-4864	10
Temazepam	1
Tetrazepam	5
Triazolam	5
Zaleplone	5
Zolpidem	1
Zopiclone	5

*Application Note LCMS 130: Detection and semi-quantitative determination of designer benzodiazepines in serum using LC-MSn*

# Example: Semi-quantification of designer benzodiazepines using Toxtyper

Sample - ID: Station Trapinski  
 Submitter JK: Method benzotyper b06a semiquant  
 Analysis Name BenzoMix I QC 5ngml\_13890.d: Acquisition Date 8/9/2016 6:11:48 PM  
 Sample Description Serum

## Base Peak Chromatogramm



## Library Search Results

Cmp Name	Cmp #	Purity	RT [min]	d RT	m/z [Da]	d m/z	Intensity	Semi-quant Conc.	ID	Comment
D4-Haloperidol	9	938	1.92	0.08	380.06	0.14	4.7 E8	48 ng/ml	MS2/MS3	
D3-Doxepin	7	951	1.85	0.07	282.96	0.14	3.7 E8	49 ng/ml	MS2/MS3	
Zolpidem	1	973	1.37	0.04	308.02	0.08	10.0 E7	7 ng/mL	MS2/MS3	
Flurazepam	5	992	1.74	0.08	387.93	0.27	8.9 E7	6 ng/mL	MS2/MS3	
Prazepam	21	996	8.04	0.07	324.94	0.16	8.4 E7	5 ng/mL	MS2/MS3	
D5-Diazepam	19	991	6.52	-0.16	289.88	0.12	7.5 E7		MS2/MS3	
Medazepam	8	912	1.86	0.12	270.84	0.16	2.9 E7	6 ng/mL	MS2/MS3	
Midazolam	6	990	1.74	0.08	325.89	0.21	2.5 E7	4 ng/mL	MS2/MS3	
Loprazolam	2	947	1.59	0.07	465.00	0.10	2.2 E7	6 ng/mL	MS2/MS3	
7-Aminoflunitrazepam	4	982	1.66	0.09	283.89	0.11	1.9 E7	5 ng/mL	MS2/MS3	
Metizolam	15	946	4.26	-0.05	328.90	0.10	1.7 E7	6 ng/mL	MS2/MS3	
Diclazepam	20	940	7.22	-0.05	319.03	-0.03	1.4 E7	5 ng/mL	MS2/MS3	
Clonazepam	10	981	3.07	0.23	353.92	0.18	1.2 E7	7 ng/mL	MS2/MS3	
Flunitrazepam	16	990	4.51	0.28	313.90	0.20	8.0 E6	5 ng/mL	MS2/MS3	
Alprazolam	13	964	4.03	-0.05	308.88	0.22	6.8 E6	5 ng/mL	MS2/MS3	
Clobazam	18	999	5.23	-0.08	300.89	0.11	5.5 E6	6 ng/mL	MS2/MS3	
Triazolam	14	991	4.25	-0.08	342.91	0.09	5.3 E6	6 ng/mL	MS2/MS3	
Flubromazepam	17	842	4.90	-0.07	332.84	0.16	2.8 E6	19 ng/mL	MS2/MS3	
Clonazepam	12	959	3.79	0.27	315.80	0.20	1.3 E6	< 5 ng/ml	MS2/MS3	
Chlordiazepoxide	3	822	1.62	0.11	299.89	0.11	7.4 E5	< 5 ng/ml	MS2/MS3	MS2 unspecific
Lorazepam	11	706	3.76	0.33	321.03	-0.03	6.0 E5	< 5 ng/ml	tentative	tentative

After a one-point calibration, qualitative and semi-quantitative results can be obtained from the same Toxtyper analysis

Application Note LCMS 130: Detection and semi-quantitative determination of designer benzodiazepines in serum using LC-MSn


# State Police Lab (LKA) Brandenburg

## “Magic Mushrooms”

Sample-ID ██████████ Station Toxyter\_2.0  
 Submitter ██████████ Method Toxyter\_Custom (2.0)  
 Analysis Name Psilos Realprobe\_RA5\_01\_46.d Acquisition Date 3/22/2017 3:17:06 PM

Sample Description

### Base Peak Chromatogram



### Library Search Results

Cmp Name	Cmp #	Purity <sup>1</sup>	RT [min]	d RT	m/z [Da]	d m/z	Intensity	ID	Comment
Psilocin	3	872	2.54	-0.06	204.99	0.14	6.8 E6	MS2/MS3	
Psilocybin	2	977	1.94	-0.07	285.04	0.06	3.8 E6	MS2	
Adenosine	1	998	1.74	0.13	268.02	0.08	8.1 E5	tentative	
Progesterone	7	877	6.54	-0.28	315.22	0.01	6.7 E5	MS2	
Methamphetamine	4	986	2.94	-0.14	150.04	0.09	6.1 E5	MS2	
3-OH-phenazepam	6	861	5.36	0.22	367.20	-0.20	5.7 E5	MS2	
Acebutolol	5	870	3.48	0.03	337.21	0.00	3.0 E5	tentative	MS2 unspecific

Identification of  
**Psilocin and Psilocybin.**  
 Methamphetamine is due  
 to contamination of the  
 packaging material.

Data courtesy of Landeskriminalamt Brandenburg  
 (State Police Lab), Eberswalde, Germany

# State Police Lab (LKA) Brandenburg

## "yellow oil sample" → Sustanon®

### Toxtyper Analysis Report




13\_01\_27.d

Station Toxtyper\_2.0

Method Toxtyper\_Custom (2.0)

Acquisition Date 3/20/2017 3:13:47 PM




Identification of  
**Testosteron**  
and its esters

#### Library Search Results

Cmp Name	Cmp #	Purity'	RT [min]	d RT	m/z [Da]	d m/z	Intensity	ID
Testosteron phenylpropionate	3	961	8.19	-0.10	421.29	-0.02	1.2 E9	MS2/MS3
Testosteron isocaproate	5	939	8.55	-0.13	387.32	-0.04	6.6 E8	MS2/MS3
Testosteron propionate	2	966	7.58	-0.13	345.23	-0.01	4.5 E8	MS2/MS3
Testosteron decanoate	4	829	8.47	-0.30	443.40	-0.04	2.3 E8	MS2/MS3
Testosteron	1	910	5.51	-0.11	289.18	-0.08	5.7 E6	MS2/MS3

Data courtesy of Landeskriminalamt Brandenburg  
(State Police Lab), Eberswalde, Germany



[www.bruker.com](http://www.bruker.com)

For research use only. Not for use in clinical diagnostic procedures.