

Resolving Industrial Solvent Mixtures on Packed and Capillary GC Columns

Nonpolar SPB-1 and polar SUPELCOWAX 10 capillary columns separate complex mixtures of most commonly used solvents with good resolution. Resolution on 0.53mm ID columns is similar to that on 0.25mm ID or 0.32mm ID columns, but 0.53mm ID columns accept up to 2000ng of each analyte, compared to about 500ng for 0.32mm ID columns. Thus, these columns are useful for analyses of samples containing both concentrated and dilute components. Alternatively, packed GC columns offer various combinations of selectivity and capacity. 80/120 Carbowax B/3% SP-1500 and 60/80 Carbowax B/1% SP-1510, respectively, separate 27 and 22 of 32 widely used solvents.

Key Words

- solvents ● industrial solvents ● water pollutants
- priority pollutants ● alcohols ● ketones ● esters
- aromatics ● chlorinated hydrocarbons

Introduction


Because many different solvents are used in industrial processes, there are many different situations requiring GC analysis. The appropriate column must be chosen from among nonpolar and polar capillary columns, or from several packed columns. The advantages of specific columns are described here.

Solvent Analyses On Capillary Columns

Both nonpolar SPB™-1 (bonded SE-30) and polar SUPELCOWAX™ 10 (bonded CARBOWAX® PEG 20M) capillary columns resolve most commonly used solvents, but with somewhat different selectivity. To determine which phase is appropriate for your sample, compare the performance of these columns with the same solvent mixture (Figure A).

Figure A also shows that complex solvent mixtures can be resolved on a 0.53mm capillary column, as well as on a 0.25mm or 0.32mm ID column. Sample resolution is almost the same, as a comparison

Figure A. SUPELCOWAX 10 and SPB-1 Capillary Columns Provide Alternative Selectivities for Solvents




794-0791,0792

Figure B. Common Solvents on an SPB-1 Column

Column: **SPB-1, 30m x 0.32mm ID, 1.0µm film**
 Cat. No.: **24045-U**
 Oven: 30°C (8 min) to 125°C at 4°C/min, hold 5 min
 Carrier: helium, 25cm/sec (set at 30°C)
 Det.: FID, 250°C
 Inj.: 0.2µL approximately equal proportions of each component, split 200:1, 250°C

- | | |
|---------------------------|----------------------------|
| 1. Methanol | 20. 1,4-Dioxane |
| 2. Methyl formate | 21. Ethyl Cellosolve |
| 3. Ethanol | U Unknown |
| 4. Acetone | 22. Methyl isobutyl ketone |
| 5. Isopropyl alcohol | 23. Toluene |
| 6. Ethyl formate | 24. Isobutyl acetate |
| 7. Methylene chloride | 25. Mesityl oxide |
| 8. n-Propyl alcohol | 26. n-Butyl acetate |
| 9. Methyl ethyl ketone | 27. Diacetone alcohol |
| 10. sec-Butyl alcohol | 28. 5-Methyl-2-hexanone |
| 11. Ethyl acetate | 29. Ethylbenzene |
| 12. Chloroform | 30. m-Xylene |
| 13. Tetrahydrofuran | 31. p-Xylene |
| 14. Methyl Cellosolve | 32. Cyclohexanone |
| 15. Isobutyl alcohol | U Unknown |
| 16. 1,1,1-Trichloroethane | 33. o-Xylene |
| 17. Isopropyl acetate | 34. Cellosolve acetate |
| 18. n-Butyl alcohol | 35. Butyl Cellosolve |
| 19. 2-Nitropropane | 36. Isoamyl acetate |



713-0795

of the chromatogram for the SPB-1 column in Figure A and the similar mixture in Figure B will reveal, but sample capacity for a 0.53mm ID column is far greater than for a narrow bore column – 2000ng per analyte, versus 500ng for a 0.32mm ID column (or 20,000ng per component for a packed column). Consequently, a 0.53mm ID column will resolve solvent mixtures with both concentrated and dilute components, such as a primary solvent containing trace contaminants. When you dilute such a sample to avoid overloading a 0.25mm or 0.32mm column with the concentrated component, the trace components can be difficult to detect or quantify. The large capacity of a 0.53mm ID column eliminates the need to dilute concentrated components. Furthermore, concentrated solvents will not damage the SPB-1 or SUPELCOWAX 10 bonded phase – you can introduce samples onto the column by split, splitless, or on-column injection.

A 0.53mm ID capillary column can be installed in instruments designed for packed columns and used with packed column injection and detection devices, as easily as a packed column. They will accept high carrier gas flow rates without significant loss of sample resolution (Figure C). When used with packed column flow rates, a 0.53mm column's efficiency is about equal to that of a

packed column. An analysis comparable to that from the packed column, however, requires significantly less time. When maximum efficiency is needed, the flow rate must be reduced toward the ideal rate of 2.5mL/minute (helium), and make-up gas must be supplied to the detector. Our injector and detector conversion kits enable you to convert a packed column system to accept 0.53mm ID capillary columns simply, inexpensively, and reversibly. For more about how to incorporate wide bore capillary columns into a packed column system, refer to our catalog.


Solvent Analyses on Packed Columns

Many industrial solvents can be separated on either 80/120 Carbowax™ B/3% SP-1500 or 60/80 Carbowax B/1% SP™-1510 packing. Although both of these packings resolve complex mixtures well, the Carbowax B/SP-1500 packing resolves more solvents from a single mixture. Under appropriate conditions, this packing separates 27 of 32 widely used solvents (Figure D). Minimal peak tailing by the low molecular weight alcohols confirms that the packing is chemically inert. The narrow, symmetrical, late eluting peaks (e.g., ethylbenzene) indicate good column

Figure C. A 0.53mm ID Capillary Column Can Be Used at Relatively High Flow Rates

Column: **SUPELCOWAX 10, 30m x 0.53mm ID, 1.0µm film**
 Cat. No.: **25301-U**
 Oven: 40°C (5 min) to 200°C at 5°C/min
 Carrier: helium, 5mL/min
 Det.: FID, 250°C
 Inj.: 1µL mixed solvents, split (50:1), 250°C

- | | |
|-------------------------------|-------------------------------|
| 1. Hexane | 26. 1,2-Dichloroethane |
| 2. 1,1-Dichloroethylene | 27. n-Butyl acetate |
| 3. Acetone | 28. 2-Hexanone (MBK) |
| 4. Methyl acetate | 29. Isobutanol |
| 5. trans-1,2-Dichloroethylene | 30. Isoamyl acetate |
| 6. Tetrahydrofuran | 31. p-Xylene |
| 7. Carbon tetrachloride | 32. m-Xylene |
| 8. 1,1,1-Trichloroethane | 33. 1-Butanol |
| 9. 1,1-Dichloroethane | 34. Methyl Cellosolve |
| 10. Ethyl acetate | 35. Amyl acetate |
| 11. Methanol | 36. o-Xylene |
| 12. Isopropanol acetate | 37. Isoamyl alcohol |
| 13. Methyl ethyl ketone | 38. Chlorobenzene |
| 14. Isopropanol | 39. Cellosolve |
| 15. Methylene chloride | 40. Styrene |
| 16. Benzene | 41. Cyclohexanone |
| 17. Propyl acetate | 42. N,N-Dimethylformamide |
| 18. Trichloroethylene | 43. Cyclohexanol |
| 19. Methyl isobutyl ketone | 44. Butyl Cellosolve |
| 20. Isobutyl acetate | 45. 2-Methylcyclohexanol |
| 21. Chloroform | 46. 1,2-Dichlorobenzene |
| 22. Tetrachloroethylene | 47. 1,1,2,2-Tetrachloroethane |
| 23. 2-Butanol | 48. o-Cresol |
| 24. Toluene | 49. p-Cresol |
| 25. 1,4-Dioxane | 50. m-Cresol |



713-0791

efficiency. At $k = 12$ (methyl ethyl ketone), the efficiency of a Carbopack B/3% SP-1500 column is 600-700 theoretical plates per foot.

A column that resolves the compounds in Figure D can be very useful to solvent analysts. For example, Cellosolves® are widely used as water-soluble solvents and as nonchlorinated replacements for chlorinated solvents, but these compounds are not resolved on all columns used for solvent analyses. A Carbopack B/3% SP-1500 column separates Cellosolves from each other and from the other components of the complex test mixture. Similarly, this column resolves *o*-, *m*-, and *p*-xylene which, like the Cellosolves, do not separate on most general purpose solvent columns.

A Carbopack B/3% SP-1500 column also separates solvents well under conditions that minimize analysis time. In Figure D, 27 of 32 solvents were resolved in 44 minutes by temperature programming the column at 4°C/minute. When the temperature programming rate was increased to 8°C/minute, analysis time was reduced to 25 minutes and 24 of 32 solvents were resolved (Table 1). Thus, you can reduce solvent analysis time significantly with only a slight sacrifice in separation.

A Carbopack B/3% SP-1500 column separates many compounds not shown in Figure D. A number of chlorinated compounds, and other materials, are no longer used as solvents, but still must be monitored by industrial hygienists. Some of these compounds, and their retention times, are listed in Table 1. This packing also is useful for analyzing volatile priority pollutants in wastewater (Table 1). The elution order for these pollutants on Carbopack B/3% SP-1500 differs from that on 60/80 Carbopack B/1% SP-1000, the primary column used for this analysis (1). By using a Carbopack B/3% SP-1500 column as a backup, you can quantify *cis*-1,3-dichloropropene, tetrachloroethylene, and tetrachloroethane. These pollutants are not isolated on the SP-1500 packing, when compared to the SP-1000 packing, but they coelute with different compounds.

The other Carbopack B packing, 60/80 Carbopack B/1% SP-1510, separates 22 of the 32 solvents in the test mixture (Figure E). Peaks are nearly symmetrical, even for the xylenes and other late-eluting compounds. At $k = 15$ (methyl ethyl ketone) column efficiency is 500-600 theoretical plates per foot. Separation on Carbopack B/1% SP-1510 is not quite as good as on Carbopack B/3% SP-1500, but the two packings differ in selectivity. Consequently, you can use a Carbopack B/1% SP-1510 column to resolve solvents not separated on Carbopack B/3% SP-1500, or to confirm results obtained on the latter column.

Either Carbopack B packing offers distinct advantages when compared to other packings used for general solvents analyses. For example, solvent resolution is better on Carbopack B/3% SP-1500 (Figure D) than on a widely used packing, Carbopack C/0.1% SP-1000 (Figure F). Overall, solvent resolution is comparable for Carbopack B/1% SP-1510 (Figure E) and the Carbopack C packing (Figure F). Both Carbopack B packings have more than twice the sample capacity of the Carbopack C material: a 10' x 1/8" column of either Carbopack B packing can accept up to 65µg of each sample component, but a comparable column of the Carbopack C packing can accept only 30µg of each component (Figure G).*

* Sample capacity was determined by progressively increasing the amount of a nonpolar test probe (4% pentane in methylene chloride) injected onto each column, and measuring peak asymmetry at 10% of peak height. The point at which peak asymmetry sharply increased was defined as the maximum sample capacity. Peaks with an asymmetry range of 0.90-1.10 were considered symmetrical.

Table 1. 80/120 Carbopack B/3% SP-1500 Column Monitors Many Solvents and Pollutants

Industrial Solvents		
Solvent	4°C/min temp. rise	t _R (min) 8°C/min temp. rise
Methanol	2.7	1.8
Methyl formate	4.2	2.7
Ethanol	5.1	2.9
Acetone	7.3	4.1
Methylene chloride	7.6	4.3
Isopropyl alcohol	8.6	4.6
Ethyl formate	9.3	5.0
Methyl acetate	9.3	5.0
n-Propyl alcohol	10.5	5.4
Tetrahydrofuran (THF)	13.5	7.1
Methyl ethyl ketone (MEK)	14.2	7.4
Methyl Cellosolve	15.4	7.8
sec-Butyl alcohol	15.8	8.0
Ethyl acetate	16.2	8.2
Isobutyl alcohol	16.6	8.3
2-Nitropropane	18.3	9.3
n-Butyl alcohol	18.6	9.8
Ethyl Cellosolve	21.9	11.0
n-Propyl acetate	24.2	12.2
sec-Butyl acetate	28.9	14.6
Cyclohexanone	29.9	15.4
Isobutyl acetate	30.6	15.4
n-Butyl acetate	32.2	16.2
Toluene	32.2	16.3
Mesityl oxide	32.2	16.3
Cellosolve acetate	36.4	18.6
Butyl Cellosolve	37.1	19.1
Isoamyl acetate	37.8	19.5
Ethylbenzene	38.5	20.1
m-Xylene	41.7	22.5
p-Xylene	42.8	23.5
o-Xylene	43.3	24.0

Volatile Wastewater Pollutants		Hazardous Solvents	
Pollutant	t _R (min) 4°C/min temp. rise	Solvent	t _R (min) 4°C/min temp. rise
Bromochloromethane	12.2	Chloroform	14.6
trans-1,2-Dichloroethylene	13.1	Pentane	15.2
1,1-Dichloroethane	13.1	1,1,1-Trichloroethane	18.3
1,2-Dichloroethane	16.0	1,4-Dioxane	18.6
Bromodichloromethane	20.4	Carbon tetrachloride	19.9
1,2-Dichloropropane	22.4	Isopropyl acetate	21.4
trans-1,3-Dichloropropene	22.4	Benzene	21.8
Chlorodibromomethane	22.4	Hexane	23.9
cis-1,3-Dichloropropene	24.3	1,1,2-Trichloroethane	24.4
Trichloroethylene	25.7	Methyl isobutyl ketone (MIBK)	28.3
1-Chloro-2-bromopropane	27.5	Methyl isoamyl ketone (MIAK)	36.4
Bromoform	31.2		
1,1,2,2-Tetrachloroethylene	33.0		
1,4-Dichlorobutane	34.4		
1,1,2,2-Tetrachloroethane	34.4		
Chlorobenzene	36.2		

Packing: **80/120 Carbopack B/3% SP-1500**
 Cat. No.: **11813** (15g bottle)
 Column: 10' x 1/8" SS
 Oven: 70°C to 235°C at 4°C/min or 8°C/min
 Carrier: nitrogen, 20mL/min
 Det.: FID
 Inj.: 0.3µL (industrial and hazardous solvents) or
 0.1µL (wastewater pollutants),
 approximately equal volumes of neat components

Solvents in Figures D-F

1. Methanol	12. Methyl Cellosolve	23. n-Butyl acetate
2. Methyl formate	13. sec-Butyl alcohol	24. Toluene
3. Ethanol	14. Ethyl acetate	25. Mesityl oxide
4. Acetone	15. Isobutyl alcohol	26. Cellosolve acetate
5. Methylene chloride	16. 2-Nitropropane	27. Butyl Cellosolve
6. Isopropyl alcohol	17. n-Butyl alcohol	28. Isoamyl acetate
7. Ethyl formate	18. Ethyl Cellosolve	29. Ethylbenzene
8. Methyl acetate	19. n-Propyl acetate	30. m-Xylene
9. n-Propyl alcohol	20. sec-Butyl acetate	31. p-Xylene
10. Tetrahydrofuran (THF)	21. Cyclohexanone	32. o-Xylene
11. Methyl ethyl ketone (MEK)	22. Isobutyl acetate	

Figure D. An 80/120 Carbopack B/3% SP-1500 Column Resolves Most Components of a 32 Solvent Mixture

Packing: **80/120 Carbopack B/3% SP-1500**
 Cat. No.: **11813-U** (15g bottle)
 Column: 10' x 1/8" SS
 Oven: 70°C to 235°C at 4°C/min
 Carrier: nitrogen, 20mL/min
 Det.: FID
 Inj.: 0.3µL, approximately equal volumes of neat components




Figure E. Selectivity for Solvents on Carbopack B/1% SP-1510 Differs from That for Carbopack B/3% SP-1500

Packing: **60/80 Carbopack B/1% SP-1510**
 Cat. No.: **11809** (15g bottle)
 Column: 10' x 1/8" SS
 Oven: 100°C to 225°C at 8°C/min
 Carrier: nitrogen, 20mL/min
 Det.: FID
 Inj.: 0.3µL, approximately equal volumes of neat components




Figure F. Industrial Solvents on a Carbopack C/0.1% SP-1000 Column

Packing: **80/100 Carbopack C/0.1% SP-1000**
 Cat. No.: **11820** (15g bottle)
 Column: 10' x 1/8" SS
 Oven: 70°C to 225°C at 4°C/min
 Carrier: nitrogen, 20mL/min
 Det.: FID
 Inj.: 0.3µL, 3% each component by volume (approx. 12µg each)





Figure G. Sample Capacity of Carbopack B Packings Is Twice That for Carbopack C




Despite the superiority of the Carbopack B packings, Carbopack C/0.1% SP-1000 is useful for separating complex mixtures of alcohols, ketones, esters, aromatics, and chlorinated hydrocarbons. Alcohol peaks do not tail. Depending on the separation needed, this packing can be used in place of the Carbopack B packings, or to confirm results obtained on them. Table 2 lists retention data for solvents on Carbopack C/0.1% SP-1000. This packing's capacity is less than the capacities of the Carbopack B packings (as discussed above), and sample volume should be restricted to 0.1-0.5µL.

SP-2100 methyl silicone is another useful phase for separating a variety of solvents, with elution order generally according to increasing boiling points. We suggest a mixed phase of 20% SP-2100 combined with 0.1% Carbowax® 1500 (to reduce tailing of alcohols) on 100/120 SUPELCOPORT™ (Figure H, Table 2). The packing is perhaps most useful for separating complex mixtures of chlorinated solvents (Figure I), or for resolving impurities from a single chlorinated solvent.

Figure H. Solvents on a 20% SP-2100/0.1% Carbowax 1500 Column


Packing: **20% SP-2100/0.1% Carbowax 1500 on 100/120 SUPELCOPORT**
 Cat. No.: **11821** (20g bottle)
 Column: 10' x 1/8" SS
 Oven: 100°C
 Carrier: nitrogen, 20mL/min
 Det.: FID
 Inj.: 0.5µL, approximately equal volumes of neat components



797-0431

Figure I. Chlorinated Solvents on a 20% SP-2100/0.1% Carbowax 1500 Column

Packing: **20% SP-2100/0.1% Carbowax 1500 on 100/120 SUPELCOPORT**
 Cat. No.: **11821** (20g bottle)
 Column: 10' x 1/8" SS
 Oven: 100°C
 Carrier: nitrogen, 20mL/min
 Det.: FID
 Inj.: 0.5µL, approximately equal volumes of neat components



797-0432

Despite the inclusion of Carbowax 1500 in this packing, methanol and ethanol tail when present in low concentrations. SP-2100/Carbowax 1500 columns also are poorly suited to separating mixtures of aliphatic solvents combined with alcohols, ketones, etc. The hydrocarbons tend to be eluted along with the other components, unless the aliphatics have either lower or higher boiling points. Furthermore, the Carbowax 1500 restricts the upper temperature limit of this packing to 175°C.


Special Purpose Packed Columns

Relative to other packings, 20% SP-2401/0.1% Carbowax 1500 on 100/120 SUPELCOPORT delays the elution of ketones. For example, methyl ethyl ketone is eluted after ethyl acetate, with good resolution (Figure J). In contrast, 20% SP-2100/0.1% Carbowax 1500 barely separates this pair (Figure H). Retention times for solvents on the SP-2401 packing at 100°C are given in Table 2. At higher temperatures, isophorone is eluted in a reasonable analysis time.

When used as a stationary phase, THEED (tetrahydroxyethylenediamine) retards elution of alcohols relative to other classes of compounds. It also allows hydrocarbons to elute rapidly, compared to other compounds with similar boiling points. Figure K shows the rapid elution of heptane and octane, esters, aromatics, and finally alcohols from 15% THEED on 100/120 Chromosorb® W AW. Retention data for other compounds are listed in Table 2. We do not recommend using a silanized support with THEED because it would be extremely difficult to coat with the THEED stationary phase.

Figure J. Solvents on a Ketone-Retarding Column

Packing: **20% SP-2401/0.1% Carbowax 1500 on 100/120 SUPELCOPORT**
 Cat. No.: **11822** (20g bottle)
 Column: 10' x 1/8" SS
 Oven: 100°C
 Carrier: nitrogen, 20mL/min
 Det.: FID
 Inj.: 0.1µL, approximately equal volumes of neat components



797-0433

Table 2. Retention Times for Common Solvents on Packed Column Phases (t_R in min)

Compound	Column				
	80/100 Carbowax C/ 0.1% SP-1000	10% SP-1000	20% SP-2100/ 0.1% Carbowax 1500	20% SP-2401/ 0.1% Carbowax 1500	15% THEED
Alcohols					
Methyl	0.6	4.0	1.1	1.1	3.0
Ethyl	0.9	4.1	1.2	1.2	3.4
Isopropyl	1.5	4.2	1.4	1.4	3.0
n-Propyl	1.8	6.2	1.8	1.7	4.9
Isobutyl	4.0	7.8	2.5	2.2	5.4
n-Butyl	5.5	10.1	3.0	2.5	7.6
Diacetone	Decomposes	3.2	8.0	14.4	Not eluted
Cellosolves					
Methyl	3.4	13.1	2.7	2.8	12.1
Ethyl	8.6	15.8	4.0	3.8	12.3
Butyl	Late	39.4	12.4	9.7	27.0
Ketones					
Acetone	1.3	3.1	1.4	2.2	1.4
MEK	3.2	3.9	2.1	3.1	1.8
MIBK	21.9	5.9	4.6	6.1	2.3
DIBK	Late	12.4	—	—	—
Mesityl oxide	Late	10.8	6.4	7.6	4.1
Cyclohexanone	—	—	11.1	18.3	Not eluted
Isophorone	—	—	43.9	—	Not eluted
Ethers					
Ethyl	2.7	2.0	—	—	—
Dioxane	4.4	7.7	—	—	—
Acetates					
Ethyl	4.5	3.5	2.2	2.5	1.3
Isopropyl	9.0	3.6	3.0	3.0	1.3
n-Propyl	Late	4.8	3.9	4.0	1.7
Isobutyl	Late	5.7	5.4	5.5	1.8
n-Butyl	Late	7.2	6.8	6.5	2.2
Cellosolve	Late	21.7	11.1	12.0	6.5
Formates					
Methyl	0.9	2.9	1.2	1.3	—
Ethyl	1.8	3.3	1.6	1.8	—
Chloroalkanes					
Carbon tetrachloride	4.9	3.6	—	—	—
1,1,1-Trichloroethane	4.4	—	3.0	1.7	1.3
Methylene chloride	1.3	4.1	1.7	1.0	1.2
Trichloroethylene	8.9	5.3	—	—	—
Chloroform	2.8	5.7	—	—	—
Tetrachloroethylene	Late	6.3	—	—	—
1,2-Dichloroethane	3.2	7.2	—	—	—
1,2,2-Trichloroethane	Late	17.8	—	—	—
Aromatics					
Benzene	9.8	4.6	3.2	2.2	1.4
Toluene	Late	6.9	5.8	3.5	1.9
Ethylbenzene	Late	10.0	—	—	—
m- & p-Xylene	Late	10.8	10.5	5.5	2.6
o-Xylene	Late	13.3	12.1	6.6	3.3
Styrene	Late	18.6	—	—	—
Aliphatics					
n-Hexane	Late	2.1	2.4	1.3	0.8
n-Heptane	Late	2.4	4.0	1.8	0.8
Isooctane	Late	2.2	—	—	—
n-Octane	Late	3.0	7.1	2.6	0.9
n-Nonane	—	—	13.0	3.9	1.1
n-Decane	Late	5.6	24.3	6.2	1.3
Miscellaneous					
Carbon disulfide	2.6	2.5	—	—	—
THF	2.4	3.9	2.7	—	1.6
2-Nitropropane	4.6	10.2	3.6	7.0	3.5
DMF	8.1	33.5	5.4	—	—

80/100 Carbowax C/0.1% SP-1000: 6' x 1/8" SS column, 20mL/min flow rate (nitrogen). 10% SP-1000 on 80/100 SUPELCOPORT: 20' x 1/8" SS column, 30mL/min flow rate. Other packings: 10' x 1/8" SS columns, 20mL/min flow rate. Oven (all columns): 100°C

Monitoring Water in Solvents on a Packed Column

Determining the amount of water in a sample by GC is difficult. With many columns, the water peak will tail severely. The degree of tailing is influenced by a variety of factors, including the nature of the column packing and the composition of the column tubing and of the instrument inlet. Diatomite supports cause tailing, even when acid washed and dimethylchlorosilane treated. However, water does not tail on porous polymers (e.g., Chromosorb 101) because these supports are polar in nature.

Water in solvents can be quantified rapidly with Chromosorb 101 packed in a short metal tube. The analysis of water in toluene in Figure L was performed on a GC equipped with a thermal conductivity detector. The Chromosorb 101 column also can be used to separate some solvent mixtures.


Metal column tubing has been blamed for water peak tailing, but this depends on both the quality of the tubing and the concentration of water in the sample. We have been able to detect as little as 0.01% water in solvents, using columns made with our chromatographic quality stainless steel tubing. Some analysts prefer using glass or Teflon® tubing to obtain a symmetric water peak. The overall efficiency of a glass column is higher than that of stainless steel, but efficiency of Teflon columns is poor.

A metal inlet in an instrument can be highly adsorptive and thus cause tailing. Avoid this by using a PureCol™ glass insert in the inlet or by injecting the sample directly onto the column.

A calibration curve for water can be obtained by injecting onto the column equal volumes of several water standards. Plot peak height versus water concentration (Figure M). The calibration mixtures should be prepared from a water-free solvent and dried with a good drying agent, such as Molecular Sieve 5A. The syringe used to inject the sample also should be water-free. For details on an internal standard procedure for water, see Hogan, *et al.* (2). References 3-5 also deal with determining water in samples.

Figure K. Solvents on an Alcohol-Retarding Column


Packing: **15% THEED on 100/120 Chromosorb WAW**
Cat. No.: **11823** (20g bottle)
Column: 10' x 1/8" SS
Oven: 80°C
Carrier: nitrogen, 20mL/min
Det.: FID
Inj.: 0.5µL, approximately equal volumes of neat components



797-0434

Figure L. Trace Water in Toluene


Packing: **60/80 Chromosorb 101**
Cat. No.: **20213** (50g bottle)
Column: 3' x 1/8" SS
Oven: 165°C
Carrier: helium, 20mL/min
Det.: TCD
Inj.: 0.1µL



797-0441

Figure M. Peak Height vs. Concentration of Water in a Solvent

Packing: **60/80 Chromosorb 101**
Cat. No.: **20213** (50g bottle)
Column: 3' x 1/8" SS
Oven: 165°C
Carrier: helium, 20mL/min
Det.: TCD
Inj.: 0.1µL



797-0442

References

1. *Federal Register*, 44, No. 233, Dec. 3 (1979).
 2. Hogan, J.M., R.A. Engel, and H.F. Stevenson, *Anal. Chem.*, 42, 249 (1970).
 3. Bennett, O.F., *Anal. Chem.*, 36, 684 (1964).
 4. Gvozdovich, T.N., G.S. Grinberg, L.V. Zuyeva, and Y.I. Yashin, *Petroleum Chemistry (USSR)*, 12, No. 2, 120 (1972) English.
 5. Hollis, O.L. and W.W. Haye, *J. Gas Chromatogr.*, 4, 235 (1966).
- References not available from Supelco.

Ordering Information:

Stock Packed Stainless Steel Columns for Solvents Analyses

General Config. (bend to fit most GCs)	Hewlett-Packard 5700, 5992-3 GC-MS (config. 5)	Hewlett-Packard 5880, 5890 5987, 6890 (config. A)	Perkin-Elmer 900, 3920 Sigma 1,2,3 (not on-column injection)	Perkin-Elmer 8300, 8400, 8500 8600, 8700, Auto System (not on-column injection)	Varian 3700 Vista series (FID)	Varian 3300/3400 Vista series (FID)
3% SP-1500 on 80/120 Carbowax B (Cat. No. 1-1813), 10' x 1/8" column						
12592	12593	12594	12595-U	13734-U	12596	13735-U
0.1% SP-1000 on 80/100 Carbowax C (Cat. No. 1-1820), 6' x 1/8" column						
12495-U	12499	12500-U	12496	13736-U	12497	13737
10% SP-1000 on 80/100 SUPELCOPORT (Cat. No. 1-1872), 10' x 1/8" column						
12537	12541	13752-U	12538	13753	12539	13754
20% SP-2100/0.1% Carbowax 1500 on 100/120 SUPELCOPORT (Cat. No. 1-1821), 10' x 1/8" column						
12718-U	12740-U	12804-U	12751	13773	12773	13774

Stock Packings for Solvents Analyses

Description	Cat. No.
80/120 Carbowax B/3% SP-1500, 15g	11813-U
60/80 Carbowax B/1% SP-1510, 15g	11809
80/100 Carbowax C/0.1% SP-1000, 15g	11820
GP 20% SP-2100/0.1% Carbowax 1500 on 100/120 SUPELCOPORT, 20g	11821
GP 20% SP-2401/0.1% Carbowax 1500 on 100/120 SUPELCOPORT, 20g	11822
15% THEED on 100/120 Chromosorb W AW, 20g	11823
10% SP-1000 on 80/100 SUPELCOPORT, 20g	11872
Chromosorb 101	
60/80, 50g	20213
80/100, 50g	20214
100/120, 50g	20215

GP – Indicates packing has been tested for specific analysis shown in this bulletin.

Nukol™ capillary columns provide symmetrical peaks and good resolution for glycols analyses. For information about these columns, please see our catalog.

Capillary Columns

Description	Cat. No.
SPB-1	
30m x 0.25mm ID, 0.25µm film	24028-U
60m x 0.25mm ID, 0.25µm film	24030-U
30m x 0.32mm ID, 0.25µm film	24044
30m x 0.32mm ID, 1.0µm film	24045-U
60m x 0.32mm ID, 0.25µm film	24046
15m x 0.53mm ID, 1.5µm film	25302-U
30m x 0.53mm ID, 1.5µm film	25303
30m x 0.53mm ID, 5.0µm film	25345-U
60m x 0.53mm ID, 5.0µm film	25349
SUPELCOWAX 10	
30m x 0.25mm ID, 0.25µm film	24079
60m x 0.25mm ID, 0.25µm film	24081
30m x 0.32mm ID, 0.25µm film	24080-U
60m x 0.32mm ID, 0.25µm film	24082
15m x 0.53mm ID, 1.0µm film	25300-U
30m x 0.53mm ID, 1.0µm film	25301-U
60m x 0.53mm ID, 2.0µm film	25376

For other column lengths and film thicknesses, see our catalog.

Trademarks

Carbowax, Nukol, PureCol, SP, SPB, SUPELCO, SUPELCOPORT, SUPELCOWAX – Sigma-Aldrich Co.

CARBOWAX, Cellosolve – Union Carbide Corp.

Chromosorb – Manville Corp.

Teflon – E. I. du Pont de Nemours & Co., Inc.

Fused silica columns manufactured under HP US Pat. No. 4,293,415.

80/120 Carbowax B/3% SP-1500 was developed in cooperation with Dr. A. DiCorcia of the University of Rome.

BULLETIN 824

For more information, or current prices, contact your nearest Supelco subsidiary listed below. To obtain further contact information, visit our website (www.sigma-aldrich.com), see the Supelco catalog, or contact Supelco, Bellefonte, PA 16823-0048 USA.

ARGENTINA · Sigma-Aldrich de Argentina, S.A. · Buenos Aires 1119 AUSTRALIA · Sigma-Aldrich Pty. Ltd. · Castle Hill NSW 2154 AUSTRIA · Sigma-Aldrich Handels GmbH · A-1110 Wien
 BELGIUM · Sigma-Aldrich N.V./S.A. · B-2880 Bornem BRAZIL · Sigma-Aldrich Quimica Brasil Ltda. · 01239-010 São Paulo, SP CANADA · Sigma-Aldrich Canada, Ltd. · 2149 Winston Park Dr., Oakville, ON L6H 6J8
 CZECH REPUBLIC · Sigma-Aldrich s.r.o. · 186 00 Praha 8 DENMARK · Sigma-Aldrich Denmark A/S · DK-2665 Vallensbaek Strand FINLAND · Sigma-Aldrich Finland/YA-Kemia Oy · FIN-00700 Helsinki
 FRANCE · Sigma-Aldrich Chimie · 38297 Saint-Quentin-Fallavier Cedex GERMANY · Sigma-Aldrich Chemie GmbH · D-82041 Deisenhofen GREECE · Sigma-Aldrich (o.m.) Ltd. · Ilioupoli 16346, Athens
 HUNGARY · Sigma-Aldrich Kft. · H-1067 Budapest INDIA · Sigma-Aldrich Co. · Bangalore 560 048 IRELAND · Sigma-Aldrich Ireland Ltd. · Dublin 24 ISRAEL · Sigma Israel Chemicals Ltd. · Rehovot 76100
 ITALY · Sigma-Aldrich s.r.l. · 20151 Milano JAPAN · Sigma-Aldrich Japan K.K. · Chuo-ku, Tokyo 103 KOREA · Sigma-Aldrich Korea · Seoul MALAYSIA · Sigma-Aldrich (M) Sdn. Bhd. · Selangor
 MEXICO · Sigma-Aldrich Química S.A. de C.V. · 50200 Toluca NETHERLANDS · Sigma-Aldrich Chemie BV · 3330 AA Zwijndrecht NORWAY · Sigma-Aldrich Norway · Torshov · N-0401 Oslo
 POLAND · Sigma-Aldrich Sp. z o.o. · 61-663 Poznań PORTUGAL · Sigma-Aldrich Quimica, S.A. · Sintra 2710 RUSSIA · Sigma-Aldrich Russia · Moscow 103062 SINGAPORE · Sigma-Aldrich Pte. Ltd.
 SOUTH AFRICA · Sigma-Aldrich (pty) Ltd. · Jet Park 1459 SPAIN · Sigma-Aldrich Quimica, S.A. · 28100 Alcobendas, Madrid SWEDEN · Sigma-Aldrich Sweden AB · 135 70 Stockholm
 SWITZERLAND · Supelco · CH-9471 Buchs UNITED KINGDOM · Sigma-Aldrich Company Ltd. · Poole, Dorset BH12 4QH
 UNITED STATES · Supelco · Supelco Park · Bellefonte, PA 16823-0048 · Phone 800-247-6628 or 814-359-3441 · Fax 800-447-3044 or 814-359-3044 · email: supelco@sial.com

H

Supelco is a member of the Sigma-Aldrich family. Supelco products are sold through Sigma-Aldrich, Inc. Sigma-Aldrich warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product for a particular use. Additional terms and conditions may apply. Please see the reverse side of the invoice or packing slip.

BKU