

Application Data Sheet

No. 10

System Gas Chromatograph

Sulfur Analyzer Nexis GC-2030PFPD2 GC-2014PFPD2

This method is for determining the sulfide compounds in gasoline using a pulsed flame photometric detector (PFPD) and capillary column. This system is composed of one split/splitless injection port, one capillary column and one PFPD. The system includes LabSolutions GC workstation software.

Analyzer Information

System Configuration:

Capillary Inlet / Capillary column / PFPD detector

Sample Information:

Sulfur compounds in light petroleum liquids, such as H₂S, COS, SO₂, mercaptans, aromatic sulfur compounds and sulfides

Methods met:

ASTM-D6228

Concentration Range:

| No. | Name of Compound | Concentration Range | |
|-----|------------------|---------------------|------------|
| | | Low Conc. | High Conc. |
| 1 | H ₂ S | 0.05ppmV | 100ppmV |
| 2 | COS | 0.05ppmV | 100ppmV |
| 3 | MeSH | 0.05ppmV | 100ppmV |
| 4 | EtSH | 0.05ppmV | 100ppmV |
| 5 | DMS | 0.05ppmV | 100ppmV |
| 6 | CS ₂ | 0.05ppmV | 100ppmV |
| 7 | PrSH | 0.05ppmV | 100ppmV |
| 8 | BuSH | 0.05ppmV | 100ppmV |

Detection limits may vary depending on the sample. Please contact us for more consultation.

System Features

- Sulfur analysis in light petroleum liquids and gasoline
- Sample lines including injection port inert in order to avoid absorption
- High selectivity for sulfur
- Equimolar, simplifies quantification of unknowns

Typical Chromatograms




Fig. 1 Chromatogram of PFPD

First Edition: November, 2017