

Application  
Data Sheet

No.21

System Gas Chromatograph

High Sensitive CO, CO<sub>2</sub>, CH<sub>4</sub> Analysis  
Nexis GC-2030CCC5  
GC-2014CCC5

This system is designed to measure a trace amount of carbon monoxide (CO), methane (CH<sub>4</sub>) and carbon dioxide (CO<sub>2</sub>) in an O<sub>2</sub> gas sample. The sample is injected automatically through a 10-port valve. First, a Porapak-N pre-column is used to cut the C<sub>2</sub> compounds. Second, Porapak functions to separate CO/CH<sub>4</sub> and CO<sub>2</sub>. CO and CH<sub>4</sub> are separated by an MS-13X column. Since a large amount of O<sub>2</sub> gas affects the lifetime of a methanizer catalyst, the O<sub>2</sub> gas needs to be removed using an additional 6-port valve. Conversely, CO<sub>2</sub> moves through the Porapak-Q. CO/CH<sub>4</sub> and CO<sub>2</sub> pass through the methanizer device and converted to methane for detection by FID. The system includes Lab Solutions GC workstation software.

**Analyzer Information**

**System Configuration:**

Three valves / four packed columns / Methanizer with FID detector

**Sample Information:**

CO, CO<sub>2</sub>, CH<sub>4</sub>

**Concentration Range:**

No.	Name of Compound	Concentration Range	
		Low Conc.	High Conc.
1	CO	1.0ppm	100ppm
2	CO <sub>2</sub>	1.0ppm	100ppm
3	CH <sub>4</sub>	1.0ppm	100ppm

Detection limits may vary depending on the sample. Please contact us for more consultation.

**System Features**

- Single channel with packed columns
- Matrix O<sub>2</sub> are removed by the third valve by using cutting technology
- Hydrocarbons and water are backflushed by the pre-column while trace CO, CO<sub>2</sub>, and CH<sub>4</sub> reach FID
- Good separate CH<sub>4</sub> and CO with MS-13X packed column
- 13 minutes analysis time

**Typical Chromatograms**




Fig. 1 Chromatogram of FID

First Edition: November, 2017