

# Application Data Sheet

## No.43

## System Gas Chromatograph

### High Speed Refinery Gas Analyzer Nexis GC-2030HSRGA1 GC-2014HSRGA1

This method is for determining the chemical composition of natural gases and similar gaseous mixtures within the composition range shown below. This test method provides data for calculating a sample's physical properties, such as its heating value and relative density, or for monitoring the concentrations of one or more of the components in a mixture. This analyzer uses a total of four valves and eight columns. The Sample is introduced into four sample loops for determination. Using a pre-column, C6+ components are back-flushed as a single peak. The valve timing then allows the hydrocarbons C3 through/to C5 to be separated individually through an Alumina capillary column and detected by FID. Finally, using MS-5A, O<sub>2</sub>, N<sub>2</sub>, CH<sub>4</sub>, and CO are separated. At the same time, CO<sub>2</sub>, C<sub>2</sub>, and H<sub>2</sub>S are separated using an Rtx-Q plot column and detected by a TCD. H<sub>2</sub> will be separated by MS-5A and, with the other components vented out, detected by another TCD using N<sub>2</sub> as carrier gas. The final analysis time is approximately six minutes. The system includes LabSolution workstation software and BTU and Specific Gravity calculation software.

#### Analyzer Information

##### System Configuration:

Four valves / eight capillary and packed columns with two TCD / one FID detectors

##### Sample Information:

He, H<sub>2</sub>, O<sub>2</sub>, N<sub>2</sub>, CO, CO<sub>2</sub>, H<sub>2</sub>S, C<sub>1</sub>-C<sub>5</sub>, C<sub>6+</sub>

##### Methods met:

ASTM-D1945, D1946, D3588, GPA-2261

##### Concentration Range:

No.	Name of Compound	Concentration Range	
		Low Conc.	High Conc.
1	He	0.01%	10.0%
2	H <sub>2</sub>	0.01%	80.0%
3	O <sub>2</sub>	0.01%	50.0%
4	N <sub>2</sub>	0.01%	50.0%
5	CH <sub>4</sub>	0.01%	80.0%
6	CO	0.01%	10.0%
7	CO <sub>2</sub>	0.01%	30.0%
8	C <sub>2</sub> H <sub>4</sub>	0.01%	10.0%
9	C <sub>2</sub> H <sub>6</sub>	0.01%	10.0%
10	C <sub>2</sub> H <sub>2</sub>	0.01%	10.0%
11	H <sub>2</sub> S	0.05%	30.0%
13	C <sub>3</sub> H <sub>8</sub>	0.01%	5.0%
14	C <sub>3</sub> H <sub>6</sub>	0.01%	5.0%
15	i-C <sub>4</sub> H <sub>10</sub>	0.01%	1.0%
16	n-C <sub>4</sub> H <sub>10</sub>	0.01%	1.0%
17	C <sub>3</sub> H <sub>4</sub>	0.01%	1.0%
18	C <sub>2</sub> H <sub>2</sub>	0.01%	1.0%
19	Other Hydrocarbons	0.01%	0.5%
20	C <sub>6</sub> plus	0.01%	0.5%

Detection limits may vary depending on the sample.  
Please contact us for more consultation.

#### System Features

- Less than 6 minutes analysis for refinery gases analysis with H<sub>2</sub>S can be carried out
- Dual TCD with FID channels for simultaneous analysis
- By using split/splitless injector, liquid hydrocarbons can be analyzed by the FID
- Good separation for H<sub>2</sub> and He, and full range capability for H<sub>2</sub>

Typical Chromatograms




Fig. 1 Chromatogram of FID-1




Fig. 2 Chromatogram of TCD-1




Fig. 3 Chromatogram of TCD-2