

Extraction of a Drugs of Abuse Panel from Whole Blood Using ISOLUTE® SLE+ Prior to UPLC-MS/MS Analysis




Figure 1. Example structures by class.

Introduction

This application note describes the extraction of 49 drugs of abuse from whole blood, prior to UPLC-MS/MS analysis. Figure 1. shows examples of these structures by class.

ISOLUTE® SLE+ Supported Liquid Extraction columns offer an efficient alternative to traditional liquid-liquid extraction (LLE) for bioanalytical sample preparation, providing high analyte recoveries, no emulsion formation, and significantly reduced sample preparation.

This application note describes an effective and efficient ISOLUTE SLE+ protocol optimized for 1 mL sample capacity column formats.

The simple sample preparation procedure delivers clean extracts, good recoveries and RSD values and LLOQ from 10 ng/mL.

Sample Preparation Procedure

Format

ISOLUTE® SLE+ 1 mL sample volume column, part number 820-0140-C

Sample Pretreatment

To 500 µL of whole blood, add 10 µL of ISTD. Allow to equilibrate and add 500 µL of 0.1% ammonium hydroxide (aq). Mix.

Sample Loading

Load 750 µL of the pre-treated whole blood onto the column and apply a pulse of vacuum or positive pressure (3–5 seconds) to initiate flow. Allow the sample to absorb for 5 minutes.

Analyte Extraction

Apply DCM/IPA (95/5, v/v, 2.5 mL) and allow to flow under gravity for 5 minutes. Collect in an appropriate glass tube containing 100 µL HCl in methanol (50 mM). This acts to stabilize free-base analytes in the solvent prior to evaporation. Apply a low vacuum or positive pressure (5–10 seconds) to elute any remaining extraction solvent, before applying the next aliquot.

Apply MTBE (2.5 mL) and allow to flow under gravity for 5 minutes. Apply a low vacuum or positive pressure (5–10 seconds) to elute any remaining extraction solvent, before applying the next aliquot.

Apply a final aliquot of MTBE (2.5 mL) and allow to flow under gravity for 5 minutes. Apply vacuum or positive pressure (5–10 seconds) to elute any remaining extraction solvent.

Post Elution and Reconstitution

Evaporate the extract to dryness in a stream of air or nitrogen using a TurboVap® LV (beginning at 8–10 psi or 1 L/min).

Reconstitute the extracts with 100 µL methanolic mobile phase (B) and vortex for 10 seconds before adding 400 µL aqueous mobile phase (A). Vortex for 10 seconds and transfer to suitable vials or a collection plate.

Analytes

2-OH-ethyl-flurazepam	Cocaine	Hydromorphone	Midazolam	PCP
6-MAM	Codeine	Ketamine	Morphine	Temazepam
7-amino-clonazepam	Diazepam	Lorazepam	Nitrazepam	THC-COOH
7-amino-flunitrazepam	Dihydrocodeine	LSD	Norbuprenorphine	Triazolam
Alprazolam	EDDP	MDA	Nordiazepam	Zaleplone
Amphetamine	Estazolam	MDEA	Norfentanyl	Zolpidem
Benzoyllecgonine	Fentanyl	MDMA	Norketamine	Zopiclone
Bromazepam	Flunitrazepam	Mephedrone	Oxazepam	α -OH-alprazolam
Buprenorphine	Flurazepam	Methadone	Oxycodone	α -OH-triazolam
Clonazepam	Hydrocodone	Methamphetamine	Oxymorphone	

UPLC Conditions

Instrument

Waters ACQUITY UPLC with 20 μ L loop

Column

Restek Raptor™ Biphenyl 2.7 μ m (100 x 2.1 mm) with Raptor Biphenyl EXP guard cartridge

Mobile Phase

A: 2 mM ammonium formate (aq), 0.1 % formic acid

B: 2 mM ammonium formate (in methanol), 0.1 % formic acid

Flow Rate:

0.4 mL min

Injection Volume

10 μ L (partial loop with overflow)

Sample Temperature

20 °C

Column Temperature

40 °C

Mass Spectrometry Conditions

Instrument

Premier XE triple quadrupole mass spectrometer equipped with an electrospray interface for mass analysis.

Desolvation Temperature

450 °C

Ion Source Temperature

150 °C

Positive ions acquired in the multiple reaction monitoring (MRM) mode:




Figure 2. Graphical representation of LC gradient conditions.

Table 1. Gradient Conditions.

Time	% A	% B	Curve
0.00	95	5	6
2.00	85	15	6
3.25	75	25	6
4.50	75	25	6
7.50	40	60	6
11.25	40	60	6
12.75	0	100	6
13.50	0	100	6
13.51	95	5	6
15.00	95	5	6

Table 2. MRM Conditions.

Compound	MRM Transition	Cone Voltage (V)	Collision Energy (eV)
Amphetamine	136.0 > 118.9	16	9
Amphetamine-D ₅	141.0 > 123.9	16	9
Methamphetamine	150.0 > 90.9	22	17
MDA	180.1 > 105.0	16	23
MDMA	194.1 > 163.0	20	13
MDEA	208.2 > 163.0	22	13
Hydromorphone	286.2 > 185.1	44	29
Morphine	286.2 > 201.0	42	25
Morphine-D ₃	289.2 > 201.0	42	25
BZE	290.1 > 168.0	30	18
BZE-D ₃	293.1 > 171.0	30	18
Oxycodone	302.2 > 198.1	34	37
Dihydrocodeine	302.2 > 199.1	42	33
Oxycodone	316.2 > 241.2	34	27
Mephedrone	178.1 > 160.0	35	12
Norfentanyl	233.1 > 84.0	25	19
7-amino-flunitrazepam	284.2 > 135.0	40	27
7-amino-clonazepam	286.2 > 121.0	40	30
Hydrocodone	300.2 > 199.1	46	33
Codeine	300.3 > 215.1	42	25
6-MAM	328.2 > 165.1	44	33
6-MAM-D ₃	331.2 > 165.1	44	33
Cocaine	304.2 > 182.0	30	20
Norketamine	224.1 > 124.9	20	23
EDDP	278.2 > 234.2	26	30
Zaleplone	306.2 > 264.2	40	22
Zopiclone	389.2 > 245.1	20	17
Norbuprenorphine	414.3 > 101.0	55	42

Compound	MRM Transition	Cone Voltage (V)	Collision Energy (eV)
Ketamine	238.1 > 124.9	25	27
Nitrazepam	282.2 > 236.1	40	25
Flunitrazepam	314.2 > 268.2	40	25
Clonazepam	316.1 > 270.1	40	25
α -OH-triazolam	359.1 > 331.1	45	26
Oxazepam	287.2 > 241.0	30	21
Estazolam	295.2 > 267.2	40	24
Temazepam	301.1 > 255.1	30	22
Zolpidem	308.2 > 235.1	45	35
Alprazolam	309.2 > 281.2	40	26
Methadone	310.2 > 265.2	26	15
Lorazepam	321.1 > 275.1	30	22
Bromazepam	316.1 > 182.1	40	30
α -OH-alprazolam	325.2 > 297.1	40	25
2-OH-ethyl-flurazepam	333.2 > 109.0	40	27
Triazolam	343.0 > 308.1	45	27
Nordiazepam	271.1 > 139.9	40	28
Diazepam	285.2 > 154.0	40	27
Diazepam-D ₅	290.2 > 154.0	40	27
Midazolam	326.2 > 291.2	45	29
Fentanyl	337.3 > 105.0	35	40
Flurazepam	388.2 > 315.1	35	23
Buprenorphine	468.3 > 468.3	55	5
PCP	244.2 > 158.9	20	15
LSD	323.8 > 222.8	30	25
THC-COOH-D ₃	348.2 > 302.2	25	20
THC-COOH	345.2 > 299.2	25	19

Results

Method performance was assessed by spiking whole blood (375 μ L) with 5 ng of the respective analytes, equating to a concentration of 13 ng/mL of each analyte when extracting

375 μ L of pre-treated whole blood. The percentage analyte recoveries for the various drug classes can be seen in Figure 3. RSDs were all below 10% as shown in Figure 4.




Figure 3. Percentage recovery of the application analytes.




Figure 4. Percentage relative standard deviation of the application analytes.

Calibration Curves

Calibration curves were constructed by spiking whole blood from 1–500 ng/mL of each analyte prior to extraction. The respective internal standards were spiked at 50 ng/mL. Quadratic effect was observed at high concentrations for a number of analytes.

However dilution and internal standards helped overcome this to achieve coefficient of determination (r^2) values greater than 0.99 for all analytes using the optimized extraction protocol. Representative curves are shown in Figures 5–8.




Figure 5. Calibration Curve for morphine.




Figure 6. Calibration Curve for amphetamine.




Figure 7. Calibration Curve for benzoylcegonine (BZE).





Figure 8. Calibration Curve for diazepam.

Table 3. Lower Limits of Quantitation (LLOQ) using optimized ISOLUTE® SLE+ extraction protocol.

Analyte	Estimated LLOQ (ng/mL) performed by Extrahera	Analyte	Estimated LLOQ (ng/mL) performed by Extrahera
Amphetamine	0.035	Ketamine	0.055
Methamphetamine	0.035	Nitrazepam	0.010
MDA	0.050	Flunitrazepam	0.025
MDMA	0.015	Clonazepam	0.025
MDEA	0.025	α -OH-triazolam	0.035
Hydromorphone	0.007	Oxazepam	0.085
Morphine	0.035	Estazolam	0.035
BZE	0.055	Temazepam	0.025
Oxymorphone	0.035	Zolpidem	0.055
Dihydrocodeine	0.010	Alprazolam	0.085
Oxycodone	0.055	Methadone	0.035
Mephedrone	0.125	Lorazepam	0.200
Norfentanyl	0.055	Bromazepam	0.055
7-amino-flunitrazepam	0.125	α -OH-alprazolam	0.200
7-amino-clonazepam	0.125	2-OH-ethyl-flurazepam	0.035
Hydrocodone	0.055	Triazolam	0.055
Codeine	0.010	Nordiazepam	0.025
6-MAM	0.015	Diazepam	0.035
Cocaine	0.055	Midazolam	0.007
Norketamine	0.035	Fentanyl	0.025
EDDP	0.025	Flurazepam	0.025
Zaleplone	0.015	Buprenorphine	0.125
Zopiclone	0.035	PCP	0.125
Norbuprenorphine	0.015	LSD	0.010
		THC-COOH	5.000

**Figure 9.** Overlaid MRM Chromatograms of application analytes at 50 ng/mL using the optimized extraction protocol.

Extract Cleanliness

An experiment was performed to evaluate the level of residual phospholipids in the final extract. Phospholipids are interfering matrix components which can mask or otherwise interfere with the quantitation of the compounds of interest in LC-MS/MS.

When the optimized method is used, ISOLUTE® SLE+ 1 mL columns show very clean total ion chromatograms (see Figure 10), compared to whole blood matrix which had been precipitated with acetonitrile only prior to LC-MS/MS analysis.




Figure 10. Total Ion Chromatograms of common phospholipid MRM transitions using the optimized extraction protocol (top), compared to blank 70/30 water/methanol (v/v) (middle) and precipitated whole blood (bottom).



Conclusions

The simple method described in this application note is suitable for extraction of a broad range of analyte classes, including THC-COOH from whole blood.

The method is easily automated using the Biotage® Extrahera™ Automation System, showing enhanced analyte recovery and reproducibility for many analytes compared to the manually processed method.

Additional Notes

1. Buprenorphine and BZE extraction recoveries are low compared to samples fortified with the analyte after extraction. However the LLOQ values in table 3 illustrate that this is not an obstacle to effective quantitation.
2. If increased sensitivity is desired on the THC-COOH quantitation, the final reconstitution volume can be modified to less than 0.5 mL.
3. Amphetamines, bath salts and ketamines can suffer loss on evaporation when drying in the more volatile free base form. To overcome this effect, 100 μ L of 50 mM HCl in MeOH is added to the collection plate/culture tubes to convert to the corresponding HCl salt forms.

Reagent Preparation

All solvents were HPLC-grade.

- » **0.1% ammonia hydroxide (aq):** Add 0.1 mL of concentrated ammonium hydroxide (28–30%) to 99.9 mL HPLC grade water.
- » **Aqueous Mobile Phase:** Weigh 126 mg and dissolve in 1 L UHPLC grade water. Add 1 mL concentrated formic acid.
- » **Methanolic Mobile Phase:** Weigh 126 mg and dissolve in 1 L UHPLC grade methanol. Add 1 mL concentrated formic acid.
- » **50 mM HCl in methanol:** Add 50 μ L concentrated hydrochloric acid to 11.95 mL HPLC grade methanol. The hydrochloric acid stock is commercially available ~12M.

Ordering Information

Part Number	Description	Quantity
820-0140-C	ISOLUTE® SLE+ 1 mL Sample Volume Columns	30
414001	Biotage® Extrahera™ Automation System	1
415040	Configuration Kit 96 Positions Dual Flow	1
414141	Extrahera clear tips	960
PPM-96	Biotage® PRESSURE+ 96 Positive Pressure Manifold (96 well)	1
C103199	TurboVap® LV Evaporator	1

Appendix

Biotage® Extrahera™ Settings

The method described in this application note was automated on the Biotage® Extrahera™, using ISOLUTE® SLE+ 1 mL columns. Performance was comparable to manually processed samples, as demonstrated in Figures 3 and 4. For the majority of analytes, recovery and RSD is improved upon compared to manual processing. This appendix contains the software settings required to configure Extrahera to run this method. An importable electronic copy of this method for Extrahera can be downloaded from www.biotage.com

Using this automated procedure, 24 samples can be processed in a total of 40 min 43 secs.

Method Name: DOA in whole blood
Sample Plate/Rack: 13 x 100 mm Test Tubes, 24
Extraction Media: ISOLUTE® SLE+ 1 mL using whole blood



< Cancel Edit SLE Method - Drugs of Abuse Whole Bloo... Save >

Method name: Drugs of Abuse Whole Blood ISOLUTE SLE+ 1 n Sample plate/rack: 13 x 100 mm Test Tube... Extraction media: SLE+ 1 mL using Whole...

Pretreatment: On
 Load: On
 Elution: On

Sample type: Whole Blood:0.1% NH4...
 Starting sample volume in plate/rack (µL): 500
 Reuse sample tips?: No

Method comment: Whole Blood:0.1% NH4OH(aq) using Columns

Settings

"Sample" Tab

Sample Type: Whole Blood
Starting Sample Volume (µL): 500
Reuse sample tips? No
Method comment:

Whole Blood:0.1% NH₄OH(aq) using Columns

Screenshot

Settings

Pre-treatment	Activated
No. of steps	1
Pause after last step	No
Dispose tips after last step	No

Solvent	
1	0.1% Ammonium Hydroxide (aq)
2	
3	
4	

	1	2	3	4
Volume (µL)	500			
Mix number of times	0			
Mix volume (µL)	0			
Wait time (min)	0			

Load	Activated
Air Push Time (s)	20
Pause after each load	No
Volume (µL)	750
Collect in position	D
Air Push Time (s)	20
Wait time (min)	5
Premix	Yes
Number of times	3

Edit SLE Method - Drugs of Abuse Whole Bloo...

Method name: **Drugs of Abuse Whole Blood ISOLUTE SLE+ 1 n** | Sample plate / rack: **13 x 100 mm Test Tube...** | Extraction media: **SLE+ 1 mL using Whole...**

Pre-treatment: On | Load: On | Elution: On

Elution (3)

Air push after last elution? Yes | Air push time (s): **20** | Dispose solvent tips after each step? No

Step	Solvent	Volume (µL)	Collect in position	Wait time (min)	Repeat (number of times)	Pause after this step?
1	DCM-IPA (95:5)	2500	B	0	1	No
2	MTBE	2500	B	0	1	No
3	MTBE	2500	B	0	1	No

Elution	Activated
No. of steps	3
Air push after last elution	Yes
Air push time (s)	20
Dispose tips after each step	No

Solvent
1 DCM-IPA (95:5)
2 MTBE
3 MTBE
4

	1	2	3	4
Volume (µL)	2500	2500	2500	
Collect in position	B	B	B	
Wait time (min)	N/A	N/A	N/A	
Repeat	1	1	1	
Pause	No	No	No	

'Advanced Settings'

0.5 bar for 300 s. Plate Dry OFF. This applies to all three elution steps

Solvent Properties

Solvent Description

1	0.1% Ammonium Hydroxide (aq)
2	DCM-IPA (95:5)
3	MTBE
4	
5	
6	
7	
8	
9	
10	



Solvent	1	2	3	4	5	6	7	8	9	10
Reservoir Type	Refillable			Non Refillable						
Capacity	N/A	N/A	N/A							
Aspiration flow rate (mL/min)	10	10	10							
Dispense flow rate (mL/min)	20	10	10							
Lower air gap flow rate (mL/min)	20	10	10							
Lower air gap volume (µL)	5	5	5							
Upper air gap flow rate (mL/min)	20	120	120							
Upper air gap volume (µL)	100	100	100							
Upper air gap dispense pause	300	300	300							
Conditioning?	Yes	Yes	Yes							
Conditioning number of times	2	2	2							
Conditioning flow rate (mL/min)	20	10	10							
Conditioning volume (%)	100	100	100							
Aspirate post dispense	Yes	Yes	Yes							
Chlorinated	No	Yes	No							
Serial dispense	No	No	No							

< Cancel Edit Sample - Whole Blood:0.1% NH4OH(aq) u... Save >

Sample	Air Gap	Aspirate
Sample name Whole Blood:0.1% NH4OH(ac)	Lower air gap flow rate (mL/min) 20	Aspirate post dispense? Yes <input type="checkbox"/>
Sample description For 1 mL SLE columns	Lower air gap volume (µL) 10	
Aspiration flow rate (mL/min) 5	Upper air gap flow rate (mL/min) 100	
Dispense flow rate (mL/min) 1	Upper air gap volume (µL) 100	
	Upper air gap dispense pause (ms) 5000	

"Sample" Screen

Sample name	Whole Blood
Sample description	For 1 mL SLE columns
Aspiration flow rate (mL/min)	5
Dispense flow rate (mL/min)	1
Lower air gap flow rate (mL/min)	20
Lower air gap volume (µL)	10
Upper air gap flow rate (mL/min)	100
Upper air gap volume (µL)	100
Upper air gap dispense pause	5000
Aspirate post dispense	Yes

< Cancel Edit Extraction Media - SLE+ 1 mL using Whole... Save >

Extraction Media	Pipetting Height
Name SLE+ 1 mL using Whole Blood	Solvent dispensation height (mm) -123.0
Manufacturer Biotage	Sample dispensation height (mm) -152.0
Part number	Aspiration height (mm) -123.0
Capacity volume (µL) 6000	Tune Pipetting Heights...
Format 24	
Comment This modified dispense height helps prevent drops clinging to the tips	

"Extraction Media" Screen

Name	ISOLUTE® SLE+ 1 mL using Whole Blood
Manufacturer	Biotage
Part number	820-0140-C
Capacity volume (µL)	6000
Format	24
Solvent dispensation height (mm)	-123.0
Sample dispensation height (mm)	-152.0
Aspiration height (mm)	-123.0
Comment	This modified dispense height helps prevent drops clinging to the tips

< Cancel Edit Sample Plate/Rack - 13 x 100 mm Test Tu... Save >

Sample Plate/Rack	Pipetting Height
Name 13 x 100 mm Test Tubes, 24	Aspiration height (mm) -197.0
Capacity volume (µL) 8000	Pre-treatment dispensation height (mm) -103.0
Format 24	Tune Pipetting Heights...

"Sample Plate/Rack" Screen

Name	13 x 100 mm Test tubes, 24
Capacity volume (µL)	8000
Format	24
Aspiration height (mm)	-197.0
Pre-treatment dispensation height (mm)	-103.0

< Cancel Edit Pipette Tip - 1000 µL Biotage tip Save >

Pipette Tip

Name
1000 µL Biotage tip

Manufacturer
Biotage

Part number
414141

Capacity (µL)
1000

Length (mm)
95

"Pipette tip" Screen	
Name	1000 µL Biotage Tip
Manufacturer	Biotage
Part number	414141
Capacity (µL)	1000
Length (mm)	95

EUROPE

Main Office: +46 18 565900
 Toll Free: +800 18 565710
 Fax: +46 18 591922
 Order Tel: +46 18 565710
 Order Fax: +46 18 565705
 order@biotage.com
 Support Tel: +46 18 56 59 11
 Support Fax: + 46 18 56 57 11
 eu-1-pointsupport@biotage.com

NORTH & LATIN AMERICA

Main Office: +1 704 654 4900
 Toll Free: +1 800 446 4752
 Fax: +1 704 654 4917
 Order Tel: +1 704 654 4900
 Order Fax: +1 434 296 8217
 ordermailbox@biotage.com
 Support Tel: +1 800 446 4752
 Outside US: +1 704 654 4900
 us-1-pointsupport@biotage.com

JAPAN

Tel: +81 3 5627 3123
 Fax: +81 3 5627 3121
 jp_order@biotage.com
 jp-1-pointsupport@biotage.com

CHINA

Tel: +86 21 2898 6655
 Fax: +86 21 2898 6153
 cn_order@biotage.com
 cn-1-pointsupport@biotage.com

KOREA

Tel: + 82 31 706 8500
 Fax:+ 82 31 706 8510
 korea_info@biotage.com
 kr-1-pointsupport@biotage.com

Distributors in other regions are listed on www.biotage.com

Part Number: AN875

© 2017 Biotage. All rights reserved. No material may be reproduced or published without the written permission of Biotage. Information in this document is subject to change without notice and does not represent any commitment from Biotage. E&OE. A list of all trademarks owned by Biotage AB is available at www.biotage.com/legal. Other product and company names mentioned herein may be trademarks or registered trademarks and/or service marks of their respective owners, and are used only for explanation and to the owners' benefit, without intent to infringe.

