

Highly Sensitive Analysis of Sulfur Compounds in Beer Using the Trap Mode of a Headspace Sampler

Yuto Nakasuji and Akara Hashimoto

User Benefits

- ◆ Using the Trap mode of HS-20 NX allows for a simple and highly sensitive analysis of concentrated headspace gas.
- ◆ By utilizing the Trap mode and SCD-2030, it becomes possible to detect trace sulfur compounds that were difficult to analyze in the Loop mode.
- ◆ This system is compatible with N₂ carrier gas, allowing analysis to be conducted without the use of helium.

Introduction

Volatile sulfur compounds, which are generated during the fermentation process of beer yeast, significantly contribute to the flavor and quality of beer. However, as these compounds exist in very small amounts, their detection requires pre-concentration techniques or the use of high-sensitivity detectors.

In a previous Application News titled "[Analysis of Volatile Sulfur Compounds in Beer Using Nexis™ SCD-2030](#)", we introduced an analysis method combining SCD (Sulfur Chemiluminescence Detector) and the Loop mode of the headspace sampler for the analysis of volatile sulfur compounds in beer.

In this Application News, we will present the results of further investigation into achieving higher sensitivity using the Trap mode of the headspace sampler.



Fig. 1 HS-20 NX + Nexis™ GC-2030 + SCD-2030

Samples and Analytical Conditions

We used two test-brewed beers as samples, differing only in the yeast (Yeast A and Yeast B).

In the procedure, 3 g of NaCl was measured and added into a vial, then 3 g of the sample was measured, added into the same vial, and sealed. The analysis was conducted by concentrating the headspace gas in the vial through five multi-injections using the Trap mode of HS-20 NX, and injection into the GC. The detailed analytical conditions are shown in Table 1.

Concentrate Headspace Gases in a Trap

The HS-20 NX can operate in loop mode or trap mode. These two modes use different methods for sampling the headspace gas. In the Loop mode, the headspace gas is collected in a measuring tube and injected into the GC. On the other hand, in the Trap mode, the headspace gas is collected in a trap tube, followed by thermal desorption and injection into the GC. The Trap mode enables high sensitivity by multiple injections from the same vial to concentrate the gas in the trap tube. The flow of analysis by the Trap mode is shown in Figure 2.

With the HS-20 NX, users can easily switch between Loop mode and Trap mode in the software, allowing for flexible usage depending on the sample and concentration.




Fig. 2 Headspace Sampling in Trap Mode

Table 1 System Configuration and Analytical Conditions

System	HS-20 NX Nexis GC-2030/SCD-2030	
HS		
Mode:	Trap (Tenax TA)	
Oven Temp.:	80 °C	
Sample Line Temp.:	100 °C	
Transfer Line Temp.:	100 °C	
Trap Cooling Temp.:	10 °C	
Trap Heating Temp.:	250 °C	
Trap Waiting Temp.:	25 °C	
Multi Injection:	5	
Vial Pressure:	80 kPa	
Dry Purge Pressure:	20 kPa	
Vial Heating Time:	35 min	
Vial Pressurization Time:	1 min	
Pressure Equilibrating Time:	0.1 min	
Loading Time:	0.5 min	
Load Equilibrating Time:	0.1 min	
Dry Purge Time:	10 min	
Injection Time:	10 min	
Needle Flush Time:	45 min	
GC		
Injection Mode:	Split	
Split Ratio:	5	
Carrier Gas:	N ₂	
Carrier Gas Control:	Const. Linear velocity (45 cm/sec)	
Column:	DB-1 (60 m × 0.32 mm I.D., 5 μm)	
Oven Program:	60 °C (3 min)_15 °C/min_240 °C (20 min)	
SCD		
Interface Temp.:	200 °C	
Electric Furnace Temp.:	850 °C	
Detector Gas:	H ₂ 100 mL	
	N ₂ 10 mL	
	O ₂ 12 mL	
	O ₃ 25 mL	

Results

Fig. 3 shows the results obtained from analyzing the beers in trap mode and loop mode. An examination of the enlarged chromatogram from 9 to 19 minutes reveals a substantial improvement in sensitivity in trap mode compared to loop mode, with multiple components that were barely detected in loop mode readily detected in trap mode.

The S/N ratios of the major peaks (peaks A to J in Fig. 3) are also shown in Table 2. The S/N results showed that trap mode improved sensitivity by a factor of 6 to 20 times compared to loop mode.

The results also showed a substantial difference in the amount of S-methylthioacetate present in the two beers (Fig. 4). S-Methylthioacetate is produced by yeast during fermentation and is known to contribute to the aroma of beer.¹⁾ The results of this analysis reveal that S-methylthioacetate potentially plays an important role in the difference in the aroma of the two beers.




Fig. 3 Chromatograms of Trap Mode and Loop Mode Results for Each Beer

Table 2 S/N Ratio of Major Peaks (A-J) in Trap Mode and Loop Mode

Peak	Yeast A		Yeast B	
	Trap	Loop	Trap	Loop
A	13.6	N.D.	18.1	N.D.
B	1896.7	287.0	104.2	14.1
C	17.6	N.D.	15.5	N.D.
D	828.5	45.0	627.2	30.9
E	34.4	N.D.	30.7	N.D.
F	79.7	9.4	78.2	8.9
G	586.8	57.5	573.7	51.9
H	76.6	10.8	42.9	6.2
I	57.3	N.D.	52.0	N.D.
J	28.3	N.D.	24.8	N.D.

Conclusion

This article examined using the HS-20 NX headspace sampler in trap mode with the SCD-2030 detection system to further improve the sensitivity of analysis of sulfur components in beer.

By analyzing in trap mode, we confirmed an improvement in sensitivity of 6 to 20 times compared to loop mode and enabled the detection of multiple compounds that were not detected in loop mode.

Analyzing two test beers brewed under identical conditions except for the brewing yeast (yeast A and yeast B) in trap mode revealed a much larger amount of S-methylthioacetate in the beer brewed with yeast A, and it identified S-methylthioacetate as the compound potentially responsible for the difference in the aroma of the two beers.




Fig. 4 Enlarged Chromatogram Showing S-Methylthioacetate

< References >

- 1) Identification and Determination of S-Methyl Thioacetate in Beer, Nippon Nogeikagaku Kaishi, Vol. 54, No. 9, 1980

< Acknowledgments >

We would like to thank Narihiro Sukuzi (CEO), Takuma Yamamiya, and Ren Takasaki of Nikenjayamochi Kadoya Honten (Ise Kadoya Brewery) and Hironori Maruyama, principal researcher at Mie Prefecture Industrial Research Institute, for their assistance with this research.

Nexis is a trademark of Shimadzu Corporation or its affiliated companies in Japan and other countries.