


Agilent 9500 ICP-QQQ with m-Lens for Ultratrace Analysis of High-Purity Reagents

Achieving low backgrounds under hot plasma conditions using ICP-QQQ preset methods



Author

Rentaro Yamashita
Agilent Technologies, Inc.

Introduction

Controlling trace-metal contaminants is essential in both semiconductor reagent analysis and metal-based materials characterization. High-purity chemicals such as nitric acid (HNO_3) and sulfuric acid (H_2SO_4), as well as high-purity metals, require analytical workflows capable of achieving sub-ppt detection limits (DLs) while maintaining excellent robustness and reproducibility. To meet these needs, the Agilent 9500 Triple Quadrupole ICP-MS (ICP-QQQ) equipped with a Dual-Cell System (DCS)¹ and the optional m-lens ion optics provides a highly stable platform, enabling low backgrounds and reliable performance under hot plasma conditions.

The m-lens is engineered to suppress background signals for elements including easily ionized elements such as Na and K that typically show elevated background equivalent concentrations (BECs) under hot plasma conditions. The 9500 ICP-QQQ with m-lens is suitable for the analysis of metal digests and other sample types that demand robust performance and low backgrounds.²

However, its benefits also extend to low-matrix, high-purity reagents, supporting sub-ppt impurity measurements of simpler matrices.³

Agilent's ICP-QQQ portfolio allows users to select instruments based on desired DLs and operational objectives. For laboratories seeking the lowest possible DLs, the Agilent 8900 Semiconductor configuration ICP-QQQ equipped with an s-lens provides exceptional sensitivity, making it suitable for laboratories operating at the extreme limits of ultratrace quantification. The 9500 ICP-QQQ with m-lens provides robust, high sensitivity measurements, ideal for routine trace level applications. This model is appropriate for labs that may not require the lowest DLs for every element but still demand high-quality, reproducible data.

In this study, we evaluated the performance of the 9500 ICP-QQQ using a method designed specifically for low-matrix, high-purity reagents. Using preset methods implemented in the instrument software for ultrapure water (UPW) and for H₂SO₄—the UPW method is also applicable to diluted HNO₃ and H₂O₂— we analyzed 1% HNO₃ and high purity H₂SO₄ separately. Sub-ppt DLs and low backgrounds were achieved for target analytes in both matrices. These results confirm that the 9500 ICP-QQQ with m-lens can be used for ultratrace impurity analysis of semiconductor-grade reagents, as well as demanding metal-matrix applications.

Experimental

Reagents, sample preparation, and calibration

UPW produced by Puric ω II (Organo Corp, Japan), Tamapure-AA-10 grade HNO₃, and Tamapure-AA-100 grade H₂SO₄ from Tama Chemicals (Japan) were used for sample preparation. Calibration standards and spike solutions were prepared using the following multi-element standards: XSTC-7, XSTC-8, and XSTC-331 (SPEX CertiPrep LLC, USA). Unless otherwise specified, all concentrations reported in this study are expressed as w/w%. Before analysis, HNO₃ samples were diluted to 1% and H₂SO₄ samples were diluted to 9.8%.

The 9500 ICP-QQQ was calibrated using a multi-point external calibration curve prepared from matrix-matched standards. The standards were prepared at 0, 5, 10, 20, and 30 ng/kg (ppt) in 1% HNO₃ and 9.8% H₂SO₄, respectively. DLs were calculated as three times the standard deviation of the blank signals, and BECs were obtained by dividing the mean blank signal by the calibration slope. A single blank solution was measured during five replicate runs, and these replicate measurements were used to calculate the DLs and BECs.

Instrumentation

The Agilent 9500 ICP-QQQ fitted with the optional m-lens and the Agilent I-AS autosampler were fully controlled using Agilent OpenLab ICP-MS software version 1.1. Designed for ultratrace metal analysis, the sample introduction system comprised a MicroFlow PFA nebulizer with I-AS probe (operated in self-aspiration mode), a temperature-controlled quartz spray chamber, and a quartz torch with a 2.5 mm inner diameter (id) injector. A platinum-tipped sampler cone with copper base and a platinum-tipped skimmer cone with nickel base for m-lens were used.

The 9500 ICP-QQQ was operated using five tune modes: no gas, Advanced Helium Mode (AHM), O₂, NH₃, and H₂. For the H₂SO₄ samples, additional specialized tunes were prepared for the determination of platinum (Pt) and zinc (Zn) to ensure optimal performance (Table 1). The use of reaction and collision gases enabled effective removal of spectral interferences originating from argon plasma species as well as matrix-related polyatomic ions. The makeup gas flow rate and the Omega lens voltage parameters were automatically optimized by the preset methods for m-lens, resulting in a stable CeO/Ce ratio that was maintained between approximately 1 and 3%.

Table 1. Agilent 9500 ICP-QQQ operating parameters.

Parameter	No Gas	No Gas (Pt)*	AHM	AHM (Zn)*	O ₂	NH ₃	H ₂
Scan Mode	Single Quad				MS/MS		
RF Power (W)	1500						
Nebulizer Gas (L/min)	0.7						
Makeup Gas (L/min)	Autotuned						
Extract 1 (V)	0						
Extract 2 (V)	-80						
Omega Lens (V)	Autotuned						
Q1 Bias (V)	-5	6	-5		0		
He Gas Flow (mL/min)	-	-	14		-	1	-
H ₂ Gas Flow (mL/min)	-	-	-	-	-	-	6
O ₂ Gas Flow (mL/min)	-	-	-	-	0.4	-	-
NH ₃ Gas Flow** (mL/min)	-	-	-	-	-	3 (30%)	-
KED (V)	5		-4.2	-3	-7		-2

* Applies to H₂SO₄ samples only. **10% NH₃ balanced with 90% He.

Results and discussion

For both acid matrices (1% HNO₃ and 9.8% H₂SO₄), all target elements exhibited excellent calibration linearity with correlation coefficients $r > 0.99$ across the respective calibration ranges. Sub-ppt DLs and BECs were achieved for most elements in both matrices, confirming the method's sensitivity for ultratrace impurity control in high-purity reagents (Table 2).

Particularly noteworthy is the performance for Zn in H₂SO₄, a measurement often considered challenging due to matrix-related backgrounds such as ³²S³²S, ³²S³⁴S, and ³⁴S³⁴S at m/z 64, 66, and 68. Using the DCS in AHM and Zn-specific tune settings, the method delivered a DL of 0.64 ppt and BEC of 0.89 ppt, demonstrating effective background suppression and robust interference control even under hot plasma conditions.

Conclusion

This study evaluated the Agilent 9500 ICP-QQQ with m-lens for ultratrace analysis of high-purity reagents using preset methods for 1% HNO₃ and 9.8% H₂SO₄. Sub- or single-figure ppt detection limits and low background equivalent concentrations were achieved for all elements, demonstrating that ultratrace measurements can be easily performed under hot-plasma conditions. The Dual-Cell System enabled effective removal of matrix-related interferences, with Advanced Helium Mode providing an effective approach for the measurement of challenging elements such as Zn in H₂SO₄. The results show that the 9500 ICP-QQQ with intuitive preset methods can deliver sensitive, low-background data, offering an accessible entry point for laboratories that are new to ultratrace analysis.

References

1. Sugiyama, N. Dual-Cell System (DCS) and Advanced Helium Mode (AHM), Agilent publication, [5994-8985EN](#)
2. Yamashita, R. Analysis of High Purity Titanium Using an Agilent 9500 ICP-QQQ, Agilent publication, [5994-9024EN](#)
3. Sakai, K.; Shimamura, Y. Ultrapure Process Chemicals Analysis by ICP-QQQ with Hot Plasma Conditions, Agilent publication, [5994-4025EN](#)

Table 2. Tune modes, measured m/z , and calculated DLs and BECs for each analyte. DLs and BECs were derived from five replicate measurements of the same blank solution ($n = 5$). Analytes with no detectable signal are reported as ND (not detected).

Analyte	Tune	Q1	Q2	1% HNO ₃		9.8% H ₂ SO ₄	
				DL (ppt)	BEC (ppt)	DL (ppt)	BEC (ppt)
Li	No gas	–	7	0.05	< DL	0.05	< DL
B	No gas	–	11	0.98	2.08	2.04	6.08
Na	NH ₃	23	23	0.16	0.17	0.19	1.42
Mg	NH ₃	24	24	0.02	< DL	0.08	< DL
Al	NH ₃	27	27	0.06	< DL	0.17	< DL
K	NH ₃	39	39	0.25	0.71	0.41	2.04
Ca	H ₂	40	40	0.04	0.06	0.16	0.20
Ti	O ₂	47	63	–	–	1.30	< DL
Ti	O ₂	48	64	0.02	< DL	–	–
V	AHM	–	51	ND	ND	–	–
V	NH ₃	51	51	–	–	0.08	< DL
Cr	NH ₃	52	52	0.06	0.06	0.14	0.24
Mn	AHM	–	55	0.04	< DL	0.12	< DL
Fe	NH ₃	56	56	0.33	0.49	0.41	1.13
Co	AHM	–	59	0.02	< DL	0.02	< DL
Ni	AHM	–	60	0.20	< DL	0.51	< DL
Cu	AHM	–	63	0.14	0.21	0.22	0.25
Zn	AHM	–	66	0.10	< DL	–	–
Zn	AHM (Zn)	–	68	–	–	0.64	0.89
As	O ₂	75	91	0.06	< DL	0.05	< DL
Sr	AHM	–	88	ND	ND	ND	ND
Mo	AHM	–	95	0.02	< DL	0.06	< DL
Cd	AHM	–	111	0.03	< DL	ND	ND
Sn	AHM	–	118	0.21	< DL	0.25	0.48
Sb	AHM	–	121	0.05	< DL	0.06	< DL
Ba	AHM	–	137	0.06	< DL	0.14	< DL
W	AHM	–	182	ND	ND	0.01	< DL
Pt	No gas	–	195	0.15	2.52	–	–
Pt	No gas (Pt)	–	195	–	–	0.13	0.34
Pb	AHM	–	208	0.03	0.03	0.10	0.26

Products used in this application

Agilent products

Product Type	Description	Part Number
Sample Introduction System	9500 ICP-MS quartz torch, 2.5 mm id for aqueous samples	M5150-67011
	9500 ICP-MS quartz connector tube, straight	M5150-67014
	9500 ICP-MS quartz spray chamber with straight exit port	M5150-67017
	MicroFlow PFA nebulizer with I-AS probe, self-aspirates at 200 µL/min	G3139-65102
Interface	ICP-MS sampler cone for 9500 ICP-MS, Pt tip with Cu base	M5150-67002
	Skimmer cone, Pt tip with Ni base for m-lens	G8400-67073
	Extraction-Omega lens assembly, m-lens, brass base	M5150-67023
Tubing Kits	Easy-fit peristaltic pump tubing, beige thermoplastic, yellow/blue, 1.52 mm id for drain	5005-0022
Bottle Kits	Waste container kit, includes a 10 L waste can, S60 StaySafe cap, fittings, and acid vapor filter	5005-0437

www.agilent.com/chem/9500icpqqq

DE-013695

This information is subject to change without notice.

© Agilent Technologies, Inc. 2026
Published in the USA, June 01, 2026
5994-9062EN