

Application Report 03

US EPA Method 8270 Semivolatiles on the Equity-5

US EPA Method 8270 describes the analysis of a wide variety of semivolatile compounds by GC/MS. In this application, the Equity-5 was used to analyze 74 compounds commonly targeted for analysis by environmental laboratories. The analysis also included 6 internal standards and 8 surrogates. The Equity column showed excellent peak shape and response for both basic and acidic compounds, including the typically low-responding compounds 2,4-dinitrophenol, pentachlorophenol, and benzidine. The low bleed of this column made it suitable for use in an MSD.

Key Words


semivolatiles, Method 8270, Equity, 28092-U

Author: K. Stenerson

Raw Data File Name:

mar02\0320004.D

Acquisition System: 5973 MSD



G001687

Conditions

Column: Equity-5, 30m x 0.25mm ID, 0.5µm
Cat. No.: 28092-U
Oven: 40°C (2 min) to 100°C @ 50°C/min to 200°C @ 10°C/min to 325°C @ 30°C/min (7.5 min)
Inj.: 280°C
MSD Interface: 325°C
Scan Range: 45-450 m/z
Flow: Pressure programmed, 20psi (0.0 min), ramp to 80psi (0.0 min), ramp to 16.5psi (3 min), ramp to 25psi (hold for remainder of run)
Injection: 1.0µL, splitless (0.61 min)
Liner: 4mm ID single taper
Sample: 50ng on-column of a 74 component semivolatile standard, 6 internal standards, and 8 surrogates

Peak IDs

1. Pyridine	13. Benzyl alcohol	25. 2,4-Dimethylphenol	57. N-Nitrosodiphenylamine
2. N-Nitrosodimethylamine	14. 1,2-Dichlorobenzene-d4 (Surr.)	26. Benzoic acid	58. Azobenzene
3. 2-Fluorophenol (Surr.)	15. 1,2-Dichlorobenzene	27. Bis(2-Chloroethoxy)methane	59. 2,4,6-Tribromophenol (Surr.)
4. Phenol-d6 (Surr.)	16. 2-Methylphenol	28. 2,4-Dichlorophenol	60. 4-Bromophenyl phenyl ether
5. Phenol	17. Bis(2-Chloroisopropyl)ether	29. 1,2,4-Trichlorobenzene	61. Hexachlorobenzene
6. Aniline	18. 4-Methylphenol	30. Naphthalene-d8 (Int. Std.)	62. Pentachlorophenol
7. Bis(2-Chloroethyl)ether	19. N-Nitroso-di-n-propylamine	31. Naphthalene	63. Phenanthrene-d10 (Int. Std.)
8. 2-Chlorophenol-d4 (Surr.)	20. Hexachloroethane	32. 4-Chloroaniline	64. Phenanthrene
9. 2-Chlorophenol	21. Nitrobenzene-d5 (Surr.)	33. Hexachlorobutadiene	65. Anthracene
10. 1,3-Dichlorobenzene	22. Nitrobenzene	34. 4-Chloro-3-methylphenol	66. Carbazole
11. 1,4-Dichlorobenzene-d4 (Int. Std.)	23. Isophorone	35. 2-Methylnaphthalene	67. Di-n-butyl phthalate
12. 1,4-Dichlorobenzene	24. 2-Nitrophenol	36. Hexachlorocyclopentadiene	68. Fluoranthene
		37. 2,4,6-Trichlorophenol	69. Benzidine
		38. 2,4,5-Trichlorophenol	70. Pyrene
		39. 2-Fluorobiphenyl (Surr.)	71. Aramite#1
		40. 2-Chloronaphthalene	72. Terphenyl-d14 (Surr.)
		41. 2-Nitroaniline	73. Aramite#2
		42. Dimethyl phthalate	74. 3,3'-Dimethylbenzidine
		43. Acenaphthylene	75. Butylbenzyl phthalate
		44. 2,6-Dinitrotoluene	76. 3,3'-Dichlorobenzidine
		45. 3-Nitroaniline	77. Bis(2-Ethylhexyl) phthalate
		46. Acenaphthene-d10 (Int. Std.)	78. Benzo(a)anthracene
		47. Acenaphthene	79. Chrysene-d12 (Int. Std.)
		48. 2,4-Dinitrophenol	80. Chrysene
		49. 4-Nitrophenol	81. Di-n-octyl phthalate
		50. Dibenzofuran	82. Benzo(b)fluoranthene
		51. 2,4-Dinitrotoluene	83. Benzo(k)fluoranthene
		52. Diethyl phthalate	84. Benzo(a)pyrene
		53. Fluorene	85. Perylene-d12 (Int. Std.)
		54. 4-Chlorophenyl phenyl ether	86. Indeno(1,2,3-cd)pyrene
		55. 4-Nitroaniline	87. Dibenz(a,h)anthracene
		56. 2-Methyl-4,6-dinitrophenol	88. Benzo(g,h,i)perylene