

Robustness of SPME Arrow Immersion Sampling: Polycyclic Aromatic Hydrocarbons in Drinking Water

Gwen Lim Sin Yee¹, Thomi Preiswerk¹, Hans-Joachim Hübschmann², Chong Chiew Mei², Günter Böhm¹
¹CTC Analytics AG, Zwingen, Switzerland, ²CTC Analytics Asia Pte Ltd, Singapore

PAL SYSTEM
 Ingenious sample handling

Abstract

Analyzing polycyclic aromatic hydrocarbons (PAHs) in drinking water becomes challenging especially for hard water (high mineral content water) and/or when sodium thiosulfate being added for treating chlorine content in drinking water during analysis. In this article, repeatability for determination of PAHs in drinking water using SPME Arrow (PDMS) immersion technique with gas chromatography mass spectrometer (GC/MS) was carried out. Single SPME Arrow had been used in this experiment to determine its durability. The repeatability (n=100) which characterized as %RSD was less than 10.93% for all analytes, ranged between 1.74 % - 10.93 %. The relative peak area of the 100th injection against the 1st injection peak area for all analytes ranged between 0.600 to 0.947; indicating SPME Arrow fiber is dutiful for at least 100 injections with immersion technique. In conclusion, fully automated SPME Arrow immersion technique with GC/MS analysis is an excellent option to analyze PAH compounds in drinking water; with higher throughput and robustness.

Keywords

Polycyclic aromatic hydrocarbons (PAHs), solid phase micro-extraction (SPME), SPME Arrow, immersion, PDMS

Introduction

Robustness of SPME and SPME Arrow used for immersion sampling is much lesser compared headspace sampling. During immersion extraction, SPME Arrow phase is fully contacted with sample matrix; e.g. drinking water and tap water which containing minerals and residual chlorine. This will affect SPME and SPME Arrow performance. In this poster, repeatability for determination of PAHs in drinking water using SPME Arrow immersion technique with GC/MS was carried out. Automated SPME Arrow workflow was optimized and one (1) SPME Arrow (PDMS) was used throughout the analysis to determine SPME Arrow durability.

Chemicals

- Pak-Mix 16 in Cyclohexane (Neochema GmbH)
- Internal Standards Mix 25 in Acetone (Dr. Ehrenstorfer GmbH)
- Benzo[a]pyrene D12 in Cyclohexane (Dr. Ehrenstorfer GmbH)
- Sodium Thiosulfate (Sigma Aldrich)
- Acetone (> 99.5% Reagent grade) (Sigma Aldrich)
- Water (HPLC grade) (J.T. Baker)

Test Solution

Switzerland tap water (drinkable, 250mg/L total minerals) spiked with

- 16 PAH compounds (100.8306 pg/mL for each PAH)
- 4 deuterated PAH compounds (136.4374 pg/mL for each deuterated PAH)


Instrumentations

- PAL RTC 120 system with 1.1 mm Ø 100 µm PDMS SPME Arrow (ARR11-P-100/20-P1)
- Shimadzu GC-2010 with Restek Rtxi-5MS Column (30 m X 0.25 mmID X 0.25 µm)
 - Inlet Temperature: 280°C
 - GC Carrier: Helium at 1.0 mL/min (Splitless)
 - GC Oven: Initial 35°C (5 minutes); 40°C/min to 150°C; 20°C/min to 250°C; 10°C/min to 305°C (22 minutes)
- GCMS-QP2010
 - Interface Temperature: 310°C
 - Source Temperature: 250°C
 - Mode: EI with SIM

Table 1: Analytes details includes compound name, retention time and quantify ion.

Compound Name	Abbrev.	Target/ISTD	Group	Retention Time, minutes	Quant. Ion, m/z
Naphthalene	N	Target	1	9.241	128
Acenaphthylene	Acy	Target	1	10.688	152
Acenaphthene-d10	Ace-d10	ISTD	1	10.832	164
Acenaphthene	Ace	Target	1	10.863	153
Fluorene	F	Target	1	11.418	166
Phenanthrene-d10	P-d10	ISTD	2	12.482	188
Phenanthrene	P	Target	2	12.508	178
Anthracene	Ant	Target	2	12.569	178
Fluoranthrene	Fl	Target	3	14.004	202
Pyrene	Pyr	Target	3	14.326	202
Benz(a)anthracene	BaAnt	Target	3	16.258	228
Chrysene-d12	Chr-d12	ISTD	3	16.275	240
Chrysene	Chr	Target	3	16.320	228
Benz(b)fluoranthene	BbFl	Target	4	18.204	252
Benz(k)fluoranthene	BkFl	Target	4	18.252	252
Benz(a)pyrene-d12	BaPyr-d12	ISTD	4	18.762	264
Benz(a)pyrene	BaPyr	Target	4	18.803	252
Indeno(1,2,3-cd)pyrene	Ipyr	Target	4	21.349	276
Dibenz(a,h)anthracene	DBahAnt	Target	4	21.436	278
Benzo(g,h,i)perylene	BghiPer	Target	4	22.043	276

Automated SPME Arrow Workflow



- 1 Adding reducing agent for water samples containing residual chlorine.
- 2 Intensive cleaning to eliminate minerals build up at SPME Arrow.

Figure 1: Automated SPME Arrow workflow modified for PAHs in water analysis.

Analysis

- 15 mL Test Solution in 20 mL vial
- Replicates of 100 samples analyzed using one (1) SPME Arrow

Results




Figure 2: TIC Chromatogram.


(A) Full TIC chromatogram.

(B) Zoom in TIC chromatogram at 17.5 minutes to 25.0 minutes.

1-N 2-Acy 3-Ace-d10 4-Ace 5-F 6-P-d10 7-P 8-Ant 9-Fl 10-Pyr 11-BaAnt 12-Chr-d12 13-Chr 14-BbFl 15-BkFl 16-BaPyr-d12 17-BaPyr 18-Ipyr 19-DBahAnt 20-BghiPer

Conclusions

SPME Arrow is robust to be used for immersion sampling technique with fully automated workflow. In PAHs in drinking water analysis, adequate cleaning procedures for SPME Arrow ensure good reproducibility results with satisfaction GC-MS response obtained for 100 sample injections.



Post-conditioning SPME Arrow at Conditioning Module.

Table 2: Reproducibility (%RSD) of relative peak area (corrected with corresponding ISTD) for each analyte (n=100).

Analyte Name/ISTD	%RSD of Relative Peak Area	Analyte Name/ISTD	%RSD of Relative Peak Area
N/Ace-d10	4.68	BaAnt/Chr-d12	2.27
Acy/Ace-d10	3.35	Chr/Chr-d12	2.15
Ace/Ace-d10	1.81	BbFl/BaPyr-d12	4.17
F/Ace-d10	2.86	BkFl/BaPyr-d12	4.45
P/P-d10	1.74	BaPyr/BaPyr-d12	2.35
Ant/P-d10	3.87	Ipyr/BaPyr-d12	9.54
Fl/P-d10	2.25	DBahAnt/BaPyr-d12	10.93
Pyr/P-d10	2.75	BghiPer/BaPyr-d12	7.60

Table 3: Ratio for Peak Area of the 100th injection against the 1st injection.

Compound	Ratio	Compound	Ratio	Compound	Ratio
N	0.832	Ant	0.854	BkFl	0.826
Acy	0.945	Fl	0.921	BaPyr-d12	0.700
Ace-d10	0.909	Pyr	0.921	BaPyr	0.724
Ace	0.839	BaAnt	0.846	Ipyr	0.678
F	0.858	Chr-d12	0.922	DBahAnt	0.600
P-d10	0.947	Chr	0.904	BghiPer	0.710
P	0.893	BbFl	0.772		

Presented at

RAFA 2019, November 05-08, Prague, Czech Republic.

References

Emanuel Coelho, Célia Ferreira and Cristina M. M. Almeida, Analysis of Polynuclear Aromatic Hydrocarbons by SPME-GC-FID in Environmental and Tap Waters, *J. Braz. Chem. Soc.*, Vol. 19, No. 6, 1084-1097, 2008.

Andreas Kremser, Maik A. Jochmann and Torsten C. Schmidt, PAL SPME Arrow—evaluation of a novel solid-phase microextraction device for freely dissolved PAHs in water, *Anal Bioanal Chem*, 408:943–952, 2016. (Open access: <http://link.springer.com/article/10.1007/s00216-015-9187-z>)

Reprints of the poster with the full list of references cited will be provided by the author upon eMail request to trpreiswerk@ctc.ch

Registered names, trademarks, etc. used, even when not marked as such, are not to be considered unprotected by law.

