

Semi-Volatile Organic Compounds Analysis using an Agilent J&W DB-5 ms Ultra Inert Column

Application Note

Environmental

Authors

Yun Zou
Agilent Technologies (Shanghai) Co.Ltd
412 YingLun Road
Waigaoqiao Free Trade Zone
Shanghai 200131 P.R.China

Chongtian Yu
Agilent Technologies Co., Ltd.
11/F Cross Tower
318 Fu Zhou Road
Shanghai 200001
China

Abstract

The analysis of basic or acidic semi-volatile environmental pollutants at low concentrations places significant demands on the analytical instrument and gas chromatographic column. This application is based on EPA Method 525.2 using Agilent J&W DB-5 ms Ultra Inert columns. The ultra inert columns show exceptionally low column bleed and outstanding column inertness. By maintaining the same selectivity as DB-5 ms columns, the DB-5 ms Ultra Inert columns can be used for the existing method without method redevelopment.



Agilent Technologies

Introduction

Analysis of low levels of semivolatiles, especially active compounds, by GC or GC/MSD continues to be challenging in environmental analysis. The United States Environmental Protection Agency (US EPA) has published numerous GC/MSD methods for organic analytes in various matrices. EPA Method 525.2 is one, which is specified for drinking water [1]. The extensive analyte list includes acids, bases and neutrals. Since many of these compounds are at trace levels, the extreme sensitivity requirements make EPA method 525.2 one of the most challenging applications. Sensitivity and linearity are two areas where laboratories are seeking improved performance. Some previous publications have discussed how to lower detection limits and improve linearity by changing instrument configuration, such as using a programmable temperature vaporizing (PTV) inlet, 6-mm drawout lens, trace ion detection (TID), and a triple-axis detector (TAD) with the MSD operating in SIM mode [2, 3]. In this application, we will discuss the impact of GC column performance on sensitivity.

Improved sensitivity or maximized signal-to-noise ratio can be achieved by either reducing the background noise or increasing the signal. Over the years, GC column manufacturers have developed low bleed columns to reduce background noise level. The signal can be increased by ensuring an inert sample flow path, which will reduce peak tailing. An inert sample flow path includes liner and seal in a split/splitless inlet, the capillary column, and any metal detector parts. Peak tailing results from an interaction of active compounds with active sites within the flow path. One of the major sources of active sites is the capillary column due to its large surface area. If the number of column active sites can be reduced, higher sensitivity should be achieved. Selecting low column bleed and a highly inert column for semivolatiles analysis is very important to improve sensitivity and to ensure accurate quantitation.

The recently introduced Agilent J&W DB-5 ms Ultra Inert GC capillary column is a new generation of 5-ms columns that exhibit the highest column inertness in the industry with exceptionally low column bleed. The improved inertness performance makes DB-5 ms Ultra Inert columns especially suitable for active compounds, trace level samples, and screening of unknown sample analyses.

Agilent J&W DB-5 ms was specified in the EPA Method 525.2 and was used in this study. In addition, the Agilent J&W DB-5 ms Ultra Inert column was used since it has the same selectivity as the Agilent J&W DB-5ms. Therefore, no method

redevelopment was needed. In addition, a direct comparison of column performance between Agilent columns and that from another GC column manufacturer was performed in this study.

Experimental

Chromatographic Conditions

The experiments were performed on an Agilent 6890N Network GC System equipped with an Agilent 5975C Series GC/MSD, and an Agilent 7683 Automatic Liquid Sampler (ALS). The split/splitless inlets were fitted with long-lifetime septa (Agilent p/n 5183-4761) and splitless deactivated dual taper direct connect liners (Agilent p/n G1544-80700). Injections were made using a 5- μ L syringe (Agilent p/n 5181-1273).

Table 1. Gas Chromatography and Mass Spectrometer Conditions

| GC Conditions | |
|--|---|
| Column: | Agilent J&W DB-5 ms Ultra Inert, 30 m \times 0.25 mm, 0.25 μ m (Agilent p/n. 122-5532UI) or Agilent J&W DB-5 ms, 30 m \times 0.25 mm, 0.25 μ m (Agilent p/n 122-5532) |
| Carrier gas: | Helium, constant flow mode, 32 cm/s at 40 $^{\circ}$ C |
| Inlet: | Pulsed splitless at 250 $^{\circ}$ C; Purge flow 30 mL/min at 1.00 min |
| Oven temperature: | 40 $^{\circ}$ C (1 min); 30 $^{\circ}$ C/min to 160 $^{\circ}$ C; 6 $^{\circ}$ C/min to 300 $^{\circ}$ C (5 min); Postrun: 320 $^{\circ}$ C hold 5 min |
| Injection size: | 1 μ L |
| MS Conditions | |
| MS: | EI, Scan |
| Solvent delay: | 4.0 min |
| MS temperature: | 300 $^{\circ}$ C (Source); 180 $^{\circ}$ C (Quad) |
| Transfer line: | 280 $^{\circ}$ C |
| Scan mode: | Mass range (45–450 amu) |
| MSD: | SIM |
| AutoSIM was used to pick ions, groups, and switching times | |
| Number of groups: | 27 |
| Compounds/group: | Varied 1 to 22 |
| Ions/group: | Varied 2 to 55 |
| Dwell times: | Varied 5 to 10 ms |
| Cycles/peak: | Minimum 10 |
| Tune: | DFTPP |

Both a multiramp temperature program and a single ramp linear temperature program were suggested in EPA method 525.2. This application used the single ramp linear temperature program, since the final oven temperature of 300 °C gave better resolution of peaks than 320 °C, which was specified in the method. Therefore, the new parameters were used for this application to test the columns, and the instrument operating parameters are listed in Table 1.

Pulsed splitless injection was used to minimize residence times of analytes in the liner, thereby reducing loss of active compounds. The inlet liner (Agilent p/n. G1544-80700) was chosen because it has shown the best performance for active compounds at low levels. It does not contain glass wool, which minimizes its contribution to the active compound degradation. This liner connects directly to the column and has a tapered top, limiting contact with metal in the inlet. Other liners can also be used, and a detailed discussion of the liners is available [4].

Full scan data acquisition was mandated by EPA Method 525. The Agilent 6890N/5975C inert GC/MSD can acquire both selected ion monitoring (SIM) and scan data in the same run. With a defined compound list, quantitation (target) ions and qualifier ions/ratios are known or can be easily determined. The use of SIM typically results in a sensitivity increase of 10 to 100 times compared to scan [5].

A source temperature of 300 °C was used instead of the typical 230 °C to 250 °C range. The higher temperature was chosen to minimize peak tailing, and therefore increase sensitivity for PAHs [6]. Lower source temperatures have historically been used to maintain performance of the active pesticides.

The system was retention time locked (RTL) to phenanthrene-d10. The primary benefit of RTL is the ability to maintain retention times after clipping or changing the column. Quantitation database and integration event times do not have to be changed. For laboratories performing SIM analyses, switching group times remain constant [7].

The column performance for EPA525.2 method was evaluated using the same instrument system and chromatographic conditions as in Table 1. Experiments were conducted on an Agilent J&W DB-5 ms column, an Agilent J&W DB-5 ms Ultra Inert column, and columns from other manufacturers respectively in both the SIM and scan modes.

Standards

EPA 525 Semivolatile Calibration Mix (1000 µg/mL, except pentachlorophenol at 4000 µg/mL), PCB Mix (500 µg/mL), EPA 525 Fortification Solution (500 µg/mL) and EPA 525 internal Standard Mix (500 µg/mL) were obtained from SUPELCO.

Pesticide samples were purchased from China National Standards Research Center at 100 µg/mL in different solvents.

Standards were made into concentrations of 0.1 to 10 ppm in dichloromethane, containing 5 ppm of ISTDs (Acenaphthylene-d10, Phenanthrene-d10 and Chrysene-d12) and SSs (1,3 Dimethyl-2-nitrobenzene, Pyrene-d10, Triphenylphosphate, Perylene-d12). Pentachlorophenol is present at four times the other analyte concentrations, as described in EPA 525.2.

Results and Discussion

Performance comparison between DB-5 ms and DB-5 ms Ultra Inert

The 10 overlaid chromatograms for the semivolatile standards are shown in Figure 1. As illustrated in Table 2, the RSD% is less than 0.098% when comparing the retention times for all 50 standards on these two columns. This same selectivity allows for direct replacement of the DB-5 ms Ultra Inert column to any existing methods developed with the DB-5 ms columns.

However, for resolution and peak height of most of semivolatile standards, slightly better performance was achieved with the DB-5 ms Ultra Inert column compared to the DB-5 ms (Figure 2). According to EPA Method 525.2, the system must provide baseline resolution of phenanthrene and anthracene; and benz[a]anthracene and chrysene should be separated by a valley whose height is less than 25% of the average peak height of these two compounds. They are closely eluting compounds with the same quantitation ions. As shown in Figure 2, the phenanthrene and anthracene pair and the benz[a]anthracene and chrysene pair were 100% baseline resolved on the DB-5 ms and DB-5 ms Ultra Inert columns, respectively. A measurement of the peak tailing was done using the US Pharmacopeia tailing factor (T_r). This is calculated using the following formula [8].

$$T_f = W_{5.0} / (T_w \times 2)$$

where: T_w = distance between peak front and retention time of peak (T_R) at 5% of peak height, units are the same as used for $W_{5.0}$.

$W_{5.0}$ = width at 5% of height

Tailing factors (T_f) of these compounds listed in Figure 2 were close to 1. The closer the T_f value is to 1, the more symmetrical a peak. The tailing factors (T_f) could be used to evaluate the inertness performance of the columns, because column

activity leads to poor peak shapes, (tailing), poor quantification, and in extreme cases complete disappearance of active analytes at low levels. Pentachlorophenol is an aggressive compound that often shows low response and tailing peaks through adsorption onto active sites in the flow path during analysis. However, compared to the DB-5 ms and competitive 5-ms columns, the improved inertness performance provided by the DB-5 ms Ultra Inert column resulted in the most symmetrical and sharpest peaks. This increased the response for pentachlorophenol, as shown in Figure 3.




Figure 1. 10 overlaid TIC for semivolatiles, 2 ng-on column. (black-DB-5 ms Ultra Inert column, red-DB-5 ms).

Table 2. Fifty Semivolatile Compound Separation using DB-5 ms and DB-5 ms Ultra Inert Columns

| Peak number | Compound | Ava. RT of DB-5 ms (n=5) | Ava. RT of DB-5 ms Ultra inert (n=5) | *RSD% of RT | Rs of DB-5 ms | Rs of DB-5 ms Ultra inert | m/z |
|-------------|---------------------------------------|--------------------------|--------------------------------------|-------------|---------------|---------------------------|---------|
| 1 | Isophorone | 5.647 | 5.642 | 0.068 | — | — | 82 |
| 2 | 1,3 Dimethyl-2-nitrobenzene(SS) | 6.047 | 6.044 | 0.061 | 8.4 | 8.41 | 134 |
| 3 | Dichlovos | 6.351 | 6.349 | 0.059 | 9.05 | 9.03 | 109 |
| 4 | Hexachlorocyclo-pentadiene | 7.129 | 7.126 | 0.065 | 21.58 | 21.8 | 237 |
| 5 | Dimethyl phthalate | 8.084 | 8.087 | 0.071 | 23.14 | 24.1 | 163 |
| 6 | 2,6-Dinitrotoluene | 8.195 | 8.196 | 0.074 | 2.51 | 2.52 | 165 |
| 7 | Acenaphthylene | 8.342 | 8.342 | 0.080 | 3.18 | 3.17 | 152 |
| 8 | Acenaphthylene-d10(IS #1) | 8.611 | 8.612 | 0.076 | 5.57 | 5.74 | 164 |
| 9 | 2-Chlorobiphenyl | 8.737 | 8.739 | 0.088 | 2.6 | 2.74 | 188 |
| 10 | 2,4-Dinitrotoluene | 8.993 | 8.996 | 0.085 | 5.44 | 5.6 | 165 |
| 11 | Diethyl phthalate | 9.572 | 9.579 | 0.098 | 12.21 | 12.35 | 149 |
| 12 | Fluorene | 9.804 | 9.808 | 0.088 | 4.65 | 4.6 | 166 |
| 13 | 2,3-Dichlorobiphenyl | 11.153 | 11.159 | 0.095 | 24.87 | 25.63 | 222/152 |
| 14 | Hexachlorobenzene | 11.218 | 11.219 | 0.090 | 1.14 | 1.12 | 284 |
| 15 | Pentachlorophenol | 11.795 | 11.798 | 0.092 | 9.78 | 10.07 | 266 |
| 16 | Phenanthrene-d10(IS #2) | 12.357 | 12.363 | 0.088 | 9.28 | 9.34 | 188 |
| 17 | Phenanthrene | 12.426 | 12.432 | 0.091 | 1.16 | 1.13 | 178 |
| 18 | Anthracene | 12.585 | 12.591 | 0.091 | 2.66 | 2.71 | 178 |
| 19 | 2,4,5-Trichlorobiphenyl | 13.133 | 13.140 | 0.089 | 9.27 | 9.61 | 256 |
| 20 | Heptachlor | 14.001 | 14.008 | 0.090 | 14.36 | 14.80 | 100 |
| 21 | Di-n-butyl phthalate | 14.587 | 14.600 | 0.095 | 9.99 | 10.14 | 149 |
| 22 | 2,2',4,4'-Tetrachlorobiphenyl | 14.733 | 14.741 | 0.085 | 2.51 | 2.42 | 292 |
| 23 | chlorpyrifos | 14.882 | 14.892 | 0.088 | 2.50 | 2.53 | 197/97 |
| 24 | 2,2',3',4,6-Pentachlorobiphenyl | 16.247 | 16.255 | 0.083 | 22.89 | 23.28 | 326 |
| 25 | Butachlor | 17.058 | 17.070 | 0.081 | 13.17 | 13.74 | 176/160 |
| 26 | Pyrene-d10(SS) | 17.153 | 17.163 | 0.078 | 1.49 | 1.41 | 212 |
| 27 | Pyrene | 17.214 | 17.223 | 0.076 | 0.87 | 0.91 | 202 |
| 28 | p,p'-DDE | 17.901 | 17.913 | 0.075 | 11.03 | 11.34 | 246 |
| 29 | 2,2',4,4',5,6'-Hexachlorobiphenyl | 18.077 | 18.088 | 0.074 | 2.79 | 2.74 | 360 |
| 30 | p,p'-DDD | 19.162 | 19.176 | 0.075 | 7.54 | 7.04 | 235/165 |
| 31 | o,p'-DDT | 19.235 | 19.250 | 0.078 | 1.17 | 1.18 | 235/165 |
| 32 | Benzyl butyl phthalate | 20.157 | 20.177 | 0.078 | 14.76 | 15.08 | 149 |
| 33 | p,p'-DDT | 20.301 | 20.316 | 0.073 | 2.39 | 2.35 | 235/165 |
| 34 | bis(2-Ethylhexyl)adipate | 20.731 | 20.755 | 0.084 | 7.17 | 7.25 | 129 |
| 35 | Triphenylphosphate (SS) | 20.851 | 20.873 | 0.082 | 1.95 | 1.95 | 326/325 |
| 36 | 2,2'3,3',4,4'6-Heptachlorobiphenyl | 21.737 | 21.754 | 0.070 | 13.80 | 13.79 | 394/396 |
| 37 | Ben[α]anthracene | 21.805 | 21.806 | 0.033 | 1.18 | 1.23 | 228 |
| 38 | Chrysene-d12 (IS #3) | 21.834 | 21.835 | 0.039 | 0.74 | 0.75 | 240 |
| 39 | 2,2',3,3',4,5,6,6'-Octachlorobiphenyl | 21.903 | 21.904 | 0.036 | 0.94 | 1.04 | 430/428 |
| 40 | Chrysene | 21.928 | 21.928 | 0.034 | 0.80 | 0.80 | 228 |
| 41 | bis(2-Ethylhexyl)phthalate | 22.605 | 22.632 | 0.081 | 9.44 | 10.04 | 149 |
| 42 | cis-Permethrine | 24.692 | 24.720 | 0.077 | 35.23 | 36.04 | 183 |
| 43 | trans-Permethrin | 24.906 | 24.936 | 0.080 | 3.58 | 3.63 | 183 |
| 44 | Benzo[b]fluoranthene | 25.704 | 25.725 | 0.064 | 12.10 | 12.29 | 252 |
| 45 | Benzo[k]fluoranthene | 25.802 | 25.824 | 0.067 | 1.30 | 1.40 | 252 |
| 46 | Benzo[a]pyrene | 26.766 | 26.789 | 0.064 | 12.69 | 13.60 | 252 |
| 47 | Perylene-d12(SS) | 26.966 | 26.991 | 0.067 | 2.60 | 2.74 | 264 |
| 48 | Indeno[1,2,3-c,d]pyrene | 30.434 | 30.470 | 0.077 | 40.75 | 40.99 | 276 |
| 49 | Dibenz[a,h]anthracene | 30.594 | 30.635 | 0.086 | 1.69 | 1.76 | 278 |
| 50 | Benzo[g,h,i]perylene | 31.373 | 31.414 | 0.088 | 7.74 | 8.07 | 276 |

*RSD% was calculated by using the retention times for each compound on the DB-5 ms and DB-5 ms Ultra Inert columns.




Figure 2. Resolution of critical pairs.




Figure 3. Extracted ion chromatograms of pentachlorophenol, m/z 266.

Performance Comparison Between DB-5 ms Ultra Inert and Competitive Column

A column performance comparison of the DB-5 ms Ultra Inert column and a 5-ms column from another manufacturer is shown in Figure 4. Elution order of some compounds were not the same for these two columns, such as peak 25 (butachlor), peak 36 (2,2'3,3',4,4'6-heptachlorobiphenyl) and peak 39 (2,2',3,3',4,5,6,6'-octachlorobiphenyl). As demonstrated in Figure 5, much lower column bleed was achieved on the DB-5 ms Ultra Inert column at high temperature compared to the competitive column. The lower bleed profile on the DB-5 ms Ultra Inert column provides greater signal-to-noise, lower detector maintenance, and potential for longer column life. In addition, indeno[1,2,3-c,d]pyrene (peak 48) and dibenz[a,h]anthracene (peak 49) are very difficult to baseline resolve, which was achievable on the DB-5 ms Ultra Inert column.

There are many active compounds such as pentachlorophenol and pesticides in these semi-volatile standards. As discussed

before, the peak symmetry could be used to evaluate the inertness performance of the columns. The T_f values for pentachlorophenol, dichlovos, heptachlor, chlorpyrifos and p,p'-DDE peaks are included in the chromatograms in Figure 3 and Figure 6. The DB-5 ms Ultra Inert column shows a significant reduction in the tailing factors. For example, the tailing factor of pentachlorophenol on the DB-5ms Ultra Inert column is 1.07 whereas the T_f is 1.28 for the competitive 5-ms column (Figure 3).

As shown in Figure 6, the competitive 5-ms column shows severe peak tailing, especially for dichlovos and chlorpyrifos, where the T_f values are 0.80 and 0.77, respectively. In addition, due to peak tailing, there is a noticeable loss of response with these two peaks on the competitive 5-ms column. In contrast, quite symmetrical peak shapes were obtained for dichlovos and chlorpyrifos with the DB-5 ms Ultra Inert column, with tailing factors being very close to 1. Furthermore, strong responses allowed for accurate quantitation especially at trace levels. The DB-5 ms Ultra Inert column demonstrates the benefits that more inert surfaces can have for the analysis




Figure 4. Comparison of Semivolatile Standards using DB-5 ms Ultra Inert column and 5-ms column from other manufacturer, 1 ng on-column.




Figure 5. Expanded view from Figure 3 shows lower bleed and better separation of late-eluting compounds on DB-5 ms Ultra inert column.




Figure 6. Extracted ion chromatograms of pesticides.

of active compounds.

Conclusions

According to the USEPA Method 525.2 criteria, the Agilent J&W DB-5 ms Ultra Inert column can provide excellent chromatography for semivolatiles compounds, including difficult to analyze active compounds at trace levels. The higher inertness and lower column bleed compared to a competitive column showed its superior performance for active compounds. In addition, the DB-5 ms Ultra Inert column has the same selectivity as the DB-5 ms column, making it very straightforward to translate any existing method from the DB-5 ms column to the DB-5 ms Ultra Inert column without method redevelopment.

References

1. USEPA Method 525.2: Methods for the Determination of Organic Compounds in Drinking Water Supplement III, National Exposure Research and Development U.S. EPA, Cincinnati, OH (August 1995).
2. M. Szelewski, "Femtogram GC/MSD Detection Limits for Environmental Semivolatiles Using a Triple-Axis Detector", Agilent Technologies, publication 5989-8342EN, www.agilent.com/chem
3. M. Szelewski, "Drinking Water Semivolatiles Analysis using the 6890N/5975B inert GC/MSD," Agilent Technologies, publication 5989-5421EN, www.agilent.com/chem.
4. M. Szelewski, B. Wilson, and P. Perkins, "Improvements in the Agilent 6890/5973 GC/MSD System for Use with USEPA Method 8270," Agilent Technologies publication 5988-3072EN, www.agilent.com/chem.
5. C. K. Meng, "Identification and Quantitation of Pesticides in the Parts-per-Trillion Range Using Retention Time Locking and GC/MS," Agilent Technologies, publication 5988-4392EN, www.agilent.com/chem.
6. M. Szelewski, "Synchronous SIM/Scan Low- Level PAH Analysis Using the Agilent Technologies 6890/5975 Inert

GC/MSD," Agilent Technologies, publication 5989-4184EN, www.agilent.com/chem.

7. K. Weiner, N. Mata, and P. Wylie, "Retention Time Locking with the G1701BA MSD Productivity ChemStation," Agilent Technologies, publication 5968-3433E, www.agilent.com/chem.
8. "Understanding Your ChemStation," Agilent Technologies manual part number G2070-91125, page 246, Edition 07/08.

For More Information

For more information on our products and services, visit our Web site at www.agilent.com/chem.

www.agilent.com/chem

Agilent shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Information, descriptions, and specifications in this publication are subject to change without notice.

© Agilent Technologies, Inc., 2009
Printed in the USA
October 13, 2009
5990-4823EN



Agilent Technologies