



Triglycerides

Determination of triglycerides on a fused silica capillary column

Application Note

Food Testing & Agriculture

Authors

Agilent Technologies, Inc.

Introduction

Gas chromatography using a custom-made Agilent CP-Sil 5 column separates 15 triglycerides in 30 minutes.




Agilent Technologies

Conditions

Technique : GC-capillary
Column : Agilent CP-Sil 5, 0.32 mm x 10 m fused silica WCOT
CP-Sil 5 (0.12 μm) (Custom-made)
+ 0.32 mm x 2 m PD - deactivated fused silica as
retention gap
Temperature : 70 $^{\circ}\text{C}$ \rightarrow 240 $^{\circ}\text{C}$, 50 $^{\circ}\text{C}/\text{min}$ \rightarrow 340 $^{\circ}\text{C}$, 5 $^{\circ}\text{C}/\text{min}$
Carrier Gas : H_2 , 20 kPa (0.2 bar, 2.9 psi), 55 cm/s
Injector : on-column
Detector : FID
T = 375 $^{\circ}\text{C}$
Sample Size : 0.2 μL

Peak identification

1. cholestane
2. C_{31}
3. C_{33}
4. C_{35}
5. C_{37}
6. C_{39}
7. C_{41}
8. C_{43}
9. C_{45}
10. C_{47}
11. C_{49}
12. C_{51}
13. C_{53}
14. C_{55}
15. C_{57}



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A00082



Agilent Technologies