



Hydrocarbons, $C_1 - C_2$

Impurity determination in ethylene

Application Note

Materials Testing & Research

Authors

Agilent Technologies, Inc.

Introduction

Gas chromatography using an Agilent CP- Al_2O_3/Na_2SO_4 column separates five C_1 to C_2 hydrocarbon impurities in ethylene in eight minutes.



Agilent Technologies


Conditions

Technique : GC-capillary
Column : Agilent CP-Al₂O₃/Na₂SO₄, 0.53 mm x 50 m fused silica PLOT Al₂O₃/Na₂SO₄ (df = 10 µm) (Part no. CP7568)
Temperature : 60 °C (3 min) → 150 °C, 10 °C/min
Carrier Gas : He, 100 kPa (1.0 bar, 14 psi)
Injector : Splitter, 1:100
T = 200 °C
Detector : FID
T = 300 °C
Sample Size : 1 mL
Concentration Range : as listed below

Courtesy : Dow Chemical Canada,
Western Canada Division, R & D Lab,
Jim Luong and Steve Craik

Peak identification

1. methane	85 ppm
2. ethane	800 ppm
3. ethylene	99% balanced
4. propane	20 ppm
5. propylene	20 ppm
6. acetylene	30 ppm



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A00579



Agilent Technologies