



Hydrocarbons, $C_1 - C_6$

Trace analysis of acetaldehyde in hydrocarbon: 'symmetry of acetaldehyde peak'

Application Note

Materials Testing & Research

Authors

Agilent Technologies, Inc.

Introduction

Gas chromatography using an Agilent PoraPLOT Q-HT column is used for the trace analysis of acetaldehyde in a hydrocarbon matrix, with separation in 12 minutes.



Agilent Technologies


Conditions

Technique : GC-wide-bore
Column : Agilent PoraPLOT Q-HT, 0.53 mm x 25 m, used silica
PLOT PoraPLOT Q-HT (df = 20 µm)
(Part no. CP7559)
Temperature : 40 °C → 260 °C, 10 °C/min
Carrier Gas : He, 60 kPa (0.6 bar, 8.7 psi)
Injector : Split
Detector : FID

PoraPLOT Q type porous polymers are known for their high inertness. Polar and non-polar compounds elute with good asymmetry, which makes trace analysis possible. Acetaldehyde elutes very well from a PoraPLOT Q-HT column without interference from hydrocarbons.

Peak identification

1. methane
2. ethane
3. propane
4. acetaldehyde
5. butane
6. pentane
7. hexane



www.agilent.com/chem

This information is subject to change without notice.

© Agilent Technologies, Inc. 2011

Printed in the USA

31 October, 2011

First published prior to 11 May, 2010

A01297



Agilent Technologies